Instrument Instructions for Use

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430, USA
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

©2016 Howmedica Osteonics Corp.

0090-9-393 Rev AA

Labeling Symbols

Attention, see Instructions for Use

Do not Reuse

STERILE R

Sterilized using Irradiation

STERILE GP

Sterilized using Hydrogen Peroxide

STERILE EO

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

EC REP

Authorized Representative in the European Community

REF

Catalog Number

LOT


Batch Code

SN

Serial Number

CONTENTS

English	1
Français	4
Deutsch	7
Italiano	10
Español	13
Svenska	16
日本語	19
Nederlands	22
Suomi	25
Dansk	28
Norsk	31
Português	34
Ελληνικά	37
Polski	40
한국어	43
中文	46

Instrument Instructions for use

IMPORTANT INFORMATION FOR DOCTORS AND OPERATING ROOM STAFF

For detailed information concerning the identification of the product (such as system association, Cat.No.) please refer to the marking on the product and/or the labelling of the package.

Ensure that you are familiar with the possible uses, compatibility and correct handling of the instrument; please remember that product systems may be subject to alterations that affect the compatibility of the instrument with other instruments or with implants.

For your information, avail yourself of the training courses and publications offered (e.g. theatre techniques).

Special comments for application

- Only use the instrument for its intended purpose.
- Always treat the instrument carefully to avoid superficial damage or alterations to the geometry.
- The design of the instrument may not be modified in any way.
- Ensure that the drilling and cutting tools to be used are sharp.
- Before the beginning of the operation, ensure that all components prepared for the operation function correctly with each other.
- In the course of the operation, repeatedly check the connections between the implant and the instruments, or between the instruments themselves, required for the precise positioning and fixing are secure.

Warnings

- Care must be taken not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing differing design parameters, varying tolerances, different materials and manufacturing specifications, HOWMEDICA OSTEONICS Instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting implant.
- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials that are visible.

Adverse Effects

- Malfunctions of instruments including, but not limited to, weld fracture, thread galling and material fatigue have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement, spinal or trauma surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.

- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impaction of the component into the preparation.

Cleaning, Maintenance and Sterilization

- This reusable instrument is **not** sterile.
- Remove all instruments from their packaging prior to cleaning and sterilization.
- All parts must be cleaned and sterilized before use. Cleaning of the reusable instrument prior to sterilization should include: presoaking, manual cleaning, ultrasonication (e.g. for complex parts), and washing/ disinfecting. Enzymatic cleaners, manual detergents and neutral cleaners should be used.
- During and after cleaning, inspect hard-to-reach areas such as lumens within the instrument to help ensure that entrapped debris and/or soil are completely removed. If debris and/or soil are evident, then reclean the instrument until there is no evidence of debris and/or soil.
- Rinsing should be performed after each cleaning step to help ensure cleaners and/or debris/soil are removed.
- Howmedica Osteonics' reusable instruments are recommended to be sterilized by moist heat sterilization.
- Instruments with articulating surfaces must be tested for movement. A moist heat compatible, medical grade lubricant should be applied to all articulating joints prior to sterilization.
- After cleaning and prior to sterilization, the reusable instrument should be double-wrapped or packaged in CSR sterilization wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique).
- Biological indicators and/or steam integrators should be used within the sterilizer load to help ensure the proper sterilization conditions of time, temperature and saturated steam have been met.
- The parameters for sterilization are:

Method:	Moist-Heat Sterilization
Cycle:	Pre-Vacuum (Pre-Vac)
Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Pressure:	2-15 PSIA
Dry-Time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Howmedica Osteonics Corporation has validated the above recommended sterilization cycle based on AAMI/ ANSI/ISO guidelines and recommendations. Other sterilization methods and cycles may also be suitable. However, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution.

- EtO sterilization and cold sterilization techniques are not recommended.

Detailed manual and automatic cleaning instructions, including times and temperatures, are provided in Stryker Orthopaedics Literature LSTP1-B available at www.stryker.com/orthopaedics/cleaning. Disassembly and

assembly instructions applicable to all Stryker Orthopaedics re-usable instruments are all also provided in this brochure.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Howmedica, Osteonics, Stryker. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ▶
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◀ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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