<p><i>Note: This portion of the document should not appear on produced Labels or IFUs.</i></p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p><b>Dimensions:</b> 8.5 inch (width) x 11 inch Booklet</p>	<p><b>Color/Material/Finish:</b> Black Graphics on White Background 20# Bond or Equivalent</p>	<p><b>Label Stock:</b> N/A</p>	
			<p><b>Description/Type:</b> Instructions For Use</p>	
			<p><b>Part Number:</b> 0206-600-700</p>	<p><b>Rev.</b> E</p>

## R<sub>x</sub> ONLY

## Introduction

This instructions for use manual contains the information intended to ensure the safe, effective, and compliant use of your product. This manual is intended for in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

## Indications For Use

The Dual Speed Cement Injection Gun is designed for use with the Stryker High Vacuum Cement Injection System. The injection gun is used to apply bone cement during total-hip and total-knee arthroplasty.

ENGLISH (EN)	
ESPAÑOL (ES)	SVENSKA (SV)
DEUTSCH (DE)	DANSK (DA)
FRANÇAIS (FR)	SUOMI (FI)
ITALIANO (IT)	PORTUGUÊS (PT)
NEDERLANDS (NL)	日本語 (JA)

## Safety Directives

### WARNINGS:

- Before using this equipment, or any component compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
- Only healthcare professionals trained and experienced in the use of this medical device should operate this equipment.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage. DO NOT use any equipment if damage is apparent.
- Upon initial receipt and before each use, clean and sterilize the equipment as indicated. See the *Processing Instructions* section.
- Before using this cement gun in a clinical environment, practice mixing bone cement and using the cement gun.
- DO NOT apply pressure in the femoral canal if the bone cement might penetrate or flow uncontrollably into the femoral cortex, arterial or venous system.
- The cartridge will not withstand the excessive pressure that is created when doughy cement is extruded. DO NOT attempt to extrude the cement if excessive force is necessary. In this situation, remove the application tip and hand pack the cement. Excessive force may cause the cartridge to break. Pieces of the broken cartridge may cause patient and/or operating room staff injury.
- Prepare the acetabulum to ensure the bone cement does not flow into the pelvic cavity.

**CAUTION:** DO NOT repair this cement gun. In case of operating difficulties, return to Stryker for repair.

## For Use With

This section identifies components intended to be used with the equipment to obtain a safe combination.

DESCRIPTION	REF
Advanced Cement Mixing (ACM) Kit	0306-XXX-000 series
High Vacuum Cement Mixing and Injection Kit	0206-553-000
180 Gram Cement Cartridge	0206-520-000

## Accessories

### ⚠ WARNINGS:

- Bone cement is intended for single use only. DO NOT reuse, reprocess, or repackage the bone cement. Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.
- DO NOT modify any equipment without the authorization of the manufacturer.
- Use only Stryker approved accessories with this product. These accessories are designed to withstand normal forces and pressures that are created during bone cement application.

**NOTE:** For a complete list of accessories, contact your Stryker sales representative. Outside the US, contact your nearest Stryker subsidiary.

The following Stryker-approved accessories are sold separately:

DESCRIPTION	REF
Umbrella Nozzle	0206-521-000

## Instructions

1. Turn the ram to disengage the drive system, then pull back the ram until its disc rests against the body of the cement gun (Figure 1).


Figure 1

2. Push the speed selector back and then either up or down to select the desired speed (Figure 2). See *Cement Delivery Settings*.


Figure 2

3. Place a filled cartridge with an application nozzle into the cartridge interface (Figure 3). Holding the cement gun in a vertical position, squeeze the trigger to fill the application nozzle until the bone cement reaches the tip. The cement gun is now ready for bone cement delivery.


Figure 3

## Cement Delivery Settings

**FILL** – The FILL setting delivers 12 to 15 mL of bone cement per squeeze of the trigger. The FILL setting is used to rapidly fill the application nozzle and distal femur.

**DISPENSE** – The DISPENSE setting delivers 2 to 3 mL of bone cement per squeeze of the trigger. The DISPENSE setting is used to inject cement at a high pressure into trabecular bone.

# Processing Instructions

## Point of Use (Post-Surgery)

### CAUTIONS:

- DO NOT use saline to wet or soak the equipment before transport to the decontamination processing area.
- If transport to the decontamination processing area is delayed, cover the equipment with a damp cloth or spray the equipment with a pretreatment foam as often as necessary to maintain moisture. The pretreatment foam will minimize the drying of soil and facilitate later decontamination processing. DO NOT allow pretreatment foam to dry on the equipment.
- Avoid mechanical damage during transport. DO NOT mix heavy devices with delicate devices.

Use absorbent wipes to remove gross soil from the equipment. Clean the equipment as soon as practical, typically within two hours, to preclude extended or repeat cleaning procedures.

## Cleaning

### WARNINGS:

- Use filtered water for diluting detergents and for rinsing the equipment. Mineral residues from hard water can stain the equipment and/or prevent effective cleaning and decontamination.
  - To clean the equipment, use specifically formulated detergents only.
  - Read, understand, and follow the indications, instructions, and safety information supplied with the detergent for correct handling and use of the product.
  - ALWAYS prepare the detergent solution according to the manufacturer's recommendations. Pay particular attention to the concentration used and the total dispersion.
  - Clean the equipment as indicated before first and every use. Use the cleaning methods as indicated in these instructions. Other cleaning methods may prevent proper sterilization of the equipment.
- CAUTIONS:**
- Prior to cleaning, separate all detachable components and remove single use cutting accessories from the cement gun.
  - Use personal protective equipment (PPE) at all times during cleaning. DO NOT disassemble the cement gun for cleaning or sterilization.
  - DO NOT allow bone cement to harden on the ram shaft.
- Poor water quality can adversely affect the life of medical devices. ALWAYS follow the water quality requirements per Association for the Advancement of Medical Instrumentation (AAMI) TIR 34.
  - ALWAYS use a detergent that is suitable for use on aluminum and stainless steel surfaces.
  - ALWAYS handle the equipment with care. DO NOT drop the equipment.
  - DO NOT use solvents, lubricants, or other chemicals, unless otherwise specified.
  - DO NOT use ultrasonic cleaning equipment.
  - DO NOT immerse or soak any equipment in liquid. Equipment may be placed briefly under running water to ensure thorough wetting and contact with liquid while actuating moving parts.
  - DO NOT use pipe cleaners or cotton swabs to clean lumens or connector receptacles.
  - ALWAYS make sure the detergent solution is completely rinsed from the interior and exterior of the equipment before drying the equipment.
  - Use of compressed air is only recommended for drying of equipment.
- Remove all gross soil from the equipment using absorbent wipes or a soft, lint-free cloth moistened with the prepared detergent solution. Device may be placed briefly under running water or flushed with a water-filled syringe to assist in the removal of gross soil. Bone cement may be scraped from the outer surface.
  - Thoroughly clean the equipment. Use suitable brushes and detergent solution to clean all surfaces. Pay particular attention to rough surfaces, crevices, and difficult-to-reach areas where soil may be shielded from brushing. Actuate all moving parts to clean hidden surfaces.
  - Thoroughly rinse the equipment with warm, running water until all traces of detergent solution are removed. Hold the equipment on an incline, distal end pointing down. Pay particular attention to rough surfaces, hinges, blind holes, and joints between mating parts. Actuate all moving parts to rinse hidden surfaces.
  - Visually inspect the equipment for any remaining soil or detergent solution. If soil or detergent solution remains, repeat the cleaning procedure using fresh detergent solution.
  - Allow the equipment to drain on absorbent wipes. If possible, orient the equipment vertically to assist in drainage.
  - Dry the equipment with a soft, lint-free cloth or medical-grade compressed air, or by heating the equipment in an oven below 110 °C [230 °F].
  - After cleaning, inspect the equipment.

## Sterilization

### WARNINGS:

- ALWAYS use a chemical indicator within every sterilization load to make sure the proper sterilization conditions of time, temperature, and saturated steam or hydrogen peroxide are achieved.
- Prior to sterilization, separate all detachable components from the cement gun.
- ALWAYS use new sterilization wrap to enclose the equipment. Do not reuse sterilization wrap.
- Sterilize the equipment as indicated before first and every use.
- Use the sterilization methods as indicated in these instructions. Using other sterilization methods may prevent proper sterilization of the equipment and/or damage the equipment.
- Follow the recommended dry times to prevent moisture from accumulating inside the equipment. Moisture may prevent proper sterilization and/or damage the equipment.
- After sterilization, allow the equipment to cool to room temperature prior to use. Failure to comply may result in a burn injury and/or damage to the equipment.

### NOTES:

- Stryker has validated several sterilization cycles for the sterilization of this equipment. However, sterilizer design and performance can affect the efficacy of the process. Healthcare facilities should verify the process they use, employing the actual equipment and operators that routinely process the equipment.
- The final responsibility for verification of sterilization techniques lies directly with the hospital. To ensure the efficacy of hospital processing, all cycles and methods should be verified for different sterilization chambers, wrapping methods, and/or various loading configurations.
- If wet trays or equipment are discovered after sterilization, a change in the product load configuration or a longer dry time may be necessary.
- The packaging material will maintain the sterility of the equipment after exposure.

## Steam Sterilization

- Enclose the equipment using a sterilization wrap that is suitable for the equipment, such as a grade 500 or higher, before sterile processing. Follow the AAMI and Association of periOperative Registered Nurses (AORN) recommended guidelines for appropriate wrapping configurations.
- Perform one of the following validated steam sterilization cycles:

Wrapping Method	Cycle	Sterilization Temperature	Minimum Exposure Time	Minimum Dry Time
Double Wrapped	Dynamic Air Removal (Pre-vacuum)	132 °C [270 °F]	4 minutes	30 minutes
		134 °C [273 °F]	3 minutes*	30 minutes
	Gravity	121 °C [250 °F]	30 minutes	30 minutes
		132 °C [270 °F]	15 minutes	30 minutes

\*Minimum exposure time may be extended to 18 minutes.

## Immediate-Use Steam Sterilization

### WARNINGS:

- DO NOT place equipment into a sterilization case or sterilization container for immediate-use steam sterilization.
- After cooling to room temperature, ALWAYS immediately use equipment sterilized by immediate-use steam sterilization. Sterility assurance cannot be maintained for unwrapped equipment.

**CAUTION:** Stryker does not recommend immediate-use steam sterilization for routine sterilization of medical devices. Immediate-use steam sterilization should only be used when individual devices require immediate sterilization and use.

Perform one of the following validated steam sterilization cycles:

Wrapping Method	Cycle	Sterilization Temperature	Minimum Exposure Time	Minimum Dry Time
Un-wrapped	Dynamic Air Removal (Pre-vacuum)	132 °C [270 °F]	4 minutes	No dry time
		134 °C [273 °F]	3 minutes	No dry time
	Gravity	132 °C [270 °F]	10 minutes	No dry time

# Storage and Handling

## Sterile Equipment

### WARNINGS:

- ALWAYS transport equipment with care to prevent damaging the sterile barrier.
- ALWAYS store processed equipment in a controlled environment and avoid extremes in temperature and moisture.
- Excessive handling of equipment will increase the likelihood of damaging the sterile barrier and may lead to contamination.

## Disposal/Recycle

Follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

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**stryker**<sup>®</sup>

## Signature History Report

Version	Signature	Role	Event Date	Vote
E	Sidow, Brian	CA_Owner	2017-07-10 14:38:03 EDT	Continue
E	Pearson, Dawn	CA_Tech Pubs-Tech Writer	2017-07-11 08:44:37 EDT	Approve
E	Khurmi, Paminder	CA_Regulatory Affairs	2017-07-11 11:02:39 EDT	Approve
E	Sidow, Brian	CA_Clinical Sciences	2017-07-11 16:58:13 EDT	Approve
E	Atkins, Dilara	CA_Marketing Manager	2017-07-13 08:26:05 EDT	Approve
E	Adar Yaron	CA_Product Design Engineer	2017-07-14 16:01:38 EDT	Approve
E	Savard, Ryan	CA_Quality Engineer	2017-07-18 14:53:40 EDT	Approve
E	Lindah, Travis	CA_Product Safety	2017-07-20 09:32:12 EDT	Approve
E	Honeck, Holly (Contractor)	Editor	2017-07-20 10:40:56 EDT	Continue

NOTE: IF THIS REPORT IS BLANK, PLEASE SEE WINDCHILL ECN FOR SIGNATURES.

## Related Information

CR Number	CR0319338
CR Class	DMR-Commercially Released
Originating Entity	INST: Kalamazoo
Origin	Product Improvement
Initial CR Approval Date	06/27/2017 16:16:50
CN Number	CN171492
CA Number	CA206588
NPD Disposition Plan	N/A