<p>NOTE: This portion of the document should not appear on produced Labels or IFUs.</p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p><b><u>Dimensions:</u></b></p> <p>4.75 inch (width) x 9 inch Booklet</p>	<p><b><u>Color/Material/Finish:</u></b></p> <p>Black Graphics on White Background 60# Offset, 20# Bond or Equivalent</p>	<p><b><u>Label Stock:</u></b></p> <p>N/A</p>	
	<p><b><u>Print Location:</u></b></p> <p>SPR</p>	<p><b><u>Suppliers/Services:</u></b></p> <p>TransPerfect MDS</p>	<p><b><u>Description/Type:</u></b></p> <p>Instructions For Use</p>	
			<p><b><u>Part Number:</u></b></p> <p>0206-701-700</p>	<p><b><u>Rev.</u></b></p> <p>S</p>

**EN**

**IMPORTANT INFORMATION: Instructions for Cement Restrictors when supplied individually or with the BioPrep Bone Preparation Kits.** For use with Stryker Cement Injection Systems.

**User/Patient Safety**

**WARNINGS:**

- Read and understand these instructions for use. The BioPrep Bone Preparation Kit contains brushes and sponges which allow the medullary canal and acetabulum to be cleaned prior to placement of bone cement and/or prosthesis. It also contains sculps which may be used to remove excess cement around the prosthesis.
- Instructions for using a cement restrictor with the Stainless Steel Inserter (REF 0206-725-000) are shown below. Familiarization with the cement restrictor and inserter prior to use is important. For further information, contact your Stryker sales representative or Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.
- Upon initial receipt and before use, inspect the package for damage and confirm the integrity of the sterile barrier. DO NOT use any equipment if damage is apparent or the sterile barrier has been compromised.
- DO NOT reuse, reprocess, or repackage a single use device. A single use device is intended for a single use only. A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the single use device resulting in operational failure. Critical product information may be lost if the single use device is repackaged. Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.
- ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of surgical waste.

**Cement Restrictor Description and Specifications**

The cement restrictor prevents flow of cement distal to the restrictor and confines the cement in a closed space. The use of the restrictor is essential for the application of low viscosity cement by means of a cement applicator. Using the restrictor also prevents cement flow into the distal part of the medullary canal and facilitates subsequent removal of cement if this becomes necessary.

Cement restrictors may be supplied pre-assembled to the inserter.

An estimate of the canal diameter can be made based on the diameter of the last intramedullary reamer used.

The cement restrictor is made of polyethylene.

Patient chart labels are provided for restrictors.

For further information, contact Stryker customer service and request a copy of document REF 0206-701-709.


**Instructions**

**NOTE:** If desired, a vent tool may be used to create a vent hole in the end of the cement restrictor. Discard any debris.

1. Remove contents from package and handle using aseptic technique.

**NOTE:** If the cement restrictor is not provided pre-assembled to the inserter, thread the restrictor onto the inserter with a clockwise motion until bottomed, approximately 3-1/2 turns.

2. The inserter is marked with rings in 1 cm increments starting at 100 mm. Determine insertion depth by placing the tip of the prosthesis at the 2 cm mark on the inserter as shown in Figure 1.

**Figure 1**

The depth markers are accurate to within +/- 0.019 cm.

3. Using the inserter, guide the cement restrictor into the medullary canal until resistance is met. Then tap the top of the inserter with a mallet or similar device to reach the predetermined depth obtained in step 2. Refer to Figure 2.

**CAUTION:** Excessive bending, twisting, or pounding may cause the inserter to break.

**Figure 2**

4. When the correct depth has been reached, unscrew the inserter from the restrictor with a counterclockwise motion, approximately 3-1/2 turns. Withdraw the inserter from the medullary canal.

**NOTE:** If necessary, the cement restrictor may be removed by rethreading the inserter to the restrictor and pulling it out of the canal. The restrictor may be adjusted as many times as necessary.

5. Continue with cement application per standard surgical protocol.

ENGLISH (EN)  
ESPAÑOL (ES)  
DEUTSCH (DE)  
FRANÇAIS (FR)  
ITALIANO (IT)  
NEDERLANDS (NL)  
SVENSKA (SV)  
DANSK (DA)  
SUOMI (FI)

PORTUGUÊS (PT)  
NORSK (NO)  
POLSKI (PL)  
ΕΛΛΗΝΙΚΑ (EL)  
日本語 (JA)  
中文 (ZH)  
한국어 (KO)  
TÜRKÇE (TR)

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**stryker**<sup>®</sup>

## Signature History Report

Version	Signature	Role	Event Date	Vote
S	Gomaa Sara	CA_Owner	2015-12-02 08:44:13 EST	Continue
S	Pearson, Dawn	Editor	2016-01-08 13:09:42 EST	Continue
S	Gomaa Sara	CA_Owner	2016-01-11 11:46:18 EST	Continue
S	Reyes Santiago, Luis	CA_Quality Engineer	2016-01-12 10:04:43 EST	Approve
S	Morgan, Carianne	CA_Product Safety	2016-01-12 14:44:02 EST	Approve
S	Khurmi, Paminder	CA_Regulatory Affairs	2016-01-13 09:18:08 EST	Approve
S	Santana Gonzalez, Michael	CA_Product Design Engineer	2016-01-18 20:52:59 EST	Approve
S	Curtis, Megan	CA_Marketing Manager	2016-02-03 11:28:26 EST	Approve
S	Pulliam Sylvia	Document Control	2016-02-03 12:14:39 EST	Complete
S	Petyko, Holly (Contractor)	Editor	2016-02-03 13:59:57 EST	Continue
S	Gomaa Sara	CA_Owner	2016-03-31 22:16:30 EDT	Approve

NOTE: IF THIS REPORT IS BLANK, PLEASE SEE WINDCHILL ECN FOR SIGNATURES.

## Related Information

CR Number	CR0264051
CR Class	DMR-Commercially Released
Originating Entity	INST: Kalamazoo
Origin	Process Improvement
Initial CR Approval Date	08/24/2015 15:34:48
CN Number	CN126853
CA Number	CA156463
NPD Disposition Plan	