<p><i>Note: This portion of the document should not appear on produced Labels or IFUs.</i></p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p>Dimensions: 4 inch (width) x 6 inch Booklet</p>	<p>Color/Material/Finish: Black Graphics on White Background 20# Bond or Equivalent</p>	<p>Label Stock: N/A</p>	
			<p>Description/Type: Instructions For Use</p>	
			<p>Part Number: 0206-725-707</p>	<p>Rev. F</p>

stryker[®] Stainless Steel Inserter

REF 0206-725-000

R_x ONLY

CE 0197

Instructions For Use

IMPORTANT INFORMATION: File in your records.

This insertion tool is designed to assist with the installation of Stryker Cement Restrictors.

User/Patient Safety

WARNINGS:

- Familiarization with this Inserter and Stryker Cement Restrictors prior to use is important. Read and understand these instructions and those supplied with the restrictors. If you require additional information, contact your Stryker sales representative or Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.
- Levering with the inserter within the intramedullary canal and applying undue pressure may cause damage to bone tissue and/or to the instrument.
- Prior to use, the inserter should be inspected for damage, especially at the threaded end. Do not use if damage is apparent.
- Upon initial receipt and before each use, clean and sterilize the equipment as indicated.

NOTE: The depth markers located along the length of the inserter are accurate to within +/- 0.013 cm.

Cleaning

WARNINGS:

- ALWAYS provide personal protective equipment (PPE) for processing personnel according to the instructions and safety data sheets (SDS) supplied with the detergent. Wear PPE at all times during processing.
- During transport to decontamination processing area, pay particular attention to sharp, cutting edges to avoid injury.
- Use filtered water for diluting detergents and for rinsing the equipment. Mineral residues from hard water can stain the equipment and/or prevent effective cleaning and decontamination. Poor water quality can adversely affect the life of medical devices. ALWAYS follow the water quality requirements per Association for the Advancement of Medical Instrumentation (AAMI) TIR 34.
- ALWAYS make sure the detergent solution is completely rinsed from the interior and exterior of the equipment before drying the equipment.

Stryker Instruments

4100 E. Milham
Kalamazoo, Michigan
(USA) 49001
1-269-323-7700
1-800-253-3210

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)

NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)
SUOMI (FI)
PORTUGUÊS (PT)
日本語 (JA)

- Read, understand, and follow the indications, instructions, and safety information supplied with the detergent for correct handling and use of the product. ALWAYS prepare the detergent solution according to the manufacturer's recommendations. Pay particular attention to the concentration used and the total dispersion.

CAUTIONS:

- DO NOT use saline to wet or soak the equipment before transport to the decontamination processing area.
- If transport to the decontamination processing area is delayed, cover the equipment with a damp cloth or spray the equipment with a pretreatment foam as often as necessary to maintain moisture. The pretreatment foam will minimize the drying of soil and facilitate later decontamination processing. DO NOT allow pretreatment foam to dry on the equipment.
- Avoid mechanical damage during transport. DO NOT mix heavy devices with delicate devices.

NOTE: Clean the equipment as soon as practical, typically within two hours, to preclude extended or repeat cleaning procedures.

1. Remove all gross soil from the equipment using absorbent wipes or a soft, lint-free cloth moistened with the prepared detergent solution.
2. Thoroughly clean the equipment. Use suitable brushes and detergent solution to clean all surfaces. Pay particular attention to the threaded end and any other rough surfaces, crevices, and difficult-to-reach areas where soil may be shielded from brushing.
3. Thoroughly rinse the equipment with warm, running water until all traces of detergent solution are removed. Hold the equipment on an incline, distal end pointing down. Pay particular attention to rough surfaces. A final rinse of the equipment using deionized or filtered water is recommended.
4. Visually inspect the equipment for any remaining soil or detergent solution. If soil or detergent solution remains, repeat the cleaning procedure using fresh detergent solution.
5. Allow the equipment to drain on absorbent wipes. If possible, orient the equipment vertically to assist in drainage.
6. Dry the equipment with a soft, lint-free cloth or medical-grade compressed air, or by heating the equipment in an oven below 110 °C [230 °F].

Sterilization**WARNINGS:**

- ALWAYS use a chemical indicator within every sterilization load to make sure the proper sterilization conditions of time, temperature, and saturated steam or hydrogen peroxide are achieved.

- ALWAYS use new sterilization wrap to enclose the equipment. Do not reuse sterilization wrap.
- Follow the recommended dry times to prevent moisture from accumulating inside the equipment. Moisture may prevent proper sterilization and/or damage the equipment.
- After cooling to room temperature, ALWAYS immediately use equipment sterilized by immediate-use steam sterilization. Sterility assurance cannot be maintained for unwrapped equipment.
- After sterilization, allow the equipment to cool to room temperature prior to use. Failure to comply may result in a burn injury and/or damage to the equipment.

CAUTION: Stryker does not recommend immediate-use steam sterilization for routine sterilization of medical devices. Immediate-use steam sterilization should only be used when individual devices require immediate sterilization and use.

NOTES:

- Stryker has validated several sterilization cycles for the sterilization of this equipment. However, sterilizer design and performance can affect the efficacy of the process. Healthcare facilities should verify the process they use, employing the actual equipment and operators that routinely process the equipment.
 - The final responsibility for verification of sterilization techniques lies directly with the hospital. To ensure the efficacy of hospital processing, all cycles and methods should be verified for different sterilization chambers, wrapping methods, and/or various loading configurations.
 - If wet trays or equipment are discovered after sterilization, a change in the product load configuration or a longer dry time may be necessary.
 - Validation is based on the AAMI protocol. International sterilization parameters are per the following standards: Australia/New Zealand per AS/NZS 4187; Netherlands per Field Standard for Loaner Instruments, Revision 03.02, April 2008; Europe and the United Kingdom per EN ISO 17664; Canada per CSA ISO 17664
1. Enclose the equipment using a sterilization wrap that is suitable for the equipment, such as a grade 500 or higher, before sterile processing. Follow the AAMI and AORN recommended guidelines for appropriate wrapping configurations.
 2. Perform one of the following validated steam sterilization cycles:

VALIDATED STEAM STERILIZATION CYCLE PARAMETERS

Wrapping Method	Cycle	Sterilization Temperature	Minimum Exposure Time	Minimum Dry Time
Wrapped	Dynamic Air Removal (Pre-vacuum)	132 °C [270 °F]	4 minutes	30 minutes
Unwrapped				No dry time
Unwrapped	Gravity	132 °C [270 °F]	10 minutes	No dry time

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Signature History Report

Version	Signature	Role	Event Date	Vote
F	Keeler, Stephanie	CA_Owner	2016-03-22 09:37:49 EDT	Continue
F	Pearson, Dawn	Editor	2016-03-22 09:55:38 EDT	Continue
F	Keeler, Stephanie	CA_Owner	2016-04-05 09:22:39 EDT	Continue
F	Jauw, James	CA_Quality Engineer	2016-04-05 09:31:40 EDT	Approve
F	Kulkarni Kedar	CA_Product Design Engineer	2016-04-05 11:38:14 EDT	Approve
F	Brown Donna	CA_Product Safety	2016-04-06 18:21:25 EDT	Approve
F	Cassar, Jason	CA_Marketing Manager	2016-04-07 08:29:39 EDT	Approve
F	Ashton, Kristi	CA_Regulatory Affairs	2016-04-11 16:41:30 EDT	Approve
F	Keeler, Stephanie	CA_Owner	2016-05-06 13:36:32 EDT	Withdraw
F	Keeler, Stephanie	CA_Owner	2016-05-06 14:11:59 EDT	Return to Editor
F	Pearson, Dawn	Editor	2016-05-10 12:24:24 EDT	Continue
F	Keeler, Stephanie	CA_Owner	2016-05-12 14:19:52 EDT	Continue
F	Jauw, James	CA_Quality Engineer	2016-05-12 14:28:02 EDT	Approve
F	Ashton, Kristi	CA_Regulatory Affairs	2016-05-13 07:47:16 EDT	Approve
F	Brown Donna	CA_Product Safety	2016-05-13 11:27:28 EDT	Approve
F	Sidow, Brian	CA_Clinical Sciences	2016-05-13 11:50:40 EDT	Approve
F	Harris, Lindsay	CA_Marketing Manager	2016-05-16 11:59:13 EDT	Approve
F	Kulkarni Kedar	CA_Product Design Engineer	2016-05-16 17:14:59 EDT	Approve
F	Campbell Desiree	Document Control	2016-05-17 08:22:55 EDT	Complete
F	Schrauben Haley	Editor	2016-06-16 11:17:29 EDT	Continue

Signature History Report

Version	Signature	Role	Event Date	Vote
F	Keeler, Stephanie	CA_Owner	2016-12-21 11:58:38 EST	Approve

NOTE: IF THIS REPORT IS BLANK, PLEASE SEE WINDCHILL ECN FOR SIGNATURES.

Related Information

CR Number	CR0229294
CR Class	DMR-Commercially Released
Originating Entity	INST: Kalamazoo
Origin	UDI
Initial CR Approval Date	01/15/2015 15:35:05
CN Number	CN142812
CA Number	CA163515
NPD Disposition Plan	