<p><i>Note: This portion of the document should not appear on produced Labels or IFUs.</i></p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p>Dimensions: 8.5 inch (width) x 11 inch Booklet</p>	<p>Color/Material/Finish: Black Graphics on White Background 20# Bond or Equivalent</p>	<p>Label Stock: N/A</p>	
			<p>Description/Type: Instructions For Use</p>	
			<p>Part Number: 0210-100-700</p>	<p>Rev. T</p>

stryker®

InterPulse® Irrigation System

Handpiece (Battery Power Pack)

REF 0210-XXX-000 Series


Handpiece (Rechargeable Power Pack)

REF 0211-XXX-000 Series

Instructions For Use

R_x ONLY

STERILE R

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)
NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)
SUOMI (FI)
PORTUGUÊS (PT)
NORSK (NO)
POLSKI (PL)
ΕΛΛΗΝΙΚΑ (EL)
日本語 (JA)
中文 (ZH)
한국어 (KO)

Introduction

This *Instructions For Use* manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, and surgical technologists.

Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

If additional information, especially safety information, or in-service training is required, contact your Stryker sales representative or call Stryker customer service.

Trademarks not the property of Stryker Corporation are the property of their respective owners.

Indications For Use

The Stryker InterPulse Irrigation System is used for wound debridement, soft tissue debridement, and cleansing of medical, clinical, or surgical sites. This includes cleansing of bone in surgical procedures, hydrodebridement of chronic wounds, debridement of loose skin from burns, and the cleansing of trauma wounds.

Contraindications

None known.

For Use With

This section describes system components that must be used with the equipment described in this manual to create a safe and effective system. The InterPulse Rechargeable Power Pack Handpiece requires the following system components:

DESCRIPTION	REF
InterPulse Rechargeable Power Pack	0211-022-000
InterPulse Power Pack Charger	0211-024-000

All InterPulse handpieces require the following system components:

DESCRIPTION	REF
InterPulse Tips or Wound Care Replacement Sets	0210-XXX-XXX Series
Stryker IV Pole (or equivalent)	0296-003-000

User/Patient Safety

WARNINGS:

- Before using any system component, or any component compatible with this system, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with the system components prior to use.
- Only trained and experienced healthcare professionals should use this equipment.
- Upon initial receipt, inspect the packaging and each component for damage. DO NOT use any equipment if damage is apparent or the sterile barrier has been compromised. See the *Inspection and Testing* section for inspection criteria.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.
- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment like this system. Install and place the system into service according to the EMC information in this manual. Portable and mobile RF communications equipment can affect the function of the system. See the *Specifications* section.

Accessories

This section describes system components that may be ordered to replace original equipment that is damaged, worn, or must be replaced. This section may also contain optional components used with the system.

WARNINGS:



- Use only Stryker-approved system components and accessories, unless otherwise specified.
- Using other electronic components and accessories may result in increased electromagnetic emissions or decreased electromagnetic immunity of the system.
- DO NOT modify any system component or accessory.
- DO NOT reuse, reprocess, or repackage a single use device. A single use device is intended for a single use only. A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the single use device resulting in operational failure. Critical product information may be lost if the single use device is repackaged. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.

The following Stryker-approved accessories are sold separately:

DESCRIPTION	REF
InterPulse Tips, sterile, disposable, 12 per package	0210-0XX-X00 Series
InterPulse Splash Shields, sterile, disposable, 12 per package	0210-034-000
InterPulse Wound Care Replacement Sets, sterile, disposable, 12 per package	0210-31X-X00 Series
IV Pole	0296-003-000

NOTE: For a complete list of accessories, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Features

A	Applied Part(s) – The tip assembly and handpiece body are the applied parts as defined by the manufacturer according to the standards listed under <i>Product Safety Certification</i> in the <i>Specifications</i> section.
B	Full Speed Lock Mechanism – The mechanism is used to hold the trigger at the high irrigation flow setting. For high flow, squeeze the trigger and push the button. To unlock, squeeze the trigger again.
C	Suction Clamp – The clamp is used to open and close the suction tubing pathway to control suction at the handpiece tip assembly.
D	Suction Tubing – The suction tubing provides a suction pathway from the handpiece tip assembly to the vacuum fluid collection canister connected to the facility vacuum source.
E	Irrigation Tubing/Power Cord Combination – The irrigation tubing provides a fluid pathway from the irrigation bag to the handpiece tip assembly. The power cord conducts electrical energy from the battery power pack to the handpiece.
F	Trigger – An activated trigger provides irrigation to the surgical/wound site.
G	Tip Locking Mechanism – The mechanism is used to engage and release the tip assembly for the handpiece.
H	Battery Power Pack Housing – The battery power pack housing stores the battery power pack and allows easy removal of the battery pack. The battery pack provides power to the handpiece. (REF 0210-XXX-000 Series)
I	Irrigation Tube Spike – The spike is used to pierce and connect the handpiece irrigation tubing to the irrigation bag.
J	Irrigation Tube Spike Cap – The protective spike cap is installed over the sharp irrigation tube spike to maintain sterility and prevent damage or injury.
K	Power Cord and Connector – The connector is installed into an InterPulse Rechargeable Power Pack. (REF 0211-XXX-000)
L	Suction Tube Adapter – The adapter is used for a large diameter port on a vacuum fluid collection canister and may be removed if necessary.

Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

SYMBOL	DEFINITION
	General warning sign
	Battery lid removal tabs (REF 0210-XXX-000 Series only)
	Type B Applied Part
	Direct current (DC)

Instructions

To Assemble the Handpiece

WARNINGS:

- For wound care applications performed outside the sterile field where the handpiece is used on a single patient repeatedly, **ALWAYS** disinfect the handpiece, and replace the tip assembly and suction tubing before each use. See the *Cleaning* section and the instructions for use supplied with the *Wound Care Tubing Replacement Sets*.
- DO NOT** reuse the handpiece for surgical applications or applications requiring a sterile field. See the instructions for use supplied with the *Wound Care Tubing Replacement Sets*.
- ALWAYS** use aseptic technique when handling the sterile handpiece, tip assembly, and irrigation tubing.
- ALWAYS** use aseptic technique when removing the irrigation tube spike cap. Failure to comply may cause the irrigation tube spike to become contaminated before insertion into the irrigation bag.
- DO NOT** allow fluid from the irrigation bag to drip onto a battery power pack. Failure to comply may cause electrical shock.
- DO NOT** submerge the battery power pack into liquid or irrigation solution.
- DO NOT** use a rechargeable power pack that gives off an odor, generates heat, becomes discolored, or leaks.

NOTES:

- Make sure a standard IV pole is equipped with an irrigation bag and a fully charged InterPulse Rechargeable Power Pack, if a power pack is required.
 - For optimal performance and maximum run time, use rechargeable power packs that have been charged within 48 hours.
 - Make sure the top of the fluid contained in the irrigation bag is less than 36 inches above the handpiece at all times after connecting the handpiece.
- If using a handpiece with a rechargeable power pack, install the handpiece power cord plug into the receptacle of the InterPulse Rechargeable Power Pack REF 0211-022-000 (figure 1). See the instructions for use supplied with the InterPulse Power Pack Charger.
 - Using aseptic technique, remove the irrigation tube spike cap from the handpiece irrigation tubing and install the irrigation tube spike into the irrigation bag aseptic peel port.

Instructions

To Assemble the Handpiece (continued)

NOTES:

- Using a regulated facility vacuum source and a fluid collection canister with the InterPulse Irrigation System is optional and depends on the InterPulse handpiece configuration and tubing assembly.
- If using suction, make sure a fluid collection canister is connected properly to a regulated facility vacuum source per local hospital protocol. Adjust the regulated facility vacuum source per local recommendations/regulations governing procedure-specific suction levels.
- To facilitate assembly, pull away and separate the irrigation tubing from the co-bonded suction tubing enough to allow easy access to the corresponding connection points.
- See the instructions for use supplied with the *Wound Care Tubing Replacement Sets* if suction tubing must be installed in the handpiece.

- Connect the handpiece suction tubing to the fluid collection canister, if required.

NOTE: The suction tube adapter may be removed to connect the handpiece suction tubing to a smaller diameter suction port on a vacuum fluid collection canister.

- Without touching the tip lock mechanism, install the appropriate tip assembly into the handpiece until an audible 'click' is heard. See the *Accessories* section for tip assembly information.


Figure 1 – To Assemble the Handpiece

To Operate Handpiece

WARNINGS:

- ALWAYS follow the current local recommendations and/or regulations governing infection control protocols for pulsatile lavage treatments. Failure to comply may result in the transmission of dangerous bacteria into the air and onto surfaces from wounds.
- ALWAYS place tubing sets away from personnel traffic areas to eliminate inadvertent tubing set kinking and/or handpiece movement.
- ALWAYS use aseptic technique when placing the handpiece in a safe and secure location when not in use.
- DO NOT place the handpiece on a patient.
- For wound care applications performed outside the sterile field where the handpiece is used on a single patient repeatedly, ALWAYS disinfect the handpiece, and replace the tip assembly and suction tubing before each use. See the *Cleaning* section and the instructions for use supplied with the *Wound Care Tubing Replacement Sets*.
- If using suction, ALWAYS follow current local recommendations/regulations governing procedure-specific suction levels to remove fluid waste safely from a surgical site.

- Squeeze the trigger enough to allow fluid to reach the end of the handpiece tip assembly before operating.
- To start irrigation, squeeze the trigger.

NOTE: If desired, push the full speed lock mechanism to hold the trigger in the full speed irrigation setting. To disengage the mechanism, squeeze the trigger.

- If using suction, slide the suction clamp to open and close the suction tubing pathway to control suction flow. If the handpiece fails to provide suction, see the *Troubleshooting* section.
- To stop irrigation and end the operation of the handpiece, release the trigger.

To Disassemble the Handpiece

WARNINGS:

- DO NOT reuse single use, disposable accessories. All disposable accessories are single use only.
- Handling biohazard waste is potentially dangerous. ALWAYS follow current local regulations governing biohazard waste to safely handle and dispose of surgical fluid waste.
- A tip assembly and suction tubing will contain surgical waste after use. ALWAYS handle these disposable accessories as potentially infectious materials after use. ALWAYS wear gloves and protective eye wear when removing and disposing of these accessories.

NOTE: Make sure all the surgical fluid waste has been removed from the tip assembly and suction tubing before disconnecting the suction tubing from the fluid collection canister.

- Press the tip locking mechanism and remove the tip assembly from the handpiece (figure 2).


Figure 2 – To Remove the Tip Assembly

- Dispose of the contaminated suction tubing and tip assembly properly.

NOTE: Before removing the irrigation tube spike from the irrigation bag, take measures to prevent the irrigation bag from leaking.

- Remove the irrigation tube spike from the irrigation bag for final handpiece disposal.
- If using a handpiece with a rechargeable power pack, remove the handpiece power cord plug from the power pack. Clean the power pack as required. See the instructions for use supplied with the InterPulse Power Pack Charger.
- Dispose of the handpiece properly. See the *Disposal/Recycle* section. If using a handpiece with a removable battery pack, see the *To Remove the Battery Pack* section.

Inspection and Testing

WARNINGS:

- DO NOT use any equipment if damage is apparent or the inspection criteria are not met.
- DO NOT disassemble or service this equipment, unless otherwise specified.

NOTE: For service, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

INTERVAL	INSPECTION CRITERIA	ACTION
Before Use	Check equipment for damage, wear, or missing components.	If damage is apparent, replace the equipment.
	Check power cord and tubing for cuts or punctures.	
	Check the handpiece and battery power pack housing for cracks.	

NOTE: If any component must be discarded, see the *Disposal/Recycle* section.

Cleaning (Between Wound Care Use - Single Patient Use Only)

WARNINGS:

- DO NOT immerse any system component in liquid.
- DO NOT sterilize any system component. See the instructions for use supplied with the InterPulse Power Pack Charger.

CAUTIONS:

- DO NOT use solvents, lubricants, or other chemicals, unless otherwise specified.
 - DO NOT allow liquids or moisture into the electrical connections of the rechargeable power pack or battery power pack.
1. Remove the tip assembly from the handpiece. Properly discard of the tip assembly.
 2. Wipe the handpiece with a soft cloth dampened with a nonabrasive, hospital disinfectant.
 3. Dry the equipment with a lint-free towel.
 4. Install new suction tubing and tip assembly to the handpiece. See the *Instructions* section.

Troubleshooting

WARNING: DO NOT disassemble or service this equipment unless otherwise specified.

NOTE: If you require service, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

PROBLEM	CAUSE	SOLUTION
The handpiece fails to provide suction.	The handpiece is not connected to a facility vacuum source or the vacuum source is adjusted incorrectly.	Make sure the handpiece is connected to a facility vacuum source and the vacuum regulator is adjusted correctly.
	The handpiece is damaged.	Replace the handpiece.
Sporadic electrical interference is experienced.	Electrical noise is present.	Turn off all the electrical equipment not in use in the room.
		Relocate the electrical equipment to maximize the distance between the equipment. Increase spatial distance.
		Plug equipment into different outlets.

NOTE: If any component must be discarded, see the *Disposal/Recycle* section.

Storage and Handling

CAUTION: ALWAYS store the equipment within the specified environmental condition values throughout its useful life. See the *Specifications* section.

To ensure the longevity, performance and safety of this equipment, use of the original packaging material is recommended when storing or transporting this equipment.

Disposal/Recycle

WARNINGS:

- ALWAYS follow the current local recommendations and/or regulations governing environmental protection and risks associated with recycling or disposing of the equipment at the end of its useful life.
- ALWAYS decontaminate the equipment and battery pack exposed to infectious material before sending them to a waste treatment facility.

To comply with European Community Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU, this device should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.

To comply with European Community Batteries Directive 2006/66/EC, this device has been designed for safe removal of the batteries by the end user using the instructions provided herein. Infected units and batteries should be decontaminated before they are sent for recycling.

To Remove the Battery Pack (figure 3)

WARNINGS:

- DO NOT remove the battery pack in the presence of a patient.
 - DO NOT cut the irrigation tubing/power cord to remove the battery power pack from the handpiece. Failure to comply may cause electrical shock, excessive heat and/or sparks and result in injury and/or fire.
1. Push the two battery lid removal tabs and remove the lid from the battery power pack housing.
 2. Disconnect the leads.
 3. Remove the battery pack.
 4. Recycle or dispose of the battery pack and equipment properly.


Figure 3 – To Remove the Battery Pack

Specifications

Model:	REF 0210-XXX-000 Series	REF 0211-XXX-000 Series
Description:	Handpiece (battery power pack)	Handpiece (rechargeable power pack)
Dimensions: (handpiece)	18.4 cm [7.25 inch] length 3.1 cm [1.23 inch] width 2.7 cm [5.0 inch] height	
Mass:	0.78 kg [1.71 lb] (handpiece with battery power pack)	0.51 kg [1.13 lb] (handpiece only)
Tubing Length:	2.82 m [9.25 ft] (irrigation) 3.05 m [10 ft] (suction)	
Electrical Power:	12 V $\overline{\text{---}}$, eight AA 1.5 V alkaline batteries, internally powered	InterPulse Rechargeable Power Pack, 7.9 V $\overline{\text{---}}$
European Conformity:	  0197 REF 0210-XXX-XXX Series (including E suffix), excluding: REF 0210-100-000S4, REF 0210-114-000S5, REF 0210-165-000, REF 0210-218-100, REF 0210-900-000, and REF 0210-918-000.	  0197 REF 0211-100-000, REF 0211-100-000E, REF 0211-110-000, REF 0211-110-000E, REF 0211-022-000, and REF 0211-024-000 only.
Ingress Protection (IP):	IPX0	
Equipment Type:	 Type B Applied Part	
Mode of Operation:	Continuous	

Product Safety Certification:

For REF 0211-1XX-000 Series with or without suffix E, REF 0211-2XX-XXX Series with and without suffix E, and REF 0211-022-000 only:

Canadian Standards Association (CSA) International

Canadian Standards Association (CSA)

CAN/CSA-C22.2 No. 601.1-M90, *Medical Electrical Equipment – Part 1: General Requirements for Safety* (Adopted IEC 60601-1 2d Ed [90])

CAN/CSA C22.2 601.1S1-94, Supplement No 1-94 to CAN/CSA C22.2 601.1-M90

CAN/CSA C22.2 601.1B-98, Amendment 2 to CAN/CSA C22.2 601.1-M90

Underwriters Laboratories (UL)

UL 60601-1 (1st Ed), *Medical Electrical Equipment – Part 1: General Requirements for Safety*

For REF 0210-0XX-XXX, REF 0210-1XX-XXX, REF 0210-2XX-XXX, REF 0210-9XX-XXX, all E and S parts included; and REF 0207-050-190:

Canadian Standards Association (CSA) International

C US

International Electrotechnical Commission (IEC)

IEC 60601-1 A1: R:2012 (3.1 Ed), *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; Corrigendum 1 (2006); Corrigendum 2 (2007); Amendment 1 (2012)

Canadian Standards Association (CSA)

CAN/CSA-C22.2 No. 60601-1:08, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* (Adopted IEC 60601-1:2005 + CORR. 1)

CAN/CSA-C22.2 No. 60601-1:08 TC 2:2011 (Corrigendum 2), *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* (Adopted IEC 60601-1:2005 + CORR. 1)

CAN/CSA-C22.2 No. 60601-1:14, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* (3.1 Ed) (Adopted IEC 60601-1:2005 Edition 3.0 + AMENDMENT 1, 2012-07, MOD)

CAN/CSA-C22.2 No. 60601-1-6:11, *Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance – Collateral Standard: Usability*

American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)

ANSI/AAMI ES60601-1:2005 (IEC 60601-1:2005, MOD), *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*





ANSI/AAMI ES60601-1:2005/C1:2009, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance – Corrigendum C1*

ANSI/AAMI ES60601-1:2005/A2:2010, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance – Amendment A2*

ANSI/AAMI ES60601-1:2005/(R)2012, AND C1:2009 AND A2:2010(R)2012 (Consolidated text — Edition 3.1), *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* (IEC 60601-1:2005+A1:2012, MOD)

European Committee for Electrotechnical Standardization (CENELEC)

CENELEC/EN 60601-1:2006/A1:2013 (3.1 ED), *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance – incorporates Amendment A12*

Environmental Conditions:**Temperature Limitation:****Operation****Humidity Limitation:****Atmospheric Pressure Limitation:****Storage and Transportation****Guidance and manufacturer's declaration - electromagnetic emissions**

The InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series is intended for use in the electromagnetic environment specified below. The customer or the user of the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration - electromagnetic immunity

The InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series is intended for use in the electromagnetic environment specified below. The customer or the user of the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = [1.2] \sqrt{P}$ <p>150 kHz to 80 MHz</p> $d = [1.2] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = [2.3] \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  <p>(Non-ionizing electromagnetic radiation)</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	

NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series is used exceeds the applicable RF compliance level above, the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Specifications (continued)

Guidance and manufacturer's declaration - electromagnetic immunity			
The InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series is intended for use in the electromagnetic environment specified below. The customer or the user of the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±2 kV, ±4 kV, ±6 kV contact ±2 kV, ±4 kV, ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV on DC power lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	Not Applicable	Not Applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Not Applicable	Not Applicable

NOTE: U_T is the alternating current mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series			
The InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the InterPulse Irrigation System, REF 0210-XXX-000E Series and REF 0211-XXX-000E Series as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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CAN/CSA-C22.2 No. 60601-1-6:11, *Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance* – Collateral Standard: Usability

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





ANSI/AAMI ES60601-1:2005/C1:2009, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* – Corrigendum C1

ANSI/AAMI ES60601-1:2005/A2:2010, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* – Amendment A2

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CENELEC/EN 60601-1:2006/A1:2013 (3.1 ED), *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* – incorporates Amendment A12

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stryker[®]

Signature History Report

Version	Signature	Role	Event Date	Vote
T	Matos Perez, Ivan	CA_Owner	2016-10-06 16:06:19 EDT	Continue
T	Matos Perez, Ivan	CA_Product Design Engineer	2016-10-06 16:08:32 EDT	Approve
T	Stewart, Adam	CA_Product Safety	2016-10-06 16:10:09 EDT	Approve
T	Ashton, Kristi	CA_Regulatory Affairs	2016-10-06 16:29:43 EDT	Approve
T	McKay, David	CA_Marketing Manager	2016-10-07 10:13:28 EDT	Approve
T	Gonzalez Glorycelis	CA_Quality Engineer	2016-10-10 13:51:26 EDT	Approve
T	Pearson, Dawn	Editor	2016-10-10 16:10:05 EDT	Continue

NOTE: IF THIS REPORT IS BLANK, PLEASE SEE WINDCHILL ECN FOR SIGNATURES.

Related Information

CR Number	CR0303398
CR Class	DMR-Commercially Released
Originating Entity	GQO: Arroyo, PR
Origin	Other
Initial CR Approval Date	09/20/2016 16:47:38
CN Number	CN154745
CA Number	CA183913
NPD Disposition Plan	