# Neptune<sup>®</sup> 2 Waste Management System Disposable Manifolds



4-Port Manifold **REF** 0702-020-000

Specimen Collection 4-Port Manifold **REF** 0702-020-001

Single Port Manifold **REF** 0702-025-000

Instructions For Use

R<sub>x</sub> ONLY

ENGLISH (EN)  
ESPAÑOL (ES)  
DEUTSCH (DE)  
FRANÇAIS (FR)  
ITALIANO (IT)  
NEDERLANDS (NL)  
SVENSKA (SV)  
DANSK (DA)  
SUOMI (FI)  
PORTUGUÊS (PT)  
NORSK (NO)  
POLSKI (PL)  
ΕΛΛΗΝΙΚΑ (EL)  
日本語 (JA)  
中文 (ZH)  
한국어 (KO)

## Indications For Use

The Neptune 2 Waste Management System is intended to be used in the operating room, pathology, surgical centers, and doctor's offices to collect and dispose of surgical fluid waste as well as collect smoke generated from electrocautery or laser devices.

## Contraindications

The Neptune 2 Waste Management System is **contraindicated against:**

- Connection directly to **chest tubes**.
- Connection to **closed wound drainage systems**.

## Contact Information

For additional information, including safety information, or in-service training, contact your Stryker sales representative or call Stryker Neptune Customer Service at 1-800-550-7836.

## For Use With

The Stryker Neptune 2 Disposable Manifold (manifold) for use with the Neptune 2 Waste Management System, specifically the rover.

## Description

The manifold is a non-sterile, disposable device that provides a fluid path from the suction lines to the rover canister. The manifold is for single patient use only. Manifolds with specimen collection capability have lids that allow the retrieval of the specimen.

## Symbol Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

SYMBOL	DEFINITION
	GENERAL WARNING SIGN

SYMBOL	DEFINITION
	SINGLE PATIENT USE ONLY

## User/Patient Safety

### WARNINGS:

- Before using any system component, or any component compatible with this system, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with the system components prior to use. Failure to comply may result in patient and/or healthcare staff injury.
- Only trained and experienced healthcare professionals should use this equipment. Failure to comply may result in patient and/or healthcare staff injury.
- Upon initial receipt, inspect each component for damage. DO NOT use any equipment if damage is apparent. Failure to comply may result in patient and/or healthcare staff injury.
- Use only Stryker-approved system components and accessories, unless otherwise specified. Failure to comply may result in patient and/or healthcare staff injury.
- DO NOT disassemble, modify, service, or repair any system component or accessory, unless otherwise specified. Call Stryker Neptune Customer Service.
- The manifold is for SINGLE PATIENT USE ONLY. DO NOT sterilize or reuse. DO NOT reuse, reprocess, or repackage a single use device. A single use device is intended for a single use only. The single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the device resulting in operational failure. Critical product information may be lost if the device is repackaged. Failure to comply may lead to infection or cross-infection and result in patient or healthcare staff injury.

- ALWAYS follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

- ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of surgical waste.
- BLOODBORNE PATHOGEN HAZARD
  - ALWAYS wear personal protective equipment (PPE) when operating or handling this equipment.
  - ALWAYS leave tubing attached to the manifold and close unused ports during disposal.
  - ALWAYS follow local regulations regarding proper handling and disposal of biohazard waste.
 Failure to comply may cause infection and result in healthcare staff injury.

## Instructions

### WARNINGS:

- DO NOT use a manifold that has been dropped or mishandled (Figure 1). Failure to comply may cause leakage and result in patient and/or healthcare staff injury.
- DO NOT use a Y-connector to add more than one suction line to each port (Figure 2). Failure to comply may cause cross-contamination and result in patient injury.

Figure 1 Handle Manifold with Care


Figure 2 One Suction Line Per Port

**NOTE:** Failure to close all the unused manifold ports will decrease the vacuum level efficiency of the rover.

## To Install the Manifold

1. Align the keyed manifold with the slotted receptacle and push the manifold into the receptacle of the desired canister of the rover (Figure 3).
2. Rotate the manifold toward the lock symbol until the manifold can no longer be rotated and is locked.


Figure 3 Install and Lock Manifold

3. Attach suction tubing to the manifold port(s) (Figure 4). Use 5/16 in. [7.9 mm] inner diameter suction tubing or 1/4 in. [6.4 mm] inner diameter suction tubing with a 5/16 in. [7.9 mm] adapter, as required.


Figure 4 Attach Tubing To Manifold

4. Close all the unused manifold ports using their corresponding caps (Figure 5).


Figure 5 Close All Unused Ports

## To Remove and Dispose of the Manifold

1. With suction active, gather the suction tubing toward the manifold port to purge the tubing of fluid waste (Figure 6). DO NOT remove any attached suction tubing from the manifold (Figure 7).


Figure 6 To Gather Suction Tubing


Figure 7 Keep Tubing Attached During Disposal

2. Turn the SUCTION LIMIT SETTING dial on the rover to zero for the canister with the manifold to be removed.
3. Push the SUCTION button on the rover control panel to stop fluid suction.
4. Rotate the manifold until it is unlocked, then pull the manifold and attached suction tubing out of the receptacle, always keeping the manifold in a horizontal orientation (Figure 8).
5. Properly dispose of the used manifold and attached suction tubing always maintaining the manifold in a horizontal orientation.


Figure 8 To Remove Manifold and Tubing

## To Use the Manifold with Specimen Collection

If the retrieval of specimen is required, use a manifold with specimen collection capability.

**CAUTION:** Repeated removal of the lid from a manifold with specimen collection capability may damage the lid seal and reduce the vacuum level efficiency.

1. Remove the manifold's lid after the surgical procedure (Figure 9).
2. Remove the filter from the manifold to collect the specimen (Figure 10).


Figure 9 Remove Manifold Lid


Figure 10 Remove Filter From Manifold

3. Properly dispose of the used manifold after collecting the specimen.

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