The Stryker logo is displayed in white lowercase letters on a dark gray rectangular background. The word "stryker" is followed by a registered trademark symbol (®).

**Endoscopic
Hardware**
User Manual

English	2
Français	9
Deutsch.....	18
Italiano	27
Português	36
Español.....	45
Nederlands	54
Dansk.....	63
Suomi.....	71
Norsk	79
Svenska	87
Polski	95
Ελληνικά	104
Türkçe.....	113
Русский	121
简体中文	130
日本語	137
한국어	145

Before Using Your Instrument

Carefully unpack the instrument and check if any damage occurred during shipment.

Read this user manual thoroughly and be familiar with its contents prior to using this equipment.

Test this equipment prior to a surgical procedure. Ensure device has not been bent, and that any locking or attachment surfaces function properly. These devices were tested at the factory before shipment.

Attempt no internal repairs or alterations. See instructions for authorized service information.

Some hardware devices have device-specific locking mechanisms. There is no guarantee that instruments selected for minimum channel width alone will be compatible in combination.

1.1 Indications

Arthroscopic Hardware:

Arthroscopic Hardware consists of manual arthroscopic instruments intended for use in surgical procedures within a joint.

Scope Accessories:

Scope accessories are intended for use with a Stryker scope.

1.2 Contraindications

There are no known contraindications.

Warnings:

- Federal Law (USA) restricts this device to sale by or on the order of a physician only. The user must be a qualified physician having complete knowledge of the use of this equipment.
- Do not contact these instruments with any powered shaver blade or place instrument in the path of a surgical laser beam. Doing so can cause damage to the instrument which may lead to patient injury.
- This device is shipped non-sterile. This device must be cleaned and sterilized prior to the first use and after every subsequent use.
- Open all instruments with features that have an open/close position (such as jaws, stopcocks, lumens, etc.) prior to cleaning, disinfection, or sterilization. If they are not in the open position, incomplete sterilization may result.
- Use only the sterilization cycles outlined in this document. Using unspecified sterilization cycles may damage the device and/or result in incomplete sterilization.
- Failure to follow cleaning instructions may result in inadequate sterilization.

- The sterilization parameters presented in this document apply only to devices sterilized outside of a sterilization tray. When using a sterilization tray, consult the instructions provided with the tray for proper sterilization parameters.
- Wear appropriate protective equipment: gloves, eye protection, etc.

1.3 Reprocessing

These reprocessing instructions are provided in accordance with ISO 17664, AAMI TIR 12, AAMI ST79, and AAMI ST81. While they have been validated by Stryker as being capable of preparing the device for re-use, it remains the responsibility of the processor to ensure that the reprocessing, as actually performed (using equipment, materials, and personnel in the reprocessing facility), achieves the desired result. This normally requires validation and routine monitoring of the process. Stryker recommends users observe these standards when reprocessing medical devices.

1.3.1 Containment and Transportation

Reprocess the device as soon as reasonably practical following use.

1.3.2 Manual Cleaning

1. Disconnect the device from all other devices.
2. If applicable, disassemble the device.
3. **Wipe**
 - a. Wipe excess soil from the device with disposable paper towels.
4. **Soak in Enzymatic Detergent**
 - a. Prepare an enzymatic detergent solution¹ with lukewarm tap water, according to manufacturer's recommendations.
 - b. Wipe the entire surface of the device using a soft clean cloth dipped in the detergent solution.
 - c. Immerse the device in the detergent solution, ensuring the solution contacts all inner and outer surfaces of the device. Use a syringe to inject any lumens with 50 ml of solution.
 - d. Soak the device in the solution for 15 minutes.

5. **Brush**

Note: Use an appropriately sized brush.

- **For the interior of the device, use a bottle brush with bristles that extend the diameter of the device lumen.**
 - **For the exterior of the device, a 2 Spectra M16 Soft-Bristle brush is recommended.**
- a. Thoroughly brush the exterior of the device, focusing on mated and rough surfaces, with a soft-bristled brush.
 - b. Brush around all movable parts while moving them through all extreme opened and closed positions.
 - c. Inject all lumens, crevices, and mated surfaces five times with 50 ml of detergent solution.
 - d. With the device immersed in the detergent, thoroughly brush the interior and exterior of the device, using an appropriately sized bottle brush

for the interior and an appropriately sized soft-bristled brush for the exterior.

6. Rinse

- a. Remove the device from the detergent solution and rinse it with water² at ambient temperature until all visible detergent residue is removed.
- b. Flush all lumens, crevices, and mated surfaces five times with 50 ml of water.
- c. After all detergent residues are removed, rinse the device for an additional 30 seconds.
- d. Drain any excess water from the device by holding it at an incline.

7. Soak in Non-Enzymatic Detergent

- a. Prepare a non-enzymatic detergent solution³ in lukewarm tap water according to manufacturer's recommendations.
- b. Immerse the device in the detergent solution, ensuring the solution contacts all inner and outer surfaces of the device. Use a syringe to inject lumens with 50 ml of solution.
- c. Soak the device in the solution for 15 minutes.

8. Brush

Note: Use an appropriately sized brush.

- **For the interior of the device, use a bottle brush with bristles that extend the diameter of the device lumen.**
- **For the exterior of the device, a 2 Spectra M16 Soft-Bristle brush is recommended.**

- a. Thoroughly brush the exterior of the device, focusing on mated and rough surfaces, with a soft-bristled brush.
- b. Brush around all movable parts while moving them through all extreme opened and closed positions.
- c. Inject all lumens, crevices, and mated surfaces five times with 50 ml of detergent solution.
- d. With the device immersed in the detergent, thoroughly brush the interior and exterior of the device, using an appropriately sized bottle brush for the interior and an appropriately sized soft-bristled brush for the exterior.

9. Rinse

- a. Remove the device from the detergent solution and rinse it with water² at ambient temperature until all visible detergent residue is removed.
- b. Flush all lumens, crevices, and mated surfaces five times with 50 ml of water.
- c. After all detergent residues are removed, rinse the device for an additional 30 seconds.
- d. Drain any excess water from the device by holding it at an incline.

10. Dry

- a. Dry the device using a clean cloth. Filtered pressurized air⁴ can be used to assist drying.

11. Inspect

- a. Visually inspect the device, including all internal surfaces, for remaining soil. Use an endoscopic camera and endoscope if necessary to see the inner surface of the lumen.
- b. If visible soil remains, repeat the manual cleaning procedure, focusing on those areas.

- ¹ Cleaning was validated using Enzol at 1 oz/gal at 35-40°C.
² Cleaning was validated using reverse osmosis/deionized (RO/
³ Recommend cleaning using Prolystica 2x Concentrate Neutral Detergent at 1/8 oz/gal at 35-40°C.
⁴ Drying was validated using pressurized air at 20 psi.

1.3.3 Automated Cleaning

1. Disconnect the device from all other devices.
2. If applicable, disassemble the device.
3. **Wipe**
 - a. Wipe excess soil from the device with disposable paper towels.
4. **Soak in Enzymatic Detergent**
 - a. Prepare an enzymatic detergent solution¹ with lukewarm tap water, according to manufacturer's recommendations.
 - b. Wipe the entire surface of the device using a soft clean cloth dipped in the detergent solution.
 - c. Immerse the device in the detergent solution, ensuring the solution contacts all inner and outer surfaces of the device. Use a syringe to inject any lumens with 50 ml of solution.
 - d. Soak the device in the solution for 15 minutes.

5. Brush

Note: Use an appropriately sized brush.

- **For the interior of the device, use a bottle brush with bristles that extend the diameter of the device lumen.**
- **For the exterior of the device, a 2 Spectra M16 Soft-Bristle brush is recommended.**

- a. Thoroughly brush the exterior of the device, focusing on mated and rough surfaces, with a soft-bristled brush.
 - b. Brush around all movable parts while moving them through all extreme opened and closed positions.
 - c. Inject all lumens, crevices, and mated surfaces five times with 50 ml of detergent solution.
 - d. With the device immersed in the detergent, thoroughly brush the interior and exterior of the device, using an appropriately sized bottle brush for the interior and an appropriately sized soft-bristled brush for the exterior.
6. **Rinse**
- a. Remove the device from the detergent solution and rinse it with water² at ambient temperature until all visible detergent residue is removed.
 - b. Flush all lumens, crevices, and mated surfaces five times with 50 ml of water.
 - c. After all detergent residues are removed, rinse the device for an additional 30 seconds.
 - d. Drain any excess water from the device by holding it at an incline.

7. Automated Wash

- a. Place the devices in an automated washer, on an incline to aid in drainage.
- b. Program the washer with the following parameters, then activate the wash:

Pre-Wash	
Recirculation Time	2 Minutes

Water Temperature	Cold tap water
Detergent Type	NA

Enzyme Wash	
Recirculation Time	2 Minutes
Water Temperature	Hot tap water
Detergent Type	Enzymatic detergent ²

Wash 1	
Recirculation Time	2 Minutes
Water Temperature	60°C (set point)
Detergent Type	Non-Enzymatic detergent ³

Rinse 1	
Recirculation Time	2 Minutes
Water Temperature	Heated 60°C
Detergent Type	NA

Dry	
Recirculation Time	7 Minutes
Temperature	115°C
Detergent Type	NA

8. Dry

- a. Dry the device using a clean cloth. Filtered pressurized air⁴ can be used to assist drying.

9. Inspect

- a. Visually inspect the device, including all internal surfaces, for remaining soil. Use an endoscopic camera and endoscope if necessary to see the inner surface of the lumen.
- b. If visible soil remains, repeat the manual cleaning procedure, focusing on those areas.

¹Cleaning was validated using Enzol at 1 oz/gal at 35-40°C.

²Cleaning was validated using reverse osmosis/deionized (RO/DI) water.

³Recommended cleaning using Prolystica 2x Concentrate Neutral Detergent at 1/8 oz/gal at 35-40°C.

⁴Drying was validated using pressurized air at 20 psi.

1.3.4 Disinfection (optional)

1. Disinfect the device in a disinfecting solution that has one of the following active ingredients:
 - ≥ 2.4% glutaraldehyde (with a minimum soak time of 45 minutes at 25°C).¹
2. Prepare the disinfecting solution according to the manufacturer's instructions.

3. Disassemble and fully immerse the device in the disinfecting solution, filling all lumens and crevices for the required time and at the appropriate temperature.
4. Rinse and flush the device in a bath of PURW. Immerse the device, shake it underwater, and allow it to soak for one minute. Repeat the rinse and flush two more times, using a fresh bath each time.
5. Dry all parts with a sterile, lint-free cloth immediately after rinsing.

¹Cidex® Activated Dialdehyde Solution is validated for disinfection efficacy.

1.3.5 Thermal Disinfection (optional)

Follow the steps below to disinfect the device using thermal disinfection.

1. Place the instruments into the washer disinfectant.
2. Program the washer using the following parameters:

Thermal Disinfection
Recirculation Time: 1 minute
Water Temperature: 90°C

1.3.6 Autoclave (Steam) Sterilization

Warning: Drying time depends on several variables, including altitude, humidity, type of wrap, preconditioning, size of chamber, mass of load, material of load, and placement in chamber. Users must verify that drying time set in their autoclave yields dry surgical equipment.

1. Open all stopcocks and valves for each device.
2. If using a sterilization tray, place the device(s) in the tray and follow instructions provided with the tray. Use only trays that are approved for sterilization with steam.
3. Double wrap each device (or tray, if being used) prior to sterilization.
4. If sterilizing the devices without a tray, use the following parameters.

	Pre-vacuum (U.S)	Pre-vacuum (EU)
Wrapping	double*	double*
Temperature	132°C (270°F)	134°C (274°F)
Sterilization Time	4 minutes	3 minutes
Dry Time	60 minutes	60 minutes

*Kimberly-Clark K600 KIMGUARD® sterilization wrap was used during sterilization validation.

Warning: Immediate-use steam sterilization is intended only for emergency situations. In the event that immediate-use reprocessing is required, the following instructions should be used.

	"Immediate-Use" Pre-vacuum (U. S)	"Immediate-Use" Pre-vacuum (EU)
Wrapping	-	-
Temperature	132°C (270°F)	134°C (274°F)
Sterilization Time	4 minutes	3 minutes
Dry Time	-	-

1.4 Maintenance and Storage

Hardware should be protected from bending and crushing. Visible bends and dents may prevent device from functioning properly. If there is any visible damage to the device, please return the unit to Stryker for repair.

Hardware should be properly cleaned and dried prior to storage. If there is any particulate matter on the device, carefully remove it prior to storage. Stainless steel components must be kept clean and dry when not in use to prevent rusting.

1.5 Expected Service Life

The expected service life of this family of products is 478 reprocessing cycles.

1.6 Authorized Service

To ensure the restoration of your instrument to original specifications, repairs should be directed to Stryker Endoscopy. Please contact your sales representative or distributor for product service outside of the United States. Please contact your local sales representative or distributor for specific warranty information outside of the United States.

Authorized Service
(United States Customers)

If service is needed either during or after the warranty period: Contact Stryker Endoscopy at (800) 624-4422 or contact your local Stryker Endoscopy sales representative.

Package all the components carefully in the original shipping container if possible.

Ship the components prepaid and insured to:

Stryker Endoscopy Customer Service

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Stryker Endoscopy Customer Service
Attention: Repair Department
5900 Optical Court
San Jose, CA 95138

Symbol Definitions / Définition des symboles / Erläuterung der Symbole / Definizione dei simboli / Explicação dos símbolos / Definiciones de los símbolos / Verklaring van symbolen / Symbolforklaring / Symbolien selitykset / Symbolforklaringer / Symboldefinitioner / Definicje symboli / Ορισμοί συμβόλων / Simge Tanımları / Описание символов / 符号定义 / 記号の定義 / 기호 정의

<div style="border: 1px solid black; padding: 2px; display: inline-block;">REF</div>	<p>Catalog number • Référence catalogue • Katalognummer • Numero di catalogo • Número de catálogo • Número de referencia • Catalogusnummer • Katalognummer • Luettelonumero • Katalognummer • Katalognummer • Numer katalogowy • Αριθμός καταλόγου • Katalog numarası • Номер по каталогу • 产品目录编号 • カタログ番号 • 카탈로그 번호</p>
	<p>Complies with the requirements of directive 93/42/EEC • Conforme aux exigences de la directive 93/42/CEE • Entspricht den Anforderungen der EU-Richtlinie 93/42/EWG • Conforme ai requisiti della Direttiva 93/42/CEE • Em conformidade com os requisitos da directiva 93/42/CEE • Cumple los requisitos de la directiva 93/42/CEE • Voldoet aan de vereisten van Richtlijn 93/42/EEG • Overholder kravene i direktiv 93/42/EØF • Täyttää direktiivin 93/42/ETY vaatimukset • Oppfyller kravene i direktiv 93/42/EØF • Uppfyller kraven i direktiv 93/42/EEG • Zgodność z wymogami dyrektywy 93/42/EWG • Συμμορφώνεται με τις απαιτήσεις της οδηγίας 93/42/ΕΟΚ • 93/42/EEC sayılı yönerge gerekliliklerine uygundur • Соответствует требованиям директивы 93/42/EEC • 符合指令 93/42/EEC 的相应要求 • 93/42/EEC の指令要件に準拠 • 지침 93/42/EEC의 요구사항을 준수합니다</p>
<div style="border: 1px solid black; padding: 2px; display: inline-block;">EC</div> <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-left: 10px;">REP</div>	<p>Authorized representative in the European Community • Représentant autorisé pour la Communauté européenne • Bevollmächtigter Vertreter in der EG • Rappresentante autorizzato nella Comunità europea • Representante autorizado na Comunidade Europeia • Representante autorizado en la Comunidad Europea • Gemachtigde vertegenwoordiger in de Europese Gemeenschap • Autoriseret repræsentant i EU • Valtuutettu edustaja Euroopan yhteisössä • Autorisert representant i EU • Auktoriserad representant i EU • Autoryzowany przedstawiciel w Unii Europejskiej • Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα • Avrupa Topluluğundaki yetkili temsilci • Уполномоченный представитель в Европейском сообществе • 欧共同体授权代表 • 欧州共同体における認定代理店 • 유럽 공동체의 공인 대리점</p>
	<p>Federal (USA) law restricts this device to sale by or on the order of a physician • Conformément à la loi fédérale américaine, ce dispositif ne peut être vendu que par un médecin ou sur prescription médicale • Laut US-Bundesgesetz darf dieses Produkt nur durch einen Arzt bzw. auf ärztliche</p>

	<p>Anordnung verkauft werden • Le leggi federali degli Stati Uniti limitano la vendita di questo dispositivo solo da parte di un medico oppure su sua prescrizione • A lei federal dos EUA limita a venda deste dispositivo a um médico ou mediante prescrição médica • La ley federal de EE. UU. limita el uso de este dispositivo a un médico o bajo su autorización expresa. • Krachtens de federale wetgeving van de Verenigde Staten mag dit hulpmiddel uitsluitend door of op voorschrift van een arts worden gekocht • Efter amerikansk lovgivning må dette udstyr kun sælges og ordineres af læger • Yhdysvaltain liittovaltion lain mukaan tämän laitteen saa myydä ainoastaan lääkäri tai lääkärin määräyksestä • I henhold til amerikansk lovgivning skal dette produktet kun selges til lege eller etter forordning fra lege • Enligt federal lag (USA) får denna produkt endast säljas av läkare eller på läkares ordination • Prawo federalne Stanów Zjednoczonych zezwala na sprzedaż tego produktu wyłącznie lekarzowi lub na jego zlecenie • Η ομοσπονδιακή νομοθεσία (των Η. Π.Α.) επιτρέπει την πώληση της διάταξης αυτής μόνο σε ιατρό ή κατόπιν εντολής ιατρού • Federal (ABD) yasalara göre bu cihaz yalnızca bir doktor tarafından veya doktor talimatıyla satılabilir • В соответствии с федеральным законодательством США это устройство может быть продано только врачом или по его указанию • 美国联邦法律限制本设备仅能由医生销售或遵照医嘱销售 • 美国連邦法により、本装置の販売は医師による（または医師の指示に基づく）注文に限定されます • 미국 연방법에 따라 본 기구는 의사가 직접 판매하거나 의사의 지시 하에서만 판매할 수 있습니다</p>
	<p>Manufacturer • Fabricant • Hersteller • Fabbricante • Fabricante • Fabricante • Fabrikant • Producent • Valmistaja • Produzent • Tillverkare • Producent • Κατασκευαστής • Üretici • Производитель • 制造商 • 製造者 • 제조사</p>
	<p>Date of Manufacturer • Date de fabrication • Herstellungsdatum • Data di produzione • Data de fabrico • Fecha de fabricación • Fabricagedatum • Fremstillingsdato • Valmistuspäivämäärä • Produksjonsdato • Tillverkningsdatum • Data produkcji • Ημερομηνία κατασκευής • Üretim Tarihi • Дата производства • 制造日期 • 製造年月日 • 제조일</p>
	<p>Non-sterile • Non stérile • Nicht steril • Non sterile • Não esterilizado • No estéril • Niet-steriel • Usteril • Steriloimaton • Usteril • Osteril • Wyrób niesterylny • Μη αποστειρωμένο • Steril Değil • Нестерильно • 非无茵 • 非滅菌 • 비멸균</p>
	<p>Consult instructions for use • Consulter le mode d'emploi • Gebrauchsanweisung beachten • Consultare le istruzioni per l'uso • Consultar as instruções de utilização • Consulte las instrucciones de uso • Raadpleeg de gebruiksaanwijzing • Se brugsanvisningen • Katso käyttöohjeet • Les bruksanvisningen • Läs bruksanvisning • Sprawdzić w instrukcji stosowania • Συμβουλευτείτε τις οδηγίες χρήσης • Kullanım talimatlarına başvurun • Обратитесь к инструкции по использованию • 请查阅《使用说明》 • 使用説明書を参照してください • 사용 지침 참조</p>

	<p>Caution • Mise en garde • Vorsicht • Attenzione • Precaução • Precaución • Let op • Forsigtig • Muistutus • Forholdsregel • Försiktighet • Przestroga • Προσοχή • Dikkat • Предостережение • 注意 • 注意 • 주의 사항</p>
<div style="border: 1px solid black; padding: 2px; display: inline-block;">LOT</div>	<p>Batch Code • Code lot • Chargennummer • Codice di lotto • Código de lote • Código de lote • Partijcode • Batchkode • Eränumero • Batchkode • Batchnummer • Kod partii • Κωδικός παρτίδας • Parti Kodu • Код партии • 批号 • バッチコード • 배 치 코드</p>
<div style="border: 1px solid black; padding: 2px; display: inline-block;">QTY</div>	<p>Quantity • Quantité • Menge • Quantità • Quantidade • Cantidad • Hoeveelheid • Antal • Määrä • Antall • Antal • Ποσότητα • Adet • Количество • 数量 • 数量 • 수량</p>
<div style="border: 1px solid black; padding: 2px; display: inline-block;">Made in Germany</div>	<p>Country of origin • Pays d'origine • Herkunftsland • Paese di origine • País de origem • País de origen • Land van oorsprong • Oprindelsesland • Valmistusmaa • Opprinnelsesland • Ursprungsland • Kraj pochodzenia • Χώρα προέλευσης • Menşei ülke • Страна происхождения • 原产国 • 原産国 • 원산지</p>

Manufactured By:
Stryker Endoscopy
5900 Optical Court
San Jose, CA 95138 USA
1-800-624-4422

U.S. Patents: www.stryker.com/patents

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks:
the Stryker Logo. All other trademarks are trademarks of their respective owners or holders.