INSTRUCTION FOR USE

## **SPY-PHI Articulating Arm**

HH1001 Main Arm / HH1002 Table Clamp / HH1003 Imaging Head Holder

EN | FR | DE | IT | PT | PT-BR | ES | NL | CZ |  
DA | FI | NO | SV | PL | EL | TR | RU | JA | KO |  
ZH-CN | ZH-TW | RO

## Table of Contents

SPY-PHI Articulating Arm Instruction for Use [EN]	3
SPY-PHI Articulating Arm Mode d'emploi [FR]	14
SPY-PHI-Gelenkarm Bedienungsanleitung [DE]	26
Braccio articolato SPY-PHI Istruzioni per l'uso [IT]	39
Braço articulado SPY-PHI Instruções de utilização [PT]	51
Braço articulado SPY-PHI Instruções de uso [PT-BR]	63
Braço articulado SPY-PHI Instrucciones de uso [ES]	74
SPY-PHI scharnierende arm Gebruiksaanwijzing [NL]	86
Kloubové rameno SPY-PHI Návod k použití [CZ]	97
SPY-PHI bevægelig arm Brugsanvisning [DA]	108
SPY-PHI-nivelvarsi Käyttöohjeet [FI]	119
SPY-PHI-leddarm Bruksanvisning [NO]	130
SPY-PHI vridbar arm Bruksanvisning [SV]	141
Ramię przegubowe SPY-PHI Instrukcja obsługi [PL]	152
Αρθρωτός βραχίονας SPY-PHI Οδηγίες χρήσης [EL]	164
SPY-PHI Mafsallı Kol Kullanım talimatları [TR]	177
Шарнирно-сочлененная консоль SPY-PHI Инструкция по применению [RU]	188
SPY-PHI Articulating Arm 使用説明書 [JA]	200
SPY-PHI 굴절암 사용 설명서 [KO]	211
SPY-PHI 关节臂 使用说明 [ZH-CN]	222
SPY-PHI 多節臂 使用說明 [ZH-TW]	232
Instrucțiuni de utilizare braț articulat SPY-PHI [RO]	242

INSTRUCTION FOR USE

**SPY-PHI Articulating Arm**

HH1001 Main Arm / HH1002 Table Clamp /  
HH1003 Imaging Head Holder

Access Surgical Innovations  
141 California Ave, Suite 102  
Half Moon Bay, CA 94019 USA

**Distributed by:**

Novadaq Technologies ULC  
8329 Eastlake Drive, Unit 101  
Burnaby, BC V5A 4W2  
Canada


**Customer and Technical Service**

*Customer Support*

Phone: (800) 624 - 4422 Toll Free (USA only)  
Fax: (800) 729 - 2917 Toll Free (USA only)  
Email: [EndoCustomerSupport@stryker.com](mailto:EndoCustomerSupport@stryker.com)

*Technical Support*

Phone: (877) 478-7953 Toll Free  
Fax: (408) 754 2037  
Email: [endots@stryker.com](mailto:endots@stryker.com)

Access Surgical Innovations, B.V.  
Kantstraat 19. NL-5076 NP Haaren  
P.O. Box NL-5060 GA Oisterwijk

## SPY-PHI Articulating Arm

### Table of Contents

- 1. Definitions**
- 2. Introduction**
- 3. Intended Use**
- 4. Contraindications**
- 5. Warnings**
- 6. Possible Adverse Effects**
- 7. Storage**
- 8. Product Description and Use**
  - 8.1. Main Arm (HH1001)**
  - 8.2. Table Clamp (HH1002)**
  - 8.3. Imaging Head Holder (HH1003)**
- 9. Inspection**
- 10. Cleaning and Sterilization**
  - 10.1. Manual Cleaning Instructions**
  - 10.2. Automated Cleaning Instructions**
  - 10.3. Sterilization Instructions**
- 11. Complaint Handling/ Reporting**

## 1. Definitions

### Explanation of Symbols

Reference number

Batch Code

Product Quantity


Compliant with European  
Medical Device Regulation

Date of Manufacture

Manufacturer

Authorized Representative in  
the European Community

Product Supplied Non-Sterile

Indicates product is a Medical Device

## 2. Introduction

The SPY-PHI Articulating Arm consists of three components: Main Arm (HH1001), Table Clamp (HH1002) and Imaging Head Holder (HH1003).

## 3. Intended Use

The SPY-PHI Articulating arm is intended as an accessory to the SPY Portable Handheld Imaging System, HH9000, for use during surgical procedures or other applications where the user wishes to position the imaging head in a fixed location while imaging, instead of holding it by hand.

**Main Arm (HH1001):** This component is used to position and hold the Imaging Head Holder (HH1003).

**Table Clamp (HH1002):** This component is used to mount the Main Arm (HH1001) to the table and hold position.

**Imaging Head Holder (HH1003):** This component is used to hold the HH9000 Imaging Head.

## 4. Contraindications

This product is not intended for use except as indicated above.

## 5. Warnings

1. CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
2. The product must only be used by trained medical personnel capable of judging and controlling danger to patients.
3. The products will be delivered non sterile. Prior to the initial use, and each following use, the products must be cleaned, disinfected, and sterilized as well as checked for visible irregularities and malfunctions according to the instructions given in this Instructions for Use.
4. Components of the SPY-PHI Articulating Arm (including table clamp, main arm and imaging head holder) are only to be used with each other to make up the SPY-PHI Articulating Arm. Components of the SPY-PHI Articulating Arm may not be used with other manufacturer's products.
5. End of life is normally determined by wear and damage due to use.
6. Use of this component for any purpose, or in any manner other than those described here may cause component damage or failure which could result in serious patient injury or death. If needed, all SPY- PHI Articulating Arm metal components or fragments thereof can be located by means of an X-Ray.

7. To maintain intended clamping capacity of the table rail clamp, do not tighten the table rail clamping knob when the main arm column is not fully installed.
8. The Rail Clamps are not insulated. If grounding of patient is NOT ACCEPTABLE, DO NOT use the system.
- 9. DO NOT FORCE BLACK KNOB OF THE MAIN ARM PAST STOP WHEN OPENING/LOOSENING.**

## 6. Possible Adverse Effects

Delayed surgical procedure.

## 7. Storage

Components should be stored in a clean and dry area. Inspect each component prior to use for functionality and damage. When necessary, dispose of products in accordance with national regulations and approved hospital practices for surgical component disposal.

## 8. Product Description and Use

### 8.1 Main Arm (HH1001)

Provides three-point positioning of the imaging head into the surgical field. The articulating segments of the Main Arm are connected to a rigid column.

To position and lock in place:

1. Hold the distal end with one hand and loosen the central knob by turning it counterclockwise.
2. Position the imaging head as needed and turn the knob clockwise to lock in place. Check to see that the Main Arm is secure enough to hold its position.

### 8.2 Table Clamp (HH1002)

Provides an attachment point to the operating table side rail. Accepts and holds the column of the main arm in position and allows for full range of rotation.

To attach to the rail:

1. Open the rail clamping knob until the L-brackets fit over the rail, then tighten to secure.
2. Insert the column of the Main Arm into the column opening, adjust the column to the desired angle, and tighten the column clamping knob when the main arm is in the desired position.

### 8.3 Imaging Head Holder (HH1003)

Connects to the distal end of the main arm to provide an attachment point for the SPY-PHI imaging head.

To connect the imaging head holder to the distal end of the main arm:

1. Insert the imaging head hex shaft into the chuck until it is fully seated.
2. Twist the collar clockwise to lock the quick connect accessory into place. Check to see that the shaft is secured, and the collar is in the locked position.

*To disassemble: Reverse the above steps.*

To position the imaging head and lock in place:

3. Rest the neck of the draped imaging head against the thin rigid band of the holder (as shown).

4. With one hand, pinch the sliding u-piece until the imaging head is firmly held in place.

5. With your other hand, turn the knob to tighten and lock.

*To loosen: Reverse the above steps.*

## 9. Inspect Before Use

All products should be inspected throughout their lifetime to ensure proper function and performance. If product does not pass the inspection, do not use.

### Main Arm (HH1001)

1. Visually inspect component for damage or cracks.
2. Check to make sure that the arm is rigid at all three joints by turning central tightening knob clockwise.
3. Insert arm column into table clamp, turn column clamping knob clockwise and ensure that it holds securely.
4. Check quick connect distal end to make sure the collar locks and unlocks accessory securely.

### Table Clamp (HH1002)

1. Visually inspect component for damage or cracks.
2. L-Bracket on Table Clamp should clamp securely to table rail.
3. Insert Main Arm into column opening, turn central tightening knob clockwise and make sure the post is held securely inside the column opening.

### Imaging Head Holder (HH1003)

1. Visually inspect instruments for damage or cracks.
2. Quick connect shaft should freely engage and disengage in and out of the quick connect distal end

## 10. Cleaning and Sterilization

### NOTE:

1. *The color of the SPY-PHI Articulating Arm components may vary due to the anodizing process or alloy used. Shading or loss of color may also occur after sterilization. This is not a defect in the products or material and will not affect the performance of your high-quality product.*
2. *Automated cleaning is not suitable for instruments with ball joints. Such instruments should only undergo a manual cleaning prior to sterilization.*

**CAUTION:** Only cleaning and sterilization procedures as indicated in this IFU should be used for cleaning and sterilization.

### Supplies and Equipment Needed:

- Ultrasonic Cleaner
- FDA Cleared CSR Wrap
- pH Neutral Cleaner such as Prolystica manufactured by Steris Corp. or equivalent
- 3-1000 MILTEX Component Cleaning Brush, Nylon Bristles or equivalent (available from Integra Miltex)
- Chemical Indicator Tape

### 10.1 Manual Cleaning Instructions

**NOTE:** *Turn tightening knob of the main arm clockwise to tighten ball joint prior to placing arm in ultrasonic cleaner (sonicator).*

1. Rinse each component individually with a steady stream of warm tap water (16.9°C to 18.9°C) until gross contaminants are removed. This process should take approximately 1-2 minutes.
2. Place each component into an ultrasonic cleaner (sonicator) containing an enzymatic, pH neutral detergent solution and warm tap water (25°C to 35°C) prepared according to the detergent manufacturer's instructions and sonicate for 10 minutes.
3. In a wash container, large enough to immerse each component, prepare a wash solution using an enzymatic, pH neutral detergent solution and warm tap water (25°C to 35°C) prepared according to the detergent manufacturer's instructions.
4. Transfer each component to a manual wash container and fully immerse in the cleaning solution prepared in Step 3.

5. While still submerged, any visible contamination and debris should be removed by scrubbing each component with a soft nylon bristle brush until visibly clean, paying particular attention to hard-to-clean areas such as crevices and joints. This process should take approximately 1-2 minutes.
6. Rinse with flowing, cold deionized water (17°C to 21°C) for 30-60 seconds until no visible soil remains.
7. Dry each component using clean, absorbent, low lint wipes to remove excess rinse water.

## 10.2 Automated Cleaning Instructions

**CAUTION:** Use only washer/disinfector machines that have been validated in accordance with ISO 15883. **AUTOMATED CLEANING IS NOT SUITABLE FOR MAIN ARM. AUTOMATED CLEANING IS ONLY SUITABLE FOR TABLE CLAMPS AND IMAGING HEAD HOLDER.**

Perform pre-cleaning to remove gross contaminants as follows:

1. Prepare a wash solution using an enzyme pH neutral detergent, in a wash container with tap water (27°C to 33°C) using the minimum concentration recommended by the detergent manufacturer.
2. Submerge and soak in wash solution for a minimum of 1 minute.
3. While still submerged, remove visible soil by scrubbing with a soft nylon bristle brush for a minimum of 4 minutes until no visible soil is observed.
4. Rinse with flowing, cold deionized water (18.5°C to 20.2°C) for a minimum of 30 seconds for each component.
5. Load components into washer/disinfector in accordance with the manufacturer's instructions.
6. Arrange components with curved surfaces and cannulations facing downward to prevent pooling of water.
7. Operate the washer/disinfector cycle according to the manufacturer's instructions.

Recommended minimal washer /disinfector parameters:

	Temperature	Time
Heated Wash	60°C (140°F)	2 Min.
Heated Tap Water Rinse	60°C (140°F)	20 Sec.
Heated Deionized Water Rinse	82°C (180°F)	1 Min.
Forced Air Drying	116°C (240°F)	9 Min.

### 10.3 Sterilization Instructions

**NOTE:**

1. For components with moving parts, lubricate joints with a steam-permeable, water soluble component lubricant prior to sterilization.
2. Components should be sterilized in the opened or unlocked position. Central knob of the articulating arms must be opened for sterilization.
3. Components should be sterilized by standard cycles using steam with the procedures listed below.

Preparation for Sterilization:

*Note: Use only FDA Cleared CSR Sterilization Wrap for component packaging.*

- Each component should be wrapped in two layers of FDA Cleared CSR Wrap using the simultaneous double-wrapping equal fold technique.
- Chemical Indicator Tape should be used to secure packaging and for labeling the contents. Note: Chemical Indicator Tape will change color or display diagonal stripes when exposed to temperatures of 121°C (250°F).

*Note: Only use wrappers validated for use in double simultaneous wrapping.*

Two single-layer wrappers or one bonded double wrapper can be used. Simultaneous double wrapping refers to when two layers of wrap material are wrapped together simultaneously using the fold technique described below:

1. Place the edge of the wrapper parallel with the edge of the table. Place the device in the center of the wrapper parallel with the edge of the wrapper.
2. Fold the edge of the wrapper over the top of the contents, covering the entire item. Fold the edge back over itself (toward the technician) to form a cuff. This will facilitate aseptic opening of the pack when used.

3. Bring the upper edge of the wrap down completely over the contents and fold it back on itself to form another cuff, overlapping the original.
4. Fold the left edge of the wrapper snugly over the pack and back on itself to form a cuff.
5. Fold the right side of the wrapper over the pack, overlapping the previous fold and folded back to form a cuff, or tuck under and secure with indicator tape.

#### Gravity

Sterilize at 121°C (250°F) for 30 minutes with a 30-minute dry time.

#### Prevacuum United States Standards

Sterilize at 132°C (270°F) for 4 minutes with a 30-minute dry time.

#### Prevacuum EU Standards

Sterilize at 134°C (273°F) for 3 minutes with a 30-minute dry time.

Autoclave temperatures should not surpass 280°F (137°C), as the handle, insulation or other non-metallic parts may be affected. (Note: The steam autoclave manufacturer may be contacted to confirm appropriate temperatures and sterilization times.)

## 11. Complaint Handling/ Reporting

Any adverse event involving the SPY-PHI Articulating Arm should be reported to Stryker and the country specific regulatory authorities immediately. To report an event to Stryker, call:

(800) 624 - 4422 Toll Free (USA only) or Email:  
EndoCustomerSupport@stryker.com

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### Distribué par:

Novadaq Technologies ULC  
8329 Eastlake Drive, Unit 101  
Burnaby, BC V5A 4W2  
Canada


### Assistance client/ Assistance technique

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