
# stryker<sup>®</sup>

## 32" 4K Surgical Display

REF 0240-031-050

CE **Rx** ONLY

# Contents

English .....	EN-1
Français.....	FR-29
Deutsch .....	DE-59
Italiano.....	IT-89
Português.....	PT-119
Español .....	ES-149
Nederlands .....	NL-179
Dansk .....	DA-209
Suomi.....	FI-237
Norsk.....	NO-265
Svenska.....	SV-293
Polski.....	PL-321
Ελληνικά .....	EL-351
Türkçe .....	TR-385
Română .....	RO-415
Русский .....	RU-445
简体中文.....	CHS-479
日本語.....	JA-507
한국어.....	KO-535

# Table of Contents

<b>Warnings and Cautions .....</b>	<b>EN-3</b>
Warnings .....	EN-3
Cautions.....	EN-4
<b>About Your Device.....</b>	<b>EN-6</b>
Intended Use and Indications for Use .....	EN-6
Intended Conditions for Use.....	EN-6
Contraindications.....	EN-6
Package Contents .....	EN-7
Device Features.....	EN-8
<b>Setup.....</b>	<b>EN-12</b>
Connections.....	EN-12
Basic Video Setup.....	EN-14
<b>Operation.....</b>	<b>EN-15</b>
On-Screen Display (OSD).....	EN-15
OSD Menus.....	EN-17
Troubleshooting.....	EN-19
<b>Cleaning and Maintenance.....</b>	<b>EN-20</b>
General Recommendations.....	EN-20
Cleaning .....	EN-20
Preventive Maintenance.....	EN-20
Inspection.....	EN-21
Electrical Safety Testing .....	EN-21
Storage.....	EN-21
Expected Service Life.....	EN-21
Disposal .....	EN-21
<b>Technical Specifications .....</b>	<b>EN-22</b>
General Description .....	EN-22
Classification and Approvals.....	EN-23
Compliance .....	EN-23
<b>Electromagnetic Compatibility.....</b>	<b>EN-24</b>
<b>Symbols and Definitions .....</b>	<b>EN-28</b>

# Warnings and Cautions

Please read this manual and follow its instructions carefully. The words warning, caution, and note carry special meanings and should be carefully reviewed:

**Warning:** Indicates measures to avoid potential serious injury to the user and the patient.

**Caution:** Indicates risks to the equipment. Failure to follow cautions may result in product damage.

**Note:** Provides special information to clarify instructions or present additional useful information.

## Warnings

To avoid potential serious injury to the user and the patient, please note the following warnings:

1. Read this manual thoroughly and be familiar with its contents prior to using this device.
2. Federal law (United States of America) restricts this device to sale by, or on the order of, a physician.
3. Carefully unpack the device and check if any damage occurred during shipment.
4. This device is non-sterile and therefore should not be placed in the sterile field.
5. Do not place the device or any other heavy object on the power cord. Damage to the cable can cause fire or electric shock.
6. To avoid electric shock, avoid removing the bezel.
7. This device should not be used adjacent to or stacked with other devices. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
8. Ensure weight compatibility with mounting apparatus. Do not mount the display on a cart arm unless it is labeled with a maximum load  $\geq 11.9$  kg/26.3 lb.
9. Test this device prior to a surgical procedure. This device was fully tested at the factory before shipment.
10. Do not attempt internal repairs or adjustments not specifically detailed in this manual. Ensure that readjustments, modifications, and/or repairs are carried out by persons authorized by Stryker Endoscopy.
11. Do not put any object into the panel. If this occurs, unplug the device and have it checked by qualified personnel before operating it any further.
12. Use appropriate caution to prevent contact with fluids if the device is being used with a power supply in patient environments.
13. The use of cables and/or other accessories with this device, other than those specified, may result in increased emissions or decreased immunity of this device.
14. This device has no means to be incorporated in an IT-network in the clinical environment.
15. The protective screen is made of tested high-resistance PMMA (polymethyl methacrylate). Nonetheless there is the possibility that it may crack if subject to strong impacts. Evaluate and prevent the risk of possible breakages of the protective screen by correctly handling and positioning the monitor in the operating room.

## Cautions

1. To achieve grounding reliability, connect the device to an AC adapter that is connected to a hospital-grade power cord, and ensure the power cord is plugged into a grounded power outlet.
2. Do not sterilize the device, as the delicate electronics cannot withstand this procedure.
3. Use only the proprietary power supply for the display. If you use an extension cable, completely secure the connection to the DC power cord of the power supply.
4. Never operate the device immediately after transportation from a cold location to a warm location.
5. To connect to an international power supply, use an attachment plug appropriate for the power outlet, as outlined in the “Technical Specifications” section of this manual.
6. Unplug the device if it is not to be used for an extended period of time. To disconnect the power cord from the display, unscrew the plug first, then pull the cord out by the plug. Never pull the cord itself.
7. Do not expose the device to moisture or apply liquid cleaners directly to the screen. Spray the cleaning solution onto a soft cloth and clean gently. For further detail, refer to the “Cleaning and Maintenance” section of this manual.
8. Allow adequate air circulation to prevent internal heat buildup. Do not place the device on surfaces (rugs, blankets, etc.) or near materials (curtains, draperies) that may block the ventilation slots. The device is cooled by natural convection and has no fan.
9. Do not touch the patient with signal input or output connectors. Equipment with SIP/SOP connectors should either comply with IEC 60601-1 and/or IEC 60601-1-1 harmonized national standards or the combination should be evaluated for safety.
10. To ensure electromagnetic compatibility, refer to the “Electromagnetic Compatibility” section of this manual. The 32” 4K Surgical Display (0240-031-050) must be installed and operated according to the EMC information provided in this manual.
11. Pay close attention to the cleaning instructions in this manual. A deviation may cause damage.
12. Do not install the device near sunlight, excessive dust, mechanical vibration, or shock.
13. Do not position the device so that it is difficult to disconnect the power cord from the supply mains.
14. Do not operate with the glass device screen facing downward.
15. Handle the device with care. Do not strike or scratch the screen.
16. Changes or modifications not expressly approved by the party responsible for compliance could void the user’s authority to operate the device.
17. **Note: This device has been tested and found to comply with the limit for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. There is no guarantee that interference will not occur in a particular installation, which can be determined by turning the device off and on. The user is encouraged to try to correct the interference by one or more of the following measures:**


- **Reorient or relocate the receiving device.**
- **Increase the separation distance between the device.**
- **Connect the device to an outlet on a circuit different from that to which the other device(s) are connected.**
- **Consult the manufacturer or field service technician for help.**

The warranty is void if any of these warnings or cautions is disregarded.

This product contains electrical waste or electronic equipment. It must not be disposed of as unsorted municipal waste and must be collected separately.

# About Your Device

The 32" 4K Surgical Display is a wide screen LED surgical display that can support a maximum resolution of 4K (4096 x 2160). The display supports the following video inputs: digital RGB (DVI), HDMI 1.4 and HDMI 2.0. It supports serial communication via the RS232 port and SDC port. It also has two USB ports that each supply power (5 V, 1 A) for accessories and peripherals (not used for data transfer).

## Intended Use and Indications for Use

The 32" 4K Surgical Display is intended for video display during surgical procedures including arthroscopy (orthopedic surgery), laparoscopy (general and gynecological surgery), thoracoscopy, endoscopy (general, gastroenterological, and ENT surgery) and general surgery. The display is a non-sterile reusable device not intended for use in the sterile field. The display is intended for use by qualified physicians and qualified operating room personnel having complete knowledge of the surgical procedures being used.


## Intended Conditions for Use

This display is intended to be used in the near-patient environment.

## Contraindications

There are no known contraindications for this device.

# Package Contents


Reference	Part Number	Package Contents
1	0240-031-050	32" 4K Surgical Display
2	–	Cable Cover
3	–	(4) M4 x 16 mm VESA screws
4	0240-031-004	Medical Power Supply Model : BPM150S24F11 (Bridgepower)
5	–	Hospital-grade AC power cord

Reference	Optional Accessories
0240-031-004	Medical Power Supply Model: BPM150S24F11 (Bridgepower)
0240-031-051	32" 4K Surgical Display Cover
0240-030-951	15-ft. (5 pin) DC extension cable Model: 1501047*** (Bridgepower) ***: blank or 001~999
0240-030-952	75-ft. (5 pin) DC extension cable Model: 1501047*** (Bridgepower) ***: blank or 001~999

# Device Features

## Front panel

**1. SDC recording status**

LED shines red when SDC is in record mode (requires SDC compatibility; contact a Stryker representative for conditions).

**2. Display Screen**

Shows video image

**3. Power switch (soft)**

Powers the display ON and OFF.

**4. Power LED\***

**Indicates current status:**

Green - Display is powered on or is in screen saver mode

Blinking Green - Display is in sleep mode

**5. PIP**

Accesses the Picture in Picture adjustment menu.

**6. Bright**

Accesses the Brightness adjustment menu.

**7. Specialty**

Accesses the Specialty selection menu.

**8. Input**

Accesses the Input selection menu.

**9. Rotary control**

Accesses the on-screen and navigates through its function.

**10. Logo Light**


**Indicates current status:**

No illumination - Off

Continuous white illumination - On


Slowly fading white illumination - Sleep mode (see Sleep Timer settings for detail)

## Rear Panel

- 1. Accessory mount (3)** Provide an access point for mounting optional accessories. (accepts 1/4 - 20 screws)
- 2. VESA mounting holes** Provide access points for mounting the display using provided screws. (100x100 mm, 200x100 mm, accepts M4 x 16 mm screws)
- 3. Power switch (hard)** Turns the input DC power On and OFF.
- 4. Handles** Aid in display positioning.  
**Caution: The handles are not intended to bear the entire weight of the display.**
- 5. Cable cover hinges** Attach the bottom of the cable cover to the display.
- 6. Velcro straps** Straps aid in cable management.
- 7. Cable cover** Covers and conceals cables.
- 8. Cable cover clips** Attach the top of the cable cover to the display.
- 9. Accessory Mount holes** Provide access points for mounting optional video accessories. (accepts M4 x 16 mm screws)

## Cable Cover


### Installing the Cable Cover

1. Align the left and right hinges of the cable cover onto the bottom rear of the display.
2. Snap on the top section of the cable cover to the aligning clips.

### Removing the Cable Cover

1. Pinch the left and right clips and pull the cable cover towards you.
2. Remove the cable cover from the left and right hinges.

## Display Handles

**Caution:** The handles are intended to aid in positioning the display, not for transporting the display. The handles should not bear the full weight of the display.

### Removing the Display Handles

1. Using a 3 mm hex key, loosen the two M4 x 25 mm screws and gently pull the handle away from the display.

### Installing the Display Handles


1. Align the handle with the screw holes on the rear of the display.
2. Using a 3 mm hex key, install the two M4 x 25 mm screws to attach the handle.

# Setup

Stryker Endoscopy considers instructional training, or inservice, an integral part of this device. Your local Stryker Endoscopy sales representative will perform at least one inservice at your convenience to help set up your device and instruct you and your staff on its operation and maintenance. To schedule an inservice, contact your local Stryker Endoscopy representative after your device has arrived.

## Connections


Video input and output signals are connected to the rear of the display, as illustrated below:

1. **Power Connector (24V)**
2. **DVI**
3. **HDMI (HDMI 1.4)**
4. **HDMI 4K (HDMI 2.0)**
5. **RS-232 communication port**
6. **SDC device control port**
7. **ACC (Accessory) POWER**


## Connecting the Power Supply

1. Connect the power supply to the 24 V input on the display.
2. Connect the AC power cord to the power supply\*.
3. Connect the AC power, using the supplied hospital-grade power cord.
4. (Optional, not shown) Connect an extension cord between the power supply and display.
5. Install cable cover.

\* Power supply information: Model Number: BPM150S24F11, Manufacturer: Bridgepower Corp.

## Basic Video Setup

1. Route the video output 1 from the camera to the SDC DVI input.
2. Route the video output 1 from the SDC DVI output to the DVI input on the display.


**Note: Stryker recommends a backup connection directly from the camera to the display.**

# Operation

Operate the display using the rotary control and the four buttons located on the front panel. A list of the display controls and their functions is provided below.

## On-Screen Display (OSD)

### Accessing the On-Screen Display

#### To use the four front-panel buttons:

1. **PIP:** Press to activate Picture Mode (Picture in Picture, Picture by Picture, Picture on Picture)
2. **Bright:** Press to activate the Brightness adjustment menu.
3. **Specialty:** Press to activate the Specialty adjustment menu.
4. **Input:** Press to activate the Input selection menu.

**Note: Press and hold any of the four front-panel buttons to exit the OSD menu.**

#### Use the Rotary Control to navigate the on-screen menus once they are activated:

- **Push** — Accesses/selects on-screen display menu.
- **Turn Right/Left** — With the on-screen display menu activated, turning increases/decreases the value of the selected parameter.
- **Push and Hold** — Exits on-screen display menu.

## Operating On-Screen Display

The device OSD helps navigate through various device menus.

1. Press the Rotary Control to activate the OSD menu.
2. Rotate the Rotary Control to move up or down through the menu.  
The parameter is highlighted when selected.
3. Press the Rotary Control to enter the next level OSD.
4. Rotate the Rotary Control to increase or decrease the value of the selected parameter, or to make a selection on different options.
5. To exit the OSD menu screen from the second or third level OSD menu, select the Exit option. To completely exit the OSD, press and hold the Rotary Control. If no keys are pressed, the OSD will automatically exit after the factory-set predetermined time (the time is customizable).

## OSD Menus

### Specialty

Menu Item	Description	Range
<b>Color</b>	Choose between color profiles for Standard, Arthro A, Arthro B, Lap A, Lap B, ENV, ENT, PACS, Norm or GYN. Color temperature RGB and gamma adjustment are available on each profile except PACS and Norm.	—
<b>Red</b>	Red balance	-128 – 127
<b>Green</b>	Green balance	-128 – 127
<b>Blue</b>	Blue balance	-128 – 127
<b>Gamma</b>	Gamma value	1.0 – 2.5, S0, S1, S2
<b>Enhancement</b>	Color adjustment to enhance sharpness	Off, Low, Mid, High

### Brightness Settings

Menu Item	Description	Range
<b>Brightness</b>	Increase or decrease the brightness.	0 – 100
<b>Contrast</b>	Increase or decrease the contrast.	0 – 100
<b>Image Sharpness</b>	Set image sharpness.	1 – 10

### Image Effect

Menu Item	Description
<b>Scale Mode</b>	Choose scale mode between Fill All, Fill To Aspect, Fill All to Aspect, One To One, V-Fill or H-Fill.
<b>Freeze Frame</b>	Enable or disable the freeze frame.
<b>Mirror</b>	Enable or disable the mirror function.
<b>PIP</b>	Enable PIP (picture in picture) function.
<b>POP</b>	Enable POP (picture on picture) function.
<b>PBP</b>	Enable PBP (picture by picture) function.

## Advanced Settings

Menu Item	Description
<b>Key Lock</b>	Key lock on: Disables all key functions (except Specialty and Inputs selection) Press and hold the rotary knob to turn off key lock.
<b>Auto Source Select</b>	Scans inputs until an active video source is detected. Auto Source Select is disabled during PIP/POP/PBP mode.
<b>Sleep Timer</b>	On: The display enters sleep mode if no active video source is detected.
	Off : The display will not enter sleep mode.
	Timer: Set the time until the display enters sleep mode: 30 or 60 minutes
<b>OSD Control</b>	Controls OSD (On Screen Display) Menu Position, Background, and Timeout
<b>Restore Factory Settings</b>	Sets to factory default
<b>Com Port</b>	Selects which port to use for communications. Select SDC for SDC device control; select RS-232 for SPI router control.
<b>Converter power</b>	ON: Enables power to converter through HDMI port
	OFF: Disables power to converter through HDMI port
<b>Time Table</b>	Enable: Default signal detection and display output
	Disable: Recommended for detecting signal in iSuite setups

## Information

Menu Item	Description
<b>User Name Entry</b>	Enter custom name to be displayed in screen saver mode
<b>Serial Number</b>	Display device serial number
<b>Runtime</b>	Display current device total run time
<b>Input</b>	Display current input format

**Note: Actual on-screen display values may vary with updated version of the firmware and user settings.**

## Troubleshooting

Before returning your display for service, consult the troubleshooting list below:

Problem	Current Status	Remedy
<b>No picture</b>	Power LED on	Using the OSD Menu, adjust the brightness and contrast to maximum or reset them to their default settings.
		Check the power supply functionality.
	Power LED off	Ensure the power switch at the front and rear of the display are set to ON.
		Check if the AC power cord is properly connected to the AC adapter and outlet.
		Check that the power supply is fully connected and functioning properly.
	Power LED blinking	Display is in sleep mode. Connect ACTIVE sources or change INPUT to wake the display.
Converters not powering on	Check that the converter power setting is ON in the OSD menu.	
	Turn the hard power switch off and on again.	
<b>Abnormal picture</b>	Oversized, undersized or missing video	Adjust settings in Scale Mode.
		Wait a few seconds after initial sync of video signals, or power cycle the display.
<b>OSD Error message</b>	Resolution not supported	Ensure the video source has a compatible output resolution.

# Cleaning and Maintenance

## Warning

To avoid electric shock and potentially fatal injury, unplug the display and power supply from the electrical outlet before cleaning.

## Caution

- Do not spray cleaning liquid directly onto the display or the power supply as product damage may result. Spray on the cloth before wiping the unit.
- Do not immerse the display or power supply in any liquid as product damage will result.
- Do not use corrosive cleaning solutions to clean the display or power supply as product damage may result.
- Do not sterilize the display or power supply as product damage may result.

## General Recommendations

- Keep the display clean to prolong its operational lifetime
- LCD panel performances may deteriorate in the long-term. Periodically check that it is correctly operating.
- Periodically check the tightness of the VESA mount screws. If not sufficiently tight, the display may detach from the arm, which may result in injury or equipment damage.

## Cleaning

**Note: Take extra care when cleaning the display screen. Excess liquid or drips that enter the bottom of the screen or the power supply receptacle may result in product damage.**

To clean the display or the power supply:

1. If the display cover is in place, remove the cover prior to cleaning.
2. Apply standard disinfectant or mild detergent to a dry sterile cloth.
3. Wipe the display or power supply.
4. Visually inspect the external surface of the device for cleanliness, focusing on hard-to-reach areas. If visible soil remains, repeat steps 1 – 3.

## Preventive Maintenance

With the display disconnected from mains perform the following periodical check:

- Check the integrity of the power cord and inspect its routing, so that it is not under the risk of being punched or cut.
- Clean the area around the power plug, dust and liquids may result in fire.
- Clean the ventilation slot of the display, dust can obstruct the air flow and cause temperature increase of the electronics.

## Inspection

Inspect the device on a continual basis for unacceptable deterioration such as (but not limited to) corrosion, discoloration, pitting, cracked seals, or abnormal noises. If a problem is observed or suspected, the device should be returned for service.

## Electrical Safety Testing

Testing of ground protective earthing from the display is not necessary. Protective earthing is assured by the certified external power supply.

## Storage

Never store the device in a non-ventilated, humid environment. This can damage the delicate electronics in the device.

## Expected Service Life

The expected service life of the device is four years.

## Disposal

This product contains electrical waste or electronic equipment. It must not be disposed of as unsorted municipal waste and must be collected separately in accordance with applicable national or institutional related policies relating to obsolete electronic equipment.

Dispose of any system accessories according to normal institutional practice relating to potentially contaminated items.

# Technical Specifications

## General Description

Relative Humidity Item		Description	
LCD Panel	<b>Description</b>	31.1" (812 mm) Diagonal	
	<b>Native Resolution</b>	4096 (H) dots x 2160 (V) lines	
	<b>Display colors</b>	1,073,741,824 colors	
	<b>Pixel Pitch</b>	0.1704 mm x 0.1704 mm	
	<b>Response Time (typical)</b>	Rise time: 11 ms Fall time: 9 ms	
	<b>Viewing Angle</b>	horizontal/vertical: 178 °	
Brightness and Contrast	<b>Brightness</b>	525 cd/m <sup>2</sup>	
	<b>Contrast</b>	1500 : 1	
Input/Output		Input	Output
		1 x DVI 1 x HDMI (HDMI 1.4) 1 x HDMI 4K (HDMI 2.0) 1 x RS-232 (SPI router control) 1 x SDC device control interface (USB)	2 x USB (5 V, 1 A power only; not used for data transfer)
Temperature	<b>Operating</b>	50 ° – 104 °F (10 ° – 40 °C)	
	<b>Transport and Storage</b>	0 ° – 140 °F (-18 ° – 60 °C)	
Relative Humidity	<b>Operating</b>	25 – 75 %	
	<b>Transport and Storage</b>	15 – 90 %	
Electrical	<b>Power Adapter</b>	Input: 100 – 240 VAC; 50 – 60 Hz; 2.5 A Output: 24 V DC; 6.25 A (150 W Max) Model Number: BPM150S24F11	
	<b>Power Consumption</b>	20 – 100.8 Watts Marked Rating: 24V DC; 4.2A	

Relative Humidity Item		Description
Current/ Voltage Rating	110V ± 10V power outlets	Select a power supply cord that is UL Listed and C.S.A Certified, type SJT or SVT, 3 – conductor, 18 AWG, terminated in a molded on hospital grade plug cap rated 110 V ± 10 V, 15 A, with a minimum length of six feet.
	220V ± 20V power outlets	Select a power supply cord that is internationally harmonized and marked “<HAR>”, 3 – conductor, 0.75 mm <sup>2</sup> minimum wire, rated 220 V ± 20 V, 10 A with a PVC insulated jacket. The cord must have a molded on plug cap rated 220 V ± 20 V, 10 A. The cord and plug cap must be suitable for medical use.
Weight with Cable Cover (approximate)		23.1 lb (10.5 kg)
Unit Dimensions (W x H x D)		756.7 x 453 x 77.2 mm 29.7 x 17.8 x 3 in
VESA Mounting Interface Dimensions		VESA 100 x 100 mm
		VESA 200 x 100 mm
Accessory Mounting Interface Dimensions		75 x 75 mm (x 2)

## Classification and Approvals

Class 1 Equipment

Medical equipment with respect to electric shock, fire, and mechanical hazards only in accordance with ANSI/AAMI ES60601-1 and CAN/CSA C22.2 No. 60601.1.

IP23: Protection against access to hazardous parts from fingers or similar objects, protection from ingress of spraying water (less than 60° from vertical)

Continuous Operation

## Compliance

FCC Regulations: FCC Part 15 Class B

FCC Identifier: QVXAMM320ES

Canadian Notice: CAN ICES-3 (B)/NMB-3(B)

CE medical device class I

**Note: Please contact your local Stryker Endoscopy sales representative for information on changes and new products.**

**Notice to the user and/or patient: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

# Electromagnetic Compatibility

Like other electrical medical equipment, the 32" 4K Surgical Display requires special precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC), the display must be installed and operated according to the EMC information provided in this manual. The display has been designed and tested to comply with IEC 60601-1-2 requirements for EMC with other devices.

## Warning

When this device is connected with other electrical equipment, leakage currents may be additive. To minimize total leakage current per patient, ensure that all systems are installed according to the requirements of IEC 60601-1-1.

## Caution

Portable and mobile RF communications equipment may affect the normal function of the display.

Do not use cables or accessories other than those provided with the display, as this may result in increased electromagnetic emissions or decreased immunity to such emissions.

If the display is used adjacent to or stacked with other equipment, observe and verify normal operation of the display in the configuration in which it will be used prior to using it in a surgical procedure. Consult the tables below for guidance in placing the display.

### Guidance and Manufacturer's Declaration: Electromagnetic Emissions

The 32" 4K Surgical Display is intended for use in the electromagnetic environment specified below. The customer or the user of the display should ensure it is used in such an environment.

Emissions test	Compliance	Electromagnetic Environment - guidance
RF emissions CISPR 11	Group 1	The display use RF energy only for their internal function; therefore, their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The display is suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This system is intended for use by health care professionals only. This system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the system or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ flicker emissions IEC 61000-3-3	Complies	

**Guidance and Manufacturer's Declaration: Electromagnetic Immunity**

The 32" 4K Surgical Display is intended for use in the electromagnetic environment specified below. The customer or the user of the display should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV (input a.c. power port) ± 1 kV (signal ports)	± 2 kV (input a.c. power port) ± 1 kV (signal ports)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	(± 0.5, ± 1, ± 2) kV; line to earth (± 0.5, ± 1) kV; line to line	(± 0.5, ± 1, ± 2) kV; line to earth (± 0.5, ± 1) kV; line to line	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips and interruptions  IEC 61000-4-11	<ul style="list-style-type: none"> <li>• 0% <math>U_r</math>, 0.5 cycle</li> <li>• 0% <math>U_r</math>, 1 cycle</li> <li>• 70% <math>U_r</math>, 25 cycles</li> <li>• 0% <math>U_r</math>, 5 Sec</li> </ul> <p>Supplementary information: If the Rated voltage range &lt;25% of the lowest rated input voltage, one rated input voltage. Otherwise, minimum and maximum rated voltage. EUT powered at one of the Nominal input frequencies. ME EQUIPMENT and ME SYSTEMS with power input voltage selection by transformer taps shall be tested at only one tap setting.</p>	<ul style="list-style-type: none"> <li>• 0% <math>U_r</math>, 0.5 cycle</li> <li>• 0% <math>U_r</math>, 1 cycle</li> <li>• 70% <math>U_r</math>, 25 cycles</li> <li>• 0% <math>U_r</math>, 5 Sec</li> </ul> <p>Supplementary information: If the Rated voltage range &lt;25% of the lowest rated input voltage, one rated input voltage. Otherwise, minimum and maximum rated voltage. EUT powered at one of the Nominal input frequencies. ME EQUIPMENT and ME SYSTEMS with power input voltage selection by transformer taps shall be tested at only one tap setting.</p>	Mains power quality should be that of a typical commercial or hospital environment. If the user of the transmitter requires continued operation during power mains interruptions, it is recommended that the Wireless Transmitter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power-frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<p><b>Note:</b> <math>U_r</math> is the AC mains voltage prior to application of the test level.</p>			

**Guidance and Manufacturer's Declaration: Electromagnetic Immunity**

The 32" 4K Surgical Display is intended for use in the electromagnetic environment specified below.  
The customer or the user of the display should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms</p> <p>10 V/m 80MHz to 2.7 GHz</p>	<p>3 Vrms</p> <p>10 V/m 80MHz to 2.7 GHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the display, including their cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended Separation Distance: <math>d = 2\sqrt{P}</math> 80 MHz to 2.7 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range(b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>(a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 32" 4K Surgical Display is used exceeds the applicable RF compliance level above, the display and transmitter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the display.</p> <p>(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Tested specifications for immunity to RF wireless communications equipment						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						
<b>Note:</b> Portable RF Communication equipment should be used no closer than 30cm to 32"4K Surgical Display. Otherwise, degradation of the performance of this equipment could result.						

# Symbols and Definitions

The following symbols appear on the product, its labeling, or the product packaging. Each symbol carries a special definition, as defined below:

Consult Instructions for Use

Attention: See Instructions for Use

Serial Number

Quantity

Authorized European Representative

Country of Origin

Catalog number

Direct Current

Date of Manufacture

Legal manufacturer

Do Not Get Device Wet

DC Power Control Switch

Maximum Stacking

This Side Up

Fragile

For Indoor Use Only

No Serviceable Parts

Efficiency Level

Japan PSE Mark Denan

Device recycling code (applicable in China)

UL, Functional Safety Recognized Component

Medical Device

China Compulsory Certificate Mark

Unique Device Identification

Federal law restricts this device to sale by or on the order of a physician

Complies with the requirements of EC directives/regulations

**IP23**

Enclosure provides protection against:

- Access to hazardous parts from fingers or similar objects
- Spraying water (less than 60° from vertical)

Medical Equipment is in accordance with AAMI ES60601-1 and CSA C22.2 No. 60601-1-08-CAN/CSA in regards to electric shock, fire hazards, and mechanical hazards

This product contains electrical waste or electronic equipment. It must not be disposed of as unsorted municipal waste and must be collected separately.

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Stryker or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: the **Stryker logo**. All other trademarks are trademarks of their respective owners or holders.
Barco, Inc.