**Connected OR Cart 230V / Carro Connected OR de 230 V /  
Connected OR 流动车 230 V / Connected OR 카트 230 V /  
Xe đẩy Connected OR 230 V**

Instructions for Use / Instrucciones de uso / 使用说明 / 사용 설명서 / Hướng dẫn sử dụng

**REF** 240-099-156

**R<sub>X</sub> ONLY**

# Contents

Overview of the Connected OR Cart 230V .....	2
Intended Use .....	4
Warnings and Cautions .....	5
Cautions .....	5
Warnings .....	5
Using the Cart .....	6
Putting the Cart into Transport Position .....	6
Moving the Cart .....	6
Powering the Cart On .....	8
Powering the Cart Off .....	8
Using the I/O Panel .....	8
Maintenance and Service .....	9
Cleaning the Cart .....	9
Adjusting the Connected OR Mount (Monitor Arm) .....	10
Resetting the Breaker Switch .....	12
Replacing the Power Cord .....	12
Replacing Monitor Cables .....	14
Disposal .....	14
Electromagnetic Compatibility .....	15
China RoHS Declaration .....	19
Symbols .....	20
Technical Specifications .....	21

# Overview of the Connected OR Cart 230V

The Connected OR Cart is a mobile, medical-grade cart designed to accommodate a combination of Stryker Advanced Imaging Modality (AIM) System, PINPOINT Endoscopic Fluorescence Imaging System, or SPY Portable Handheld Imaging (SPY-PHI) System components. The cart is intended to be used during or between surgical or clinical procedures. The Connected OR Cart provides adjustable shelves, an articulated monitor arm, lockable steering wheels, handles for transporting and positioning, isolated electrical power, a central power switch, an external I/O panel for network or video connections, and holders for accessories such as CO<sub>2</sub> tanks, imaging heads, and foot switches.


Figure 1 – Connected OR Cart overview, right side

- |  |                                 |
|--|---------------------------------|
| 1. Display monitor   | 6. SPY-PHI cable hook           |
| 2. Monitor stow and positioning handle                                 | 7. Cart transporting handle     |
| 3. Connected OR Mount monitor arm                                      | 8. Cart positioning handle      |
| 4. Monitor stow bracket (for use when the monitor is in stow position) | 9. Adjustable shelves           |
| 5. SPY-PHI imager holder   | 10. CO <sub>2</sub> tank holder |


Figure 2 – Connected OR Cart overview, left side

- |                                   |                       |
|-----------------------------------|-----------------------|
| 1. Monitor positioning handles    | 5. Foot switch holder |
| 2. Connected OR Mount monitor arm | 6. IFU holder         |
| 3. Cart power button              | 7. Rear wheel         |
| 4. Accessory basket               | 8. Front wheel        |


Figure 3 – Connected OR Cart overview, rear

- |                             |  |
|-----------------------------|--|
| 1. I/O panel                | 4. Rear door latch                       |
| 2. Rear door                | 5. Power cord hooks                      |
| 3. Transport position label | 6. Breaker switch location (inside cart) |

### Intended Use

The Cart is intended to store and transport components of minimally invasive and open-field imaging systems such as the Stryker Advanced Imaging Modality (AIM) System, PINPOINT Endoscopic Fluorescence Imaging System, or SPY Portable Handheld Imaging (SPY-PHI) System) during intraoperative use. The Cart also provides power, video, and data connections to the components.

The Cart should be operated only by qualified personnel such as surgeons and clinic support staff.

There are no known contraindications for the Connected OR Cart.

## Warnings and Cautions

Please read and follow these Instructions for Use carefully. The words warning, caution, and note carry special meaning and should be carefully reviewed:

Warning	Indicates risks to the safety of the patient or user. Failure to follow warnings may result in injury to the patient or user.
Caution	Indicates risks to the equipment. Failure to follow cautions may result in product damage.
Note	Clarifies the instructions or presents additional useful information.

### Cautions

To avoid potential damage to this device, please note the following cautions:

1. Only Stryker-approved components and accessories may be used with the Connected OR Cart.
2. Ensure the cart is in transport position before moving it.
3. Ensure that modifications and/or repairs are carried out by persons authorized by Stryker.

### Warnings

To avoid potential serious injury to the user and the patient and/or damage to this device, please note the following general warnings.

1. Read these Instructions for Use thoroughly, especially the warnings, and be familiar with its contents before using the cart.
2. Federal law restricts this device to sale by or on the order of a physician.
3. To avoid the risk of electric shock, only connect equipment explicitly specified and approved by Stryker to the multiple-socket outlet (power outlets) on the Connected OR Cart.
4. To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth.
5. Do not position the cart so that it is difficult to disconnect the power cord from supply mains.
6. Do not connect additional multiple-socket outlets or extension cords to the system.
7. The Connected OR Mount (monitor arm) may have pinch points: keep hands clear.
8. Failure to properly adjust the Connected OR Mount could result in equipment damage or injury.
9. No modification of this equipment is permitted. Do not modify this equipment without authorization of the manufacturer. If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
10. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

## Using the Cart


The Connected OR Cart can be moved to the desired location and positioned by one person. The Connected OR Mount articulated monitor arm enables the monitor to be positioned in the optimal viewing position. Once in position, the cart can be powered on using the central power button and connected to additional devices using the I/O panel.

### Putting the Cart into Transport Position

Before moving the cart, ensure it is in transport position.

To put the cart into transport position:


1. Rotate the monitor 90 degrees to the right side of the cart.
2. Lower the monitor on the arm.
3. Hook the monitor stow handle in the stow bracket.

### Moving the Cart


To move the cart:

1. Disconnect the cart from supply mains and secure the power cord to the cord hooks on the rear door. Ensure any network or video cables are disconnected from the external I/O panel on the rear of the cart.
2. Ensure that any external components, such as imaging heads, are disconnected.
3. Release the brakes on all four wheels by lifting the foot lever on each wheel.
4. Push the cart from the rear using the transporting handles (rear-facing handles).
5. Move the cart into position, taking care not to block access to supply mains.
6. When the cart is in its final position, lock all four wheels by pushing the foot levers down.
7. Release the monitor from the stow bracket by pushing the monitor down and away from the cart.
8. Using the monitor positioning handles, position the monitor for optimal viewing.

The monitor arm can be moved to the other side of the cart by releasing the rotation lock. Turn the rotation lock knob on the lateral section of the arm counterclockwise to release the arm.

## Powering the Cart On

To power the Connected OR Cart on:

1. Connect the power cord to supply mains with protective earth.
2. On the front of the cart, in the top-left corner, push the power button.

The power button illuminates, and power is supplied to any consoles or accessories installed on the cart.

## Powering the Cart Off

To power the Connected OR Cart off:

1. On the front of the cart, in the top-left corner, push the power button.

The power button illumination turns off and power to the consoles installed on the cart is disconnected.

2. Disconnect the power cord from supply mains.

## Using the I/O Panel

An I/O panel on the rear of the cart provides the following connections:

- HDMI
- USB
- DVI
- Ethernet (RJ45)


Figure 4 – I/O Panel

## Maintenance and Service

The following regular maintenance may be performed as needed. Component part lists, circuit diagrams, descriptions, and calibration instructions are available upon request.

**Note: Secure the Connected OR Mount in the monitor stow bracket before removing the monitor.**

### Cleaning the Cart

It is recommended that the cart be cleaned after each use.


To clean the Connected OR Cart:

1. Disconnect the cart from supply mains.
2. Wipe the cart with a damp cloth or sponge and any of the following:
  - mild soap solution
  - mild detergent solution
  - 70% ethyl or isopropyl alcohol
  - multipurpose germicidal disinfectant

## Adjusting the Connected OR Mount (Monitor Arm)

**Warning:** Failure to properly adjust the Connected OR Mount could result in equipment damage or injury.

**Caution:** Do not overtighten tensioning screws. Damage to the bushings may occur.

### 1. Adjusting Counterbalance

If the monitor drifts up or down on its own, the counterbalance may need to be adjusted.

To adjust the counterbalance:

1. With the monitor in viewing position, insert a 7/32" hex key at the adjustment point shown in 1 above.
2. Turn the hex key clockwise or counterclockwise until the monitor remains where it is positioned without drifting up or down.

### 2. Adjusting Monitor Swing Tension

To adjust monitor swing tension at the monitor-end of the arm:

1. Insert a 3/32" hex key at either of the points shown in 2 above.
2. Turn clockwise to add tension or counterclockwise to remove tension.

### **3. Adjusting Radial Arm Swing Tension**

To adjust radial arm swing tension:

1. Remove the radial arm cover. Using a 3/32" hex key, remove the two screws securing the top cover.
2. Insert a 5/32" hex key in the point shown in 3 above.
3. Turn the tensioning screw clockwise to add tension or counterclockwise to remove tension.
4. Replace the arm cover.

### **4. Adjusting Monitor Tilt Tension**

If the monitor droops, adjust the tilt tension:

1. Tilt the monitor upward and insert a 7/32" hex key in the adjustment point shown.
2. Turn the tensioning screw clockwise or counterclockwise until the monitor no longer droops when positioned.

### **5. Adjusting Upper Arm Pivot Tension**

To adjust the upper arm pivot tension:

1. Insert a 3/32" hex key at the point shown in 7 above.
2. Turn the tensioning screw clockwise to add tension or counterclockwise to remove tension.

## Resetting the Breaker Switch

To reset the breaker switch:

1. Open the rear door of the cart.


Figure 5 – Location of the breaker switch



2. Locate the breaker switch and reset it.

## Replacing the Power Cord


The Connected OR Cart power cord is replaceable.

To replace the power cord:


1. Power the cart off and disconnect the power cord from supply mains.
2. Open the rear door of the cart.
3. Detach the power cord plate by removing the two screws using a 3 mm hex key.

4. Disconnect the power cord from the cart outlet.
5. Pull the power cord assembly through the opening.

6. Push the new power cord assembly through the plate opening in the cart and connect it to the cart outlet.
7. Attach the power cord plate to the cart using the screws and 3 mm hex key.




8. Close the rear door of the cart.


## Replacing Monitor Cables

To replace a monitor cable (monitor power cable, HDMI cable, DVI cable):

1. Disconnect the cart from supply mains.
2. Disconnect the cable to be replaced from inside the cart and from the monitor.
3. Using a 3/8" nut driver, remove the cable housing from the underside of the top of the cart.
4. Pull the disconnected cable(s) up through the cable port in the top of the cart, removing any rubber grommets as necessary.
5. Remove the Cable Manager from the arm.

6. Using a 3/32" hex key, remove the top cover on the radial section of the arm (as shown).

7. Remove the cable to be replaced from the arm.
8. Install the replacement cable in the arm as shown and reinstall the top cover.
9. Route the replacement cable down through the top cable port and reconnect it inside the cart.
10. Reinstall the cable housing to the underside of the top of the cart.
11. Reconnect the cables to the monitor.

## Disposal

All Connected OR Cart components shall be disposed of in compliance with local, regional and national regulations. Do not dispose any part in general waste.

## Electromagnetic Compatibility

Like other electrical medical equipment, the Connected OR Cart model 240-099-155 and 240-099-156, requires precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC), the Connected OR Cart must be installed and operated according to the EMC information provided in this manual.

This equipment is for use in a professional healthcare environment. Not for use in the RF shielded room of a medical electrical system for magnetic resonance imaging, where the intensity of EM disturbances is high.

This equipment is not likely susceptible to interference from HF surgical instruments in the Special Environment of being in close proximity to an active HF surgical instrument. In the case that HF surgical interference is observed, adjust the separation distance of the equipment.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The Connected OR Cart, model 240-099-155 and 240-099-156, is intended for use in the electromagnetic environment specified below. The customer or the user of the Connected OR Cart should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function; therefore, its RF emissions are very low and are not likely to cause any interference in nearby equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	
		<b>Warning:</b> This system is intended for use by health care professionals only. This system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the system or shielding the location.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The Connected OR Cart, model 240–099–155 and 240-099-156, is intended for use in the electromagnetic environment specified below. The customer or the user of the Connected OR Cart should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment — Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. If disruption of video occurs, the device may require restart.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If disruption of video occurs, the device may require restart.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% $U_t$ 0.5 cycle 0% $U_t$ 1 cycle 70% $U_t$ 25 cycles 0% $U_t$ 5 Sec	0% $U_t$ 0.5 cycle 0% $U_t$ 1 cycle 70% $U_t$ 25 cycles 0% $U_t$ 5 Sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Connected OR Cart, model 240–099–155, requires continued operation during power mains interruptions, it is recommended that the Connected OR Cart be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: $U_t$ is the a.c. mains voltage prior to application of the test level.			
Data extracted from IEC 60601-1-2:2014, Table 4 and Table 5			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The Connected OR Cart, model 240-099-155 and 240-099-156, is intended for use in the electromagnetic environment specified below. The customer or the user for the Connected OR Cart should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Portable and mobile RF communications equipment should be used no closer to any part of the Connected OR Cart, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.			
Conducted RF IEC 61000-4-6	3 Vrms 0.15 MHz to 80 MHz  6 Vrms in ISM bands between 0.15 MHz and 80 MHz	3 Vrms 0.15 MHz - 80 MHz  6 Vrms in ISM bands between 0.15 MHz and 80 MHz	Recommended separation distance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	$d = \sqrt{P}$ 80 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and d is the recommended separation distance in meters (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey (a) should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol:
			
(a) Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Connected OR Cart is used exceeds the applicable RF compliance level above, the Connected OR Cart should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Connected OR Cart. (b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			
Data extracted from IEC 60601-1-2:2014, Table 5			

Tested specifications for immunity to RF wireless communications equipment						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						
Note: Portable RF Communication equipment should be used no closer than 30 cm to the Connected OR Cart. Otherwise, degradation of the performance of this equipment could result.						
Data extracted from IEC 60601-1-2:2014, Table 9						

# China RoHS Declaration

This section provides a hazardous material declaration based on the requirements of the SJ/T11364 standard.


O: Indicates that the hazardous substances contained in all homogeneous materials of this part are below the limit requirements of GB/T26572.

X: Indicates that the hazardous substance contained in at least one homogeneous material used in this component is higher than limitation requirements of GB/T26572. The suggested Environmental Friendly Use Period for all components is equal to the expected lifetime of the Connected OR Cart (5 years).

Component	Pb	Hg	Cd	Cr(VI)	PBBs	PBBDEs	EFUP
Metal Parts	O	O	O	O	O	O	5
Plastic Parts	O	O	O	O	O	O	5
Power Cords	O	O	O	O	O	O	5
Cables/ Connectors	X	O	O	O	O	O	5
Circuit Breaker	X	O	X	O	O	O	5
Relays	O	O	X	O	O	O	5
Switches/ Distributors	O	O	O	O	O	O	5
Mechanical Hardware	O	O	O	O	O	O	5
Fasteners	O	O	O	O	O	O	5
Transformer	O	O	O	O	O	O	5
Packaging Materials	O	O	O	O	O	O	5
Labeling/ Literature	O	O	O	O	O	O	5

# Symbols

## Device and Package Labeling

	Warning/Caution: See Instructions for Use		Legal manufacturer
	Catalogue number		Date of manufacture
	Serial number		Federal law restricts this device to sale by or on the order of a physician
	Consult Instructions for Use		Made in United States
	Mass of mobile medical equipment		This product contains electrical waste or electronic equipment. It must not be disposed of as unsorted municipal waste and must be collected separately.
	Quantity		Circuit breaker
	Ambient temperature range		Relative humidity range
	Atmospheric pressure		Radio frequency hazard
	Fragile; handle with care		This side up
	SGS Q-mark with contract number		Environmentally Friendly Use Period (in years): 5

# Technical Specifications

The following table lists technical specifications for the Connected OR Cart.

Mass	Empty cart: 328 lbs / 149 kg Loaded cart: 570 lbs / 258 kg
Power Consumption	230 V AC, 50/60 Hz, 1700 VA input
Maximum Permitted Electrical Load through Multiple Socket Outlets	7 A
Electrical Safety Classification	Class 1 Equipment No applied part Continuous Operation
Environmental Specifications	Operating: Temperature: +10° C to +30° C Relative humidity: 15% to 93% Atmospheric Pressure: 70 kPa to 106 kPa
	Storage and transport: Temperature: -18° C to +60° C Relative humidity: 15% to 90% Atmospheric Pressure: 50 kPa to 106 kPa
Video and Network Outputs	HDMI USB DVI RJ45

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: The Stryker Logo. All other trademarks are trademarks of their respective owners or holders.