
<p><i>Note: This portion of the document should not appear on produced Labels or IFUs.</i></p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p>Dimensions: 8.25 inch (width) x 5.5 inch Booklet</p>	<p>Color/Material/Finish: Black Graphics on White Background 20# Bond or Equivalent</p>	<p>Label Stock: N/A</p>	
			<p>Description/Type: Instructions For Use</p>	
			<p>Part Number: 0306-105-702</p>	<p>Rev. C</p>

R_x ONLY

CE 0197

STERILE R

10 Gauge **REF** 0306-105-000
11 Gauge **REF** 0306-115-000
13 Gauge **REF** 0306-135-000
9 inch 10 Gauge **REF** 0306-195-000

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)
NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)
SUOMI (FI)
PORTUGUÊS (PT)
NORSK (NO)

Indications For Use

Stryker's Bone Biopsy Kit is used with the ACM® Percutaneous Delivery System and the PCD® Precision System. The kit requires the use of a four-faceted tip introduction needle.

The Bone Biopsy Kit is used to remove a sample of bone tissue for diagnostic purposes.

User/Patient Safety

WARNINGS:

- Only trained and experienced healthcare professionals should use this equipment. Before using any system component, or any component compatible with this system, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with the system components prior to use.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.
- Upon initial receipt and before use, inspect the package for damage and confirm the integrity of the sterile barrier. DO NOT use if package is damaged.
- Upon initial receipt and before use, inspect each component for damage. DO NOT use if any component is damaged.
- Before using an access cannula/stylet assembly, ALWAYS make sure the diameter of the access cannula is compatible with the selected biopsy needle.
- The components of this system should only be manipulated while under fluoroscopic guidance.

- DO NOT reuse, reprocess, or repackage this device. This device is intended for a single use only. This device may not withstand chemical, chemical vapor, or high-temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the device resulting in operational failure. Critical product information may be lost if the device is repackaged. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.
- Follow the current local regulations governing the handling and disposal of sharps and potential biohazard waste.

Kit Descriptions

Each kit contains a coaxial biopsy needle, a syringe, and an obturator.


DESCRIPTION	REF
10 Gauge Bone Biopsy Kit	0306-105-000
11 Gauge Bone Biopsy Kit	0306-115-000
13 Gauge Bone Biopsy Kit	0306-135-000
9 Inch 10 Gauge Bone Biopsy Kit	0306-195-000

Required Component

(1) Introduction Needle (stylet with access cannula)

NOTE: Introduction needles are available in multiple sizes to accommodate various application needs. Introduction needles may be obtained from the Stryker PCD Precision System, the ACM Percutaneous Delivery System or ordered separately. Contact your Stryker Instruments sales representative for a complete list of introduction needles.

Features

A	Coaxial Biopsy Needle
B	Obturator
C	Syringe
D	Plunger

Instructions

To Prepare the Patient

NOTE: If the introduction needle's stylet is removed from its cannula, it must be reinserted properly to ensure that it is secure before use. The silver marker on the cap of the stylet must align with the arrow on the cannula's handle and both parts must lock together to be secure.

- Prepare the patient per standard medical procedures.
- Select an introduction needle. Ensure the size of the coaxial biopsy needle will be compatible with the introduction needle cannula.

To Insert the Access Cannula

1. Hold the introduction needle firmly; stabilize the needle by placing your palm against the cannula's handle and the index finger against its shaft (figure 1).


Figure 1 – Insert the Access Cannula/Styler Assembly

2. Under fluoroscopic guidance, insert the access cannula/styler assembly into the surgical site.
3. Twist the access cannula/styler assembly or use a gentle hammering technique to advance the access cannula/styler assembly into the cortical bone of the vertebral body.

NOTE: Generally, decreased resistance indicates entry into the marrow cavity or fracture site.

4. Just after entering the marrow cavity, immediately stop advancing the introduction needle.
5. Remove the styler by twisting its cap and pulling the styler straight out of the cannula (figure 2).


Figure 2 – Remove the Styler

To Insert the Biopsy Needle

1. Under fluoroscopic guidance, insert the biopsy needle into the access cannula (figure 3).


Figure 3 – Insert the Biopsy Needle

2. Turn the biopsy needle's handle clockwise to lock the biopsy needle in the cannula.
3. Twist the biopsy needle or use a gentle hammering technique to advance the biopsy needle into the marrow cavity. The biopsy needle will capture a bone tissue sample.
4. When the biopsy needle is properly positioned and a sample of the desired size is captured, stop advancing the biopsy needle.

To Remove the Sample

1. Attach the kit syringe to the Luer connector on the biopsy needle (figure 4).


Figure 4 – Attach the Syringe

2. Pull the plunger to create a vacuum and retain the sample in the biopsy needle.
3. Turn the biopsy needle handle to unlock the biopsy needle from the cannula; turn the biopsy needle in the cannula to break off the sample.
4. Twist and remove the biopsy needle from the access cannula.
5. Remove the syringe from the biopsy needle.
6. Use the obturator to remove the sample from the biopsy needle (figure 5).


Figure 5 – Use the Obturator to Remove the Sample

7. Twist and remove the access cannula from the surgical site.
8. Properly dispose of each component used in surgery.

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