stryker®

VertePort® Cement Cannulas and Access Cannulas

REF	0306-080-000	0306-191-000
	0306-100-000	0306-230-000
	0306-101-000	0306-310-000
	0306-110-000	0306-330-000
	0306-111-000	0306-410-000
	0306-114-000	0306-500-000
	0306-130-000	0306-510-000
	0306-131-000	0306-511-000
	0306-140-000	0306-530-000
	0306-190-000	

Instructions For Use

R_x ONLY CE 0197

STERILE R

REF 0306-080-000
REF 0306-114-000
REF 0306-140-000
REF 0306-310-000
REF 0306-330-000
REF 0306-511-000
REF 0306-530-000


STERILE EO

REF 0306-100-000
REF 0306-101-000
REF 0306-110-000
REF 0306-111-000
REF 0306-130-000
REF 0306-131-000
REF 0306-190-000
REF 0306-191-000
REF 0306-230-000
REF 0306-410-000
REF 0306-500-000
REF 0306-510-000

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)
NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)

SUOMI (FI)
PORTUGUÊS (PT)
NORSK (NO)
POLSKI (PL)
ΕΛΛΗΝΙΚΑ (EL)
日本語 (JA)
中文 (ZH)
한국어 (KO)

Introduction

This *Instructions For Use* manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, especially safety information, or in-service training, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Indications For Use

Stryker VertePort Cement Cannulas and Access Cannulas are used to provide and maintain access to the vertebral body for the injection of bone cement/augmentation material.

Contraindications

None known.

For Use With

VertePort Cement Cannulas and Access Cannulas are intended for use with the following Stryker products:

DESCRIPTION	REF
VertaPlex® HV Radiopaque Bone Cement	0406-622-000

VertaPlex Radiopaque Bone Cement	0406-4XX-XXX series
SpinePlex® Radiopaque Bone Cement	0406-2XX-XXX series
AutoPlex® System Disposable Cement Mixer Kit	060X-XXX-XXX series
PCD® Precision System Kit	050X-XXX-XXX series
Cortoss Bone Augmentation Material (Resin)	2101-XXXX series
Cortoss Delivery Gun	2110-0008
Cortoss Mix-Tips, Luer (QTY 3)	2110-0031
Aliquot® Side-Port Syringe, 1cc	2110-0513
Aliquot Flexible Extension, 2 inch, (QTY 3)	2110-0507

Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

SYMBOL	DEFINITION
	General warning sign
	Indexing mark; indicates the direction the handle will steer. Located on top of bevel stylet only.
8G	Indicates 8 gauge cannula size - light grey
10G	Indicates 10 gauge cannula size - dark blue
11G	Indicates 11 gauge cannula size - green
13G	Indicates 13 gauge cannula size - dark grey


User/Patient Safety

WARNINGS:

- Before using any system component, or any component compatible with this system, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with system components prior to use.
- Only trained and experienced healthcare professionals should use this equipment.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.
- Upon initial receipt, inspect each component for damage. DO NOT use any component if damage is apparent.
- Sterile only if package is unopened and undamaged.
- Single use only. DO NOT reuse or resterilize.
- DO NOT reuse, reprocess, or repackage this device. This device is intended for a single use only. This device may not withstand chemical, chemical vapor, or high-temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the device resulting in operational failure. Critical product information may be lost if the device is repackaged. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.
- Use only Stryker-approved system components and accessories.
- DO NOT modify any system component or accessory.

Features

A	Stylet – A solid rod with a sharp point used to penetrate tissue and the vertebral body when installed in the access cannula.
B	Access Cannula – A cannula designed to provide and maintain access to the vertebral body.
C	Cement Obturator – A solid rod with a plastic hub used to manually inject any bone cement/augmentation material contained within the cement cannula.
D	Cement Cannula – A cannula designed to slide within the previously placed access cannula for precise delivery of bone cement/augmentation material under fluoroscopic visualization.

Instructions

To Insert the Access Cannula

WARNINGS:

- Incorrect placement of the access cannula/stylet assembly may result in rupture of the aorta and/or nerve damage.
- The marks on the length of the access cannula are for reference only and are not intended to replace the use of fluoroscopic guidance.
- The marks on the length of the access cannula indicate the distance from the distal end of the access cannula and do not correspond to the distance from the sharp, distal tip of the stylet.
- Follow the current local protocol governing the preparation of the patient for percutaneous bone cement/augmentation material delivery.

NOTE: If desired, incise the surgical site to facilitate insertion of the access cannula/stylet assembly.

1. Remove the protective sheath from the access cannula/stylet assembly.
2. Under fluoroscopic guidance, insert the access cannula/stylet assembly into the surgical site (figure 1).


Figure 1 – Insert the Access Cannula/ Stylet Assembly

3. Twist the access cannula/stylet assembly or use a gentle hammering technique to advance the access cannula/stylet assembly through the cortical bone of the vertebral body.

NOTE: Generally, decreased resistance indicates entry into the cancellous bone of the vertebral body.

4. Unlock and remove the stylet from the access cannula (figure 2).


Figure 2 – Remove the Stylet from the Access Cannula

To Inject Polymethylmethacrylate (PMMA) Bone Cement

WARNINGS:

- Read, understand, and follow the indications for use, instructions, and WARNING information supplied with the bone cement and the mixing/delivery system.
 - ALWAYS prepare the bone cement according to the manufacturer's recommendations.
 - Monitor the amount and location of injected bone cement with extreme care.
1. Connect the mixing/delivery system to the cement cannula.
 2. Fill the cement cannula with bone cement.
 3. Disconnect the mixing/delivery system from the cement cannula.

NOTE: As an alternative, leave the mixing/delivery system connected, and inject the bone cement directly from the system. See the instructions for use supplied with the system.

- Under fluoroscopic guidance, insert the cement cannula into the access cannula. Place the tip of the cement cannula in the desired injection location inside the vertebral body.
- Under fluoroscopic guidance, inject the bone cement using the cement obturator (figure 3).


Figure 3 – Inject Bone Cement Using the Cement Obturator

- When the bone cement has been fully injected, remove the cement cannula from the access cannula.

To Inject Cortoss Bone Augmentation Material

WARNINGS:

- Read, understand, and follow the indications for use, instructions, and WARNING information supplied with the augmentation material and the mixing/delivery system.
- ALWAYS prepare the augmentation material according to the manufacturer's recommendations.
- DO NOT connect a Cortoss mixing/delivery system directly to the access cannula for injection. Direct connection of the Cortoss system does not allow for controlled injection of the bone augmentation material.

- Monitor the amount and location of injected augmentation material with extreme care.

- Connect the Cortoss Mix Tip to the cement cannula (figure 4).


Figure 4 – Connect the Cortoss Mix Tip to the Cement Cannula

- Fill the cement cannula with augmentation material.
- Disconnect the mix tip from the cement cannula.
- Under fluoroscopic guidance, insert the cement cannula into the access cannula. Place the tip of the cement cannula in the desired injection location inside the vertebral body.
- Under fluoroscopic guidance, inject the augmentation material using the cement obturator (figure 5).


Figure 5 – Inject Augmentation Material Using the Cement Obturator


- When the augmentation material has been fully injected, remove the cement cannula from the access cannula.

To Remove the Access Cannula

WARNINGS:

- ALWAYS remove the access cannula from the vertebral body immediately after bone cement/augmentation material injection. If the bone cement/augmentation material is allowed to harden, the access cannula may be difficult to remove and/or a bone cement/augmentation material spike may be created.
 - Follow the current local regulations governing the safe handling and disposal of sharps and potential biohazard waste.
1. When access to the surgical site is no longer required, install the stylet into the access cannula.

2. Twist and then remove the access cannula/stylet assembly from the surgical site before the bone cement/augmentation material hardens.
3. Properly dispose of each component used in surgery.

Specifications

8 Gauge

MODEL	REF
iVAS® 8 Gauge Access Cannula	0306-080-000
8 Gauge VertePort Cement Cannula	0306-310-000

10 Gauge

MODEL	REF
10 G Match-Ground Four Facet Tip Introduction Needle 5 inch (12.7 cm)	0306-100-000
10 G Match-Ground Bevel Tip Introduction Needle 5 inch (12.7 cm)	0306-101-000
10 G Match-Ground Four Facet Tip Introduction Needle 9 inch (22.86 cm)	0306-190-000
10 G Match-Ground Bevel Tip Introduction Needle 9 inch (22.86 cm)	0306-191-000
VertePort 10 G System	0306-400-000
VertePort Cement Cannula for Kit 10 G	0306-410-000
iVAS 10 Gauge Access Cannula	0306-530-000

11 Gauge

MODEL	REF
11 G Match-Ground Four Facet Tip Introduction Needle 5 inch (12.7 cm)	0306-110-000
11 G Match-Ground Bevel Tip Introduction Needle 5 inch (12.7 cm)	0306-111-000
11 G Match-Ground Four Facet Tip Introduction Needle 4 inch (10.16 cm) with Bevel Stylet	0306-114-000
iVAS 11 Gauge Access Cannula	0306-330-000
VertePort Cement Cannula for Kit 11 G	0306-510-000
11 Gauge VertePort Cement Cannula	0306-511-000

13 Gauge

MODEL	REF
13 G Match-Ground Four Facet Tip Introduction Needle 5 inch (12.7 cm)	0306-130-000
13 G Match-Ground Bevel Tip Introduction Needle 5 inch (12.7 cm)	0306-131-000
iVAS 13 Gauge Access Cannula	0306-230-000

14 Gauge

MODEL	REF
14 G Match-Ground Four Facet Tip Introduction Needle 4.4 inch (11.18 cm) with Bevel Stylet	0306-140-000

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