## Surgical Helmets

Flyte<sup>®</sup> Helmet

REF 0408-600-000

Flyte Helmet with Fiber-optic Headlight


REF 0408-635-000

Flyte Helmet with Battery-powered Headlight

REF 0408-645-000

Instructions For Use

**R<sub>x</sub> ONLY**

ENGLISH (EN)  
ESPAÑOL (ES)  
DEUTSCH (DE)  
FRANÇAIS (FR)  
ITALIANO (IT)  
NEDERLANDS (NL)  
SVENSKA (SV)  
DANSK (DA)  
SUOMI (FI)  
PORTUGUÊS (PT)  
NORSK (NO)  
POLSKI (PL)  
ΕΛΛΗΝΙΚΑ (EL)

## Introduction

This *Instructions For Use* manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians.

Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

If additional information, especially safety information, or in-service training is required, contact your Stryker sales representative or call Stryker customer service.

**NOTE:** The user and/or patient should report any serious product-related incident to both the manufacturer and the Competent Authority of the European Member State where the user and/or patient is established.

Trademarks not the property of Stryker Corporation are the property of their respective owners.

## Compliance Statements

The CE marking affixed to the equipment indicates compliance with the following European Community Regulations:

- Medical Device Regulation (EU) 2017/745
- Regulation (EU) 2016/425 - Personal Protective Equipment (PPE)

This PPE has been EC-type examined by the following notified body:

BSI Group the Netherlands B.B.  
John M. Keynesplein 9,  
1066 EP Amsterdam, The Netherlands, BSI 2797

The lenses were developed and tested in accordance with the requirements of BS EN 166: 2002.

## Marking Definitions — S1663F

CODE	DEFINITION
S	Code of Manufacturer, Stryker
166	European Norm Standard - Personal Eye Protection Specification
3	Liquid Droplets
F	Mechanical Strength, Low Energy Impact

## Indications For Use

The Stryker Flyte Helmets are surgical devices that pull in and circulate fresh air throughout the helmet to aid user comfort.

## Contraindications

None known.

## Limitations of Use

### WARNINGS:

- The surgical helmet is not intended to be used as a respiratory protection device.
- The materials contained within this device are not known allergens. However, if an allergic reaction occurs, contact a physician.

## For Use With

Rechargeable Flyte Power Pack and the Flyte Power Cord. See *Accessories* section.

## Description

In addition to pulling in and circulating fresh air throughout the surgical helmet, each helmet has microphone and speaker capability. Stryker-approved dictation recorders, MP3 players or other compatible communication devices may be connected to the helmet via a communication cord.

## User/Patient Safety

### WARNINGS:

- Before using any component, or any component compatible with this equipment, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the components prior to use.
- Only trained and experienced healthcare professionals should use this equipment.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage. DO NOT use any equipment if damage is apparent or the inspection criteria are not met. See the *Inspection and Testing* section for inspection criteria.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.
- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment like the Flyte helmet. Install and place the helmet into service according to the EMC information in this manual. Portable and mobile radio frequency (RF) equipment, such as wireless phones, can affect the function of the helmet. See the *Specifications* section.
- DO NOT stack or place equipment adjacent to the product. If such a configuration is necessary, observe the configuration to ensure that electromagnetic interference does not degrade performance.
- DO NOT use the product in a magnetic resonance imaging (MRI) environment. Using the product in an MRI environment could affect the function of the system.

## Accessories

This section describes components that may be ordered to replace original equipment that is damaged, worn, or must be replaced. This section may also contain optional components used with the surgical helmet.

### WARNINGS:


- Use only Stryker-approved components and accessories, including earphones and dictation equipment, unless otherwise specified.
- Using other electronic components and accessories may result in increased electromagnetic emissions or decreased electromagnetic immunity of the equipment.
- DO NOT modify any component or accessory, unless otherwise specified.
- DO NOT reuse, reprocess, or repackage a single use device. A single use device is intended for a single use only. A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the single use device resulting in operational failure. Critical product information may be lost if the single use device is re-packaged. Failure to comply may lead to infection or cross-infection and result in patient and/or health care staff injury.

The following Stryker-approved accessories are sold separately:

DESCRIPTION	REF
Flyte Lithium-Ion Extended Life Power Pack	0408-660-000
Flyte Comfort Pads – A set of five disposable pads are installed in the inner helmet to provide additional comfort; eight sets per case	0408-210-000
Flyte Helmet Power Cord	0408-600-300
Flyte Helmet Dictation Power Cord – When using a Stryker-approved device such as a dictation recorder or MP3 player, replace the power cord with a dictation power cord.	0408-605-300
SONY Earphone-Right	0400-615-600
SONY Earphone-Left	0400-615-610
PANASONIC Earphone-Right	0400-615-620
PANASONIC Earphone-Left	0400-615-630

**NOTE:** For a complete list of accessories, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

# Features

<b>A</b>	<b>Front Flexible Air Duct</b> – This vent channels air to the front of the helmet.
<b>B</b>	<b>Height Adjustment Strap</b> – This strap is adjusted to provide a comfortable fit. Pull to extend or push to retract the length of the strap for fit and comfort.
<b>C</b>	<b>Headlight Adjustment Handle</b> – Available on the fiber-optic and battery-powered headlight helmet models only. The handle is used to direct the headlight beam.
<b>D</b>	<b>Hook</b> – This hook aids in positioning the hood or toga lens on the helmet.
<b>E</b>	<b>Headlight</b> – Available on the fiber-optic and battery-powered headlight helmet models only.
<b>F</b>	<b>Light Shroud</b> – Available on the fiber-optic and battery-powered headlight helmet models only. The shroud prevents glare from the shield and may be trimmed or removed as desired.
<b>G</b>	<b>Microphone Indicator Light</b> – Illuminates green while the microphone is on.
<b>H</b>	<b>Power Loss Indicator Light</b> – Flashes on and off red when the power pack has a reduced amount of power remaining.
<b>I</b>	<b>Variable Control Buttons</b> – Press these buttons to select a fan speed level. There are six speed levels available. For helmets with a battery-powered headlight, press both buttons simultaneously until the light turns on or off, as desired.
<b>J</b>	<b>Microphone</b> – A built-in microphone provides voice transmission to compatible devices such as a dictation recorder.
<b>K</b>	<b>Microphone Button</b> – Press this button to toggle the microphone on and off (mute).
<b>L</b>	<b>Attachment Points (two magnets)</b> – These magnets secure the disposable hood or toga lens to the helmet.
<b>M</b>	<b>Comfort Pads</b> – These removable pads provide a comfortable fit.

<b>N</b>	<b>Adjustable Air Nozzle Tip</b> – This vent channels air to the user’s back. Rotate the vent to direct flow toward or away from the neck for comfort.
<b>O</b>	<b>Headband Adjustment Knob</b> – Customize the fit by turning the knob clockwise to tighten and counterclockwise to loosen the headband.
<b>P</b>	<b>Fan Intake</b> – The fan pulls in air and circulates the air throughout the helmet.
<b>Q</b>	<b>Fiber-optic Cable</b> – Available on the fiber-optic headlight helmet model only. This cable routes light from the fiber-optic light source to the headlight.
<b>R</b>	<b>Audio Jack</b> – When using a dictation device, use this jack to connect an earphone.
<b>S</b>	<b>Power Cable</b> – This built-in cable is compatible with the power cord and holster combination, and the dictation power cord.
<b>T</b>	<b>Power Cord with Holster</b> – This cord is supplied with all helmet models. The holster holds standard and extended life power packs. The holster swivels to hang from your waist band for either left or right handed access.
<b>U</b>	<b>Power Cable Plug</b> – This plug is used to connect to the appropriate power cord. See the <i>Accessories</i> section.
<b>V</b>	<b>ACMI Fiber-optic Cable Plug with Cap</b> – Available on the fiber-optic headlight helmet model only. This plug is used to connect to an ACMI compatible fiber-optic light source.

## Symbol Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

Increase Air Flow	Decrease Air Flow	Microphone Mute	General Warning Sign	Refer to Instruction Manual/ booklet	Direct Current (DC)

## Instructions

### WARNINGS:

- Perform recommended inspection and testing as indicated in the *Inspection and Testing* section.
- ALWAYS operate the equipment within the specified environmental condition values. See the *Specifications* section.
- This equipment is suitable to use in a professional healthcare facility environment.
- ALWAYS wear the surgical helmet with a surgical toga or surgical hood to obtain the required personal protection.
- DO NOT cut or puncture the fiber-optic cable.
- ALWAYS turn the light source power off before connecting or disconnecting the fiber-optic cable to or from the source.
- ALWAYS hold the strain relief near the end of the fiber-optic cable when connecting or disconnecting the cable from the light source.

### CAUTIONS:

- DO NOT scratch the headlight lens.
- DO NOT stretch or kink the fiber-optic cable. Failure to comply will damage the cable and reduce the light output permanently.
- If using a dictation recorder or MP3 player, ALWAYS make sure the device has standard 3.5 mm port(s). Failure to comply may result in helmet or accessory damage.

**NOTE:** Helmet position will vary from person to person depending on head size and personal comfort preferences.

### To Don a Surgical Helmet

**WARNING:** DO NOT connect the power cable plug to the power pack before donning the helmet.

**NOTE:** ALWAYS perform the following steps in sequence. Failure to comply may cause the red power loss indicator to flash. See the *Troubleshooting* section.

- Apply the comfort pads inside the helmet as required (figure 1).
- Connect the helmet power cable to the appropriate power cord (figures 2 and 2a). See the *Accessories* section.

### NOTES:

- To disconnect the power cord, pull the connector sleeve back before separating the connector plug and receptacle (figure 2).
  - If using a dictation power cord, connect the pink-band plug into the microphone receptacle and the green-band plug into the speaker receptacle (figure 2a). Only one connector plug may be used based on the recorder or player receptacles available.
  - If using a dictation recorder or MP3 player, connect an earphone to the helmet audio jack (figure 2b as required). See the *Accessories* section.
- Loosen the headband adjustment knob (figure 3).

**CAUTION:** DO NOT pull to extend or push to retract the front flexible air duct. Failure to comply may result in product damage.

- Put the helmet on your head while extending the height adjustment strap for a comfortable fit (figure 4).
- Turn the headband knob to secure the helmet (figure 5).

### NOTES:

- Make sure the helmet fits securely on top of your head.
- Make sure the face mask ties are accessible to allow easy removal of the face mask when donning the hood or toga.

**WARNING:** Keep your fingers away from the fan during operation.

- Slide a fully charged power pack into the holster by aligning the grooves on the power pack with the holster. The power pack will snap into place (figure 6).

### NOTES:

- Make sure the indicator lights momentarily flash and the fan turns on automatically. See the *Troubleshooting* section if required.
  - If using a helmet with a battery-powered headlight, install an extended life power pack into the holster.
- Hang the holster from your waistband and swivel the holster to the desired orientation (figure 7).
  - Press the variable control buttons to adjust the fan speed for optimal air flow. See the *Troubleshooting* section if required (figure 8).
  - Rotate the adjustable air nozzle tip for optimal air flow (figure 9).

### NOTES:

- If using a dictation recorder or MP3 player, see the instructions for use supplied with the device for volume adjustment, then place the earphone in your ear (figure 2b).
- You are now ready to scrub for surgery and don a hood or toga. See the instructions for use supplied with the hood or toga.
- If the face shield lens begins to fog during use, see the *Troubleshooting* section.

### To Use a Surgical Helmet with a Headlight

#### To Use a Fiber-optic Headlight

##### CAUTIONS:

- ALWAYS remove the protective cap from the plug of the ACMI fiber-optic cable before installing the plug into the fiber-optic light source receptacle.
- ALWAYS put the protective cap onto the plug of the ACMI fiber-optic cable after use and before storage.

##### NOTES:

- If using a fiber-optic headlight, turn the light beam on and off and adjust the brightness at the light source. See the instructions for use supplied with the light source.
  - A clamp may be used to fasten the fiber-optic cable to the gown to support the weight of the cable.
- Remove the protective cap on the male plug of the ACMI fiber-optic cable and install the plug into any compatible ACMI fiber-optic light source receptacle (figure A).
  - Adjust the headlight beam direction with the headlight adjustment handle (figure C).

#### To Use a Battery-powered Headlight

**NOTE:** Using the battery-powered headlight will decrease the power pack capacity by half. See the instructions for use supplied with the power pack for capacity values.

- Make sure a power pack is installed.
- Press both variable air flow control buttons simultaneously until the light turns on or off, as desired (figure B).
- Adjust the headlight beam direction with the headlight adjustment handle (figure C).

### To Doff a Surgical Helmet

- Flex the tab back to release the power pack, then lift the power pack from the holster (figure D).
- Loosen the headband knob and remove the helmet from your head.
- Remove and discard the comfort pads from the helmet. Clean the helmet and cord. See the *Cleaning* section. To clean the power pack, see the instructions for use supplied with the Flyte charger.

## Inspection and Testing

### WARNINGS:

- Only individuals trained and experienced in the maintenance of reusable medical devices should inspect and test this equipment.
- DO NOT use any equipment if damage is apparent or the inspection criteria are not met.
- DO NOT disassemble or service this equipment, unless otherwise specified.

**NOTE:** For service, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.


INTERVAL	INSPECTION CRITERIA	ACTION
Before Use	Inspect equipment for damage, wear, or missing components.	If damage is apparent, replace the equipment.
	Inspect power cord and fiber-optic cable for cuts or punctures.	
	Inspect headlight lens for scratches.	
	Inspect the power pack housing for cracks.	
	Check fan motor for excessive noise.	

**NOTE:** If any component must be discarded, see the *Disposal/Recycle* section.


# Instructions

## Figures


1

3


2

4


2a


5

2b


6


# Instructions

## Figures (continued)

7


A

8

B


9

C

D

## Cleaning

### WARNINGS:

- ALWAYS remove the power pack, any type power cord, earphone, and comfort pads from the helmet before cleaning.
- DO NOT sterilize any component.

### CAUTIONS:

- DO NOT use solvents, lubricants, or other chemicals, unless otherwise specified.
- DO NOT allow liquids or moisture into the fan or electrical connections of the helmet.
- DO NOT allow debris to enter the helmet fan.

### NOTES:

- The cord may be the power cord or the dictation power cord.
- To clean the power pack, see the instructions for use supplied with the Flyte charger.

### To Clean the Surgical Helmet and Cord

1. Remove all the accessories from the helmet. Discard the comfort pads.
2. Wipe the helmet, cord, and power pack holster with a soft cloth dampened with a nonabrasive, hospital disinfectant.
3. Dry the equipment with a lint-free towel or medical-grade compressed air.
4. Clean the headlight lens with lens tissue. See the instructions for use supplied with the lens tissue.
5. See the *Storage and Handling* section.

## Storage and Handling

**CAUTION:** ALWAYS store the equipment within the specified environmental condition values throughout its useful life. See the *Specifications* section.

To ensure the longevity, performance and safety of this equipment, use of the original packaging material is recommended when storing or transporting this equipment.

## Disposal/Recycle

**WARNING:** ALWAYS follow the current local recommendations and/or regulations governing environmental protection and risks associated with recycling or disposing of the equipment at the end of its useful life.

To comply with European Community Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU, this device should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.

## Troubleshooting

**WARNING:** DO NOT service this equipment, unless otherwise specified. If you require service, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

PROBLEM	CAUSE	ACTION
The face shield lens is fogging.	Not enough airflow is circulating throughout the helmet.	Press the variable control buttons to increase the fan speed if possible.
	The power pack requires replacement.	Replace the power pack with a fully charged power pack.
	The helmet is damaged.	Replace the helmet.
The indicator lights fail to flash and/or the fan fails to turn on during start up.	The power pack requires replacement.	Replace the power pack with a fully charged power pack.
	The power cord with holster is damaged.	Replace the power cord with holster.
	The helmet is damaged.	Replace the helmet.
The helmet will not allow airflow adjustment.	The helmet is damaged.	Replace the helmet.
The fan fails to provide enough or any airflow.	The fan speed is not adjusted correctly.	Press the variable control buttons to increase the fan speed.
	The helmet is damaged.	Replace the helmet.
The green microphone indicator flashes rapidly for several seconds.	The power pack may be dislodged or the fan may be experiencing physical interference.	Make sure the power pack is secure in the holster. Remove any fan obstruction.
The red power loss indicator light flashes on and off.	The steps within the <i>To Don a Surgical Helmet</i> section were performed out of sequence.	See the <i>To Don a Surgical Helmet</i> section and perform the steps in sequence.
	The power pack requires replacement.	Replace the power pack with a fully charged power pack to prevent the loss of air circulation in the helmet.
Sporadic electrical interference is experienced.	Electrical noise is present.	Turn off all the electrical equipment not in use in the room.
		Relocate the electrical equipment to maximize the distance between the equipment. Increase spatial distance.
		Plug equipment into different outlets.

**NOTE:** If any component must be discarded, see the *Disposal/Recycle* section.

## Specifications

<b>Model:</b>	REF 0408-600-000	REF 0408-635-000	REF 0408-645-000
<b>Description:</b>	Flyte Helmet	Flyte Helmet with Fiber-optic Headlight	Flyte Helmet with Battery-powered Headlight
<b>European Conformity:</b>		Not Applicable	
<b>Mass:</b>	0.47 kg [16.6 oz]	0.68 kg [24.0 oz]	0.52 kg [18.3 oz]
<b>Fiber-optic Cable Length:</b>	Not Applicable	3.05 m [10 ft]	Not Applicable
<b>Fan Speed:</b>	2,250 - 3,500 rpm		
<b>Input Voltage:</b>	6 - 8.4 V $\equiv$		
<b>Ingress Protection (IP):</b>	IPX0		
<b>Mode of Operation:</b>	Continuous		

**Product Safety Certification:**

CSA International

### International Electrotechnical Commission (IEC)

IEC 60601-1:2012, Ed. 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); Amendment 1 (2012)

IEC 60601-1-6:2010+A1:2013, Ed. 3.1, *Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance - Usability*

IEC 62366:2007+A1:2014, Ed. 1.1, *Medical Devices - Part 1: Application of Usability Engineering to Medical Devices*

### Canadian Standards Association (CSA)

CAN/CSA-C22.2 No. 60601-1:2014, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; (IEC 60601-1:2005+A1:2012, MOD)

### American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)

ANSI/AAMI ES60601-1:2005/(R) 2012, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; Consolidated Reprint (2009)/(R) 2012; Amendment 2 (2010)/(R) 2012; Amendment 1 (2012)

### European Committee for Electrotechnical Standardization (CENELEC)

EN 60601-1:2006+A12:2014, Ed. 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment A11 (2011); IEC Amendment 1 (2013); IEC Corrigendum 3 (2014); CENELEC Amendment A12 (2014)

**Environmental Conditions:**

#### Operation

#### Storage and Transportation



**Temperature Limitation**

**Humidity Limitation**

**Atmospheric Pressure Limitation**

## Specifications

Guidance and manufacturer's declaration - electromagnetic emissions		
Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 are intended for use in the electromagnetic environment specified below. The customer or the user of Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  <b>NOTE:</b> The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer adequate protection to radio frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Not Applicable	Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 are suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:   <b>WARNING:</b> This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. Mitigation measures may be necessary, such as reorienting or relocating the Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 or shielding the location.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration - electromagnetic immunity			
Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 are intended for use in the electromagnetic environment specified below. The customer or the user of Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V 0.15 MHz – 80 MHz  6 V in ISM bands between 0.15 MHz and 80 MHz  80% AM at 1 kHz	3 V 0.15 MHz – 80 MHz  6 V in ISM bands between 0.15 MHz and 80 MHz  80% AM at 1 kHz	IEC 60601-1-2 3rd Edition:  Portable and mobile RF equipment should be used no closer to any part of the Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz  Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m).  IEC 60601-1-2 4th Edition:   <b>WARNING:</b> Portable RF equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Flyte Helmet REF 0408-600-000, 0408-635-000, or 0408-645-000, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> .  Interference may occur in the vicinity of equipment marked with the following symbol:    (Non-ionizing electromagnetic radiation)
Radiated RF IEC 61000-4-3	<sup>b</sup> 3 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz  27 V/m 385 MHz, pulse modulation 18 Hz, Maximum power = 1.8 W  28 V/m 450 MHz, FM ± 5 kHz deviation, 1 kHz sine, Maximum power = 2 W  9 V/m 710, 745, 780, 5240, 5500, 5785 MHz, pulse modulation 217 Hz, Maximum power = 0.2 W  28 V/m 810, 870, 930 MHz, pulse modulation 18 Hz, Maximum power = 2 W  28 V/m 1720, 1845, 1970, 2450 MHz, pulse modulation 217 Hz, Maximum power = 2 W	<sup>b</sup> 3 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz  27 V/m 385 MHz, pulse modulation 18 Hz, Maximum power = 1.8 W  28 V/m 450 MHz, FM ± 5 kHz deviation, 1 kHz sine, Maximum power = 2 W  9 V/m 710, 745, 780, 5240, 5500, 5785 MHz, pulse modulation 217 Hz, Maximum power = 0.2 W  28 V/m 810, 870, 930 MHz, pulse modulation 18 Hz, Maximum power = 2 W  28 V/m 1720, 1845, 1970, 2450 MHz, pulse modulation 217 Hz, Maximum power = 2 W	

NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Flyte Helmet REF 0408-600-000, 0408-635-000, and 0408-645-000 are used exceeds the applicable RF compliance level above, the Flyte Helmet REF 0408-600-000, 0408-635-000, and 0408-645-000 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Flyte Helmet REF 0408-600-000, 0408-635-000, and 0408-645-000.

<sup>b</sup>Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## Specifications (continued)

Guidance and manufacturer's declaration - electromagnetic immunity			
Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 are intended for use in the electromagnetic environment specified below. The customer or the user of Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 should assure that it is used in such an environment. <b>NOTE:</b> The values provided in the table below have changed due to 60601-1-2 4th Edition requirements.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV at 100 kHz repetition frequency for power supply lines ± 1 kV at 100 kHz repetition frequency for input/output lines	Not Applicable	Not Applicable
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Not Applicable	Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 0% $U_T$ (100% dip in $U_T$ ) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% $U_T$ (100% dip in $U_T$ ) for 1 cycle at 0° 70% $U_T$ (30% dip in $U_T$ ) for 25 and 30 cycles at 0° <5% $U_T$ (>95% dip in $U_T$ ) for 5 s 0% $U_T$ (100% dip in $U_T$ ) for 250/300 cycles	Not Applicable	Not Applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m, 30 A/m at 50 and 60 Hz	3 A/m, 30 A/m at 50 and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE:  $U_T$  is the a.c. mains voltage prior to application of the test level.

IEC 60601-1-2 3rd Edition: Recommended separation distances between portable and mobile RF equipment and the Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000			
Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 are intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF equipment (transmitters) and Flyte Helmet REF 0408-600-000, 0408-635-000 and 0408-645-000 as recommended below, according to the maximum output power of the equipment.			
Rated maximum output power of transmitter  W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.74
1	1.2	1.2	2.3
10	3.8	3.8	7.4
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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### International Electrotechnical Commission (IEC)

IEC 60601-1:2012, Ed: 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); Amendment 1 (2012)  
 IEC 60601-1-6:2010+A1:2013, Ed: 3.1, *Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance - Usability*  
 IEC 62366:2007+A1:2014, Ed: 1.1, *Medical Devices - Part 1: Application of Usability Engineering to Medical Devices*

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
EN 60601-1:2006+A12:2014, Ed: 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment 2 (2010)/(R) 2012; Amendment 1 (2012); IEC Amendment 1 (2013); IEC Corrigendum 3 (2014); CENELEC Amendment A12 (2014)

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
### International Electrotechnical Commission (IEC)

IEC 60601-1:2012, Ed: 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); Amendment 1 (2012)

IEC 60601-1-6:2010+A1:2013, Ed: 3.1, *Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance - Usability*

IEC 62366:2007+A1:2014, Ed: 1.1, *Medical Devices - Part 1: Application of Usability Engineering to Medical Devices*

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



EN 60601-1:2006+A12:2014, Ed: 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment A11 (2011); IEC Amendment 1 (2013); IEC Corrigendum 3 (2014); CENELEC Amendment A12 (2014)

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