<p><i>Note: This portion of the document should not appear on produced Labels or IFUs.</i></p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p>Dimensions: 8.5 inch (width) x 11 inch Booklet</p>	<p>Color/Material/Finish: Black Graphics on White Background 20# Bond or Equivalent</p>	<p>Label Stock: N/A</p>	
			<p>Description/Type: Instructions For Use</p>	
			<p>Part Number: 0506-583-700</p>	<p>Rev. M</p>

PCD Precision System with 10 G Needle **REF** 0506-482-000

PCD Precision System with 11 G Needle **REF** 0506-483-000

PCD Precision System with 13 G Needle **REF** 0506-485-000

PCD Precision System without Needle **REF** 0506-486-000

PCD Precision System with 90 Degree Extension Tube and VertaPlex® HV Bone Cement **REF** 0507-586-000

Instructions For Use

R_x ONLY  0197 

STERILE R

Not made with natural rubber latex

Introduction

This instructions for use manual contains information intended to ensure the safe and effective use of your product. This manual is intended for in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Definitions

SYMBOL	DEFINITION
	General warning sign
	Unlock direction
PUSH	Push
ROTATE	Rotate
	Stop
	Transfer Inject

Indications for Use

Stryker's PCD System kits provide components for mixing bone cement and delivering it percutaneously.

Contraindications

None known

User/Patient Safety

WARNINGS:

- Only trained and experienced healthcare professionals should use this equipment. Before using any system component or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with system components prior to use.
- Upon initial receipt and before use, inspect the package for damage and confirm the integrity of the sterile barrier. DO NOT use any equipment if damage is apparent or the sterile barrier has been compromised.
- DO NOT use the equipment after the expiration date printed on the package. The equipment may not be safe or effective after the expiration date.
- DO NOT reuse, reprocess, or repackage a device that is intended for single use only. A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a contamination risk and compromise structural integrity resulting in operational failure. Critical product information may be lost during repackaging. Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.
- Use only Stryker-approved components and accessories, unless otherwise specified. DO NOT modify any component or accessory.
- Monitor patients carefully for any change in blood pressure during and immediately following the application of bone cement.

- Treating multiple vertebral levels may increase the risk of sudden drop in blood pressure, particularly if more than three vertebral levels are treated in a single operation.
- Consider other conventional therapies before using a percutaneous vertebroplasty or kyphoplasty procedure.
- Assess the risk/benefit for patients with fractures of the vertebral body or walls of the pedicles due to disease or trauma. Have immediate surgical support available for urgent decompressive surgery, if required.
- Follow the current local regulations governing biohazard waste to safely handle and dispose of sharps.
- ALWAYS follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.
- Stryker bone cement is supplied with all of the necessary components for safe and effective use. DO NOT add a radiopacifier (e.g. barium), antibiotics, or other drugs or materials to the bone cement. The addition of such substances may alter the bone cement properties and result in patient injury.

Instructions

Prepare and Organize Equipment

Vacuum Pump

NOTES:

- You may use either Stryker's High Vacuum Pump REF 0206-000-000 or Automatic High Vacuum Pump REF 0206-500-000.
- If using the High Vacuum Pump, it can be sterilized and placed within the sterile field.
- If using the Automatic High Vacuum Pump, connect the nitrogen source to the vacuum pump and set the pump outside the sterile field.

Bone Cement

NOTES:

- Cement mixing and working times are approximate and may vary due to room temperature and user technique.

ENGLISH (EN)	SUOMI (FI)
ESPAÑOL (ES)	PORTUGUÊS (PT)
DEUTSCH (DE)	NORSK (NO)
FRANÇAIS (FR)	POLSKI (PL)
ITALIANO (IT)	ΕΛΛΗΝΙΚΑ (EL)
NEDERLANDS (NL)	日本語 (JA)
SVENSKA (SV)	中文 (ZH)
DANSK (DA)	한국어 (KO)


Prepare and Organize Equipment (continued)

- To extend working time of cement, you may remove cement monomer from the outer packaging and refrigerate to 35 °F (1.7 °C) before use. DO NOT refrigerate the powder. Maintain an ambient temperature in the operating room between 60-68 °F (15.6-20 °C).
- Mixing capacity is 10 ml of liquid monomer and 20 g of powder.

Mixing and Delivery Components

WARNING: Reduce nitrogen source pressure regulator to zero before connecting the vacuum tube to the Automatic High Vacuum Pump.

- Use proper sterile technique to remove the mixing and delivery components from the packages.
- Attach the vacuum tube from package to vacuum pump.

- Attach the extension tube to the delivery cartridge in the orientation shown. Make sure the connection is secure.

Mixing the Cement


NOTE: Cement mixing and working times are approximate and may vary due to room temperature and user technique.

- Set the mixing unit on the table.
- Place the funnel on the mixing unit.
- Pour the bone cement powder into the mixing unit and then add the liquid monomer.


- Remove the funnel.
- Screw the mixing blade assembly on the mixing unit, turning it clockwise to secure.


NOTE: To ensure proper mixing, tighten the mixing blade assembly until it can no longer be turned.

- Attach the vacuum hose to the vent.


- Activate the vacuum source until the gauge reads 20-22 inHg.


NOTE: For optimal mixing, hold the unit at an angle during the next step.

- Hold the mixing unit in one hand while simultaneously pumping and twisting the handle. To ensure thorough blending, continually twist while pumping the entire stroke length. Straight up and down movements will only pack down the cement. Observe the mixture as blending occurs. Blending takes approximately 30 seconds.


- Detach the vacuum tube from the unit.
- Press and rotate the blade release (PUSH/ROTATE) button. Withdraw the plunger shaft from the mixing unit.


NOTE: The mixing blade will remain in the mixing unit.

- Turn the delivery cartridge with extension tube counterclockwise onto the mixing unit.

- Transfer the cement from the mixing unit to the delivery cartridge by twisting the outer body of the unit clockwise (Transfer/Inject direction).


- As you transfer cement from the mixing unit, observe that cement moves down the nozzle and into the tube. At a certain point, you will hear a click. Keep twisting the outer body of the mixing unit until the entire length of tube is primed with cement.

NOTE: If twisting becomes difficult, slow down. This will ease cement transfer.

NOTE: Avoid an air gap in the end of the tubing prior to needle attachment.

- Attach the extension tube to the introduction needle with the Luer connector. Make sure the fitting is secure.


Cement Injection

WARNING: The physician must monitor the amount and placement of injected cement with extreme care. During application, variations in incremental and total cement delivery may occur because cement is a compressible material. A buildup of pressure can cause cement to suddenly flow and overflow the fracture site which may result in injury to the patient.

NOTE: Once cement has reached the needle, note the starting location of the obturator piston when you begin to deliver cement as you can estimate the volume of delivered cement by using the 1 cc markers on the cement cartridge.

- Continue to turn the outer body of the mixing unit with care while applying cement to the fracture.

Disassembly

- When cement application is complete, detach the Luer connector from the introduction needle.

WARNING: The needle must be immediately removed. A needle left in the surgical site after cement hardens could become very difficult to remove and cause injury to the patient.

- Twist the needle to remove it from the surgical site.
- Discard each component from the kit.

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