
Consolidated Operating Room Equipment  
Sterilization Cases

REF 5400-276-000  
5400-277-000  
5400-278-000

Instructions For Use

R<sub>x</sub> ONLY

ENGLISH (EN)  
ESPAÑOL (ES)  
DEUTSCH (DE)  
FRANÇAIS (FR)  
ITALIANO (IT)  
NEDERLANDS (NL)  
SVENSKA (SV)  
DANSK (DA)  
SUOMI (FI)  
PORTUGUÊS (PT)  
NORSK (NO)  
POLSKI (PL)  
ΕΛΛΗΝΙΚΑ (EL)  
日本語 (JA)  
中文 (ZH)  
한국어 (KO)

## Introduction

This instructions for use manual contains information intended to ensure the safe, effective and compliant use of your product. Keep and consult this reference manual during the life of the product.

This manual is intended for in-service trainers, biomedical equipment technicians, and central supply/sterile processing technicians.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

## Indications For Use

The Stryker Sterilization Case is designed to contain an insert tray and protect reusable medical devices during steam sterilization, storage and transportation. The Sterilization Case must be wrapped in a validated sterilization wrap prior to sterilization to maintain sterility. The Stryker Insert Tray is designed to provide a convenient holder for reusable medical devices during automated washing and steam sterilization. This Insert Tray may also be placed in a rigid container that does not require a validated wrap.

## Contraindications

None known.

## Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.


| SYMBOL  | DEFINITION           |
|---|----------------------|
|  | General warning sign |

## For Use With

The trays are intended for use with CORE handpieces, attachments, and accessories. Each tray will accommodate a variety of handpiece, attachment, and accessory combinations.


## Features

### Sterilization Case, Small (REF 5400-276-000)


|   |                |
|---|----------------|
| A | Lid            |
| B | Insert Tray    |
| C | Mat            |
| D | Base           |
| E | Latches (two)  |
| F | Handles (four) |

**Sterilization Case, Medium (REF 5400-277-000)**

|          |                |
|----------|----------------|
| <b>A</b> | Lid            |
| <b>B</b> | Mats (two)     |
| <b>C</b> | Insert Tray    |
| <b>D</b> | Base           |
| <b>E</b> | Latches (two)  |
| <b>F</b> | Handles (four) |


**Sterilization Case, Large (REF 5400-278-000)**


|          |                      |
|----------|----------------------|
| <b>A</b> | Lid                  |
| <b>B</b> | Insert Trays (three) |
| <b>C</b> | Base                 |
| <b>D</b> | Handles (two)        |
| <b>E</b> | Latches (four)       |

## To Remove the Lid from the Base

1. Push down on one of the latch release buttons.


2. Slide the latch away from the center of the case to unlock the latch.

3. Repeat steps to unlock each of the remaining latches.
4. Remove the lid from the base.

## To Secure the Lid to the Base

1. Place the mats inside the insert tray(s), if applicable, and stack the insert tray(s) inside of the base.
2. Put the lid on the base.
3. Slide each latch toward the center of the case until it securely locks into place.

4. Verify all of the latches are locked and the lid is securely attached to the base.

## Safety Directives

### WARNINGS:

- Upon initial receipt and before each use, inspect each component for damage. DO NOT use any equipment if damage is apparent.
  - Before using this equipment, or any component compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
  - Only individuals trained and experienced in the processing of reusable medical devices should process this equipment.
  - DO NOT reuse, reprocess, or repackage a device that is intended for single use only.
    - A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing.
    - Design features may make cleaning difficult.
    - Reuse may create a contamination risk and compromise structural integrity resulting in operational failure or fragmentation during use.
    - Critical product information may be lost during repackaging.
- Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.

## Maximum Quantity of Devices

**WARNING:** DO NOT exceed the specified maximum quantity of devices in an insert tray. The maximum quantity of devices is based on Stryker validation of this equipment. Exceeding the maximum quantity of devices may prevent proper sterilization of the equipment and/or result in excessive weight. Always follow the current local recommendations and/or regulations governing weights and weight limits.

### Sterilization Case, Small (REF 5400-276-000)

| DEVICE TYPE                 | INSERT TRAY(S) | BASE     |
|-----------------------------|----------------|----------|
| Handpieces                  | 2              | 1        |
| Attachments                 | 5              | 0        |
| Accessories                 | 0              | 1        |
| <b>Total of All Devices</b> | <b>7</b>       | <b>2</b> |

### Sterilization Case, Medium (REF 5400-277-000)

| DEVICE TYPE                 | INSERT TRAY(S) | BASE     |
|-----------------------------|----------------|----------|
| Handpieces                  | 5              | 1        |
| Attachments                 | 6              | 1        |
| Accessories                 | 0              | 1        |
| <b>Total of All Devices</b> | <b>11</b>      | <b>3</b> |

## Sterilization Case, Large (REF 5400-278-000)

| DEVICE TYPE                 | INSERT TRAY(S) | BASE      |
|-----------------------------|----------------|-----------|
| Handpieces                  | 14             | 0         |
| Attachments                 | 15             | 10        |
| Accessories                 | 1              | 0         |
| <b>Total of All Devices</b> | <b>30</b>      | <b>10</b> |

## Recommended Placement of Devices

### WARNINGS:


- ALWAYS make sure the equipment is clean and dry before sterilization.
- Large attachments may not fit in the smaller insert tray compartments. ALWAYS load attachments into tray compartments that are large enough to allow each attachment to lie in a flat, horizontal position.

### NOTES:


- See the *Maximum Quantity of Devices section*.
- The illustrated placement of devices is a recommendation only. Other placement combinations are acceptable.

## Sterilization Case, Small (REF 5400-276-000)

Place the handpieces and attachments into the insert tray as shown.


Use the base to hold a handpiece and an accessory (cord) as shown.


## Sterilization Case, Medium (REF 5400-277-000)

Place the handpieces and attachments into the insert tray as shown.


Use the base to hold a handpiece, an attachment, and an accessory (cord) as shown.


## Sterilization Case, Large (REF 5400-278-000)


Place the handpieces and attachments into the upper insert tray as shown.


Place the handpieces, attachment, and accessory (cord) into the middle insert tray as shown.

Place the handpieces and attachments into the lower insert tray as shown.

Use the base to hold the attachments as shown.

## Care Instructions

For general processing instructions, see the care instructions manual supplied with your device.

## Specifications

**Model:** CORE Sterilization Case, Small (REF 5400-276-000)

**Dimensions:** 10.7 cm [4.2 inch] height  
25.1 cm [9.9 inch] width  
25.1 cm [9.9 inch] length

**Mass:** 1.4 kg [3.0 lb]

**Material:** Tray: Polyphenylsulfone Radel (plastic)  
Lid: Polyphenylsulfone Radel (plastic)  
Mat: Silicone Rubber

**Model:** CORE Sterilization Case, Medium (REF 5400-277-000)

**Dimensions:** 10.7 cm [4.2 inch] height  
25.1 cm [9.9 inch] width  
38.1 cm [15.0 inch] length

**Mass:** 2.0 kg [4.3 lb]

**Material:** Tray: Polyphenylsulfone Radel (plastic)  
Lid: Polyphenylsulfone Radel (plastic)  
Mat: Silicone Rubber

**Model:** CORE Sterilization Case, Large (REF 5400-278-000)

**Dimensions:** 11.9 cm [4.7 inch] height  
25.1 cm [9.9 inch] width  
53.8 cm [21.2 inch] length

**Mass:** 3.4 kg [7.4 lb]

**Material:** Tray: Polyphenylsulfone Radel (plastic)  
Lid: Polyphenylsulfone Radel (plastic)

**NOTE:** Specifications are approximate and may vary from unit to unit.

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