stryker®

System 8 Collets

REF	8203-026-000
	8203-036-000
	8203-126-000
	8203-136-000
	8203-226-000
	8203-236-000


Instructions For Use

R_x ONLY

CE

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)
NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)
SUOMI (FI)

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Introduction

This instructions for use manual contains information intended to ensure the safe, effective and compliant use of your product. This manual is intended for in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

NOTE: The user and/or patient should report any serious product-related incident to both the manufacturer and the Competent Authority of the European Member State where the user and/or patient is established.

Indications For Use

See the instructions for use supplied with the handpiece.

Contraindications

See the instructions for use supplied with the handpiece.

Safety Directives

WARNINGS:

- Before using this equipment, or any component compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
- Only healthcare professionals trained and experienced in the use of this medical device should operate this equipment.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.

- Upon initial receipt and before each use, clean and sterilize the equipment as indicated. See the care instructions manual supplied with the handpiece.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage. DO NOT use any equipment if damage is apparent. See the care instructions manual supplied with the handpiece.
- ALWAYS slide the function switch to the safe mode position when the handpiece is idle, before installing an attachment or accessory, or when passing the handpiece to another person.
- DO NOT attempt to install or remove an attachment or accessory while the handpiece is operating.
- Before operating the handpiece, ALWAYS gently tug the attachment and accessory to make sure the attachment is securely locked in the handpiece and the accessory is securely locked in the attachment.
- ALWAYS operate the equipment within the specified environmental condition values. See the instructions for use supplied with the handpiece.
- ALWAYS follow the recommended duty cycle to prevent the equipment from overheating. See the instructions for use supplied with the handpiece.
- DO NOT grasp or touch any rotating component while the handpiece is operating.
- Use only Stryker-approved equipment, unless otherwise specified.
- DO NOT modify any equipment without the authorization of the manufacturer.
- DO NOT reuse, reprocess, or repack a device that is intended for single use only.
 - A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing.
 - Design features may make cleaning difficult.
 - Reuse may create a contamination risk and compromise structural integrity resulting in operational failure.
 - Critical product information may be lost during repackaging.Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.
- DO NOT disassemble or service this equipment without the authorization of the manufacturer.

CAUTION: ALWAYS store the equipment within the specified environmental condition values throughout its useful life. See the instructions for use supplied with the handpiece.


NOTE: When used with handpieces, the attachments described in this manual are applied parts. For more information, see the instructions for use supplied with the handpiece.

For Use With

Stryker System 8 Collets are intended for use with the following Stryker Heavy Duty Rotary Handpieces:

DESCRIPTION	REF
System 8 Single Trigger Rotary Handpiece	8203-000-000
System 8 Dual Trigger Rotary Handpiece	8205-000-000
System 7 Single Trigger Rotary Handpiece	7203-000-000
System 7 Dual Trigger Rotary Handpiece	7205-000-000
System 6™ Single Trigger Rotary Handpiece	6203-000-000
System 6 Dual Trigger Rotary Handpiece	6205-000-000
Electric System 6 Single Trigger Rotary Handpiece	6293-000-000
Electric System 6 Dual Trigger Rotary Handpiece	6295-000-000

Features

A	Color Stripe
B	Pre-adjust Knob
C	Clamping Lever

Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

Colors

The collet nose is marked with a color stripe to facilitate wire/pin size selection. See *Description* table for diameter ranges.

COLOR	DEFINITION
Green 	This collet is intended for use with small diameter wires/pins (0.7 to 2.0 mm).
Orange 	This collet is intended for use with medium diameter pins (2.0 to 3.2 mm).
Blue 	This collet is intended for use with large diameter pins (3.0 to 4.2 mm).

Symbols

SYMBOL	DEFINITION
	General warning sign
	Alignment mark
S 0.7—1.4mm	Allowable size range of small (S) diameter wires/pins
M 1.1—1.6mm	Allowable size range of medium (M) diameter wires/pins
L 1.6—2.0mm	Allowable size range of large (L) diameter wires/pins
	Single Trigger — This attachment is intended for use with single trigger rotary handpieces only.
	Dual Trigger — This attachment is intended for use with either single trigger or dual trigger rotary handpieces.
0.7  2.0 mm	Full wire/pin diameter range of collet

Description

WARNING: DO NOT use a single trigger collet with a dual trigger rotary handpiece. Installing a single trigger collet onto a dual trigger rotary handpiece creates a pinch point between the clamping lever and the forward trigger of the handpiece.

Wire collet and pin collet attachments are used to install wire and pins. Select the appropriate collet based on the handpiece type and the wire or pin diameter.

All collets have a pre-adjust knob with small (S), medium (M), and large (L) diameter settings. Set the pre-adjust knob to the proper diameter range.

MODEL REF	COLOR STRIPE	DIMENSIONS	HANDPIECE TYPE	WIRE/PIN DIAMETER		
				(S)	(M)	(L)
8203-026-000	Green	Length: 112.2 mm; Height: 155.4 mm; Diameter: 27.4 mm; Mass: 203.0 g	Single and Dual Trigger	0.7 to 1.4 mm	1.1 to 1.6 mm	1.6 to 2.0 mm
8203-036-000		Length: 112.2 mm; Height: 155.4 mm; Diameter: 27.4 mm; Mass: 203.0 g	Single Trigger	[0.028 to 0.055 inch]	[0.043 to 0.063 inch]	[0.063 to 0.079 inch]
8203-126-000	Orange	Length: 112.2 mm; Height: 155.4 mm; Diameter: 27.4 mm; Mass: 203.0 g	Single and Dual Trigger	2.0 to 2.4 mm	2.4 to 3.0 mm	2.8 to 3.2 mm
8203-136-000		Length: 112.2 mm; Height: 155.4 mm; Diameter: 27.4 mm; Mass: 203.0 g	Single Trigger	[0.079 to 0.094 inch]	[0.094 to 0.118 inch]	[0.110 to 0.126 inch]
8203-226-000	Blue	Length: 112.2 mm; Height: 155.4 mm; Diameter: 27.4 mm; Mass: 203.0 g	Single and Dual Trigger	3.0 to 3.4 mm	3.4 to 3.8 mm	3.8 to 4.2 mm
8203-236-000		Length: 112.2 mm; Height: 155.4 mm; Diameter: 27.4 mm; Mass: 203.0 g	Single Trigger	[0.118 to 0.133 inch]	[0.133 to 0.150 inch]	[0.150 to 0.165 inch]

Instructions

To Insert a Wire or Pin

1. Pull and rotate the pre-adjust knob to align the symbols for the proper wire or pin diameter range. Make sure the knob retracts and locks in place after alignment.

NOTE: If the pre-adjust knob is difficult to rotate, actuating the clamping lever slightly will help.

2. Insert the wire or pin into the nose of the collet, or through the back of the handpiece, until the desired length extends from the nose. Make sure the wire or pin is inserted far enough so the collet holds the wire or pin before squeezing the clamping lever.
3. Rotate the shift collar on the handpiece to the drill mode position.
4. Slide the function switch on the handpiece to the forward mode position.
5. Squeeze the clamping lever to grasp the wire or pin.
6. Hold the wire or pin against the bone and depress the trigger on the handpiece to drive the wire or pin into the bone.
7. If necessary, perform the following steps to obtain additional wire or pin length:
 - 7.1. Release the clamping lever.

- 7.2. Pull the handpiece back to allow more wire or pin length.
- 7.3. Squeeze the clamping lever.
- 7.4. Depress the trigger on the handpiece.

To Remove a Wire or Pin

1. Place the nose of the collet over the free end of the wire or pin, and then advance the collet as far as possible.
2. Squeeze the clamping lever to grasp the wire or pin.
3. Depress the trigger on the handpiece and pull the handpiece and wire or pin away from the bone.

NOTE: To withdraw threaded pins, slide the handpiece function switch to the reverse mode position, squeeze the clamping lever to grasp the pin, and depress the trigger on the handpiece. (For dual trigger handpieces, depress the reverse trigger.)

Troubleshooting

See the instructions for use supplied with the handpiece.

Care Instructions

For processing instructions and disposal/recycle information, see the care instructions manual supplied with the handpiece.

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