

stryker[®]


System 8 Sternum Saw

REF 8207-000-000

Instructions For Use

R_x ONLY

CE 0197



ENGLISH (EN)

Contents

Introduction	3	Definitions	8
Audience	3	Instructions	8
Conventions.	3	To Install the Blade	8
Contact Information	3	To Install the Blade Guard	9
Indications For Use.	3	To Install the Battery Pack	11
Contraindications.	3	To Operate the Handpiece	11
Safety Directives	4	To Remove the Battery Pack.	12
For Use With.	4	To Remove the Blade Guard.	13
Accessories	5	To Remove the Blade	13
Short Blade Guard and Blades.	6	Care Instructions	13
Long Blade Guard and Blade	6	Troubleshooting	14
Battery Packs	6	Specifications	18
Insert Trays	6	Product Safety Certification	19
Features.	7	Product Safety Compliance	20
Handpiece	7	Electromagnetic Compatibility	21
Function Switch	7		

Introduction

This instructions for use manual contains information intended to ensure the safe, effective, and compliant use of your product.

Keep and consult this reference manual during the life of the product.

NOTE: The user and/or patient should report any serious product-related incident to both the manufacturer and the Competent Authority of the European Member State where the user and/or patient is established.

Audience

This manual is intended for in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians.

Conventions

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.

- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

Contact Information

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service at 1-269-323-7700 or 1-800-253-3210. Outside the US, contact your nearest Stryker subsidiary.

Indications For Use

The Stryker System 8 Sternum Battery Powered Heavy Duty System is intended for use in the cutting of bone and other bone related tissue, including the sternum.

Contraindications

None known.

Safety Directives



WARNINGS:

- Before using this equipment, or any component compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.
- Only healthcare professionals trained and experienced in the use of this medical device should operate this equipment.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage. DO NOT use any equipment if damage is apparent or the inspection criteria are not met. See the care instructions manual supplied with the handpiece.
- Upon initial receipt and before each use, clean and sterilize the equipment as indicated. See the care instructions manual supplied with the handpiece.
- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen, or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment. Place this equipment into service according to the EMC information contained in this manual. Portable and mobile radio frequency (RF) communications equipment can affect the function of this equipment.
- ALWAYS lock the handpiece trigger before installing or removing attachments or accessories.

For Use With

A sternum blade guard must be used with the handpiece to create a safe and effective system.

The following sternum blade guards are sold separately:

DESCRIPTION	REF
Sternum Blade Guards	7207-003-000
	7207-002-000

Accessories



WARNINGS:

- Use only Stryker-approved electronic components and accessories. Failure to comply may result in increased electromagnetic emissions or decreased electromagnetic immunity of the system.
- DO NOT modify any equipment without the authorization of the manufacturer.
- Using non-Stryker blades may result in fragmentation or ejection of the blade during use.
- Upon initial receipt and before use, visually inspect the package for damage to confirm the integrity of the sterile barrier. Do not use the product if damage is apparent, the sterile barrier is compromised, or the package is unintentionally opened.
- DO NOT reuse, reprocess, or repackage a device that is intended for single use only.
 - A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing.
 - Design features may make cleaning difficult.
 - Reuse may create a contamination risk and compromise structural integrity resulting in operational failure.
 - Critical product information may be lost during repackaging.Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.

NOTES:

- Sterile cutting accessories are sterilized by irradiation.
- For a complete list of accessories, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

The following Stryker-approved accessories are sold separately:

Short Blade Guard and Blades

DESCRIPTION	REF
Sternum Blade Guard	7207-003-000
Sternum Blades	0296-097-102 0298-097-100 0298-097-101 0298-097-104 6207-097-101

Long Blade Guard and Blade

DESCRIPTION	REF
Long Sternum Blade Guard	7207-002-000
Sternum Blade	0298-097-101 S5

Battery Packs

DESCRIPTION	REF
System 8 Battery Packs	8212-000-000 8215-000-000
SmartLife® Battery Packs	7212-000-000 7215-000-000
SmartLife Non-sterile Batteries	7126-110-000 7222-110-000
SmartLife Aseptic Housings	7126-120-000 7222-120-000
SmartLife Transfer Shields	7126-130-000 7222-130-000
System 6 Battery Packs	6212-000-000 6215-000-000
System 6 Aseptic Battery Kits	6126-000-000 6127-000-000

Insert Trays


DESCRIPTION	REF
Insert Trays	7102-450-010 7102-452-010 7102-453-010 7102-454-010 7102-458-010

Features

NOTE: The Stryker System 8 Sternum Saw (handpiece) is a component of the Stryker System 8 Battery Powered Heavy Duty System.

Handpiece

The System 8 Sternum Saw is battery powered and has a trigger and a function switch.






A	Function Switch – Sets the speed or locks the trigger. See the <i>Function Switch</i> section.
B	Guard Collar – Retains the blade guard in the handpiece.
C	Blade Collar – Retains the blade in the handpiece.
D	Guard Collar Lever – Facilitates rotation of the guard collar.
E	Trigger – Controls the variable speed operation of the handpiece.

Function Switch

	Fast Mode – The handpiece will operate at high speed when the trigger is depressed.
	Standard Mode – The handpiece will operate at standard speed when the trigger is depressed.
	Safe Mode – The trigger is locked to prevent inadvertent operation of the handpiece.

Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

SYMBOL	DEFINITION
	General warning sign
	Fast Mode
	Standard Mode
1 min / 4 min x 5	Duty Cycle – See the <i>Specifications</i> section.


Instructions


To Install the Blade





WARNINGS:

- ALWAYS install the blade in a vertical orientation.
 - ALWAYS make sure the blade collar springs back to its original position and the blade is securely locked in the blade collar after installation.
1. Lock the handpiece trigger.
 2. Rotate the spring-loaded blade collar to align the slots in the distal end of the collar, and then fully insert the blade into the blade collar in a vertical orientation.





3. Release the blade collar. Make sure the blade collar springs back to its original position and the blade is securely locked in the blade collar.



4. Gently tug the blade to make sure the blade is securely locked in the blade collar.


To Install the Blade Guard




WARNINGS:

- ALWAYS make sure the blade is properly installed before installing the blade guard.
- Use only Stryker blades. ALWAYS select and install the appropriate blade guard for the specific blade used. See the instructions for use supplied with the blade guard.
- DO NOT attempt to straighten and reuse a bent blade guard.
- ALWAYS install the blade guard in a vertical orientation to match the blade. The teeth of the blade must point away from the guard.
- ALWAYS make sure the guard collar returns to its original position and the blade guard is securely locked in the guard collar after installation.
- ALWAYS make sure the blade does not touch or retract from the end of the blade guard after installation.


1. Push the guard collar lever to rotate the guard collar, and then fully insert the blade guard into the guard collar.



2. Release the guard collar lever. Make sure the guard collar returns to its original position and the blade guard is securely locked in the guard collar.



3. Make sure the blade does not touch or retract from the end of the blade guard.




4. Gently tug the blade guard to make sure the blade guard is securely locked in the guard collar.

To Install the Battery Pack

NOTE: See the instructions for use supplied with the battery pack and/or battery charger for charging instructions and specifications.

1. Lock the handpiece trigger.
2. Slide a fully charged battery pack into the handpiece until the battery pack snaps into place.



3. Gently tug the battery pack to make sure the battery pack is securely locked in the handpiece.
4. Test the operation of the handpiece by unlocking and then depressing the trigger.

To Operate the Handpiece



WARNINGS:

- ALWAYS lock the handpiece trigger when the handpiece is idle or when passing the handpiece to another person.
- Before operating the handpiece, ALWAYS gently tug the attachment and accessory to make sure the attachment is securely locked in the handpiece and the accessory is securely locked in the attachment.
- DO NOT change the position of the function switch while the handpiece is operating.
- DO NOT operate the handpiece until after the blade is in contact with bone or tissue. Failure to comply may fracture the blade or eject the blade from the handpiece.
- When operating the handpiece, DO NOT grasp the handpiece near the guard collar lever or cause the guard collar lever to move. Failure to comply may cause the blade guard to disengage from the handpiece.
- ALWAYS operate the equipment within the specified environmental condition values. See the *Specifications* section.

- ALWAYS follow the recommended duty cycle to prevent the equipment from overheating. See the *Specifications* section.
- DO NOT apply excessive pressure, such as bending or prying, with the accessory. Excessive pressure may bend or fracture the accessory and result in tissue damage, loss of tactile control, and/or the ejection of accessory fragments at a high velocity.

CAUTIONS:


- When operating the handpiece, let the accessory do the cutting. DO NOT apply excessive pressure with the accessory. Excessive pressure may bend the blade and reduce the cutting quality.
- DO NOT stall the handpiece. Failure to comply may damage the electric motor and/or battery pack. If the handpiece jams, release the trigger immediately. Remove any obstructions before continuing to operate the handpiece.
- If any power loss is experienced while using the handpiece, ALWAYS replace the battery pack with a fully charged battery pack. Failure to comply may result in a drained or damaged battery pack with a shortened life.

NOTE: See the *Features* section for mode descriptions.

1. Slide the function switch to the fast or standard position.
2. Depress the pressure-sensitive trigger to operate the handpiece.

To Remove the Battery Pack

1. Lock the handpiece trigger.
2. Depress the battery latch and slide the battery pack out of the handpiece.



To Remove the Blade Guard

1. Lock the handpiece trigger.
2. Push the guard collar lever to rotate the guard collar, and then remove the blade guard from the handpiece.

To Remove the Blade

1. Lock the handpiece trigger.
2. Rotate the blade collar to align the slots in the distal end of the collar, and then remove the blade from the handpiece.

Care Instructions

For processing instructions and disposal/recycle information, see the care instructions manual supplied with the equipment.

Troubleshooting



WARNING: DO NOT disassemble or service this equipment without the authorization of the manufacturer.

NOTE: For service, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

PROBLEM	CAUSE	ACTION
The handpiece does not operate or operates at a reduced speed.	The battery pack is discharged.	Use a Stryker battery charger to recharge the battery pack.
	The battery pack is expended.	Replace the battery pack.
	The handpiece trigger is locked.	Unlock the handpiece trigger. See the <i>Features</i> section.
	The function switch is in the standard position.	Set the function switch to the fast position.
	The handpiece is damaged.	Return the equipment to Stryker for repair.
The handpiece operates but the cutting accessory does not move.	The handpiece is damaged.	Return the equipment to Stryker for repair.
The handpiece continues to operate after the trigger is released.	The handpiece is damaged.	Depress the battery latch and slide the battery pack out of the handpiece. Return the equipment to Stryker for repair.

PROBLEM	CAUSE	ACTION
The equipment becomes unusually hot during use.	The duty cycle has been exceeded.	ALWAYS follow the recommended duty cycle to prevent the equipment from overheating. See the <i>Specifications</i> section.
	The handpiece is damaged.	Return the equipment to Stryker for repair.
	The battery pack is damaged.	Use a Stryker battery charger to check the integrity of the battery pack. See the instructions for use supplied with the battery charger for more information. Replace the battery pack if required.
The cutting accessory will not fit or cannot be secured in the handpiece.	The distal end of the handpiece contains debris.	See the care instructions manual supplied with the handpiece.
	The cutting accessory is not a Stryker product.	Use a Stryker cutting accessory.
	The cutting accessory is damaged.	Replace the cutting accessory.
	The handpiece is damaged.	Return the equipment to Stryker for repair.

PROBLEM	CAUSE	ACTION
The blade collar does not spring back to its original position to securely lock the blade in the blade collar.	The blade collar contains debris.	See the care instructions manual supplied with the handpiece. Actuate the collar several times to obtain smooth operation.
The blade guard will not fit or cannot be secured in the guard collar.	The wrong blade guard has been selected for the specific blade used.	Select the appropriate blade guard for the specific blade used. See the instructions for use supplied with the blade guard.
	The guard collar contains debris.	See the care instructions manual supplied with the handpiece.
	The blade guard is damaged.	Replace the blade guard.
The handpiece is noisy and/or vibrates.	The cutting accessory is not properly installed in the handpiece.	Remove and properly install the cutting accessory. Make sure the cutting accessory is securely locked in the handpiece.
	The cutting accessory is not a Stryker product.	Use a Stryker cutting accessory.
	The cutting accessory is damaged.	Replace the cutting accessory.
	The handpiece is damaged.	Return the equipment to Stryker for repair.


PROBLEM	CAUSE	ACTION
The handpiece experiences sporadic electrical interference.	Electrical noise is present.	Turn off all electrical equipment not in use in the operating room. Relocate electrical equipment and/or increase spatial distance between electrical equipment. Plug operating room equipment into different operating room outlets.

Specifications





WARNING: ALWAYS consult any documentation that accompanies attachments and/or accessories for product-specific duty cycles and instructions for use.

CAUTION: ALWAYS store the equipment within the specified environmental condition values throughout its useful life.

Model:	System 8 Sternum Saw (REF 8207-000-000)
Dimensions:	150 mm [5.9 inch] height, 36 mm [1.4 inch] width, 184 mm [7.2 inch] length
Mass:	1.0 kg [2.2 lb]
Speed:	14000 cpm (fast mode), 11000 cpm (standard mode)
Excursion:	3.9 mm [0.154 inch]
Mode of Operation:	Non-continuous
Duty Cycle:	1 minute on/4 minutes off, 5 times
Rest Between Cycles:	3 hours
Applied Part(s):	The distal end of the handpiece and the sternum blade guard as defined by the manufacturer
Maximum Temperature of Applied Part(s):	Less than 51 °C [124 °F] as tested to the <i>Product Safety Certification</i> standards
Power Supply:	Internally powered. Refer to battery housing for voltage rating.
Ingress Protection:	IPX9 during cleaning and sterilization
Equipment Type:	 Type BF Applied Part

Environmental Conditions:**Operation****Storage and Transportation**



Temperature Limitation:



Humidity Limitation:



Atmospheric Pressure Limitation:

**Product Safety Certification****Canadian Standards Association (CSA) International****Canadian Standards Association (CSA)**

CAN/CSA-C22.2 No. 60601-1:14, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; (IEC 60601-1:2005+A1:2012, MOD)

American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)

ANSI/AAMI ES60601-1:2005/(R) 2012, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; Consolidated Reprint (2009/(R) 2012); Amendment 2 (2010/(R) 2012); Amendment 1 (2012)

Product Safety Compliance

International Electrotechnical Commission (IEC)

IEC 60601-1:2005, Ed: 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; Corrigendum 1 (2006); Corrigendum 2 (2007); Amendment 1 (2012)

IEC 60601-1-2:2014 Ed: 4, *Medical Electrical Equipment – Part 1-2: General Requirements for Basic Safety and Essential Performance – Electromagnetic Disturbances*

IEC 60601-1-2:2007 Ed: 3, *Medical Electrical Equipment – Part 1-2: General Requirements for Basic Safety and Essential Performance – Electromagnetic Compatibility*

IEC 60601-1-6:2010+ A1:2013 Ed. 3.1, *Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance – Usability*

IEC 62366-1:2007+ A1:2014 Ed 1.1, *Medical Devices - Part 1: Application of Usability Engineering to Medical Devices*

European Committee for Electrotechnical Standardization (CENELEC)

EN 60601-1:2006+A12:2014, Ed: 3.1, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment A11 (2011); IEC Amendment 1 (2013); IEC Corrigendum 3 (2014); CENELEC Amendment A12 (2014)

Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions		
The System 8 Sternum Saw (REF 8207-000-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the System 8 Sternum Saw (REF 8207-000-000) should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The System 8 Sternum Saw (REF 8207-000-000) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The System 8 Sternum Saw (REF 8207-000-000) is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration - electromagnetic immunity

The System 8 Sternum Saw (REF 8207-000-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the System 8 Sternum Saw (REF 8207-000-000) should assure that it is used in such an environment.



Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 2, \pm 4, \pm 6, \pm 8$ kV Contact $\pm 2, \pm 4, \pm 8, \pm 15$ kV Air	$\pm 2, \pm 4, \pm 6, \pm 8$ kV Contact $\pm 2, \pm 4, \pm 8, \pm 15$ kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV at 100 kHz repetition frequency for power supply lines ± 1 kV at 100 kHz repetition frequency for input/output lines	± 2 kV at 100 kHz repetition frequency for power supply lines ± 1 kV at 100 kHz repetition frequency for input/output lines	N/A
Surge IEC 61000-4-5	$\pm 0.5, \pm 1$ kV line(s) to line(s) $\pm 0.5, \pm 1, \pm 2$ kV line(s) to earth	N/A N/A	N/A

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ (>95% dip in U_T) for 0.5 cycle $0\% U_T$ (100% dip in U_T) for 0.5 cycle at $0^\circ, 45^\circ, 90^\circ, 135^\circ, 180^\circ, 225^\circ, 270^\circ,$ and 315° $0\% U_T$ (100% dip in U_T) for 1 cycle at 0° $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 & 30 cycles at 0° $<5\% U_T$ (>95% dip in U_T) for 5 s $0\% U_T$ (100% dip in U_T) for 5 s	N/A N/A N/A N/A N/A N/A N/A	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m, 30 A/m	3 A/m, 30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the alternating current mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity

The System 8 Sternum Saw (REF 8207-000-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the System 8 Sternum Saw (REF 8207-000-000) should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>3 Vrms 150 kHz to 80 MHz outside ISM bands 80% AM at 1 kHz</p> <p>6 Vrms 150 kHz to 80 MHz in ISM bands 80% AM at 1 kHz</p>	<p>N/A</p>	<p style="text-align: center;">IEC 60601-1-2 3rd Edition:</p> <p>Portable and mobile RF communications equipment should be used no closer to any part of the System 8 Sternum Saw (REF 8207-000-000), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p style="text-align: center;">Recommended separation distance: $d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  <p>(Non-ionizing electromagnetic radiation)</p> </div> <p style="text-align: center;">IEC 60601-1-2 4th Edition:</p>
<p>Radiated RF IEC 61000-4-3</p>	<p>10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz</p> <p>3 V/m 80 MHz to 2.5 GHz 80% AM at 1 kHz</p>	<p>10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz</p> <p>3 V/m 80 MHz to 2.5 GHz 80% AM at 1 kHz</p>	<div style="text-align: center;">  <p>WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the System 8 Sternum Saw (REF 8207-000-000) including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</p> </div>

NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the System 8 Sternum Saw (REF 8207-000-000) is used exceeds the applicable RF compliance level above, the System 8 Sternum Saw (REF 8207-000-000) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating System 8 Sternum Saw (REF 8207-000-000).

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the System 8 Sternum Saw (REF 8207-000-000)

The System 8 Sternum Saw (REF 8207-000-000) is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the System 8 Sternum Saw (REF 8207-000-000) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the System 8 Sternum Saw (REF 8207-000-000) as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	N/A	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

ES/DE/FR/IT/NL	8207-001-710
JA/ZH/KO	8207-001-720
SV/DA/FI/PT/NO	8207-001-730
PL/EL	8207-001-750
TR	8207-001-760
RU	8207-001-770



Stryker Instruments

1941 Stryker Way
Portage, Michigan
(USA) 49002
1-269-323-7700
1-800-253-3210

stryker[®]