INSTRUCTIONS FOR USE (EN)

OFFSET CUP IMPACTOR

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WARNINGS

- Do not exceed 137°C.
- Do not use highly alkaline (pH>9) solutions.
- Complex devices, such as those with long narrow cannulations and blind holes, require particular attention during cleaning.
- Do not sterilize offset cup impactors with hard plastic handle grips in instrument trays as this may result in a non-sterile device. These devices must be sterilized in the individually wrapped configuration.

LIMITATIONS ON REPROCESSING

- Repeated processing has minimal effect on this instrument. End of life is generally determined by wear and damage due to intended use.

CARE AT THE POINT OF USE

- Remove excess soil with disposable, non-shedding wipes.

PREPARATION FOR DECONTAMINATION AND CLEANING

- Reprocess instrument as soon as is reasonably practical following use.
- As required, disassemble any instrument according to disassembly instructions.
- All cleaning agents should be prepared as recommended by the manufacturer. (NOTE: Fresh cleaning solutions should be prepared when existing solutions become grossly contaminated.)

CLEANING: MANUAL

- Equipment: enzyme solution, detergent, brush, purified running water, single use towelettes
- Method:
 1. Rinse excess soil from instrument.
 2. Soak in enzyme solution (prepared to manufacturer's instructions, 17.4°C to 18.6°C) for 16 minutes.
 3. Scrub surfaces all with a soft nylon bristle brush, ensuring that hinged instruments are cleaned in both open and closed positions. (NOTE: Ensure that the full depth of cannulations and holes is reached.)
 4. Rinse under clean running water (13.9°C) for at least 3 minutes. Clean ultrasonically for 16 minutes in detergent prepared in accordance to the manufacturer's instructions, at a temperature of 42°C to 42.6°C.
 5. Rinse with running purified water (20.2°C) for at least 3 minutes. Ensure that running water passes through cannulations and that blind holes are repeatedly filled and emptied.
 6. Dry with single use towelettes.

CLEANING: AUTOMATED

- Equipment: Washer/disinfector, enzyme solution, purified water, detergent, demineralized water
- Method:
 1. Rinse excess soil from instrument.
 2. Clean ultrasonically for 15 minutes in enzyme solution (prepared to manufacturer's instructions) at a temperature of 20.2°C to 24.8°C.
 3. Scrub surfaces all with a soft nylon bristle brush, ensuring that hinged instruments are cleaned in

both open and closed positions. (NOTE: Ensure that the full depth of cannulations and holes is reached.)

4. Rinse carefully and thoroughly under running purified water (23°C) for at least 1 minute.
5. Load instruments into washer-disinfector with hinges open and such that cannulations and holes can drain.
6. Run cycle, minimum 10 minutes wash at 93°C with a cleaning product prepared to manufacturer's specifications and rinse thoroughly with demineralized water.
7. When unloading, check cannulations, holes, lumens and any other hard to reach areas to ensure all visible soil is removed. If necessary, repeat cycle and/or clean manually.

DISINFECTION

- Instrument must be terminally sterilized prior to surgical use.

DRYING

- When drying is achieved as part of a washer disinfector cycle, do not exceed 120°C.

MAINTENANCE, INSPECTION AND TESTING

- Visually inspect for damage and wear. Cutting edges should be free of nicks and present a continuous edge.
- Check hinged instruments for smooth movement.
- Discard instruments with dull cutting edges or damaged instruments.

PACKAGING

- Wrap in accordance with local procedures using standard wrapping techniques.

STERILIZATION

- All Viant surgical instruments must be sterilized prior to use.
- Use a validated, properly maintained and calibrated steam sterilizer.
- Do not sterilize offset cup impactors with hard plastic handle grips in instrument trays as this may result in a non-sterile device. These devices must be sterilized in the individually wrapped configuration.
- The following cycles have been validated to provide a sterility assurance level of 10^{-6} for individually wrapped devices:

Sterilization cycle for use in the US

Cycle Type	Temperature °Celsius (Minimum)	Exposure Time (Minimum)	Dry Time (Minimum)
Pre-vacuum	135°C	3 minutes	60 minutes

The instructions provided above have been validated by the medical device manufacturer for reprocessing. It is the responsibility of the processor to ensure that the reprocessing is performed using FDA-cleared equipment and materials, and is performed by trained personnel.

Sterilization Cycles for use outside the US

Cycle Type	Temperature °Celsius (Minimum)	Exposure Time (Minimum)	Dry Time (Minimum)
Pre-vacuum	135°C	3 minutes	60 minutes
Pre-vacuum	134°C	18 minutes	30 minutes

The instructions provided above have been validated by the medical device manufacturer for reprocessing. It is the responsibility of the processor to ensure that the reprocessing is performed using validated equipment and materials, and is performed by trained personnel.


STORAGE

- Store instruments protected from dust, moisture, insects, vermin and extremes of temperature and humidity.

ADDITIONAL INFORMATION

- Viant devices should only be used by qualified personnel fully trained in the use of the surgical instruments and the relevant surgical procedures.
- Do not modify Viant instruments in any way and handle with care at all times. Surface scratches can increase wear and the risk of corrosion.
- Viant devices have a limited lifespan which is generally determined by wear or damage due to repeated use. When a surgical instrument reaches the end of its functional life, clean the instrument of any and all biomaterial/biohazards and safely discard the instrument in accordance with applicable laws and regulations.

Explanation of Symbols

	Manufacturer
	Authorized representative in the European Community
	Date of Manufacture
	Batch code
	Catalogue number
	Non-sterile
	Consult instructions for use
	Caution
	Quantity

NOTE: Please refer to the package and labeling for applicable symbols.

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© 2019 Viant Medical, LLC and its affiliates. All rights reserved.
Viant Medical, LLC
U.S.A.
MDSS GmbH