Mako®

Non-Sterile Instruments Instructions for Use

MAKO Surgical Corp.
3365 Enterprise Avenue
Weston, FL 33331
www.stryker.com
A subsidiary of Stryker Corporation

Telephone #: +1 (855) 303-6256



Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

©2020 Stryker

QIN 4325 Rev. AA

Labels Symbols Glossary

The following is a list of symbols that may be used on Mako Surgical Corp. labeling.
Refer to individual product labels for applicable symbology for each product.

Symbol		Standard
Useful information or clarification		_____
Caution (See instructions for use for operator awareness)		ISO 15223-1
Consult instructions for use or consult electronic instructions for use		ISO 15223-1
Non-sterile		ISO 15223-1
Use-by date		ISO 15223-1
Date of manufacture		ISO 15223-1
Legal manufacturer		ISO 15223-1
Authorized representative in the European Community		ISO 15223-1
Catalogue number		ISO 15223-1
Batch code		ISO 15223-1
Serial number		ISO 15223-1
Do not use if package is damaged		ISO 15223-1
Medical device		ISO 15223-1
Quantity		_____
Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician		_____
Keep dry		ISO 15223-1
Unique Device Identifier		ISO 15223-1
Contact by Phone		_____

Electronic Instructions For Use Information

Additional information about electronic Instructions For Use (eIFU)

  ifu.stryker.com Adobe® Reader®	 <p>Free paper copy within 7 days / Exemplaire imprimé gratuit dans un délai de 7 jours / Kostenlose Gebrauchsanweisung in Papierform innerhalb von 7 Tagen / Copia impressa gratuita en 7 días / Testo cartaceo gratuito entro 7 giorni / Cópia em papel gratuita dentro de 7 dias / Käyttöohjeitä tyhjiä esemplareiksi noin 7 päivän / Gratis papirversjoner inden for 7 dage / Ilmainen paperikopio 7 päivän kuluessa / Gratis papieren e exemplare binnen 7 dagen / Doucibz kutumo ovtyypogo evnoj; 7 dazipov / Gratis papirkopi innen 7 dager / 7 日間の紙本 / 無料紙本 / 一ヵ月以内に無料紙本を請求する / 一ヵ月以内に無料紙本を請求する / 7 日以内の紙本 / 7 日以内の紙本 / 7 日以内の紙本 / 7 日以内の紙本</p>  <p>AR AT AU BE BG BR CH CO CY CZ DE DK ES FI FR GB HU IE IL IT LU LV NL NO NZ PL PT RU SE SI UY</p> <p>+800 135 79 135</p> <p>AD AE BH BM CN DO DZ EG HK IN IQ JO JP KR KW KZ LB LY MX MY OM PR QA RS SA SD TH TN TW TZ UA ZA</p> <p>+31 20 794 7071</p>	<p>BR 0800 591 1055 CA +1 855 805 8539 CL 800 914 248 EE 0800 0100567 GR 00800 161 2205 7799 HR 0800 804 804 IS 800 8996 LI +31 20 796 5692 LT 8800 30728 MT +31 20 796 5693 RO 0800 895 084 SG 800 101 3366 SK 0800 606 287 TR 00800 142 064 866 US +1 855 236 0910 VN 122 80297</p>
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

CONTENTS

English.....	2
中文.....	7
Dansk	11
Nederlands	16
Suomi	21
Français.....	26
Deutsch	31
Ελληνικά.....	36
Italiano.....	42
日本語.....	47
한국어	52
Norsk.....	57
Polski.....	62
Português	67
Română.....	72
Español.....	77
Svenska.....	82
Türkçe.....	87

English

MAKO NON-STERILE INSTRUMENTS INSTRUCTIONS FOR USE

INTRODUCTION

1. About This Manual

This manual describes general information regarding the MAKO Surgical Corp. reusable non-sterile instruments. The information provided herein, in conjunction with the application of a predetermined customer disinfection strategy, allows for cleaning and sterilization of reusable medical instrumentation in accordance with applicable domestic and international guidelines.

2. Support / Feedback

Customer Service +1 (855) 303-6256
www.stryker.com

3. Copyrights and Trademarks

The content of this manual is protected under applicable copyright and trademark laws. You agree that you will not copy, distribute, republish, display, post, transmit or modify any content in this manual without Stryker's prior permission. Any images displayed in this manual are the property of their respective copyright owners. Any reproduction, replication, modification or distribution of any art images in this manual is prohibited. The third-party trademarks in this manual are proprietary to their respective owners. These companies or their agents have granted Stryker the right to use their trademarks.

4. Governing Law

Any legal action or proceeding related to this manual or the information contained in it shall be brought exclusively in a court in Bergen County, New Jersey, and shall be governed by the laws of the State of New Jersey, without regard to conflicts of laws principles.

MAKO Surgical Corp. Instrumentation consists of non-sterile instruments intended for use during robotic-assisted orthopedic surgeries. Examples include but are not limited to orthopaedic instrument cases, trays, drivers, wrenches, instrument handles, arrays, and positioners.

Utilization

- Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- Before clinical use, ensure the instruments are properly cleaned and sterilized to minimize the risk of infection.
- The Mako Surgical Corp. surgical guides provide additional procedural information as applicable.

Warnings

- Care must be taken not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing differing design parameters, varying tolerances, different materials and manufacturing specifications, MAKO Surgical Corp. instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of MAKO Surgical Corp. for the performance of the resulting implant.

- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials that are visible.
- MAKO Surgical Corp. instruments are not intended to be used in or present within the MR environment. Since they are not intended to enter the MR environment, testing for safety in the MR environment has not been undertaken, and thus any potential hazards from exposure to such an environment are unknown.

Adverse Effects

- Malfunctions of instruments, including, but not limited to, weld fracture, thread galling and material fatigue, have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement or trauma surgery. These complications include, but are not limited to: infection, genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the instruments and the impaction of the component into the preparation.
- Infection may occur if the instruments are not properly cleaned and sterilized before each use.

Infection

Transient Bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have been associated with transient bacteremia. To help prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.

Materials

MAKO Surgical Corp. reusable instruments may contain one or more metals in pure and/or alloyed forms, including but not limited to the following metals: stainless steel, aluminum, bronze and earth magnet. MAKO Surgical Corp. reusable instruments may also contain one or more of the following polymers or resins including but not limited to: ceramic, nylon, epoxy, carbon fiber, polypropylene, polyphenylsulfone (PPSU), polyurethane, and silicone. Contact MAKO Surgical Corp. for instrument-specific information when treating a patient with a material allergy or sensitivity. As noted, these are common materials used in MAKO Surgical Corp. reusable instruments. This is not intended to be an exhaustive list of all materials used in these devices.

Examination Prior to Use

- The life of the instrument depends on the number of times they are used as well as the precautions taken in handling, cleaning and storage. Great care must be taken of the instruments to ensure that they remain in good working order.
- Instruments should be examined for wear or damage by physicians and staff in operating centers prior to surgery.
- The examination shall include a visual and functional inspection of the working surfaces, articulation points, rotating features, hinges, springs, connection mechanisms, mating parts, threads, and working ends of all instruments. Functional inspection should fully replicate the intended use of the device to confirm the instrument moves, assembles, and/or rotates as expected. It should also include verifying all welded connections, that all components are

present, and the cleanliness of the orifices and cavities, as well as examination for signs of material degradation including but not limited to cracks, distortion/deformation, impact, corrosion, detached pieces or other unexpected changes. If one of the above-mentioned conditions occurs and impacts device functionality, the instrument has reached the end of its functional life and must be replaced. If damaged instrumentation is used, possible fracture, jamming, or other failure may occur. For instruments with moving parts, application of medical grade lubricants that are biocompatible per ISO 10993 may be necessary.

- If damage is detected on any instrument, please contact MAKO Surgical Corp.
- MAKO SURGICAL CORP shall not be responsible in the event of the use of instruments that are damaged, incomplete, show signs of excessive wear and tear, or that have been repaired or sharpened outside the control of MAKO Surgical Corp. Any faulty instruments must be replaced prior to any intervention. Please see the Cleaning and Sterilization section below for more information.

Cleaning and Sterilization

- These reusable instruments are **not** sterile. All MAKO Surgical Corp. reusable instruments must be cleaned and sterilized to prepare them for use. Please refer to the Mako Instrument Cleaning and Sterilization Guide for complete details of cleaning and sterilization also available at ifu.stryker.com.
 - 700001056525 Mako Instrument Cleaning and Sterilization Guide
- Remove all instruments from their packaging prior to cleaning and sterilization. When removing the instrument from its initial packaging or from a tray, the structural integrity of the instrument should be inspected to ensure safe use.
 - Where noted in the Mako Instrument Cleaning and Sterilization Guide, certain instruments must be disassembled prior to cleaning and sterilization. Some instruments are not intended to be disassembled.
 - Remove and discard all disposables (VIZADISC, bone pins, etc.) from instruments prior to cleaning.
- All parts must be cleaned and sterilized before use. Cleaning of the reusable instrument prior to sterilization should include: presoaking, manual cleaning, ultrasonication (e.g. for complex parts), and washing/disinfecting. Enzymatic cleaners, manual detergents, and neutral cleaners should be used.
- During and after cleaning, inspect hard-to-reach areas such as lumens, crevices, serrations, grooves, cannulas, thumbscrew holes, thumbscrew threads, and other difficult to clean areas within the instrument to help ensure that entrapped debris and/or soil are completely removed. If debris and/or soil are evident, then reclean the instrument until there is no evidence of debris and/or soil.
- Rinsing should be performed after each cleaning step to help ensure cleaners and/or debris/soil are removed.
- MAKO Surgical Corp. reusable instruments are recommended to be sterilized by moist heat sterilization.
- Instruments with articulating surfaces must be tested for movement. A moist heat compatible, medical grade lubricant should be applied to all articulating joints prior to sterilization.
- After cleaning and prior to sterilization, the reusable instrument should be double-wrapped or packaged in blue Central Sterilization Room (CSR) sterilization wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique).

- Biological indicators and/or steam integrators should be used within the sterilizer load to help ensure the proper sterilization conditions of time, temperature, and saturated steam have been met.
- The parameters for PKA, TKA and Leg Positioner sterilization are:
 - Method: Moist-Heat Sterilization
 - Cycle: Pre-Vacuum (Pre-Vac)
 - Temperature: 270° F (132° C)
 - Exposure Time: 4 minutes (minimum)
 - Pressure: 2-15 PSIA
 - Dry-Time: 30 minutes (minimum, in chamber)
 - Cool-Time: appropriate time based on relative humidity and environment of at least 45 minutes at room temperature
- The parameters for THA sterilization are:
 - Method: Moist-Heat Sterilization
 - Cycle: Pre-Vacuum (Pre-Vac)
 - Temperature: 270° F - 273° F (132°C - 134° C)
 - Exposure Time: 4 minutes (minimum)
 - Pressure: 2-15 PSIA
 - Dry-Time: 45 minutes (minimum, in chamber)
 - Cool-Time: appropriate time based on relative humidity and environment of at least 45 minutes at room temperature

MAKO Surgical Corp. has validated the above recommended sterilization cycle based on AAMI/ANSI/ISO guidelines and recommendations. Other sterilization methods and cycles may also be suitable. However, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution.

- ETO sterilization and cold sterilization techniques are not recommended.
- Select MAKO Surgical Corp. trays and instruments may be sterilized with the following parameters **outside of the United States**:
 - Method: Moist-heat sterilization
 - Cycle: Pre-vacuum
 - Temperature: 273° F- 278.6 °F (134°C - 137°C)
 - Exposure Time: 3 minutes (minimum)
 - Dry- time: 30 minutes (minimum, in chamber)
 - Cool-Time: appropriate time based on relative humidity and environment of at least 45 minutes at room temperature

MAKO Surgical Corp. has validated the above recommended OUS sterilization cycle based on ISO guidelines and recommendations. Please refer to the Mako Instrument Cleaning and Sterilization Guide for a list of instruments/trays that may be sterilized using this cycle. Additional validated cycle options are also available in the e-ifu available at ifu.stryker.com. The use of a washer-disinfector meeting the requirements of the ISO 15883 series is recommended for any and all devices that are reusable produced by Stryker.

In case of serious incident, please notify the Manufacturer, and Competent Authority in your region.

Transport & Storage Information

Store the devices in standard hospital environmental conditions unless specific requirements are defined and described on the product label. Instrument trays are not intended as sterilization trays to maintain the product sterility. As stated above, blue wrap and select rigid containers may be used to maintain sterility in transit post hospital sterilization.

Safe Disposal

If a device is being returned for evaluation, please contact your local Mako representative for shipping/handling information. If the device is not being returned to MAKO Surgical Corp., instruments are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state and local regulations.

This IFU is available electronically at ifu.stryker.com. Detailed manual and automatic cleaning instructions, including times and temperatures, are provided in the Mako Instrument Cleaning and Sterilization Guide at ifu.stryker.com. Disassembly and assembly instructions applicable to all MAKO Surgical Corp. reusable instruments are all also provided in this user guide.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. The products depicted are CE marked according to the Medical Device Directive 93/42/EEC and European Medical Device Regulation 2017/745. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Mako, Mako Surgical Corp., Stryker, and VIZADISC. All other trademarks are trademarks of their respective owners or holders.

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