Non-Sterile Instruments

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430 USA
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

©2020 Howmedica Osteonics Corp.


QIN 4383 Rev AA

Refer to product label for CE mark status and Legal Manufacturer.
The CE mark is only valid if also found on the product label.

Labeling Symbols Glossary

The following is a list of symbols that may be used on Stryker medical device labeling. Refer to individual product labels for applicable symbology for each product.

Symbol	Standard
Caution (See instructions for use for operator awareness)	 ISO 15223-1
Consult instructions for use or consult electronic instructions for use	 ISO 15223-1
Do not re-use	 ISO 15223-1
Sterilized using irradiation	 ISO 15223-1
Sterilized using hydrogen peroxide	 ----
Sterilized using ethylene oxide	 ISO 15223-1
Sterilized using Aseptic processing techniques (Aseptic fill)	 ISO 15223-1
Non-sterile	 ISO 15223-1
Do not resterilize	 ISO 15223-1
Single sterile barrier system	 ISO 15223-1
Double sterile barrier system	 ISO 15223-1
Use-by date	 ISO 15223-1

Symbol	Standard
Date of manufacture 	ISO 15223-1
Legal manufacturer 	ISO 15223-1
Authorized representative in the European Community 	ISO 15223-1
Catalogue number 	ISO 15223-1
Batch code 	ISO 15223-1
Serial number 	ISO 15223-1
MR Safe 	ASTM F2503
MR Conditional 	ASTM F2503
MR Unsafe 	ASTM F2503
Do not use if package is damaged 	ISO 15223-1
Medical device 	ISO 15223-1
Quantity 	----
Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician 	----

Symbol	Standard
Contains hazardous substances	 ISO 15223-1
Temperature limit	 ISO 15223-1
Keep dry	 ISO 15223-1
Keep away from sunlight	 ISO 15223-1
Patient identification	 ISO 15223-1
Date (of Implantation)	 ISO 15223-1
Health care center or doctor	 ISO 15223-1
Patient information website	 ISO 15223-1
Unique Device Identifier	 ISO 15223-1
Contact by Phone	 N/A

stryker

ifu.stryker.com

Adobe® Reader®

Free paper copy within 7 days/Version papier gratuite dans les 7 jours/
Kostenlose Papierversion in 7 Tagen/Copia en papel gratuita dentro de unos
7 días/Copia cartacea gratuita entro 7 giorni/Copia impressa gratuita dentro
de 7 dias/Gratis kopia på papper inom 7 dagar /Gratis papirkopi i løbet af
7 dage/Ilmainen kopio paperilla 7 päivän kuluessa/Binnen 7 dagen gratis op
papier afgedrukt exemplaar/Δωρεάν έντυπο αντίγραφο εντός 7 ημερών/
Gratis papirkopi innen 7 dager/紙のコピーを7日間以内無料提供/Bezplatna
kopia papierowa w ciągu 7 dni/7일 이내 무료 사본/可在 7 天内免费索取纸
质副本/Copie gratuite pe hârtie în termen de 7 zile/Nemokama popierinė
kopija per 7 dienas/Безплатно копие на хартиен носител до 7 дни/
Kopie na papieru zdarma během 7 dnů/7 gün içinde ücretsiz basılı kopya

**AR AT AU BE BG BR
CH CO CY CZ DE DK
ES FI FR GB HU IE IL
IT LU LV NL NO NZ
PL PT RU SE SI UY**

+800 135 79 135

**AD AE BH BM CN DO
DZ EG HK IN IQ JO JP
KR KW KZ LB LY MX
MY OM PR QA RS SA
SD TH TN TW TZ UA
ZA**

+31 20 794 7071

BR 0800 591 1055
CA + 855 805 8539
CL 800 914 248
EE 0800 0100567
GR 00800 161 2205 7799
HR 8000 804 804
IS 800 8996
LI +31 20 796 5692
LT 8800 30728
MT +31 20 796 5693
RO 0800 895 084
SG 800 101 3366
SK 0800 606 287
TR 00800 142 064 866
US + 855 236 0910
VN 122 80297

EN: Non-Sterile Instruments	1
BG: Нестерилни инструменти	6
CH: 非无菌器械	12
CS: Nesterilní nástroje	16
DA: Usterile instrumenter	21
DU: Niet-steriele instrumenten	26
FI: Steriloimattomat instrumentit	32
FR: Instruments non stériles	37
GE: Unsteriles Instrumentarium	43
GK: Μη αποστειρωμένα εργαλεία	49
IT: Strumenti non sterili	55
JA: 非滅菌器具	60
KO: 비멸균 기구	65
LT: Nesterilūs instrumentai	70
NO: Usterile instrumenter	75
PL: Narzędzia niesterylne	80
PT: Instrumentos não esterilizados	86
RO: Instrumente nesterile	92
SP: Instrumentos no estériles	97
SV: Icke-sterila instrument	103
TR: Steril Olmayan Aletler	108

English

NON-STERILE INSTRUMENTS

HOWMEDICA OSTEONICS Non-invasive Instrumentation consists of manual surgical instruments intended for use during hip, knee and osteosyntheses orthopedic surgeries. Examples include but are not limited to orthopaedic instrument cases, trays, modular instrument handles and hammers.

Utilization

- Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- Radiographic templates are available to help assist in the preoperative prediction of component size and style. Radiographic templates are non-patient contacting – please consult the template for the magnification factor.
- The Howmedica Osteonics Corp. Surgical Protocols provide additional procedural information as applicable.

Warnings

- Care must be taken not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing differing design parameters, varying tolerances, different materials and manufacturing specifications, HOWMEDICA OSTEONICS instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting implant.
- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials that are visible.
- Unless otherwise specified on the instrument, markings on instruments with a measuring function have a linear accuracy of ± 0.01 inches (± 0.25 mm) or $\pm 5^\circ$ for angular measurements.
- Stryker Orthopaedics instruments are not intended to be used in or present within the MR environment. Since they are not intended to enter the MR environment, testing for safety in the MR environment has not been undertaken, and thus any potential hazards from exposure to such an environment is unknown.

Adverse Effects

- Malfunctions of instruments, including, but not limited to, weld fracture, thread galling and material fatigue, have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement or trauma surgery. These complications include but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.

- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impaction of the component into the preparation.
- Infection. Ensure the instruments are properly cleaned and sterilized to minimize the risk of infection.

Infection

Transient Bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have been associated with transient bacteremia. To help prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.

Materials

Stryker Orthopaedics reusable instruments may contain one or more metals in pure and/or alloyed forms, including but not limited to the following metals: stainless steel, cobalt chromium, nickel, aluminum, and titanium. Stryker reusable instruments may also contain one or more of the following polymers including but not limited to: polypropylene, polyphenylsulfone (PPSU), polyetherimide, polyoxymethylene, silicone, polyamide, polytetrafluoroethylene (PTFE), polyethylene, Santoprene, and polycarbonate. Contact Stryker for instrument-specific information when treating a patient with a material allergy or sensitivity. As noted, these are common materials used in Stryker reusable instruments. This is not intended to be an exhaustive list of all materials used in these devices.

Examination Prior to Use

- The life of the instrument depends on the number of times they are used as well as the precautions taken in handling, cleaning and storage. Great care must be taken of the instruments to ensure that they remain in good working order.
- Instruments should be examined for wear or damage by physicians and staff in operating centers prior to surgery.
- The examination shall include a visual and functional inspection of the working surfaces, articulation points, rotating features, hinges, springs, connection mechanisms, mating parts, threads, and working ends of all instruments. Functional inspection should fully replicate the intended use of the device to confirm the instrument moves, assembles, and/or rotates as expected. It should also include verifying all welded connections, that all components are present, and the cleanliness of the orifices and cavities, as well as examination for signs of material degradation including but not limited to cracks, distortion/deformation, impact, corrosion, detached pieces or other unexpected changes. If one of the above-mentioned conditions occurs and impacts device functionality, the instrument has reached the end of its functional life and must be replaced. If damaged instrumentation is used, possible fracture, jamming, or other failure may occur. For instruments with moving parts, application of medical grade lubricants that are biocompatible per ISO 10993 may be necessary.
- Stryker Orthopaedics shall not be responsible in the event of the use of instruments that are damaged, incomplete, show signs of excessive wear and tear, or that have been repaired or sharpened outside the control of Stryker Joint Replacement. Any faulty instruments must be replaced prior to any surgical procedure.

Cleaning and Sterilization

- These reusable instruments are **not** sterile. All Stryker Orthopaedics reusable instruments must be cleaned and sterilized to prepare them for use. Please refer to Cleaning, Sterilization, Inspection and Maintenance of Reusable Medical Devices - LSTPI-B, available at ifu.stryker.com.
- Remove all instruments from their packaging prior to cleaning and sterilization. When removing the instrument from its initial packaging or from a tray, the structural integrity of the instrument should be inspected to ensure safe use.
- Where noted in LSTPI-B, certain instruments must be disassembled prior to cleaning and sterilization.
- All parts must be cleaned and sterilized before use. Cleaning of the reusable instrument prior to sterilization should include: presoaking, manual cleaning, ultrasonication (e.g. for complex parts), and washing/disinfecting. Enzymatic cleaners, manual detergents, and neutral cleaners should be used.
- During and after cleaning, inspect hard-to-reach areas such as lumens within the instrument to help ensure that entrapped debris and/or soil are completely removed. If debris and/or soil are evident, then reclean the instrument until there is no evidence of debris and/or soil.
- Rinsing should be performed after each cleaning step to help ensure cleaners and/or debris/soil are removed.
- Howmedica Osteonics' reusable instruments are recommended to be sterilized by moist heat sterilization.
- Instruments with articulating surfaces must be tested for movement. A moist heat compatible, medical grade lubricant should be applied to all articulating joints prior to sterilization.
- After cleaning and prior to sterilization, the reusable instrument should be double-wrapped or packaged in CSR sterilization wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique).
- Select Stryker Orthopaedics trays and instruments are compatible with Aesculap SterilContainer System (JN441 or JN442 perforated bottom with JK48X series of aluminum 2000 lids). Please refer to LSTPI-B for a list of compatible parts, and detailed instructions for use.
- Biological indicators and/or steam integrators should be used within the sterilizer load to help ensure the proper sterilization conditions of time, temperature, and saturated steam have been met.
- The parameters for sterilization are:

Method:	Moist-Heat Sterilization
Cycle:	Pre-Vacuum (Pre-Vac)
Temperature:	270° F (132° C)
Exposure Time:	4 minutes (minimum)
Pressure:	2-15 PSIA
Dry-Time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Howmedica Osteonics Corp. has validated the above recommended sterilization cycle based on AAMI/ANSI/ISO guidelines and recommendations. Other sterilization methods and cycles may also be suitable. However, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution.

- EtO sterilization and cold sterilization techniques are not recommended.

Select Stryker Orthopaedics trays and instruments may be sterilized with the following parameters **outside of the United States:**

Method:	Moist-heat sterilization
Cycle:	Pre-vacuum
Temperature:	273° F - 279 °F (134°C - 137°C)
Exposure Time:	3 minutes (minimum)
Dry-time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Howmedica Osteonics Corp. has validated the above recommended OUS sterilization cycle based on ISO guidelines and recommendations. Please refer to e-IFU LSTPI-B for a list of instruments/trays that may be sterilized using this cycle. Additional validated cycle options are also available in this e-IFU available at ifu.stryker.com.

In case of serious incident, please notify the Manufacturer, and Competent Authority in your region.

Transport & Storage Information

Store the devices in standard hospital environmental conditions unless specific requirements are defined and described on the product label. Instrument trays are not intended as sterilization trays to maintain the product sterility. CSR sterilization wraps and select rigid containers may be used to maintain sterility in transit post hospital sterilization.

Safe Disposal

If a device is being returned for evaluation, please contact your local Stryker representative for shipping/handling information. If the device is not being returned to Stryker, implant components are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state and local regulations.

This IFU is available electronically at ifu.stryker.com. Detailed manual and automatic cleaning instructions, including times and temperatures, are provided in Stryker Orthopaedics Literature LSTPI-B available at ifu.stryker.com. Disassembly and assembly instructions applicable to all Stryker Orthopaedics re-usable instruments are all also provided in this brochure.

Refer to Aesculap Instructions for Use for Care & Handling of Aesculap SterilContainer Systems (JN441 or JN442).

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Howmedica, Osteonics, Stryker, Stryker Orthopaedics. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ►
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◀ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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Howmedica Osteonics Corp. ürün etiketlerinde kullanılan kısaltmaların bir listesi aşağıdaki tabloda verilmiştir: