stryker[®]

STRYKER SYSTEM G™


Power Module

REF 7309-120-000

Instructions For Use

R_x ONLY

ENGLISH (EN)
ESPAÑOL (ES)
ITALIANO (IT)
PORTUGUÊS (PT)
TÜRKÇE (TR)
РУССКИЙ (RU)
中文 (ZH)
한국어 (KO)

Contents

Introduction	2
Audience	2
Conventions	2
Contact Information	3
Indications For Use	3
Contraindications	3
For Use With	3
Safety Directives	4
Features	5
Definitions	5
Care Instructions	6
Instructions	6
Storage and Handling	7
Disposal/Recycle	8
Product Safety Compliance	10
Specifications	11

Introduction

This instructions for use manual contains information intended to ensure the safe, effective, and compliant use of your product.

Audience

This manual is intended for in-service trainers, physicians, nurses, and surgical technologists. Keep and consult this reference manual during the life of the product.

Conventions

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

Contact Information

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service at 1-269-323-7700 or 1-800-253-3210. Outside the US, contact your nearest Stryker subsidiary.

NOTE: The user and/or patient should report any serious product-related incident to both the manufacturer and the Competent Authority of the European Member State where the user and/or patient is established.

Indications For Use

The STRYKER SYSTEM G Power Module is intended to power the STRYKER SYSTEM G handpieces. The power module is a non-sterile, rechargeable unit, which is intended to be charged using the STRYKER SYSTEM G Charger.

Contraindications

None known.

For Use With

This section identifies components intended to be used with the equipment to obtain a safe combination.

Description	REF
STRYKER SYSTEM G Rotary Drill	7305-001-000
STRYKER SYSTEM G Reciprocating Saw	7306-001-000
STRYKER SYSTEM G Sagittal Saw	7308-001-000
STRYKER SYSTEM G Transfer Shield	7304-120-000
STRYKER SYSTEM G Charger, 2-Bay	7310-120-000

NOTE: For a complete list of components or ordering information, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.


Safety Directives

WARNINGS:

- Before using this equipment, read and understand the instructions for use. Pay particular attention to safety information.
- Only trained and experienced healthcare professionals should use this equipment.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker does not recommend surgical procedure or technique.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage. DO NOT use any equipment if damage is apparent. See the care instructions manual supplied with the handpiece for inspection criteria.
- Before use, ALWAYS allow the equipment to reach the operation temperature range specified in the instructions for use supplied with the handpiece.
- To prevent the equipment from overheating, ALWAYS follow the recommended duty cycle specified in the instructions for use supplied with the handpiece.
- DO NOT allow the battery contacts to touch metal objects.

Features

A	Handle – Used to insert and remove the power module from the handpiece.
B	Battery Contacts – Interlock with the battery charger to complete an electric circuit.

Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

Symbol	Definition
	General Warning Sign
	Caution
DO NOT STERILIZE	Do not sterilize
	Consult instructions for use
	Universal Recycling Symbol

Symbol	Definition
1 V	1 Volt
1 A	1 Ampere
1 VA	1 Volt Ampere
1 Hz	1 Hertz
AC	Alternating Current
REF	Catalog Number

Care Instructions

WARNINGS:

- Upon initial receipt and before each use, clean the equipment as indicated.
- DO NOT sterilize the power module.

For processing instructions, see the care instructions manual supplied with the handpiece.

Instructions

WARNINGS:

- The power module is not supplied with a full charge. Always charge the power module before first use.
- ALWAYS use a STRYKER SYSTEM G Charger to charge power modules.

CAUTIONS:

- ALWAYS install the power module onto the charger to activate it before first use. A new power module will not function until it has been installed onto the charger.
- DO NOT store the power module installed in a handpiece. Discharge will occur even though the handpiece is not operating.

Fully charge the power module before first and every use.

NOTES:

- For charging instructions, see the instructions for use supplied with the battery charger.
- For maximum capacity, store the power module on the battery charger to maintain a full charge.
- Stryker recommends using power modules that have been charged within the last 5 days.
- To install and remove the power module from a handpiece, see the instructions for use supplied with the handpiece.
- Typical power module run time is 5 minutes or longer, depending on surgical procedure and technique.

Storage and Handling

CAUTION: ALWAYS transport and store the equipment within the specified environmental condition values:

Temperature Limitation	Humidity Limitation	Atmospheric Pressure Limitation
 -20 °C to 40 °C	 10 % to 75 %	 50 kPa to 106 kPa

Disposal/Recycle

WARNINGS:

- Power modules contain Lithium-ion material. Risk of fire, burns, and explosion. DO NOT short circuit, open, crush or incinerate. Consult instructions for use.
- Follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.

NOTE: Contact Stryker for the recycling passport for this product.

The Batteries Directive 2006/66/EC introduces new requirements from September 2008 on removability of batteries from waste equipment in EU Member States. To comply with this Directive, this device has been designed for safe removal of the batteries at end-of-life by a waste treatment facility. Infected units should be decontaminated before they are sent for recycling. In the case that it is not possible to decontaminate the unit for recycling, the hospital should not attempt to remove the batteries from waste equipment. Continued disposal of small amounts of portable batteries to landfill and incineration is allowed under the Batteries Directive 2006/66/EC and Member State regulations.

Li-ION

Power modules contain Lithium-ion material.

China RoHS Standard SJ/T 11364

To comply with China RoHS Standard SJ/T 11364, this device has been marked with the environmentally friendly use period (EFUP) number, measured in years. The device contains at least one of the listed hazardous substances above threshold.

China RoHS Disclosure Report

REF 7309-120-000

Part Name	Hazardous Substances					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium (Cr [VI])	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)
Electronic Circuit Boards	X	O	O	O	O	O

This table is prepared in accordance with the provisions of SJ/T11364.

O: Indicates that said hazardous substance contained in all of the homogenous materials for this part is below the limit requirement of GB/T26572.

X: Indicates that said hazardous substance contained in at least one of the homogenous materials used for this part is above the limit requirement of GB/T26572.

Material: In accordance with the European REACH regulation and other environmental regulatory requirements, components within the power module contain Lead Cas. No. 7439-92-1. This declaration is made in good faith and is either based on a technical evaluation, supplier data and/or laboratory testing.

Product Safety Compliance

International Electrotechnical Commission (IEC)

IEC 60601-1 A1: R:2012 (3.1 Edition): *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); Amendment 1 (2012)

IEC 60601-1-6: 2013: *USABILITY – Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability* – Edition 3.1 Consolidated Reprint.

IEC 62366: 2014: *USABILITY – Medical devices – Application of usability engineering to medical devices* - Edition 1.1 Consolidated Reprint.

IEC 60601-1-2: 2014: *EMC (Version 4.0) – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests* – Edition 4.0

IEC 62133-2: 2017: *Secondary Cells and Batteries Containing Alkaline or other Non-Acid Electrolytes – Safety Requirements for Portable Sealed Secondary Cells and for Batteries Made from Them for use in Portable Applications*

European Committee for Electrotechnical Standardization (CENELEC)

EN 60601-1:2006/A1:2013 (3.1 Edition): *Medical electrical equipment Part 1: General requirements for basic safety and essential performance* – Incorporates Amendment A12: 2014

Specifications

Model:	STRYKER SYSTEM G Power Module
Catalog Number:	7309-120-000
Dimensions:	Length: 98 mm Width: 69 mm Height: 153 mm
Mass:	0.425 kg
Maximum Speed:	21000 revolutions per minute (nominal)
Electrical Input:	Internally Powered, 12.6 V === (Direct current) maximum
Capacity:	21.6 Wh
Lithium-ion Code:	31NR 19/65
Ingress Protection (IP):	IPX0
European Conformity:	 0197

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IEC 60601-1-6: 2013: *USABILITY – Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability – Edition 3.1 Consolidated Reprint.*

IEC 62366: 2014: *USABILITY – Medical devices – Application of usability engineering to medical devices - Edition 1.1 Consolidated Reprint.*

IEC 60601-1-2: 2014: *EMC (Version 4.0) – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests – Edition 4.0*

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EN 60601-1:2006/A1:2013 (3.1 Edition): *Medical electrical equipment Part 1: General requirements for basic safety and essential performance – Incorporates Amendment A12: 2014*

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