stryker®


STRYKER SYSTEM G™

Charger, 2-Bay

REF 7310-120-000

Instructions For Use

R_x ONLY

ENGLISH (EN)
ESPAÑOL (ES)
ITALIANO (IT)
PORTUGUÊS (PT)
POLSKI (PL)
ΕΛΛΗΝΙΚΑ (EL)
TÜRKÇE (TR)
РУССКИЙ (RU)
中文 (ZH)
한국어 (KO)

Contents

Introduction	3
Audience	3
Conventions	3
Contact Information	3
Indications For Use	4
Contraindications	4
Safety Directives	4
For Use With	5
Features	6
Charger	6
Information Screen	7
Indicator Lights	8
Definitions	8
Instructions	9
To Connect the Charger	10
To Charge a Power Module	11
To Reset a Circuit Breaker	15
To Disconnect the Charger	15
Cleaning	16
Inspection	16
Storage and Handling	17
Service	17
Disposal/Recycle	18
Specifications	19
Product Safety Compliance	20
Electromagnetic Compatibility	21

Introduction

This instructions for use manual contains information intended to ensure the safe, effective, and compliant use of your product.

Audience

This manual is intended for in-service trainers, physicians, nurses, and surgical technologists. Keep and consult this reference manual during the life of the product.

Conventions

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

Contact Information

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service at 1-269-323-7700 or 1-800-253-3210. Outside the US, contact your nearest Stryker subsidiary.

NOTE: The user and/or patient should report any serious product-related incident to both the manufacturer and the Competent Authority of the European Member State where the user and/or patient is established.

Indications For Use

The STRYKER SYSTEM G Charger, 2-Bay is a two-station, modular battery charger intended to charge STRYKER SYSTEM G handpiece power modules.

Contraindications

None known.

Safety Directives

WARNINGS:

- Before using this equipment, read and understand the instructions for use. Pay particular attention to safety information.
- Only trained and experienced healthcare professionals should use this equipment.

For Use With

This section identifies components intended to be used with the equipment to obtain a safe combination.

WARNINGS:

- Use only Stryker-approved equipment, unless otherwise specified.
- Use of unapproved electrical equipment may result in increased electromagnetic emissions or decreased electromagnetic immunity of the system.


NOTES:

- Components are sold separately.
- For a complete list of components or ordering information, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

Description	REF
STRYKER SYSTEM G Power Module	7309-120-000
Power Cord, Plug Type E/F	0039-308-000
Power Cord, Plug Type G	0039-309-000
Power Cord, Plug Type B	0039-310-000
Power Cord, Plug Type N	0039-311-000
Power Cord, Plug Type I	0039-312-000
Power Cord, Plug Type D	0039-313-000
Power Cord, Plug Type M	0039-314-000


Features

Charger

A	Information Screen – Provides status information for the power module.
B	Indicator Lights – Provide status information for the power module.
C	Charger Contacts – Interlock with the power module to complete an electric circuit.
D	Power Cord Receptacle – Allows connection of the power cord.
E	Circuit Breakers – Allow recovery from momentary electrical faults.
F	Power Cord – Connects the charger to a facility power receptacle.

Information Screen

A	Counter – Indicates the number of charge cycles or the charge level for the power module.
B	Cycles – Indicates that the counter is displaying the number of charge cycles for the power module.
C	Percentage – Indicates that the counter is displaying the charge level for the power module.
D	Charging – Indicates that the power module is charging.
E	<p>Module Health – For power modules with 10 or more charge cycles, indicates the health of the power module after charging:</p> <p> Excellent – The power module is at or near the beginning of its useful life.</p> <p> Good – The power module is in the middle of its useful life.</p> <p> Poor – The power module is at or near the end of its useful life. Replacement is recommended.</p>
F	Replace – Indicates that the power module is at or near the end of its useful life. Replacement is recommended.

Indicator Lights

Light	Definition
 Charging	Charging (Amber) – Indicates that the power module is charging.
 Ready	Ready (Green) – <ul style="list-style-type: none"> ▪ When no power modules are installed, indicates that the charger is connected to facility power and ready for use. ▪ When a power module is installed, indicates that the power module is fully charged and ready for use.
 Replace	Replace (Amber) – Indicates that the power module is at or near the end of its useful life. Replacement is recommended.

Definitions




The symbols located on the equipment and/or labeling are defined in this section and/or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

Symbol	Definition
	General Warning Sign
	Caution
DO NOT STERILIZE	Do not sterilize
	Consult instructions for use
	Universal Recycling Symbol
1 V	1 Volt
1 A	1 Ampere
1 VA	1 Volt Ampere
1 Hz	1 Hertz
AC	Alternating Current
REF	Catalog Number

Instructions

WARNINGS:

- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen, or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using this equipment. Place the equipment into service according to the EMC information contained in this manual.
- Portable and mobile radio frequency (RF) communications equipment can affect the function of this equipment.
- This equipment should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, ALWAYS verify normal operation in the configuration in which it will be used.
- To avoid the risk of electric shock, ALWAYS connect this equipment to a hospital-grade, facility power receptacle with protective earth.
- DO NOT operate the charger using a voltage inconsistent with the rating on the back of the unit.
- DO NOT operate the charger with a damaged power cord or plug.
- DO NOT modify the power cord.
- ALWAYS position the equipment so that the power cord may be easily disconnected as required.
- DO NOT place the charger within the sterile field.
- ALWAYS operate the equipment within the specified environmental condition values:

Temperature Limitation	Humidity Limitation	Atmospheric Pressure Limitation
		

CAUTION: ALWAYS place the power cord where it will not be stepped on, tripped over, or otherwise subjected to damage or stress.

To Connect the Charger




1. Before operating the equipment, verify the following:

- You are using Stryker-approved equipment.
- The equipment is not damaged.
- The equipment is clean.

NOTE: If you experience difficulty during operation, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

2. Use the power cord to connect the charger to a hospital-grade, facility power receptacle. Upon connection, the charger will indicate that it is ready for use:


- The ready lights will turn on.
- “000” will display on both screens.

To Charge a Power Module

NOTES:

- Fully charge power modules before first and every use.
 - The charging sequence may take between 5 and 45 minutes, based on the existing charge level in the power module.
 - For maximum capacity, store power modules on the battery charger to maintain a full charge.
1. Install a clean, dry power module onto the charger. Make sure the power module is securely snapped into place on the charger contacts.

2. Upon installation, the charger will communicate with the power module:


- The charging light will turn on.
- “000” will flash.

- Charging
- Ready
- Replace

3. After 3 seconds, the charging sequence will begin:

- The charging light will flash.
- The charging symbol will flash.
- The charge level for the power module will display.

- Charging
- Ready
- Replace


NOTE: If the power module is damaged or at the end of its useful life, the charger will not charge the power module:

- The replace light will turn on.
- The replace symbol will display.
- The number of charge cycles for the power module will display.


- Charging
- Ready
- Replace

4. After 60 seconds, the charge level will alternate with the number of charge cycles for the power module.


NOTE: If the power module is in poor health, the charging sequence will include additional information:

- The replace light will turn on.
- The replace symbol will display.


5. When the power module is fully charged, the charging sequence will end:

- 100% charge level will display.
- The charge level will continue to alternate with the number of charge cycles for the power module.
- For power modules with 10 or more charge cycles, the module health symbol will indicate the health of the power module. See the *Features* section.

NOTE: If the power module is in poor health, the end of the charging sequence will include additional information:


- The replace light will remain on.
- The replace symbol will continue to display.

To Reset a Circuit Breaker

⚠ WARNING: DO NOT use the charger if either of the circuit breakers continues to trip and turn off power. Disconnect the power cord from the facility power receptacle and contact Stryker.


In the event of a momentary electrical fault, a circuit breaker may trip and turn off power. Push the toggle switch up to reset the circuit breaker and continue operation.

To Disconnect the Charger

CAUTION: To reduce the risk of damage to the power cord, ALWAYS pull the power cord plug, not the power cord, when disconnecting the charger.

Disconnect the power cord from the charger and the facility power receptacle.

Cleaning

 **WARNING:** Before cleaning, ALWAYS disconnect the power cord from the facility power receptacle and the charger to reduce the risk of electric shock.

CAUTIONS:

- DO NOT immerse the charger in liquid.
 - DO NOT use solvents, lubricants, or other chemicals to clean the charger or charger modules unless otherwise directed.
 - DO NOT allow water to collect on top of the charger.
 - DO NOT allow water to enter the cooling vents on the back of the charger.
 - DO NOT sterilize the charger.
1. Disconnect the power cord from the facility power receptacle and the charger.
 2. Gently wipe the external surfaces of the charger and power cord with a soft, lint-free cloth moistened with a prepared cleaning solution or non-abrasive, hospital disinfectant and water.
 3. Immediately dry the charger with a soft, lint-free cloth or compressed air, < 140 kPa.
 4. Inspect the charger and power cord for damage. See the *Inspection* section.

Inspection

 **WARNINGS:**



- Perform recommended inspection as indicated in these instructions.
- DO NOT use any equipment if damage is apparent or the inspection criteria are not met.
- Before inspection, ALWAYS disconnect the power cord from the charger to reduce the risk of electric shock.

Interval	Activity	Criteria
Upon initial receipt and before each use	Inspect the equipment	No damage or signs of wear
		No loose or missing components
		No cuts in the power cord
		No bent contacts
		No cracks in the housing

NOTE: If the equipment fails to meet the inspection criteria, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Storage and Handling

CAUTION: ALWAYS transport and store the equipment within the specified environmental condition values:

Temperature Limitation	Humidity Limitation	Atmospheric Pressure Limitation
		

Service

⚠ WARNING: DO NOT disassemble, modify, service, or repair this equipment without authorization of the manufacturer.

NOTES:

- The useful life of this equipment is dependent upon many factors including, but not limited to, the method and duration of each use, and the handling of the equipment between uses.
- Routine and careful inspection of the equipment is the best method for determining the serviceable life span of the equipment.
- Maintenance documentation for this equipment is available upon request to Stryker-authorized service personnel only.
- For service information, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

Disposal/Recycle

WARNING: ALWAYS follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.

NOTE: Contact Stryker for the recycling passport for this product.

China RoHS Standard SJ/T 11364

To comply with China RoHS Standard SJ/T 11364, this device has been marked with the environmentally-friendly use period (EFUP) number, measured in years. The device contains at least one of the listed hazardous substances above threshold.

China RoHS Disclosure Report

REF 7310-120-000

Part Name	Hazardous Substances					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium (Cr [VI])	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)
Electronic Circuit Boards	X	O	X	O	O	O

This table is prepared in accordance with the provisions of SJ/T11364.

O: Indicates that said hazardous substance contained in all of the homogenous materials for this part is below the limit requirement of GB/T26572.

X: Indicates that said hazardous substance contained in at least one of the homogenous materials used for this part is above the limit requirement of GB/T26572.

Material: In accordance with the European REACH regulation and other environmental regulatory requirements, components within the charger contain Lead Cas. No. 7439-92-1, Cadmium Cas. No. 7440-43-9, N,N-dimethylformamide Cas. No. 68-12-2, 2-(2H-Benzotriazol-2-yl)-4, 6-ditertpentylphenol Cas. No. 25973-55-1. This declaration is made in good faith and is either based on a technical evaluation, supplier data and/or laboratory testing.

Specifications

Model:	STRYKER SYSTEM G Charger, 2-Bay
Catalog Number:	7310-120-000
Dimensions:	
Height:	130 mm
Width:	245 mm
Length:	394 mm
Mass:	1.6 kg
Mode of Operation:	Continuous Operation
Electrical:	
Input:	100-240 V AC / 135-175 VA (50/60 Hz)
Output:	12.6 V --- (Direct current) / 3 A maximum
Circuit Breakers:	1 A maximum
Ingress Protection (IP):	IPX0
Equipment Type:	Class I
Ground Type:	 Protective Earth Ground
Means of Isolation from Supply Mains:	Disconnect the power cord from the charger and the facility power receptacle.
European Conformity:	

Product Safety Compliance

International Electrotechnical Commission (IEC)

IEC 60601-1 A1: R:2012 (3.1 Edition): *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); Amendment 1 (2012)

IEC 60601-1-6: 2013: *USABILITY – Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability* – Edition 3.1 Consolidated Reprint.

IEC 62366: 2014: *USABILITY – Medical devices – Application of usability engineering to medical devices* - Edition 1.1 Consolidated Reprint.

IEC 60601-1-2: 2014: *EMC (Version 4.0) – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests* – Edition 4.0

European Committee for Electrotechnical Standardization (CENELEC)

EN 60601-1:2006/A1:2013 (3.1 Edition): *Medical electrical equipment Part 1: General requirements for basic safety and essential performance* – Incorporates Amendment A12: 2014

Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions		
<p>The STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) should assure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	<p>The STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:</p> <p> WARNING: This equipment/system is intended for use in a professional healthcare facility environment. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) or shielding the location.</p>
Harmonic emissions IEC 61000-3-2	Class A	

Guidance and manufacturer's declaration - electromagnetic immunity			
NOTE: The STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 2 kV, ± 4 kV, ± 6 kV, ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 2 kV, ± 4 kV, ± 6 kV, ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 20%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV at 100 kHz repetition frequency for power supply lines ± 1 kV at 100 kHz repetition frequency for input/output lines	± 2 kV at 100 kHz repetition frequency for power supply lines ± 1 kV at 100 kHz repetition frequency for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 , ± 1 kV line(s) to line(s) ± 0.5 , ± 1 , ± 2 kV line(s) to earth	± 0.5 , ± 1 kV line(s) to line(s) ± 0.5 , ± 1 , ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $0\% U_T$ (100% dip in U_T) for 0.5 cycle at 0° , 45° , 90° , 135° , 180° , 225° , 270° , and 315° $0\% U_T$ (100% dip in U_T) for 1 cycle at 0° $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 & 30 cycles at 0° $<5\% U_T$ ($>95\%$ dip in U_T) for 5 s $0\% U_T$ (100% dip in U_T) for 5 s	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $0\% U_T$ (100% dip in U_T) for 0.5 cycle at 0° , 45° , 90° , 135° , 180° , 225° , 270° , and 315° $0\% U_T$ (100% dip in U_T) for 1 cycle at 0° $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 & 30 cycles at 0° $<5\% U_T$ ($>95\%$ dip in U_T) for 5 s $0\% U_T$ (100% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) requires continued operation during power mains interruptions, it is recommended that the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m at 50 and 60 Hz	30 A/m at 50 and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the alternating current mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity			
<p>NOTE: The STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) should assure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	<p>3 Vrms 150 kHz to 80 MHz outside ISM bands 80% AM at 1 kHz</p> <p>6 Vrms 150 kHz to 80 MHz in ISM bands 80% AM at 1 kHz</p>	<p>3 Vrms 150 kHz to 80 MHz outside ISM bands 80% AM at 1 kHz</p> <p>6 Vrms 150 kHz to 80 MHz in ISM bands 80% AM at 1 kHz</p>	<p>IEC 60601-1-2 3rd Edition:</p> <p>Portable and mobile RF equipment should be used no closer to any part of the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1.2\sqrt{P}$ <p>$d=1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>IEC 60601-1-2 4th Edition:</p>
Radiated RF IEC 61000-4-3	<p>27 V/m 385 MHz, PM at 18 Hz, Max. Power = 1.8W</p> <p>28 V/m 450 MHz, FM \pm 5 kHz deviation, 1 kHz sine, Max. Power = 2W</p> <p>9 V/m 710, 745, 780, 5240, 5500, 5785 MHz, PM 217 Hz, Max. Power = 0.2W</p> <p>28 V/m 810, 870, 930 MHz, PM 18 Hz, Max. Power = 2W</p> <p>28 V/m 1720, 1845, 1970, 2450 MHz, PM 217 Hz, Max. Power = 2W</p>	<p>27 V/m 385 MHz, PM at 18 Hz, Max. Power = 1.8W</p> <p>28 V/m 450 MHz, FM \pm 5 kHz deviation, 1 kHz sine, Max. Power = 2W</p> <p>9 V/m 710, 745, 780, 5240, 5500, 5785 MHz, PM 217 Hz, Max. Power = 0.2W</p> <p>28 V/m 810, 870, 930 MHz, PM 18 Hz, Max. Power = 2W</p> <p>28 V/m 1720, 1845, 1970, 2450 MHz, PM 217 Hz, Max. Power = 2W</p>	<p> WARNING: Portable RF equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> <p> (Non-ionizing electromagnetic radiation)</p>

NOTE 1: At 80 MHz and 800MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) is used exceeds the applicable RF compliance level above, the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000).

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

IEC 60601-1-2 3rd Edition: Recommended separation distances between portable and mobile RF equipment and the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000)			
The STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF equipment (transmitters) and the STRYKER SYSTEM G Charger, 2-Bay (REF 7310-120-000) as recommended below, according to the maximum output power of the equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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