QIN 4332 Rev. AC

**stryker**  
**Howmedica**  
**OSTEONICS**

## Sterile Instrumentation

Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ 07430, USA  
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

**CE 2797**

©2021/08 Howmedica Osteonics Corp.

QIN 4332 Rev. AC

Refer to product label for CE mark status and Legal Manufacturer.  
The CE mark is only valid if also found on the product label.

## Labeling Symbols

Attention, See Instructions for Use

Do not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

## CONTENTS

|                  |    |
|------------------|----|
| English .....    | 1  |
| Français .....   | 6  |
| Deutsch.....     | 12 |
| Italiano.....    | 18 |
| Español.....     | 24 |
| Português .....  | 30 |
| Svenska .....    | 36 |
| Suomi.....       | 41 |
| Dansk.....       | 46 |
| Nederlands ..... | 51 |
| Norsk.....       | 57 |
| Polski.....      | 62 |
| Ελληνικά.....    | 68 |
| 简体中文.....        | 74 |
| 한글 .....         | 79 |
| 日本語.....         | 85 |

Howmedica Osteonics Corp.'s Sterile Instrumentation consists of manual surgical instruments intended for use during hip or knee arthroplasty, or trauma surgeries.

### Utilization

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- Howmedica Osteonics Corp.'s trial components should be used for size determination, preparation evaluation, trial reduction and range of motion evaluation, thus preserving the integrity of the actual implants and their sterile packaging.
- Radiographic templates are available to assist in the preoperative prediction of component size and style.
- Howmedica Osteonics Corp.'s Surgical Protocols provide additional procedural information.

### Warnings

- Use caution when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing different design parameters, varying tolerances, different materials and manufacturing specifications, Howmedica Osteonics Corp.'s Sterile Instrumentation should not be used to implant any other manufacturer's components with the exception of the Pelvic Alignment Level device which can be used to implant other manufacturers' hip stems. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting implant.
- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials.

- Stryker Orthopaedics instruments are not intended to be used in or present within the MR environment. Since they are not intended to enter the MR environment, testing for safety in the MR environment has not been undertaken, and thus any potential hazards from exposure to such an environment is unknown.
- Patient post-operative pain. Inherent to all joint replacement is the risk that a patient will develop post-operative pain; pain is a commonly reported symptom regardless of the device implanted. The clinical literature reveals numerous potential causes of pain not directly related to the implant performance including, but not limited to, prior history of trauma and natural disease progression.
- For patients who present with pain following implantation of any orthopedic implant system, the physician should consider all potential causes of the symptoms identified in the clinical literature, including infection, soft tissue impingement, and possible adverse local tissue reactions associated with wear debris, metal ions or corrosion. Accurate diagnosis of the source of pain and directed, timely intervention is essential to ensuring effective treatment of pain.
- This is a single-use device and should never be reused. Reuse of a single use device may result in a myriad of risks including, but not limited to:
  1. Contaminants leading to infection.
  2. Material fragments, debris, corrosion byproducts or unintended foreign objects leading to inflammatory response.
  3. Biologic Contaminants (non-pathological) leading to inflammation.

Additionally, although the device may appear undamaged, previous use may have created nonvisible damage that could result in loss of device functionality such as:

1. Fractured device.
2. Assembly issues.

## **Adverse Effects**

- Malfunctions of instruments, including, but not limited to: weld fracture, thread galling and material fatigue, have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement or trauma surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impaction of the component into the preparation.

## **Infection**

Transient Bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have been associated with transient bacteremia. To help prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.

## **Sterilization**

- This instrument has been sterilized by gamma radiation.
- The packaging of all sterile products should be inspected for flaws in the sterile barrier before opening. In the presence of such a flaw, the product must be assumed nonsterile.
- Care should be taken to prevent contamination of the instrument. In the event of contamination, this product must be discarded.

- If the package is opened, but the product is not used, the instrument **must not** be resterilized and must be discarded or returned to the supplier.
- Sterility cannot be assured for resterilized devices, as cleaning and re-sterilization have not been verified.
- Device should not be used after the expiry date displayed on the label as packaging had not been validated beyond this date.

## **Transport & Storage Information**

The device is individually packed in protective packaging that is labelled according to its contents. Store and transport the device in the original protective packaging. Do not remove the device from the packaging until it is planned to be used. Store the device in standard hospital environmental conditions unless specific requirements are defined and described on the product label.

## **Safe Disposal**

If a device is being returned for evaluation, please contact your local Stryker representative for shipping/handling information. If the device is not being returned to Stryker, instruments are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state, and local regulations.

**CAUTION:** Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

**Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Howmedica, Osteonics, Stryker. All other trademarks are trademarks of their respective owners or holders.**

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

| Term           | Abbreviation | Term           | Abbreviation |
|----------------|--------------|----------------|--------------|
| Alpha Code     | ALPH CDE     | Neck           | NK           |
| Angle          | ANG          | Offset         | OFFST        |
| Degree         | DEG or °     | Outer Diameter | OD           |
| Diameter       | DIA          | Right          | RT ►         |
| Extra Deep     | XDP          | Screw Holes    | SCR HLS      |
| Extra Large    | XLGE         | Side           | SDE          |
| Extra Small    | XSM          | Size           | SZE          |
| Head           | HD           | Small          | SM           |
| Height         | HT           | Standard       | STD          |
| Inner Diameter | ID           | Taper          | TPR          |
| Insert         | INSR         | Thickness      | THKNS        |
| Large          | LGE          | Type           | TYP          |
| Left           | ◀ LFT        | With           | W/           |
| Length         | LNTH         | Without        | W/O          |
| Medium         | MED          |                |              |

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