Non-Sterile Instruments

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430, USA
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

CE 2797

© 2022/02 Howmedica Osteonics Corp.
QIN 4333 Rev. AE

Refer to product label for CE Mark Status and Legal Manufacturer.
The CE mark is only valid if also found on the product label.

Labeling Symbols

Attention, See Instructions for Use

Do Not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide


Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

Non-Sterile Instruments	1
Instruments non stériles	5
Nicht sterile Instrumente	10
Strumenti non sterili	15
Instrumentos no estériles.....	20
Instrumentos não esterilizados	25
Osterila instrument.....	30
Steriloimattomat instrumentit	34
Usterile instrumenter	38
Niet-steriele instrumenten	42
Ikke-sterile instrumenter.....	47
Oprządowanie niesterylne	51
Μη αποστειρωμένα εργαλεία	56
非消毒器械.....	61
비멸균 기구	66
非滅菌性手術器具.....	71

English

NON-STERILE INSTRUMENTS

HOWMEDICA OSTEONICS Instrumentation consists of manual surgical instruments intended for use during hip or knee arthroplasty, or trauma surgeries.

Utilization

- Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- HOWMEDICA OSTEONICS Trial Components should be used for size determination, preparation evaluation, trial reduction and range of motion evaluation, thus helping to preserve the integrity of the actual implants and their sterile packaging.
- Radiographic templates are available to help assist in the preoperative prediction of component size and style.
- The Howmedica Osteonics Corp. Surgical Protocols provide additional procedural information.

Warnings

- Use caution when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing differing design parameters, varying tolerances, different materials and manufacturing specifications, HOWMEDICA OSTEONICS instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting implant.
- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials that are visible.
- Unless otherwise specified on the instrument, markings on instruments with a measuring function have a linear accuracy of ± 0.01 inches (± 0.25 mm) or $\pm 5^\circ$ for angular measurements.
- Stryker Orthopaedics instruments are not intended to be used in or present within the MR environment. Since they are not intended to enter the MR environment, testing for safety in the MR environment has not been undertaken, and thus any potential hazards from exposure to such an environment is unknown.
- Patient post-operative pain. Inherent to all joint replacement is the risk that a patient will develop post-operative pain; pain is a commonly reported symptom regardless of the device implanted. The clinical literature reveals numerous potential causes of pain not directly related to the implant performance including, but not limited to, prior history of trauma and natural disease progression.
- For patients who present with pain following implantation of any orthopedic implant system, the physician should consider all potential causes of the symptoms identified in the clinical literature, including infection, soft tissue impingement, and possible adverse local tissue reactions associated with wear debris, metal ions or corrosion. Accurate diagnosis of the source of pain and directed, timely intervention is essential to ensuring effective treatment of pain.

Adverse Effects

- Malfunctions of instruments, including, but not limited to, weld fracture, thread galling and material fatigue, have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impaction of the component into the preparation.
- Infection. Ensure the instruments are properly cleaned and sterilized to minimize the risk of infection.

Infection

Transient Bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have been associated with transient bacteremia. To help prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.

Examination Prior to Use

- The life of the instruments depends on the number of times they are used as well as the precautions taken in handling, cleaning and storage. Great care must be taken of the instruments to ensure that they remain in good working order.
- Instruments should be examined for wear or damage by physicians and staff in operating centers prior to surgery.
- The examination shall include a visual and functional inspection of the working surfaces, articulation points, rotating features, hinges, springs, connection mechanisms, mating parts, threads, and working ends of all instruments. Functional inspection should fully replicate the intended use of the device to confirm the instrument moves, assembles, and/or rotates as expected. It should also include verifying all welded connections, that all components are present, and the cleanliness of the orifices and cavities, as well as examination for signs of material degradation including but not limited to cracks, distortion/deformation, impact, corrosion, detached pieces or other unexpected changes. If one of the above-mentioned conditions occurs and impacts device functionality, the instrument has reached the end of its functional life and must be replaced. If damaged instrumentation is used, possible fracture, jamming, or other failure may occur. For instruments with moving parts, application of medical grade lubricants that are biocompatible per ISO 10993 may be necessary.
- Stryker Orthopaedics shall not be responsible in the event of the use of instruments that are damaged, incomplete, show signs of excessive wear and tear, or that have been repaired or sharpened outside the control of Stryker. Any faulty instruments must be replaced prior to any surgical procedure. Please refer to e-IFU "Instructions for Inspection and Maintenance of Reusable Medical Devices" – SLI0001 for a comprehensive guide on instrument inspection and determination of an instrument's end of life. This e-IFU is available at ifu.stryker.com.

Cleaning and Sterilization

- These reusable instruments are **not** sterile. All Stryker Orthopaedics reusable instruments must be cleaned and sterilized to prepare them for use. Please refer to Cleaning, Sterilization, Inspection and Maintenance of Reusable Medical Devices - LSTPI-B, available at ifu.stryker.com.
- Remove all instruments from their packaging prior to cleaning and sterilization.
- All parts must be cleaned and sterilized before use. Cleaning of the reusable instrument prior to sterilization should include: presoaking, manual cleaning, ultrasonication (e.g. for complex parts), and washing/ disinfecting. Enzymatic cleaners, manual detergents, and neutral cleaners should be used.
- During and after cleaning, inspect hard-to-reach areas such as lumens within the instrument to help ensure that entrapped debris and/or soil are completely removed. If debris and/or soil are evident, then reclean the instrument until there is no evidence of debris and/or soil.
- Rinsing should be performed after each cleaning step to help ensure cleaners and/or debris/soil are removed.
- Howmedica Osteonics reusable instruments are recommended to be sterilized by moist heat sterilization.
- Instruments with articulating surfaces must be tested for movement. A moist heat compatible, medical grade lubricant should be applied to all articulating joints prior to sterilization.
- After cleaning and prior to sterilization, the reusable instrument should be double-wrapped or packaged in CSR sterilization wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique).
- Select Stryker Orthopaedics trays and instruments are compatible with Aesculap SterilContainer System (JN441 or JN442 perforated bottom with JK48X series of aluminum 2000 lids). Please refer to LSTPI-B for a list of compatible parts, and detailed instructions for use.

- Biological indicators and/or steam integrators should be used within the sterilizer load to help ensure the proper sterilization conditions of time, temperature, and saturated steam have been met.
- The parameters for sterilization are:

Method:	Moist-Heat Sterilization
Cycle:	Pre-Vacuum (Pre-Vac)
Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Pressure:	2-15 PSIA
Dry-Time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Howmedica Osteonics Corp. has validated the above recommended sterilization cycle based on AAMI/ANSI/ISO guidelines and recommendations. Other sterilization methods and cycles may also be suitable. However, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution.

- EtO sterilization and cold sterilization techniques are not recommended for these products.

Select Stryker Orthopaedics trays and instruments may be sterilized with the following parameters **outside of the United States**:

Method:	Moist-Heat Sterilization
Cycle:	Pre-Vacuum
Temperature:	273° F - 279° F (134° C - 137° C)
Exposure Time:	3 minutes (minimum)
Dry-time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Howmedica Osteonics Corp. has validated the above recommended OUS sterilization cycle based on ISO guidelines and recommendations. Please refer to LSTPI-B for a list of instruments/trays that may be sterilized using this cycle.

Detailed manual and automatic cleaning instructions, including times and temperatures, are provided in Stryker Orthopaedics e-IFU LSTPI-B available at ifu.stryker.com. Disassembly and assembly instructions applicable to all Stryker Orthopaedics re-usable instruments are all also provided in this e-IFU.

Refer to Aesculap Instructions for Use for Care & Handling of Aesculap SterilContainer Systems (JN441 or JN442).

Transport & Storage Information

Store the devices in standard hospital environmental conditions unless specific requirements are defined and described on the product label. Instrument trays are not intended as sterilization trays to maintain the product sterility. CSR sterilization wraps or pouches and select rigid containers may be used to maintain sterility in transit post hospital sterilization.

Safe Disposal

If a device is being returned for evaluation, please contact your local Stryker representative for shipping/handling information. If the device is not being returned to Stryker, instruments are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state, and local regulations.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Howmedica, Osteonics, Stryker. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ►
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◄ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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