The Stryker logo is displayed in white lowercase letters on a grey rectangular background. The word "stryker" is followed by a registered trademark symbol (®).

stryker[®]

**Endoscopic
Hardware
User Manual**

| | |
|------------------|-----|
| English | 3 |
| Français | 11 |
| Deutsch..... | 21 |
| Italiano | 31 |
| Português | 41 |
| Español..... | 51 |
| Nederlands | 61 |
| Dansk..... | 71 |
| Suomi..... | 79 |
| Norsk | 87 |
| Svenska | 95 |
| Polski | 103 |
| Ελληνικά | 113 |
| Türkçe..... | 123 |
| Русский | 131 |
| 简体中文 | 141 |
| 日本語 | 149 |
| 한국어 | 157 |

Before Using Your Instrument

Carefully unpack the instrument and check if any damage occurred during shipment.

Read this user manual thoroughly and be familiar with its contents prior to using this equipment.

Test this equipment prior to a surgical procedure. Ensure device has not been bent, and that any locking or attachment surfaces function properly. These devices were tested at the factory before shipment.

Attempt no internal repairs or alterations. See instructions for authorized service information.

Some hardware devices have device-specific locking mechanisms. There is no guarantee that instruments selected for minimum channel width alone will be compatible in combination.

1.1 Indications

Arthroscopic Hardware:

Arthroscopic Hardware consists of manual arthroscopic instruments intended for use in surgical procedures within a joint.

Scope Accessories:

Scope accessories are intended for use with a Stryker scope.

1.2 Contraindications

There are no known contraindications.

Warnings:

- Federal Law (USA) restricts this device to sale by or on the order of a physician only. The user must be a qualified physician having complete knowledge of the use of this equipment.
- Do not contact these instruments with any powered shaver blade or place instrument in the path of a surgical laser beam. Doing so can cause damage to the instrument which may lead to patient injury.
- This device is shipped non-sterile. This device must be cleaned and sterilized prior to the first use and after every subsequent use.
- Open all instruments with features that have an open/close position (such as jaws, stopcocks, lumens, etc.) prior to cleaning, disinfection, or sterilization. If they are not in the open position, incomplete sterilization may result.
- Use only the sterilization cycles outlined in this document. Using unspecified sterilization cycles may damage the device and/or result in incomplete sterilization.
- Failure to follow cleaning instructions may result in inadequate sterilization.

- The sterilization parameters presented in this document apply only to devices sterilized outside of a sterilization tray. When using a sterilization tray, consult the instructions provided with the tray for proper sterilization parameters.
- Wear appropriate protective equipment: gloves, eye protection, etc.

1.3 Reprocessing

These reprocessing instructions are provided in accordance with ISO 17664, AAMI TIR 12, AAMI ST79, and AAMI ST81. While they have been validated by Stryker as being capable of preparing the device for re-use, it remains the responsibility of the processor to ensure that the reprocessing, as actually performed (using equipment, materials, and personnel in the reprocessing facility), achieves the desired result. This normally requires validation and routine monitoring of the process. Stryker recommends users observe these standards when reprocessing medical devices.

1.3.1 Containment and Transportation

Reprocess the device as soon as reasonably practical following use.

1.3.2 Manual Cleaning

1. Disconnect the device from all other devices.
2. If applicable, disassemble the device.
3. **Wipe**
 - a. Wipe excess soil from the device with disposable paper towels.
4. **Soak in Enzymatic Detergent**
 - a. Prepare an enzymatic detergent solution¹ with lukewarm tap water, according to manufacturer's recommendations.
 - b. Wipe the entire surface of the device using a soft clean cloth dipped in the detergent solution.
 - c. Immerse the device in the detergent solution, ensuring the solution contacts all inner and outer surfaces of the device. Use a syringe to inject any lumens with 50 ml of solution.
 - d. Soak the device in the solution for 15 minutes.

5. **Brush**

Note: Use an appropriately sized brush.

- **For the interior of the device, use a bottle brush with bristles that extend the diameter of the device lumen.**
 - **For the exterior of the device, a 2 Spectra M16 Soft-Bristle brush is recommended.**
- a. Thoroughly brush the exterior of the device, focusing on mated and rough surfaces, with a soft-bristled brush.
 - b. Brush around all movable parts while moving them through all extreme opened and closed positions.
 - c. Inject all lumens, crevices, and mated surfaces five times with 50 ml of detergent solution.
 - d. With the device immersed in the detergent, thoroughly brush the interior and exterior of the device, using an appropriately sized bottle brush

for the interior and an appropriately sized soft-bristled brush for the exterior.

6. Rinse

- a. Remove the device from the detergent solution and rinse it with water² at ambient temperature until all visible detergent residue is removed.
- b. Flush all lumens, crevices, and mated surfaces five times with 50 ml of water.
- c. After all detergent residues are removed, rinse the device for an additional 30 seconds.
- d. Drain any excess water from the device by holding it at an incline.

7. Soak in Non-Enzymatic Detergent

- a. Prepare a non-enzymatic detergent solution³ in lukewarm tap water according to manufacturer's recommendations.
- b. Immerse the device in the detergent solution, ensuring the solution contacts all inner and outer surfaces of the device. Use a syringe to inject lumens with 50 ml of solution.
- c. Soak the device in the solution for 15 minutes.

8. Brush

Note: Use an appropriately sized brush.

- **For the interior of the device, use a bottle brush with bristles that extend the diameter of the device lumen.**
- **For the exterior of the device, a 2 Spectra M16 Soft-Bristle brush is recommended.**

- a. Thoroughly brush the exterior of the device, focusing on mated and rough surfaces, with a soft-bristled brush.
- b. Brush around all movable parts while moving them through all extreme opened and closed positions.
- c. Inject all lumens, crevices, and mated surfaces five times with 50 ml of detergent solution.
- d. With the device immersed in the detergent, thoroughly brush the interior and exterior of the device, using an appropriately sized bottle brush for the interior and an appropriately sized soft-bristled brush for the exterior.

9. Rinse

- a. Remove the device from the detergent solution and rinse it with water² at ambient temperature until all visible detergent residue is removed.
- b. Flush all lumens, crevices, and mated surfaces five times with 50 ml of water.
- c. After all detergent residues are removed, rinse the device for an additional 30 seconds.
- d. Drain any excess water from the device by holding it at an incline.

10. Dry

- a. Dry the device using a clean cloth. Filtered pressurized air⁴ can be used to assist drying.

11. Inspect

- a. Visually inspect the device, including all internal surfaces, for remaining soil. Use an endoscopic camera and endoscope if necessary to see the inner surface of the lumen.
- b. If visible soil remains, repeat the manual cleaning procedure, focusing on those areas.

- ¹ Cleaning was validated using Enzol at 1 oz/gal at 35-40°C.
² Cleaning was validated using reverse osmosis/deionized (RO/DI) water.
³ Recommend cleaning using Prolystica 2x Concentrate Neutral Detergent at 1/8 oz/gal at 35-40°C.
⁴ Drying was validated using pressurized air at 20 psi.

1.3.3 Automated Cleaning

1. Disconnect the device from all other devices.
2. If applicable, disassemble the device.
3. **Wipe**
 - a. Wipe excess soil from the device with disposable paper towels.
4. **Soak in Enzymatic Detergent**
 - a. Prepare an enzymatic detergent solution¹ with lukewarm tap water, according to manufacturer's recommendations.
 - b. Wipe the entire surface of the device using a soft clean cloth dipped in the detergent solution.
 - c. Immerse the device in the detergent solution, ensuring the solution contacts all inner and outer surfaces of the device. Use a syringe to inject any lumens with 50 ml of solution.
 - d. Soak the device in the solution for 15 minutes.
5. **Brush**

Note: Use an appropriately sized brush.

- **For the interior of the device, use a bottle brush with bristles that extend the diameter of the device lumen.**
- **For the exterior of the device, a 2 Spectra M16 Soft-Bristle brush is recommended.**

- a. Thoroughly brush the exterior of the device, focusing on mated and rough surfaces, with a soft-bristled brush.
 - b. Brush around all movable parts while moving them through all extreme opened and closed positions.
 - c. Inject all lumens, crevices, and mated surfaces five times with 50 ml of detergent solution.
 - d. With the device immersed in the detergent, thoroughly brush the interior and exterior of the device, using an appropriately sized bottle brush for the interior and an appropriately sized soft-bristled brush for the exterior.
6. **Rinse**
 - a. Remove the device from the detergent solution and rinse it with water² at ambient temperature until all visible detergent residue is removed.
 - b. Flush all lumens, crevices, and mated surfaces five times with 50 ml of water.
 - c. After all detergent residues are removed, rinse the device for an additional 30 seconds.
 - d. Drain any excess water from the device by holding it at an incline.
7. **Automated Wash**
 - a. Place the devices in an automated washer, on an incline to aid in drainage.
 - b. Program the washer with the following parameters, then activate the wash:

| | |
|--------------------|-----------|
| Pre-Wash | |
| Recirculation Time | 2 Minutes |

| | |
|-------------------|----------------|
| Water Temperature | Cold tap water |
| Detergent Type | NA |

| | |
|--------------------|----------------------------------|
| Enzyme Wash | |
| Recirculation Time | 2 Minutes |
| Water Temperature | Hot tap water |
| Detergent Type | Enzymatic detergent ¹ |

| | |
|--------------------|--------------------------------------|
| Wash 1 | |
| Recirculation Time | 2 Minutes |
| Water Temperature | 60°C (set point) |
| Detergent Type | Non-Enzymatic detergent ³ |

| | |
|--------------------|-------------|
| Rinse 1 | |
| Recirculation Time | 2 Minutes |
| Water Temperature | Heated 60°C |
| Detergent Type | NA |

| | |
|--------------------|-----------|
| Dry | |
| Recirculation Time | 7 Minutes |
| Temperature | 115°C |
| Detergent Type | NA |

8. Dry

- a. Dry the device using a clean cloth. Filtered pressurized air⁴ can be used to assist drying.

9. Inspect

- a. Visually inspect the device, including all internal surfaces, for remaining soil. Use an endoscopic camera and endoscope if necessary to see the inner surface of the lumen.
- b. If visible soil remains, repeat the manual cleaning procedure, focusing on those areas.

¹Cleaning was validated using Enzol at 1 oz/gal at 35-40°C.

²Cleaning was validated using reverse osmosis/deionized (RO/DI) water.

³Recommended cleaning using Prolystica 2x Concentrate Neutral Detergent at 1/8 oz/gal at 35-40°C.

⁴Drying was validated using pressurized air at 20 psi.

1.3.4 Disinfection (optional)

1. Disinfect the device in a disinfecting solution that has one of the following active ingredients:
 - ≥ 2.4% glutaraldehyde (with a minimum soak time of 45 minutes at 25°C).¹
2. Prepare the disinfecting solution according to the manufacturer's instructions.

3. Disassemble and fully immerse the device in the disinfecting solution, filling all lumens and crevices for the required time and at the appropriate temperature.
4. Rinse and flush the device in a bath of PURW. Immerse the device, shake it underwater, and allow it to soak for one minute. Repeat the rinse and flush two more times, using a fresh bath each time.
5. Dry all parts with a sterile, lint-free cloth immediately after rinsing.

¹Cidex® Activated Dialdehyde Solution is validated for disinfection efficacy.

1.3.5 Thermal Disinfection (optional)

Follow the steps below to disinfect the device using thermal disinfection.

1. Place the instruments into the washer disinfectant.
2. Program the washer using the following parameters:

| |
|------------------------------|
| Thermal Disinfection |
| Recirculation Time: 1 minute |
| Water Temperature: 90°C |

1.3.6 Autoclave (Steam) Sterilization

Warning: Drying time depends on several variables, including altitude, humidity, type of wrap, preconditioning, size of chamber, mass of load, material of load, and placement in chamber. Users must verify that drying time set in their autoclave yields dry surgical equipment.

1. Open all stopcocks and valves for each device.
2. If using a sterilization tray, place the device(s) in the tray and follow instructions provided with the tray. Use only trays that are approved for sterilization with steam.
3. Double wrap each device (or tray, if being used) prior to sterilization.
4. If sterilizing the devices without a tray, use the following parameters.

| | Pre-vacuum (U.S) | Pre-vacuum (EU) |
|-------------------------------|---------------------|--------------------|
| Wrapping | double* | double* |
| Temperature | 132°C (270°F) | 134°C (274°F) |
| Sterilization Time | 4 minutes | 3 minutes |
| Dry Time | 60 minutes | 60 minutes |

*Kimberly-Clark K600 KIMGUARD® sterilization wrap was used during sterilization validation.

Warning: Immediate-use steam sterilization is intended only for emergency situations. In the event that immediate-use reprocessing is required, the following instructions should be used.

| | "Immediate-Use" Pre-vacuum (U. S) | "Immediate-Use" Pre-vacuum (EU) |
|-------------------------------|-----------------------------------------|---------------------------------------|
| Wrapping | - | - |
| Temperature | 132°C (270°F) | 134°C (274°F) |
| Sterilization Time | 4 minutes | 3 minutes |
| Dry Time | - | - |

1.4 Maintenance and Storage

Hardware should be protected from bending and crushing. Visible bends and dents may prevent device from functioning properly. If there is any visible damage to the device, please return the unit to Stryker for repair.

Hardware should be properly cleaned and dried prior to storage. If there is any particulate matter on the device, carefully remove it prior to storage. Stainless steel components must be kept clean and dry when not in use to prevent rusting.

1.5 Expected Service Life

Inspect instruments and devices before each use. If a problem relating to functionality or appearance (such as but not limited to corrosion, discoloration, pitting or protrusions, binding, excessive wear) is observed or suspected, discontinue use immediately and contact your Stryker representative.

1.6 Authorized Service

To ensure the restoration of your instrument to original specifications, repairs should be directed to Stryker Endoscopy. Please contact your sales representative or distributor for product service outside of the United States. Please contact your local sales representative or distributor for specific warranty information outside of the United States.

Authorized Service

(United States Customers)

If service is needed either during or after the warranty period: Contact Stryker Endoscopy at (800) 624-4422 or contact your local Stryker Endoscopy sales representative.

Package all the components carefully in the original shipping container if possible.

Ship the components prepaid and insured to:

Stryker Endoscopy Customer Service

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks:
the Stryker Logo. All other trademarks are trademarks of their respective owners or holders.