stryker[®]

STRYKER SYSTEM G™

Rotary Drill


REF 7305-001-000

Rotary Drill Attachments

REF 7307-XXX-XXX Series
4103-XXX-XXX Series

Instructions For Use

R_x ONLY

ENGLISH (EN)
ESPAÑOL (ES)
ITALIANO (IT)
PORTUGUÊS (PT)
POLSKI (PL)
ΕΛΛΗΝΙΚΑ (EL)
TÜRKÇE (TR)
РУССКИЙ (RU)
中文 (ZH)
한국어 (KO)

Contents

Introduction	3	To Install a Power Module	8
Audience	3	Sterile Assistant	8
Conventions	3	Circulating Assistant	8
Contact Information	3	Sterile Assistant	9
Intellectual Property	3	To Install an Attachment	9
Indications For Use	3	To Install and Remove a Cutting or Fixation Tool	10
Contraindications	3	Wire and Pin Collets	10
Safety Directives	3	Keyed Chucks	11
For Use With	4	Keyless Chucks	11
Power and Care Components	4	Drill and Reamer Attachments	12
Attachments and Cutting and Fixation Tools	4	To Operate a Handpiece	12
Wire and Pin Collets	4	To Remove an Attachment	13
Keyed Chucks	4	To Remove a Power Module	13
Keyless Chucks	4	Specifications	14
Drill Attachments	4	Handpiece	14
Reamer Attachments	5	Attachments	14
Features	5	Wire and Pin Collets	14
Handpiece	5	Keyed Chucks	15
Safety Switch	6	Keyless Chucks	15
Attachments	6	Drill Attachments	15
Wire and Pin Collets	6	Reamer Attachments	16
Keyed Chucks	6	Transfer Shield	16
Keyless Chucks	6	Product Safety Compliance	16
Drill and Reamer Attachments	7	Electromagnetic Compatibility	17
Definitions	7		
Instructions	7		
To Process the Equipment	7		
To Charge a Power Module	7		

Introduction

This instructions for use manual contains information intended to ensure the safe, effective, and compliant use of your product.

Audience

This manual is intended for in-service trainers, physicians, nurses, and surgical technologists. Keep and consult this reference manual during the life of the product.

Conventions

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

Contact Information

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service at 1-269-323-7700 or 1-800-253-3210. Outside the US, contact your nearest Stryker subsidiary.

NOTE: The user and/or patient should report any serious product-related incident to both the manufacturer and the Competent Authority of the European Member State where the user and/or patient is established.

Intellectual Property

Trademarks not the property of Stryker Corporation are the property of their respective owners.

Indications For Use

The STRYKER SYSTEM G Rotary Drill (handpiece) is a surgical battery-powered instrument intended for use during a variety of orthopedic and trauma procedures to provide power to operate various accessories or attachments, which may include a pin/wire driver, a driver bit, and/or a drill bit.

This handpiece is designed for surgical procedures where hard tissue and/or bone must be reamed, drilled and/or fixated with screws, including, but not limited to, the hand wrist, elbow, shoulder, foot, ankle, knee, and hip.

Contraindications

None known.

Safety Directives

WARNINGS:

- Before using this equipment, read and understand the instructions for use. Pay particular attention to safety information.
- Only trained and experienced healthcare professionals should use this equipment.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker does not recommend surgical procedure or technique.

For Use With

This section identifies components intended to be used with the equipment to obtain a safe combination.

WARNINGS:

- Use only Stryker-approved equipment, unless otherwise specified.
- Use of unapproved electrical equipment may result in increased electromagnetic emissions or decreased electromagnetic immunity of the system.

NOTE: For a complete list of components or ordering information, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

Power and Care Components

Description	REF
STRYKER SYSTEM G Transfer Shield	7304-120-000
STRYKER SYSTEM G Power Module	7309-120-000
STRYKER SYSTEM G Charger, 2-Bay	7310-120-000
STRYKER SYSTEM G Sterilization Cases	7302-120-000
	7302-130-000

Attachments and Cutting and Fixation Tools

 **WARNING:** All cutting and fixation tools are intended for single use only.

Wire and Pin Collets

NOTE: Wire and pin collets REF 4103-063-000, 4103-126-000S3, and 4103-126-000 are not compatible with STRYKER SYSTEM G.

Attachment	REF	Tools
STRYKER SYSTEM G Wire Collet	7307-001-000	Kirschner wires and Steinmann pins from 0.7 to 1.8 mm in diameter
STRYKER SYSTEM G Pin Collet, Small	7307-002-000	Kirschner wires and Steinmann pins from 2.0 to 3.2 mm in diameter
STRYKER SYSTEM G Pin Collet, Large	7307-003-000	Kirschner wires and Steinmann pins from 3.0 to 4.2 mm in diameter

Keyed Chucks

NOTE: Each chuck is supplied with a chuck key (REF 1331-001-009)

Attachment	REF	Tools
1:1 1/4 inch (6.4 mm) Keyed Chuck	4103-131-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter
High Torque 1/4 inch (6.35 mm) Drill with Chuck 2:1 with Key	4103-180-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter
3.25:1 1/4 inch (6.4 mm) Keyed Chuck	4103-231-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter

Keyless Chucks

NOTE: The Keyless Chuck Adaptor (REF 4103-082-131) is intended for use with standard Trinkle attachments.

Attachment	REF	Tools
Keyless Chuck Adaptor	4103-082-131	Wires, pins, drill bits, and tools up to 6.4 mm in diameter
Keyless Drill Attachment	4103-133-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter

Drill Attachments


Attachment	REF	Tools
1:1 AO Small Drill	4103-110-000	Drill bits, taps, automatic screwdrivers, and tools with AO snap-lock shanks
1:1 Hudson	4103-113-000	Drill bits, reamers, taps, automatic screwdrivers, and tools with Hudson fittings
1:1 Hudson/Modified Trinkle	4103-135-000	Drill bits, reamers, automatic screwdrivers, and tools with Hudson or modified Trinkle fittings
1:1 Trinkle	4103-160-000	Drill bits, reamers, automatic screwdrivers, and tools with standard Trinkle fittings

Reamer Attachments

Attachment	REF	Tools
3.25:1 AO Large Reamer	4103-210-000	Drill bits, flexible reamers, automatic screwdrivers, and tools with AO snap-lock shanks
3.25:1 Hudson	4103-213-000	Drill bits, reamers, taps, automatic screwdrivers, and tools with Hudson fittings
3.25:1 Hudson/ Modified Trinkle	4103-235-000	Drill bits, reamers, automatic screwdrivers, and tools with Hudson or modified Trinkle fittings
3.25:1 Trinkle	4103-260-000	Drill bits, reamers, automatic screwdrivers, and tools with standard Trinkle fittings
4:1 Hudson	4103-413-000	Drill bits, reamers, taps, automatic screwdrivers, and tools with Hudson fittings
4:1 Hudson/ Modified Trinkle	4103-435-000	Drill bits, reamers, automatic screwdrivers, and tools with Hudson or modified Trinkle fittings




Features

Handpiece


A	Safety Switch – Locks one or both of the triggers.
B	Release Lever – Push to release an attachment from the handpiece.
C	Reverse Trigger – Controls the variable speed operation of the handpiece in a counterclockwise direction.
D	Forward Trigger – Controls the variable speed operation of the handpiece in a clockwise direction.
E	Door – Opens to allow installation of a power module.
F	Latch – Secures the door.

Safety Switch

	<p>Forward Mode – Only the forward trigger is functional. The reverse trigger is locked to prevent inadvertent operation of the handpiece in a counterclockwise direction.</p>
	<p>Forward/Reverse Mode – Both triggers are functional.</p>
	<p>Safe Mode – Both triggers are locked to prevent inadvertent operation of the handpiece.</p>


Attachments

Wire and Pin Collets


<p>A</p>	<p>Clamping Lever – Secures and releases the wire or pin from the attachment.</p>
<p>B</p>	<p>Pre-adjust Knob – Allows the collet to receive wires and pins of various diameters.</p>
<p>C</p>	<p>J-Notch – Secures the attachment in the handpiece.</p>

Keyed Chucks


<p>A</p>	<p>Jaws – Secure the cutting tool in the attachment.</p>
<p>B</p>	<p>Chuck Key – A tool used to open and close the jaws.</p>
<p>C</p>	<p>J-Notch – Secures the attachment in the handpiece.</p>

Keyless Chucks

<p>A</p>	<p>Lock Switch – Locks the collar to prevent inadvertent operation of the attachment (REF 4103-133-000 only).</p>
<p>B</p>	<p>Chuck – Rotates to open and close the jaws.</p>
<p>C</p>	<p>Jaws – Secure the cutting tool in the attachment.</p>
<p>D</p>	<p>J-Notch – Secures the attachment in the handpiece.</p>

Drill and Reamer Attachments

A	Retaining Collar – Slide to allow installation of a cutting tool.
B	J-Notch – Secures the attachment in the handpiece.

Definitions

The symbols located on the equipment and/or labeling are defined in this section and/or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

Symbol	Definition
	General Warning Sign
	Forward/Reverse Trigger (Action is in the direction of the arrow.)
	Forward/Reverse Trigger (Action is in the direction of the arrow.)
F/R 1 min/ 4 min X 3	Duty Cycle (See the <i>Specifications</i> section.)
	Lock
	Unlock
	Direction of Rotation
	Alignment Mark
USE ONLY STRYKER POWER MODULE	Use only Stryker power module.
STERILIZE WITH DOOR FULLY OPEN	Sterilize with door fully open.
S	Small Diameter
L	Large Diameter
	Run
	Load
SERIAL	Serial Number

Instructions

To Process the Equipment

WARNING: Upon initial receipt and before each use, process the equipment as indicated in the care instructions manual supplied with the equipment.

To Charge a Power Module

Fully charge the power module before first and every use.


NOTES:

- For charging instructions, see the instructions for use supplied with the battery charger.
- For additional information, including safety information, see the instructions for use supplied with the power module.


To Install a Power Module

Sterile Assistant


1. Set the safety switch to the safe mode position.

2. Turn the latch to the unlock position and open the door.

3. Place the sterile transfer shield over the open compartment.


4. Present the compartment to the Circulating Assistant.


Circulating Assistant

⚠ WARNING: DO NOT contaminate the handpiece.
The transfer shield guards the sterile handpiece from contamination during power module installation.

1. Hold the power module by the handle.
2. Insert the power module through the transfer shield into the open compartment.


3. Push the handle down.

⚠ WARNING: ALWAYS sterilize the contaminated transfer shield before using it again.

4. Remove the contaminated transfer shield.


Sterile Assistant


⚠ WARNINGS:

- When closing the door, hold the door by the latch. DO NOT touch the edges of the door or the internal surface of the door that faces the power module.
- ALWAYS verify the door is securely closed to prevent exposing the power module to the sterile field.

1. Holding the door by the latch, use the door to push the power module handle down until it locks into position against the power module.


2. Close the door and turn the latch to the lock position. Verify the mark on the latch aligns with the lock symbol and the door is securely closed.

NOTE: Although one power module is usually sufficient for one operation, Stryker recommends preparing a second handpiece with power module to ensure quick intraoperative exchange of handpieces under sterile conditions.

To Install an Attachment

1. Align the J-notch with the handpiece and insert the attachment until it snaps securely into place.

2. Gently pull on the attachment to verify it is secure.

To Install and Remove a Cutting or Fixation Tool

Wire and Pin Collets

1. Select the appropriate collet for the desired wire or pin diameter.

NOTE: All collets have a pre-adjust knob with small and large diameter settings.

Wire Collet (REF 7307-001-000)

Setting	Wire or Pin Diameter
Small (S)	0.7 to 1.2 mm
Large (L)	1.3 to 1.8 mm


Pin Collet, Small (REF 7307-002-000)

Setting	Wire or Pin Diameter
Small (S)	2.0 to 2.6 mm
Large (L)	2.7 to 3.2 mm


Pin Collet, Large (REF 7307-003-000)

Setting	Wire or Pin Diameter
Small (S)	3.0 to 3.6 mm
Large (L)	3.7 to 4.2 mm


2. Pull and rotate the pre-adjust knob to align the symbols for the desired wire or pin diameter. After alignment, snap the knob back into place.

3. Insert the wire or pin through the front of the collet or the back of the handpiece until the desired length extends from the collet.


4. Pull the clamping lever back to grasp the wire or pin.


5. To remove the wire or pin, or to adjust the wire length, release the clamping lever.

Keyed Chucks


1. Use the chuck key to open the jaws.

2. Insert the tool into the center of the jaws.


3. Use the chuck key to close the jaws.
4. Visually inspect the tool to make sure it is centered in the jaws.


5. Gently pull on the tool to verify it is secure.
6. To remove the tool, use the chuck key to open the jaws again.

Keyless Chucks

1. Hold the attachment and manually rotate the chuck to open the jaws.

2. Insert the tool into the center of the jaws.


3. Hold the attachment and manually rotate the chuck to close the jaws.
4. Visually inspect the tool to make sure it is centered in the jaws.

WARNING: DO NOT operate the handpiece to close the jaws on the cutting tool.


5. Hold the attachment and manually rotate the chuck to close the jaws.
6. Visually inspect the tool to make sure it is centered in the jaws.


5. Gently pull on the tool to verify it is secure.
6. To remove the tool, hold the attachment and manually rotate the chuck to open the jaws again.

Drill and Reamer Attachments


1. Slide the retaining collar back.

2. Hold the retaining collar in position and insert the tool into the attachment.

3. Release the retaining collar.

4. Gently pull on the tool to verify it is secure.
5. To remove the tool, slide the retaining collar back again.

To Operate a Handpiece

⚠ WARNINGS:

- DO NOT open the handpiece housing during surgery as the inside of the device is not sterile.
- ALWAYS follow the specified duty cycle to prevent the equipment from overheating. See the *Specifications* section.
- ALWAYS consult attachment and tool labeling for product-specific duty cycles and instructions for use.
- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen, or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using this equipment. Place the equipment into service according to the EMC information contained in this manual.
- Portable and mobile radio frequency (RF) communications equipment can affect the function of this equipment.
- This equipment should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, ALWAYS verify normal operation in the configuration in which it will be used.
- ALWAYS operate the equipment within the specified environmental condition values:

Temperature Limitation	Humidity Limitation	Atmospheric Pressure Limitation

1. Before operating the equipment, verify the following:
 - You are using Stryker-approved equipment.
 - The equipment is not damaged.
 - The equipment is clean and sterile.
 - You have properly assembled the equipment.
 - You have installed a fully-charged power module.

NOTE: If you experience difficulty during operation, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

- Set the safety switch to the forward or forward/reverse mode position.

NOTE: Set the safety switch back to the safe mode position when the handpiece is idle or when passing the handpiece to another person.


CAUTION: DO NOT apply excessive pressure or restrict the motion of the cutting or fixation tool during operation. An excessive load on the system may cause the power module to temporarily shut down. If this happens, release the trigger to resume operation or allow the power module to cool for 20 to 40 minutes.

- Pull the triggers to control the variable speed operation of the handpiece.

NOTE: Typical power module run time is 5 minutes or longer, depending on surgical procedure and technique.

To Remove an Attachment


- Push the release lever down and pull the attachment from the handpiece.

To Remove a Power Module

NOTE: Keep the power module as clean as possible because it cannot be sterilized.

- Wipe all debris from the external surfaces of the handpiece.
- Turn the latch to the unlock position and open the door.
- Use the handle to remove the power module from the handpiece.

Specifications

Handpiece

Model:	STRYKER SYSTEM G Rotary Drill
Catalog Number:	7305-001-000
Dimensions:	
Length:	168.15 mm
Width:	83.82 mm
Height:	226.10 mm
Mass:	1.043 kg
Nominal Speed:	
Drill:	825 revolutions per minute
Ream:	250 revolutions per minute
Mode of Operation:	Non-continuous
Duty Cycle:	1 minute on/4 minutes off, 3 times
Rest Between Cycles:	4 hours
Applied Part:	The handpiece, attachment, and cutting or fixation tool as defined by the manufacturer
Maximum Temperature of Applied Part:	Less than 43 °C, as tested to the <i>Product Safety Compliance</i> standards at a maximum ambient temperature of 30 °C.
Power Supply:	Internally Powered
Ingress Protection (IP):	IPX0
Equipment Type:	 Type BF Applied Part
European Conformity:	 0197

Attachments

Wire and Pin Collets

Model:	STRYKER SYSTEM G Wire Collet	STRYKER SYSTEM G Pin Collet, Small	STRYKER SYSTEM G Pin Collet, Large
Catalog Number:	7307-001-000	7307-002-000	7307-003-000
Dimensions:			
Length:	88 mm	88 mm	88 mm
Diameter:	26 mm	26 mm	26 mm
Mass:	0.27 kg	0.27 kg	0.27 kg
Gear Ratio:	1:1	1:1	1:1
Ingress Protection (IP):	IPX0		
European Conformity:			

Keyed Chucks

Model:	1:1 1/4 inch (6.4 mm) Keyed Chuck	High Torque 1/4 inch (6.35 mm) Drill with Chuck 2:1 with Key	3.25:1 1/4 inch (6.4 mm) Keyed Chuck
Catalog Number:	4103-131-000	4103-180-000	4103-231-000
Dimensions:			
Length:	105.90 mm	111.25 mm	118.62 mm
Diameter:	28.60 mm	33.80 mm	28.60 mm
Mass:	0.313 kg	0.198 kg	0.182 kg
Gear Ratio:	1:1	2:1	3.25:1
Ingress Protection (IP):	IPX0		
European Conformity:	CE		

Keyless Chucks

Model:	Keyless Chuck Adaptor	Keyless Drill Attachment
Catalog Number:	4103-082-131	4103-133-000
Dimensions:		
Length:	87.12 mm	126.50 mm
Diameter:	30.07 mm	31.95 mm
Mass:	0.185 kg	0.157 kg
Gear Ratio:	1:1	1:1
Ingress Protection (IP):	IPX0	
European Conformity:	CE	

Drill Attachments

Model:	1:1 AO Small Drill	1:1 Hudson	1:1 Hudson/Modified Trinkle	1:1 Trinkle
Catalog Number:	4103-110-000	4103-113-000	4103-135-000	4103-160-000
Dimensions:				
Length:	89.50 mm	84.33 mm	78.25 mm	79.85 mm
Diameter:	26.04 mm	26.04 mm	26.04 mm	26.04 mm
Mass:	0.185 kg	0.183 kg	0.157 kg	0.154 kg
Gear Ratio:	1:1	1:1	1:1	1:1
Ingress Protection (IP):	IPX0			
European Conformity:	CE	N/A	CE	CE

Reamer Attachments

Model:	3.25:1 AO Large Reamer	3.25:1 Hudson	3.25:1 Hudson/Modified Trinkle	3.25:1 Trinkle	4:1 Hudson	4:1 Hudson/Modified Trinkle
Catalog Number:	4103-210-000	4103-213-000	4103-235-000	4103-260-000	4103-413-000	4103-435-000
Dimensions:						
Length:	95.45 mm	100.46 mm	94.00 mm	97.30 mm	100.46 mm	94.00 mm
Diameter:	27.95 mm	27.95 mm	27.95 mm	27.70 mm	32.40 mm	32.40 mm
Mass:	0.198 kg	0.256 kg	0.182 kg	0.180 kg	0.225 kg	0.208 kg
Gear Ratio:	3.25:1	3.25:1	3.25:1	3.25:1	4:1	4:1
Ingress Protection (IP):	IPX0					
European Conformity:						

Transfer Shield

Model:	STRYKER SYSTEM G Transfer Shield
Catalog Number:	7304-120-000
Dimensions:	
Length:	150.9 mm
Width:	103.5 mm
Height:	37.0 mm
Mass:	77 g
European Conformity:	

Product Safety Compliance

International Electrotechnical Commission (IEC)

IEC 60601-1 A1: R:2012 (3.1 Edition): *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); Amendment 1 (2012)

IEC 60601-1-6: 2013: *USABILITY – Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability* – Edition 3.1 Consolidated Reprint.

IEC 62366: 2014: *USABILITY – Medical devices – Application of usability engineering to medical devices* - Edition 1.1 Consolidated Reprint.

IEC 60601-1-2: 2014: *EMC (Version 4.0) – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests* – Edition 4.0

European Committee for Electrotechnical Standardization (CENELEC)

EN 60601-1:2006/A1:2013 (3.1 Edition): *Medical electrical equipment Part 1: General requirements for basic safety and essential performance* – Incorporates Amendment A12: 2014

Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions		
The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	<p>The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:</p> <p>⚠ WARNING: This equipment/system is intended for use in a professional healthcare facility environment. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) or shielding the location.</p>
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration - electromagnetic immunity			
NOTE: The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±2 kV, ±4 kV, ±6 kV, ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±2 kV, ±4 kV, ±6 kV, ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 20%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV at 100 kHz repetition frequency for power supply lines ±1 kV at 100 kHz repetition frequency for input/output lines	N/A N/A	N/A
Surge IEC 61000-4-5	±0.5, ±1 kV line(s) to line(s) ±0.5, ±1, ±2 kV line(s) to earth	N/A N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 0% U_T (100% dip in U_T) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% U_T (100% dip in U_T) for 1 cycle at 0° 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 & 30 cycles at 0° <5% U_T (>95% dip in U_T) for 5 s 0% U_T (100% dip in U_T) for 5 s	N/A N/A N/A N/A N/A	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m at 50 and 60 Hz	30 A/m at 50 and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the alternating current mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity			
<p>NOTE: The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should assure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands 80% AM at 1 kHz 6 Vrms 150 kHz to 80 MHz in ISM bands 80% AM at 1 kHz	N/A	<p>IEC 60601-1-2 3rd Edition:</p> <p>Portable and mobile RF equipment should be used no closer to any part of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1.2\sqrt{P}$ <p>$d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>IEC 60601-1-2 4th Edition:</p> <p>⚠ WARNING: Portable RF equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> <p style="text-align: center;"> (Non-ionizing electromagnetic radiation)</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	3 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz	

NOTE 1: At 80 MHz and 800MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is used exceeds the applicable RF compliance level above, the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000).

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

IEC 60601-1-2 3rd Edition: Recommended separation distances between portable and mobile RF equipment and the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000)			
<p>The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF equipment (transmitters) and the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) as recommended below, according to the maximum output power of the equipment.</p>			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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IEC 60601-1-6: 2013: *USABILITY – Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability – Edition 3.1 Consolidated Reprint.*

IEC 62366: 2014: *USABILITY – Medical devices – Application of usability engineering to medical devices - Edition 1.1 Consolidated Reprint.*

IEC 60601-1-2: 2014: *EMC (Version 4.0) – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests – Edition 4.0*

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EN 60601-1:2006/A1:2013 (3.1 Edition): *Medical electrical equipment Part 1: General requirements for basic safety and essential performance – Incorporates Amendment A12: 2014*

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