Non-Sterile Instruments

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430 USA
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

C **€** 2797

©2021/08 Howmedica Osteonics Corp.

QIN 4330 Rev AA 08/2021

Refer to product label for CE mark status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

stryker

ifu.stryker.com

Adobe® Reader®

Free paper copy within 7 days / Version papier gratuite dans les 7 jours / Kostenlose Papierversion in 7 Tagen / Cópia en papel gratuita dentro de unos 7 dias / Cópia cartacea gratuita entro 7 giorni / Cópia impressa gratuita dentro de 7 dias / Gratis kopia på papper inom 7 dagar / Gratis papirkopi i løbet af 7 dage/Ilmainen kopio paperilla 7 päivän kuluessa/Binnen 7 dagen gratis op papier afgedrukt exemplaar / Δωρεάν έντυπο αντίγραφο εντός 7 ημερών/Gratis papirkopi innen 7 dager / 紙のコピーを7日間以内無料提供 / 7일 이내 무료 사본 / Bezpłatna kopia papierowa w ciągu 7 dni / 7일 이내 무료 사본 /可在7天内免费索取纸质副本 / Copie gratuită pe hârtie în termen de 7 zile / Немоката ропиеринэ копија per 7 dienas / Безплатно копие на хартиен носител до 7 дни / Kopie na papíru zdarma během 7 dnů / 7gün içinde ücretsiz basılı kopya

**AR AT AU BE BG BR CH CO CY CZ DE DK
ES FI FR GB HU IE IL IT LU LV NL NO NZ
PL PT RU SE SI UY**

+800 135 79 135

**AD AE BH BM CN DO DZ EG HK IN IQ JO
JP KR KW KZ LB LY MX MY OM PR QA
RS SA SD TH TN TW TZ UA ZA**

+31 20 794 7071

BR 0800 591 1055
CA +1 855 805 8539
CL 800 914 248
EE 0800 0100567
GR 00800 161 2205 7799
HR 0800 804 804
IS 800 8996
LI +31 20 796 5692
LT 8800 30728
MT +31 20 796 5693
RO 0800 895 084
SG 800 101 3366
SK 0800 606 287
TR 00800 142 064 866
US +1 855 236 0910
VN 122 80297

CONTENTS


English	4
Български	9
中文	14
Čeština	18
Dansk	22
Nederlands	26
Suomi	31
Français	35
Deutsch	40
Ελληνικά	44
Italiano	49
日本語	54
한국어	58
Lietuvių kalba	62
Norsk	66
Polski	70
Português	75
Română	80
Español	85
Svenska	90
Türkçe	94


NON-STERILE INSTRUMENTS

Labeling Symbols Glossary

The following is a list of symbols that may be used on Stryker medical device labeling. Refer to individual product labels for applicable symbology for each product.

Symbol		Standard
Caution (See instructions for use for operator awareness)		ISO 15223-1
Consult instructions for use or consult electronic instructions for use	 <small>lit.stryker.com</small>	ISO 15223-1
Do not re-use		ISO 15223-1
Sterilized using irradiation	STERILE R	ISO 15223-1
Sterilized using hydrogen peroxide	STERILE GP	N/A
Sterilized using ethylene oxide	STERILE EO	ISO 15223-1
Sterilized using Aseptic processing techniques (Aseptic fill)	STERILE A	ISO 15223-1
Non-sterile		ISO 15223-1
Do not resterilize		ISO 15223-1
Single sterile barrier system		ISO 15223-1
Double sterile barrier system		ISO 15223-1
Use-by date		ISO 15223-1
Date of manufacture		ISO 15223-1
Legal manufacturer		ISO 15223-1
Authorized representative in the European Community	EC REP	ISO 15223-1
Catalogue number	REF	ISO 15223-1
Batch code	LOT	ISO 15223-1
Serial number	SN	ISO 15223-1
MR Safe	MR	ASTM F2503
MR Conditional		ASTM F2503
MR Unsafe		ASTM F2503
Do not use if package is damaged		ISO 15223-1
Medical device	MD	ISO 15223-1
Quantity	QTY	N/A
Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician	R_x Only	N/A
Contains hazardous substances		ISO 15223-1
Temperature limit		ISO 15223-1

Symbol		Standard
Keep dry		ISO 15223-1
Keep away from sunlight		ISO 15223-1
Patient identification		ISO 15223-1
Date (of Implantation)		ISO 15223-1
Health care center or doctor		ISO 15223-1
Patient information website		ISO 15223-1
Unique Device Identifier		ISO 15223-1
Contact by Phone		N/A

HOWMEDICA OSTEONICS Corp. (herein referred to as Stryker Orthopaedics) invasive instrumentation consists of manual surgical instruments intended for use during hip, knee and osteosyntheses orthopaedic surgeries. Generic, hand-held orthopaedic instrumentation examples include but are not limited to retractors, inserters, impactors, in addition to positioning, trialing, measuring and attached to power instruments.

Utilization

- Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- Stryker Orthopaedics Trial Components should be used for size determination, preparation evaluation, trial reduction and range of motion evaluation, thus helping to preserve the integrity of the actual implants and their sterile packaging.
- The Stryker Orthopaedics Surgical Protocols provide additional procedural information as applicable.

Clinical Benefit

- The instruments provided are intended to be used in orthopaedic implant surgeries. The clinical benefit of the instruments is the facilitation of implantation and/or explantation of the orthopaedic implants.

Warnings

- Care should be taken when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing differing design parameters, varying tolerances, different materials and manufacturing specifications, Stryker Orthopaedics instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of Stryker Orthopaedics for the performance of the resulting implant.
- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials that are visible.
- Unless otherwise specified on the instrument, markings on instruments with a measuring function have a linear accuracy of ± 0.01 inches (± 0.25 mm) or $\pm 5^\circ$ for angular measurements.
- Stryker Orthopaedics instruments are not intended to be used in or present within the MR environment. Since they are not intended to enter the MR environment, testing for safety in the MR environment has not been undertaken, and thus any potential hazards from exposure to such an environment is unknown.
- Stryker Orthopaedics non-sterile instruments [or: one or more components of these instruments] may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight:

- Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless-steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects. Refer to the instrument label to determine if the instrument contains hazardous substances.

Adverse Effects

- Malfunctions of instruments, including, but not limited to, weld fracture, thread galling and material fatigue, have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement or trauma surgery. These complications include but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impaction of the component into the preparation.
- Infection. Ensure the instruments are properly cleaned and sterilized to minimize the risk of infection.

Infection

Transient Bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have been associated with transient bacteremia. To help prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.

Materials

Stryker Orthopaedics reusable instruments may contain one or more metals in pure and/or alloyed forms, including but not limited to the following metals: stainless steel, cobalt chromium, nickel, aluminum, and titanium. Stryker Orthopaedics reusable instruments may also contain one or more of the following polymers including but not limited to: polypropylene, polyphenylsulfone (PPSU), polyetherimide, polyoxymethylene, silicone, polyamide, polytetrafluoroethylene (PTFE), polyethylene, Santoprene, and polycarbonate. Contact Stryker Orthopaedics for instrument-specific information when treating a patient with a material allergy or sensitivity. As noted, these are common materials used in Stryker reusable instruments. This is not intended to be an exhaustive list of all materials used in these devices.

Examination Prior to Use

- The life of the instruments depends on the number of times they are used as well as the precautions taken in handling, cleaning and storage. Great care must be taken of the instruments to ensure that they remain in good working order.
- Instruments should be examined for wear or damage by physicians and staff in operating centers prior to surgery.
- The examination shall include a visual and functional inspection of the working surfaces, articulation points, rotating features, hinges, springs, connection mechanisms, mating parts, threads, and working ends of all instruments. Functional inspection should fully replicate the intended use of the device to confirm the instrument moves, assembles, and/or rotates as expected. It should also include verifying all welded connections, that all components are present, and the cleanliness of the orifices and cavities, as well as examination for signs of material degradation including but not limited to cracks, distortion/deformation, impact, corrosion, detached pieces or other unexpected changes. If one of the above-mentioned conditions occurs and impacts device functionality, the instrument has reached the end of its functional life and must be replaced. If damaged instrumentation is used, possible fracture, jamming, or other failure may occur. For instruments with moving parts, application of medical grade lubricants that are biocompatible per ISO 10993 may be necessary.
- Stryker Orthopaedics shall not be responsible in the event of the use of instruments that are damaged, incomplete, show signs of excessive wear and tear, or that have been repaired or sharpened outside the control of Stryker. Any faulty instruments must be replaced prior to any surgical procedure. Please refer to e-IFU "**Instructions for Inspection and Maintenance of Reusable Medical Devices**" – **SLI0001** for a comprehensive guide on instrument inspection and determination of an instrument's end of life. **This e-IFU is available at ifu.stryker.com.**

Cleaning and Sterilization

- These reusable instruments are **not** sterile. All Stryker Orthopaedics reusable instruments must be cleaned and sterilized to prepare them for use. Please refer to Cleaning, Sterilization, Inspection and Maintenance of Reusable Medical Devices - LSTPI-B, available at ifu.stryker.com.

- Remove all instruments from their packaging prior to cleaning and sterilization. When removing the instrument from its initial packaging or from a tray, the structural integrity of the instrument should be inspected to ensure safe use.
- Where noted in LSTPI-B, certain instruments must be disassembled prior to cleaning and sterilization.
- All parts must be cleaned and sterilized before use. Cleaning of the reusable instrument prior to sterilization should include: presoaking, manual cleaning, ultrasonication (e.g. for complex parts), and washing/disinfecting. Enzymatic cleaners, manual detergents, and neutral cleaners should be used.
- During and after cleaning, inspect hard-to-reach areas such as lumens within the instrument to help ensure that entrapped debris and/or soil are completely removed. If debris and/or soil are evident, then reclean the instrument until there is no evidence of debris and/or soil.
- Rinsing should be performed after each cleaning step to help ensure cleaners and/or debris/soil are removed.
- Stryker Orthopaedics reusable instruments are recommended to be sterilized by moist heat sterilization.
- Instruments with articulating surfaces must be tested for movement. A moist heat compatible, medical grade lubricant should be applied to all articulating joints prior to sterilization.
- After cleaning and prior to sterilization, the reusable instrument should be double-wrapped or packaged in CSR sterilization wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique).
- Select Stryker Orthopaedics trays and instruments are compatible with Aesculap SterilContainer System (JN441 or JN442 perforated bottom with JK48X series of aluminum 2000 lids). Please refer to LSTPI-B for a list of compatible parts, and detailed instructions for use.
- Biological indicators and/or steam integrators should be used within the sterilizer load to help ensure the proper sterilization conditions of time, temperature, and saturated steam have been met.
- The parameters for sterilization are:

Method:	Moist-Heat Sterilization
Cycle:	Pre-Vacuum (Pre-Vac)
Temperature:	270° F (132° C)
Exposure Time:	4 minutes (minimum)
Pressure:	2-15 PSIA
Dry-Time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Stryker Orthopaedics has validated the above recommended sterilization cycle based on AAMI/ANSI/ISO guidelines and recommendations. Other sterilization methods and cycles may also be suitable. However, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution.

- EtO sterilization and cold sterilization techniques are not recommended.

Select Stryker Orthopaedics trays and instruments may be sterilized with the following parameters **outside of the United States:**

Method:	Moist-heat sterilization
Cycle:	Pre-vacuum
Temperature:	273° F- 279 °F (134°C - 137°C)
Exposure Time:	3 minutes (minimum)
Dry-time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Stryker Orthopaedics has validated the above recommended OUS sterilization cycle based on ISO guidelines and recommendations. Please refer to e-IFU LSTPI-B for a list of instruments/trays that may be sterilized using this cycle. Additional validated cycle options are also available in this e-IFU available at ifu.stryker.com.

In case of serious incident, please notify the Manufacturer, and Competent Authority in your region.

Transport & Storage Information

Store the devices in standard hospital environmental conditions unless specific requirements are defined and described on the product label. Instrument trays are not intended as sterilization trays to maintain the product sterility. CSR sterilization wraps or pouches and select rigid containers may be used to maintain sterility in transit post hospital sterilization.

Safe Disposal

If a device is being returned for evaluation, please contact your local Stryker representative for shipping/handling information. If the device is not being returned to Stryker, instruments are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state and local regulations.

This IFU is available electronically at ifu.stryker.com. Detailed manual and automatic cleaning instructions, including times and temperatures, are provided in Stryker Orthopaedics e-IFU LSTPI-B available at ifu.stryker.com. Disassembly and assembly instructions applicable to all Stryker Orthopaedics re-usable instruments are all also provided in this e-IFU.

Refer to Aesculap Instructions for Use for Care & Handling of Aesculap SterilContainer Systems (JN441 or JN442).

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Howmedica, Osteonics, Stryker, Stryker Orthopaedics. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The following table contains a list of abbreviations that are used on Stryker Orthopaedics product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ►
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◀ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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