<p><i>Note: This portion of the document should not appear on produced Labels or IFUs.</i></p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p>Dimensions: 6 inch (width) x 6 inch Booklet</p>	<p>Color/Material/Finish: Black Graphics on White Background 20# Bond or Equivalent</p>	<p>Label Stock: N/A</p>	
			<p>Description/Type: Instructions For Use</p>	
			<p>Part Number: 4103-126-702</p>	<p>Rev. F</p>

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
Introduction

This instructions for use manual contains information intended to ensure the safe, effective and compliant use of your product. Keep and consult this reference manual during the life of the product.

This manual is intended for in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Indications for Use

See the instructions for use supplied with the handpiece.

Contraindications

None known.

For Use With

DESCRIPTION	REF
System 4 Rotary Handpiece	4103-000-000
System 5 Rotary Handpiece	4203-000-000
Dual Trigger Rotary Handpiece	4205-000-000

User/Patient Safety

WARNINGS:

- Before using this equipment, or any component compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
- Only healthcare professionals trained and experienced in the use of this medical device should operate this equipment.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage. DO NOT use any component if damage is apparent.
- Upon initial receipt and before each use, clean and sterilize the equipment as indicated. See the care instructions manual supplied with the handpiece.
- DO NOT disassemble or service this equipment without the authorization of the manufacturer.
- Check all moving parts for free movement. Be alert for unusual sounds or vibrations and note the operating speed.

Accessories

WARNINGS:

- Use only Stryker-approved components and accessories unless otherwise specified. DO NOT modify any component or accessory.
- DO NOT reuse, reprocess, or repackage a device that is intended for single use only. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.


See the instructions for use supplied with the handpiece.

Definitions

The symbols located on the equipment and/or labeling are defined in this section and/or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

S	Small Diameter Pins	L	Large Diameter Pins
0.079-0.102 inch [2.0-2.6 mm]		0.106-0.126 inch [2.7-3.2 mm]	

Features

A	Pre-adjust Knob – Used to prepare the collet to accept the correct pin diameter size.
B	Knob Pin Holes (2) – Used to receive the knob pin.
C	Knob Pin – Used to insert into either the small or large pre-adjust knob pin hole.
D	Nose – Allows the insertion of the pin into the collet.
E	Clamping Lever – Used to clamp or grasp the pin before insertion or removal.

Instructions

WARNINGS:

- ALWAYS place the handpiece in the SAFE mode while the handpiece is idle, before installing or removing any accessory or attachment, or when passing the handpiece to another person.
- DO NOT attempt to insert or remove any cutting accessory or attachment while the handpiece is operating.

To Install the Attachment

See the instructions for use supplied with the handpiece.

To Insert the Pin

1. Set the pre-adjust knob for a small (S) or large (L) pin diameter range by pulling and rotating the knob until S aligns with S or L aligns with L, respectively.

NOTE: Make sure the knob retracts and locks in place when the symbols align.

2. Insert the pin into the nose of the collet or into the back cap of the handpiece until the desired length pin extends from the nose.

NOTE: Make sure the pin extends far enough through the collet to allow the collet to hold the pin before squeezing the clamping lever.

3. Remove the handpiece from the SAFE mode to allow rotary drill operation.
4. Place the handpiece trigger in the FORWARD orientation as required.
5. Squeeze the clamping lever against the handpiece handle to grasp pin.
6. Hold the pin against the bone. Press the FORWARD handpiece trigger to drive the pin into the bone.
7. To obtain additional pin length perform the following steps.
 - 7.1. Release the clamping lever of the collet.
 - 7.2. Pull the handpiece back from the bone to allow more pin length between the handpiece and bone.
 - 7.3. Squeeze the clamping lever to grasp the pin.
 - 7.4. Press the FORWARD handpiece trigger to drive the pin into the bone.

To Remove the Pin

1. Place the handpiece trigger in the REVERSE orientation as required.
2. Squeeze the clamping lever to grasp the pin.
3. Pulling the handpiece and pin away from the bone, press the REVERSE handpiece trigger to remove the pin from the bone.

To Remove the Attachment

See the instructions for use supplied with the handpiece.

Care Instructions

For processing instructions and disposal/recycle information, see the care instructions manual supplied with the handpiece.

Troubleshooting

See the instructions for use supplied with the handpiece.

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4103-126-702 Rev-F

www.stryker.com

Print Date: Sep 21, 2016 02:12:36 PM
0000000958, Rev. F Effective Date: Sep 21, 2016 11:14:29 AM

Signature History Report

Version	Signature	Role	Event Date	Vote
F	Mehrotra Abhishek	CA_Owner	2016-06-22 09:04:08 EDT	Continue
F	Pearson, Dawn	Editor	2016-07-20 11:09:08 EDT	Continue
F	Mehrotra Abhishek	CA_Owner	2016-07-20 12:12:10 EDT	Continue
F	Unger, Lauren	CA_Marketing Manager	2016-07-20 12:44:23 EDT	Approve
F	Sidow, Brian	CA_Clinical Sciences	2016-07-20 12:57:37 EDT	Approve
F	Brown Donna	CA_Product Safety	2016-07-20 18:10:15 EDT	Approve
F	Raghu Varun	CA_Product Design Engineer	2016-07-21 08:49:01 EDT	Approve
F	Modi Karankumar	CA_Regulatory Affairs	2016-07-21 14:29:07 EDT	Approve
F	Wippich, Caitlin	CA_Quality Engineer	2016-07-21 17:19:03 EDT	Approve
F	Snell, Tami	Document Control	2016-07-22 08:19:56 EDT	Complete
F	Schrauben Haley	Editor	2016-08-30 11:11:44 EDT	Continue
F	Mehrotra Abhishek	CA_Owner	2016-09-20 15:44:26 EDT	Approve

NOTE: IF THIS REPORT IS BLANK, PLEASE SEE WINDCHILL ECN FOR SIGNATURES.

Related Information

CR Number	CR0212150
CR Class	DMR-Commercially Released
Originating Entity	INST: Kalamazoo
Origin	Standards Revision
Initial CR Approval Date	01/16/2015 15:25:17
CN Number	CN140178
CA Number	CA165528
NPD Disposition Plan	