<p><i>Note: This portion of the document should not appear on produced Labels or IFUs.</i></p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p>Dimensions: 5 inch (width) x 8.5 inch Booklet; folded to 2.5 inch (width) x 8.5 inch</p>	<p>Color/Material/Finish: Black Graphics on White Background 20# Bond or Equivalent</p>	<p>Label Stock: N/A</p>	
			<p>Description/Type: Instructions For Use</p>	
			<p>Part Number: 4103-133-700</p>	<p>Rev. N</p>

stryker[®]

Keyless Drill Attachment


REF 4103-133-000

Instructions For Use

R_x ONLY

CE

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)
NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)
SUOMI (FI)
PORTUGUÊS (PT)
NORSK (NO)
POLSKI (PL)
ΕΛΛΗΝΙΚΑ (EL)
日本語 (JA)
中文 (ZH)
한국어 (KO)

For Use With

For use with System 4 Rotary Handpiece REF 4103-000-000 or System 5 Rotary Handpieces REF 4203-000-000 and REF 4205-000-000.

User/Patient Safety

WARNINGS:

- Before using this equipment, or any components compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
- Upon initial receipt and before each use, operate the equipment and inspect each component for damage. DO NOT use any equipment if damage is apparent.
- Upon initial receipt and before each use, clean and sterilize the equipment as indicated.
- DO NOT disassemble or service this equipment without the authorization of the manufacturer.
- Check all moving parts for free movement. Be alert for unusual sounds or vibrations, and note the operating speed.
- Use only Stryker approved accessories. Other accessories may cause tissue damage to patient and/or injury to healthcare staff, delay surgery if not functioning optimally, or cause damage to the handpiece. DO NOT modify any accessory to fit the handpiece.
- All cutting accessories are single use only. DO NOT reuse.
- ALWAYS place the handpiece in the SAFE mode while the handpiece is idle, before installing or removing any accessory, or when passing the handpiece to another person.


NOTE: For additional information, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Instructions

WARNING: You must manually close the chuck jaws to secure the accessory. Operating the handpiece to automatically close the jaws on the accessory may result in injury to the user.

1. Open the chuck jaws by holding the attachment with the knurled ring in one hand and turning the collar with the other hand.
2. Insert the cutting accessory and make certain it is centered in the jaws.
3. Close the jaws by holding the collar and turning the knurled ring until tight.
4. Tug the cutting accessory to ensure it is secure.
5. Visually verify that the cutting accessory is centered in the jaws.
6. Install the attachment in the handpiece.

Troubleshooting

Problem: The accessory wobbles.

Cause: The accessory extends too far from the distal end of the attachment, is the wrong size, or is not properly centered in the attachment.

Action: Reinsert the accessory. If wobble occurs, the bit may be worn or damaged. Try a different accessory. If wobble persists, return the attachment for repair.

Care Instructions

For processing instructions and disposal/recycle information, see the care instructions manual supplied with the handpiece.

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0000000996, Rev. N Effective Date: Sep 21, 2016 10:33:30 AM

Signature History Report

Version	Signature	Role	Event Date	Vote
N	Mehrotra Abhishek	CA_Owner	2016-06-21 11:13:09 EDT	Continue
N	Pearson, Dawn	Editor	2016-07-20 11:19:34 EDT	Continue
N	Mehrotra Abhishek	CA_Owner	2016-07-20 12:17:49 EDT	Continue
N	Unger, Lauren	CA_Marketing Manager	2016-07-20 12:41:14 EDT	Approve
N	Sidow, Brian	CA_Clinical Sciences	2016-07-20 12:59:40 EDT	Approve
N	Brown Donna	CA_Product Safety	2016-07-20 15:05:13 EDT	Approve
N	Raghu Varun	CA_Product Design Engineer	2016-07-21 08:45:07 EDT	Approve
N	Modi Karankumar	CA_Regulatory Affairs	2016-07-21 11:19:58 EDT	Approve
N	Wippich, Caitlin	CA_Quality Engineer	2016-07-21 17:06:22 EDT	Approve
N	Pulliam Sylvia	Document Control	2016-07-22 07:48:49 EDT	Complete
N	Schrauben Haley	Editor	2016-08-30 11:14:08 EDT	Continue
N	Mehrotra Abhishek	CA_Owner	2016-09-20 15:53:02 EDT	Approve

NOTE: IF THIS REPORT IS BLANK, PLEASE SEE WINDCHILL ECN FOR SIGNATURES.

Related Information

CR Number	CR0212150
CR Class	DMR-Commercially Released
Originating Entity	INST: Kalamazoo
Origin	Standards Revision
Initial CR Approval Date	01/16/2015 15:25:17
CN Number	CN140180
CA Number	CA165562
NPD Disposition Plan	