<p>NOTE: This portion of the document should not appear on produced Labels or IFUs.</p> <p>Stryker Instruments (269) 323-7700 (800) 253-3210</p>	<p><u>Dimensions:</u></p> <p>8.5 inch (width) x 11 inch Booklet</p>	<p><u>Color/Material/Finish:</u></p> <p>Black Graphics on White Background 20# Bond or equivalent</p>	<p><u>Label Stock:</u></p> <p>N/A</p>	
	<p><u>Print Location:</u></p> <p>Print Center</p>	<p><u>Suppliers/Services:</u></p> <p>TransPerfect MDS</p>	<p><u>Description/Type:</u></p> <p>Instructions For Use</p>	
			<p><u>Part Number:</u></p> <p>5450-800-700</p>	<p><u>Rev.</u></p> <p>A</p>

stryker[®]


Sonopet[®] Sterilization and Assembly Trays

REF 5450-800-278
5450-800-279

Instructions For Use

R_x ONLY

CE

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)
NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)
SUOMI (FI)
PORTUGUÊS (PT)
NORSK (NO)
POLSKI (PL)
ΕΛΛΗΝΙΚΑ (EL)

Introduction

This instructions for use manual contains information intended to ensure the safe, effective and compliant use of your product. Keep and consult this reference manual during the life of the product.

Audience

This manual is intended for in-service trainers, biomedical equipment technicians, and central supply/sterile processing technicians.

Conventions

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

Contact Information

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

Indications For Use

The Stryker Sonopet Sterilization and Assembly Tray is designed to protect reusable medical devices during steam sterilization, storage and transportation. The Sterilization and Assembly Tray must be wrapped in a validated sterilization wrap prior to sterilization to maintain sterility or may also be placed in a rigid container that does not require a validated wrap.

Contraindications

None.

Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

SYMBOL	DEFINITION
	General warning sign
COIL CORD HERE	The designated location for placement of the coiled handpiece cord inside the tray.
	DO NOT grasp the handpiece by hand to install or remove the ultrasonic tip.

NOTE: Handpiece assembly instructions are displayed on the sterilization tray lid. See the instructions for use supplied with the handpiece.

For Use With

The tray is intended for use with Stryker Sonopet Ultrasonic handpieces, attachments, and accessories.

Description

The sterilization tray, otherwise known as a sterilization case, is designed to hold medical instruments for sterilization, storage, transportation, aseptic presentation of contents, and return of the contaminated items to the decontamination area. The typical sterilization tray consists of a bottom or base, a lid that is secured to the base by means of a latching mechanism, and carrying handles. To permit sterilant penetration into the tray, there are perforations in the base and lid. The tray may include an open basket(s) or insert tray(s), placed inside, to hold the handpieces and accessories. The sterilization tray must be completely enclosed in an validated sterilization wrap to maintain sterility.

Safety Directives


WARNINGS:

- Before using this equipment, or any component compatible with this equipment, read and understand the instructions for use. Pay particular attention to safety information. Become familiar with the equipment before use.
 - Only individuals trained and experienced in the processing of reusable medical devices should use this equipment.
 - Upon initial receipt and before each use, inspect each component for damage. DO NOT use any equipment if damage is apparent or the inspection criteria are not met. See the care instructions manual supplied with the handpiece.
 - DO NOT reuse, reprocess, or repackage a device that is intended for single use only.
 - A single use device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing.
 - Design features may make cleaning difficult.
 - Reuse may create a contamination risk and compromise structural integrity resulting in operational failure or fragmentation during use.
 - Critical product information may be lost during repackaging.
- Failure to comply may lead to infection or cross infection and result in patient and/or healthcare staff injury.

Care Instructions

For general processing instructions, see the care instructions manual supplied with the handpiece.


Features

A	Lid
B	Base
C	Handles (two)
D	Latches (two)
E	Assembly Bracket (7 mm or 10 mm)

To Remove the Lid from the Base


1. Open the latches.

2. Remove the lid from the base.

To Secure the Lid to the Base

1. Put the lid on the base and close both of the latches.

2. Verify that the latches are closed and the lid is securely attached to the base.

Maximum Quantity of Devices

WARNING: DO NOT exceed the specified maximum quantity of devices in an insert tray. The maximum quantity of devices is based on Stryker validation of this equipment. Exceeding the maximum quantity of devices may prevent proper sterilization of the equipment and/or result in excessive weight. Always follow the current local recommendations and/or regulations governing weights and weight limits.

DEVICE TYPE	MAXIMUM QUANTITY
Handpiece	1
Attachment	1
Accessories	2
Total of All Devices	4


Recommended Placement of Devices

WARNING: ALWAYS make sure the surgical instruments are clean and dry before loading the instruments into the sterilization tray. See the instructions for use supplied with the surgical instruments for instrument cleaning recommendations.


CAUTIONS:

- ALWAYS install the sterilization cap onto the handpiece plug properly to protect the plug during sterilization. See the instructions for use supplied with the handpiece.
- DO NOT drop the handpiece during preparation for sterilization.
- DO NOT strike the handpiece against a metal or surgical object during preparation for sterilization.

NOTE: See the *Maximum Quantity of Devices* section.

Place the handpiece, attachment, and accessories into the base as shown.

Specifications

Model:	Sonopet Sterilization and Assembly Tray REF 5450-800-278; 7 mm (0.28 inch) bracket REF 5450-800-279; 10 mm (0.39 inch) bracket
Size:	5.1 cm (2 inch) height 22.9 cm (9 inch) width 24.9 cm (9.8 inch) length
Weight:	1.4 kg (3.1 lb)
Material:	Lid and Base: Radel (plastic) Bracket, Handles, and Latches: Stainless Steel Insert: Silicon Rubber

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stryker[®]

Signature History Report

Version	Signature	Role	Event Date	Vote
A	Keeler, Stephanie (Contractor)	CA_Owner	2015-08-04 13:43:21 EDT	Continue
A	Pearson, Dawn	Editor	2015-08-11 13:40:27 EDT	Continue
A	Keeler, Stephanie (Contractor)	CA_Owner	2015-10-23 16:19:24 EDT	Continue
A	Thompson Kristen	CA_Clinical Sciences	2015-10-23 16:51:31 EDT	Approve
A	Morgan, Carianne	CA_Product Safety	2015-10-26 10:19:00 EDT	Approve
A	Ekola, David	CA_Quality Engineer	2015-10-26 13:25:23 EDT	Approve
A	Helgeson, Julia	CA_Regulatory Affairs	2015-10-27 09:33:29 EDT	Approve
A	Van Loon, Mackenzie	CA_Product Design Engineer	2015-10-28 10:09:44 EDT	Approve
A	Zemlick, Larry	CA_Marketing Manager	2015-10-29 15:54:43 EDT	Approve
A	Pulliam Sylvia	Document Control	2015-10-29 16:18:35 EDT	Complete
A	Petyko, Holly (Contractor)	Editor	2015-11-20 13:53:09 EST	Continue
A	Keeler, Stephanie (Contractor)	CA_Owner	2016-01-28 18:51:02 EST	Approve

NOTE: IF THIS REPORT IS BLANK, PLEASE SEE WINDCHILL ECN FOR SIGNATURES.

Related Information

CR Number	CR0229294
CR Class	DMR-Commercially Released
Originating Entity	INST: Kalamazoo
Origin	UDI
Initial CR Approval Date	01/15/2015 15:35:05
CN Number	CN115831

Related Information

CA Number	CA133940
NPD Disposition Plan	