# STRYKER SYSTEM G™

**Sterilization Case, 1 Hand Piece**


REF 7302-120-000

**Sterilization Case, 2 Hand Pieces**

REF 7302-130-000

Instructions For Use

**R<sub>x</sub> ONLY**

ENGLISH (EN)  
ESPAÑOL (ES)  
ITALIANO (IT)  
PORTUGUÊS (PT)  
TÜRKÇE (TR)  
РУССКИЙ (RU)  
中文 (ZH)  
한국어 (KO)

## Contents

Introduction . . . . . 2

    Audience . . . . . 2

    Conventions . . . . . 2

    Contact Information . . . . . 2

Indications For Use . . . . . 2

Contraindications . . . . . 2

For Use With . . . . . 2

Safety Directives . . . . . 2

Features . . . . . 3

Definitions . . . . . 3

Instructions . . . . . 3

Care Instructions . . . . . 3

    Maximum Quantity of Devices . . . . . 3

    Recommended Placement of Devices . . . . . 4

Specifications . . . . . 4

References . . . . . 4

## Introduction

This instructions for use manual contains information intended to ensure the safe, effective, and compliant use of your product.

## Audience

This manual is intended for in-service trainers, biomedical equipment technicians, and central supply/sterile processing technicians. Keep and consult this reference manual during the life of the product.

## Conventions

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

## Contact Information

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service at 1-269-323-7700 or 1-800-253-3210. Outside the US, contact your nearest Stryker subsidiary.

## Indications For Use

STRYKER SYSTEM G Sterilization Cases are designed to contain and protect reusable medical devices during sterilization, storage and transportation.

## Contraindications

None known.

## For Use With

This section identifies components intended to be used with the equipment to obtain a safe combination.

 **WARNING:** Use only Stryker-approved equipment, unless otherwise specified.

### NOTES:

- Components are sold separately.
- For a complete list of components or ordering information, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

STRYKER SYSTEM G Sterilization Cases are intended to be used with STRYKER SYSTEM G handpieces, attachments, and transfer shields.


## Safety Directives

 **WARNINGS:**

- Before using any system component, or any component compatible with this system, read and understand the instructions. Pay particular attention to safety information. Become familiar with the system components prior to use.
- Only individuals trained and experienced in the processing of reusable medical devices should use this equipment.
- See the care instructions manual supplied with your handpiece for information on the safe and effective use of this equipment.

## Features

**NOTE:** The Sterilization Case, 1 Hand Piece is shown for reference.

A	Handle
B	Latch
C	Bracket
D	Lid
E	Base
F	Perforations (Holes and Slots)


## Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.


Symbol	Definition
	General Warning Sign
	Consult instructions for use

## Instructions

To remove the lid from the base, pull the two latches up and away from the base. Unhook the latches from the lid and remove the lid.

To secure the lid to the base, hook the two latches onto the lid. Push the latches down and toward the base.

## Care Instructions

For general processing instructions, see the care instructions manual supplied with your handpiece.

## Maximum Quantity of Devices

**⚠ WARNING:** DO NOT exceed the specified maximum quantity of devices in a sterilization case. The maximum quantity of devices is based on Stryker validation of this equipment. Exceeding the maximum quantity of devices may prevent proper sterilization of the equipment and/or result in excessive weight. Always follow the current local recommendations and/or regulations governing weights and weight limits.

### Sterilization Case, 1 Hand Piece

Device Type	Maximum Quantity
Handpiece	1
Attachment	5
Accessory (Chuck Key)	2
Transfer Shield	1
<b>Total of All Devices</b>	<b>9</b>

### Sterilization Case, 2 Hand Pieces

Device Type	Maximum Quantity
Handpiece	2
Attachment	5
Accessory (Chuck Key)	2
Transfer Shield	2
<b>Total of All Devices</b>	<b>11</b>


## Recommended Placement of Devices

**NOTES:**

- See the *Maximum Quantity of Devices* section.
- The illustrated placement of devices is a recommendation only. Other placement combinations are acceptable.


### Sterilization Case, 2 Hand Pieces

**CAUTION:** ALWAYS load the handpieces with the doors in the open position.

### Sterilization Case, 1 Hand Piece

**CAUTION:** ALWAYS load the handpiece with the door in the open position.

## Specifications

<b>Model:</b>	STRYKER SYSTEM G Sterilization Case, 1 Hand Piece	STRYKER SYSTEM G Sterilization Case, 2 Hand Pieces
<b>Catalog Number:</b>	7302-120-000	7302-130-000
<b>Dimensions:</b>		
Height:	140 mm	140 mm
Length:	255 mm	255 mm
Width:	315 mm	492 mm
<b>Mass:</b>	2.27 kg	3.38 kg
<b>Material:</b>	Aluminum Stainless Steel Nylon	
<b>European Conformity:</b>	CE	

## References

This document is compliant with AAMI and International Organization for Standardization (ISO) standards related to the sterile processing of reusable medical devices.

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