**stryker**<sup>®</sup>

# STRYKER SYSTEM G™

## Rotary Drill


REF 7305-001-000

## Rotary Drill Attachments

REF 7307-XXX-XXX Series  
4103-XXX-XXX Series

Instructions For Use

**R<sub>x</sub> ONLY**

ENGLISH (EN)  
ESPAÑOL (ES)  
ITALIANO (IT)  
PORTUGUÊS (PT)  
TÜRKÇE (TR)  
РУССКИЙ (RU)  
中文 (ZH)  
한국어 (KO)

# Contents

Introduction . . . . .	3	To Install a Power Module . . . . .	8
Audience . . . . .	3	Sterile Assistant . . . . .	8
Conventions . . . . .	3	Circulating Assistant . . . . .	8
Contact Information . . . . .	3	Sterile Assistant . . . . .	9
Intellectual Property . . . . .	3	To Install an Attachment . . . . .	9
Indications For Use . . . . .	3	To Install and Remove a Cutting or Fixation Tool . . .10	
Contraindications . . . . .	3	Wire and Pin Collets . . . . .	.10
Safety Directives . . . . .	3	Keyed Chucks . . . . .	11
For Use With . . . . .	4	Keyless Chucks . . . . .	11
Power and Care Components . . . . .	4	Drill and Reamer Attachments . . . . .	.12
Attachments and Cutting and Fixation Tools . . . . .	4	To Operate a Handpiece . . . . .	.12
Wire and Pin Collets . . . . .	4	To Remove an Attachment . . . . .	.13
Keyed Chucks . . . . .	4	To Remove a Power Module . . . . .	.13
Keyless Chucks . . . . .	4	Specifications . . . . .	.14
Drill Attachments . . . . .	4	Handpiece . . . . .	.14
Reamer Attachments . . . . .	5	Attachments . . . . .	.14
Features . . . . .	5	Wire and Pin Collets . . . . .	.14
Handpiece . . . . .	5	Keyed Chucks . . . . .	.15
Safety Switch . . . . .	6	Keyless Chucks . . . . .	.15
Attachments . . . . .	6	Drill Attachments . . . . .	.15
Wire and Pin Collets . . . . .	6	Reamer Attachments . . . . .	.16
Keyed Chucks . . . . .	6	Transfer Shield . . . . .	.16
Keyless Chucks . . . . .	6	Product Safety Compliance . . . . .	.16
Drill and Reamer Attachments . . . . .	7	Electromagnetic Compatibility . . . . .	.17
Definitions . . . . .	7		
Instructions . . . . .	7		
To Process the Equipment . . . . .	7		
To Charge a Power Module . . . . .	7		

## Introduction

This instructions for use manual contains information intended to ensure the safe, effective, and compliant use of your product.

## Audience

This manual is intended for in-service trainers, physicians, nurses, and surgical technologists. Keep and consult this reference manual during the life of the product.

## Conventions

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

## Contact Information

For additional information, including safety information, in-service training, or current literature, contact your Stryker sales representative or call Stryker customer service at 1-269-323-7700 or 1-800-253-3210. Outside the US, contact your nearest Stryker subsidiary.

## Intellectual Property

Trademarks not the property of Stryker Corporation are the property of their respective owners.

## Indications For Use

The STRYKER SYSTEM G Rotary Drill (handpiece) is a surgical battery-powered instrument intended for use during a variety of orthopedic and trauma procedures to provide power to operate various accessories or attachments, which may include a pin/wire driver, a driver bit, and/or a drill bit.

This handpiece is designed for surgical procedures where hard tissue and/or bone must be reamed, drilled and/or fixated with screws, including, but not limited to, the hand wrist, elbow, shoulder, foot, ankle, knee, and hip.

## Contraindications

None known.

## Safety Directives

### **WARNINGS:**

- Before using this equipment, read and understand the instructions for use. Pay particular attention to safety information.
- Only trained and experienced healthcare professionals should use this equipment.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker does not recommend surgical procedure or technique.

## For Use With

This section identifies components intended to be used with the equipment to obtain a safe combination.

### WARNINGS:

- Use only Stryker-approved equipment, unless otherwise specified.
- Use of unapproved electrical equipment may result in increased electromagnetic emissions or decreased electromagnetic immunity of the system.

**NOTE:** For a complete list of components or ordering information, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

## Power and Care Components

Description	REF
STRYKER SYSTEM G Transfer Shield	7304-120-000
STRYKER SYSTEM G Power Module	7309-120-000
STRYKER SYSTEM G Charger, 2-Bay	7310-120-000
STRYKER SYSTEM G Sterilization Cases	7302-120-000
	7302-130-000

## Attachments and Cutting and Fixation Tools

 **WARNING:** All cutting and fixation tools are intended for single use only.

### Wire and Pin Collets

**NOTE:** Wire and pin collets REF 4103-063-000, 4103-126-000S3, and 4103-126-000 are not compatible with STRYKER SYSTEM G.

Attachment	REF	Tools
STRYKER SYSTEM G Wire Collet	7307-001-000	Kirschner wires and Steinmann pins from 0.7 to 1.8 mm in diameter
STRYKER SYSTEM G Pin Collet, Small	7307-002-000	Kirschner wires and Steinmann pins from 2.0 to 3.2 mm in diameter
STRYKER SYSTEM G Pin Collet, Large	7307-003-000	Kirschner wires and Steinmann pins from 3.0 to 4.2 mm in diameter

## Keyed Chucks

**NOTE:** Each chuck is supplied with a chuck key (REF 1331-001-009)

Attachment	REF	Tools
1:1 1/4 inch (6.4 mm) Keyed Chuck	4103-131-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter
High Torque 1/4 inch (6.35 mm) Drill with Chuck 2:1 with Key	4103-180-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter
3.25:1 1/4 inch (6.4 mm) Keyed Chuck	4103-231-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter

## Keyless Chucks

**NOTE:** The Keyless Chuck Adaptor (REF 4103-082-131) is intended for use with standard Trinkle attachments.

Attachment	REF	Tools
Keyless Chuck Adaptor	4103-082-131	Wires, pins, drill bits, and tools up to 6.4 mm in diameter
Keyless Drill Attachment	4103-133-000	Wires, pins, drill bits, and tools up to 6.4 mm in diameter

## Drill Attachments


Attachment	REF	Tools
1:1 AO Small Drill	4103-110-000	Drill bits, taps, automatic screwdrivers, and tools with AO snap-lock shanks
1:1 Hudson	4103-113-000	Drill bits, reamers, taps, automatic screwdrivers, and tools with Hudson fittings
1:1 Hudson/Modified Trinkle	4103-135-000	Drill bits, reamers, automatic screwdrivers, and tools with Hudson or modified Trinkle fittings
1:1 Trinkle	4103-160-000	Drill bits, reamers, automatic screwdrivers, and tools with standard Trinkle fittings

## Reamer Attachments

Attachment	REF	Tools
3.25:1 AO Large Reamer	4103-210-000	Drill bits, flexible reamers, automatic screwdrivers, and tools with AO snap-lock shanks
3.25:1 Hudson	4103-213-000	Drill bits, reamers, taps, automatic screwdrivers, and tools with Hudson fittings
3.25:1 Hudson/ Modified Trinkle	4103-235-000	Drill bits, reamers, automatic screwdrivers, and tools with Hudson or modified Trinkle fittings
3.25:1 Trinkle	4103-260-000	Drill bits, reamers, automatic screwdrivers, and tools with standard Trinkle fittings
4:1 Hudson	4103-413-000	Drill bits, reamers, taps, automatic screwdrivers, and tools with Hudson fittings
4:1 Hudson/ Modified Trinkle	4103-435-000	Drill bits, reamers, automatic screwdrivers, and tools with Hudson or modified Trinkle fittings




## Features

### Handpiece


<b>A</b>	<b>Safety Switch</b> – Locks one or both of the triggers.
<b>B</b>	<b>Release Lever</b> – Push to release an attachment from the handpiece.
<b>C</b>	<b>Reverse Trigger</b> – Controls the variable speed operation of the handpiece in a counterclockwise direction.
<b>D</b>	<b>Forward Trigger</b> – Controls the variable speed operation of the handpiece in a clockwise direction.
<b>E</b>	<b>Door</b> – Opens to allow installation of a power module.
<b>F</b>	<b>Latch</b> – Secures the door.

## Safety Switch

	<p><b>Forward Mode</b> – Only the forward trigger is functional. The reverse trigger is locked to prevent inadvertent operation of the handpiece in a counterclockwise direction.</p>
	<p><b>Forward/Reverse Mode</b> – Both triggers are functional.</p>
	<p><b>Safe Mode</b> – Both triggers are locked to prevent inadvertent operation of the handpiece.</p>


## Keyed Chucks

<b>A</b>	<b>Jaws</b> – Secure the cutting tool in the attachment.
<b>B</b>	<b>Chuck Key</b> – A tool used to open and close the jaws.
<b>C</b>	<b>J-Notch</b> – Secures the attachment in the handpiece.


## Attachments

### Wire and Pin Collets


<b>A</b>	<b>Clamping Lever</b> – Secures and releases the wire or pin from the attachment.
<b>B</b>	<b>Pre-adjust Knob</b> – Allows the collet to receive wires and pins of various diameters.
<b>C</b>	<b>J-Notch</b> – Secures the attachment in the handpiece.

### Keyless Chucks

<b>A</b>	<b>Lock Switch</b> – Locks the collar to prevent inadvertent operation of the attachment (REF 4103-133-000 only).
<b>B</b>	<b>Chuck</b> – Rotates to open and close the jaws.
<b>C</b>	<b>Jaws</b> – Secure the cutting tool in the attachment.
<b>D</b>	<b>J-Notch</b> – Secures the attachment in the handpiece.

## Drill and Reamer Attachments

<b>A</b>	<b>Retaining Collar</b> – Slide to allow installation of a cutting tool.
<b>B</b>	<b>J-Notch</b> – Secures the attachment in the handpiece.

## Definitions

The symbols located on the equipment and/or labeling are defined in this section and/or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

Symbol	Definition
	General Warning Sign
	Forward/Reverse Trigger (Action is in the direction of the arrow.)
	Forward/Reverse Trigger (Action is in the direction of the arrow.)
<b>F/R 1 min/ 4 min X 3</b>	Duty Cycle (See the <i>Specifications</i> section.)
	Lock
	Unlock
	Direction of Rotation
	Alignment Mark
<b>USE ONLY STRYKER POWER MODULE</b>	Use only Stryker power module.
<b>STERILIZE WITH DOOR FULLY OPEN</b>	Sterilize with door fully open.
<b>S</b>	Small Diameter
<b>L</b>	Large Diameter
	Run
	Load
<b>SERIAL</b>	Serial Number

## Instructions

### To Process the Equipment

**WARNING:** Upon initial receipt and before each use, process the equipment as indicated in the care instructions manual supplied with the equipment.

### To Charge a Power Module

Fully charge the power module before first and every use.


#### NOTES:

- For charging instructions, see the instructions for use supplied with the battery charger.
- For additional information, including safety information, see the instructions for use supplied with the power module.


## To Install a Power Module

### Sterile Assistant


1. Set the safety switch to the safe mode position.

2. Turn the latch to the unlock position and open the door.

3. Place the sterile transfer shield over the open compartment.


4. Present the compartment to the Circulating Assistant.


### Circulating Assistant

**⚠ WARNING: DO NOT** contaminate the handpiece.  
The transfer shield guards the sterile handpiece from contamination during power module installation.

1. Hold the power module by the handle.
2. Insert the power module through the transfer shield into the open compartment.


3. Push the handle down.

**⚠ WARNING:** ALWAYS sterilize the contaminated transfer shield before using it again.

4. Remove the contaminated transfer shield.


### Sterile Assistant


#### ⚠ WARNINGS:

- When closing the door, hold the door by the latch. DO NOT touch the edges of the door or the internal surface of the door that faces the power module.
- ALWAYS verify the door is securely closed to prevent exposing the power module to the sterile field.

1. Holding the door by the latch, use the door to push the power module handle down until it locks into position against the power module.


2. Close the door and turn the latch to the lock position. Verify the mark on the latch aligns with the lock symbol and the door is securely closed.

**NOTE:** Although one power module is usually sufficient for one operation, Stryker recommends preparing a second handpiece with power module to ensure quick intraoperative exchange of handpieces under sterile conditions.

### To Install an Attachment

1. Align the J-notch with the handpiece and insert the attachment until it snaps securely into place.

2. Gently pull on the attachment to verify it is secure.

## To Install and Remove a Cutting or Fixation Tool

### Wire and Pin Collets

1. Select the appropriate collet for the desired wire or pin diameter.

**NOTE:** All collets have a pre-adjust knob with small and large diameter settings.

#### Wire Collet (REF 7307-001-000)

Setting	Wire or Pin Diameter
Small (S)	0.7 to 1.2 mm
Large (L)	1.3 to 1.8 mm


#### Pin Collet, Small (REF 7307-002-000)

Setting	Wire or Pin Diameter
Small (S)	2.0 to 2.6 mm
Large (L)	2.7 to 3.2 mm


#### Pin Collet, Large (REF 7307-003-000)

Setting	Wire or Pin Diameter
Small (S)	3.0 to 3.6 mm
Large (L)	3.7 to 4.2 mm


2. Pull and rotate the pre-adjust knob to align the symbols for the desired wire or pin diameter. After alignment, snap the knob back into place.

3. Insert the wire or pin through the front of the collet or the back of the handpiece until the desired length extends from the collet.


4. Pull the clamping lever back to grasp the wire or pin.


5. To remove the wire or pin, or to adjust the wire length, release the clamping lever.

## Keyed Chucks


1. Use the chuck key to open the jaws.

2. Insert the tool into the center of the jaws.


3. Use the chuck key to close the jaws.
4. Visually inspect the tool to make sure it is centered in the jaws.


5. Gently pull on the tool to verify it is secure.
6. To remove the tool, use the chuck key to open the jaws again.

## Keyless Chucks

1. Hold the attachment and manually rotate the chuck to open the jaws.

2. Insert the tool into the center of the jaws.


3. Hold the attachment and manually rotate the chuck to close the jaws.
4. Visually inspect the tool to make sure it is centered in the jaws.

**WARNING:** DO NOT operate the handpiece to close the jaws on the cutting tool.


5. Hold the attachment and manually rotate the chuck to close the jaws.
6. Visually inspect the tool to make sure it is centered in the jaws.


5. Gently pull on the tool to verify it is secure.
6. To remove the tool, hold the attachment and manually rotate the chuck to open the jaws again.

## Drill and Reamer Attachments


1. Slide the retaining collar back.

2. Hold the retaining collar in position and insert the tool into the attachment.

3. Release the retaining collar.

4. Gently pull on the tool to verify it is secure.
5. To remove the tool, slide the retaining collar back again.

## To Operate a Handpiece

### ⚠ WARNINGS:

- DO NOT open the handpiece housing during surgery as the inside of the device is not sterile.
- ALWAYS follow the specified duty cycle to prevent the equipment from overheating. See the *Specifications* section.
- ALWAYS consult attachment and tool labeling for product-specific duty cycles and instructions for use.
- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen, or nitrous oxide.
- Take special precautions regarding electromagnetic compatibility (EMC) when using this equipment. Place the equipment into service according to the EMC information contained in this manual.
- Portable and mobile radio frequency (RF) communications equipment can affect the function of this equipment.
- This equipment should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, ALWAYS verify normal operation in the configuration in which it will be used.
- ALWAYS operate the equipment within the specified environmental condition values:

Temperature Limitation	Humidity Limitation	Atmospheric Pressure Limitation

1. Before operating the equipment, verify the following:
  - You are using Stryker-approved equipment.
  - The equipment is not damaged.
  - The equipment is clean and sterile.
  - You have properly assembled the equipment.
  - You have installed a fully-charged power module.

**NOTE:** If you experience difficulty during operation, contact your Stryker sales representative or call Stryker customer service. See the *Contact Information* section. Outside the US, contact your nearest Stryker subsidiary.

2. Set the safety switch to the forward or forward/reverse mode position.

**NOTE:** Set the safety switch back to the safe mode position when the handpiece is idle or when passing the handpiece to another person.


**CAUTION:** DO NOT apply excessive pressure or restrict the motion of the cutting or fixation tool during operation. An excessive load on the system may cause the power module to temporarily shut down. If this happens, release the trigger to resume operation or allow the power module to cool for 20 to 40 minutes.

3. Pull the triggers to control the variable speed operation of the handpiece.

**NOTE:** Typical power module run time is 5 minutes or longer, depending on surgical procedure and technique.

## To Remove an Attachment


1. Push the release lever down and pull the attachment from the handpiece.

## To Remove a Power Module

**NOTE:** Keep the power module as clean as possible because it cannot be sterilized.

1. Wipe all debris from the external surfaces of the handpiece.
2. Turn the latch to the unlock position and open the door.
3. Use the handle to remove the power module from the handpiece.

## Specifications

### Handpiece

<b>Model:</b>	STRYKER SYSTEM G Rotary Drill
<b>Catalog Number:</b>	7305-001-000
<b>Dimensions:</b>	
Length:	168.15 mm
Width:	83.82 mm
Height:	226.10 mm
<b>Mass:</b>	1.043 kg
<b>Nominal Speed:</b>	
Drill:	825 revolutions per minute
Ream:	250 revolutions per minute
<b>Mode of Operation:</b>	Non-continuous
Duty Cycle:	1 minute on/4 minutes off, 3 times
Rest Between Cycles:	4 hours
<b>Applied Part:</b>	The handpiece, attachment, and cutting or fixation tool as defined by the manufacturer
<b>Maximum Temperature of Applied Part:</b>	Less than 43 °C, as tested to the <i>Product Safety Compliance</i> standards at a maximum ambient temperature of 30 °C.
<b>Power Supply:</b>	Internally Powered
<b>Ingress Protection (IP):</b>	IPX0
<b>Equipment Type:</b>	 Type BF Applied Part
<b>European Conformity:</b>	 0197

## Attachments

### Wire and Pin Collets

<b>Model:</b>	STRYKER SYSTEM G Wire Collet	STRYKER SYSTEM G Pin Collet, Small	STRYKER SYSTEM G Pin Collet, Large
<b>Catalog Number:</b>	7307-001-000	7307-002-000	7307-003-000
<b>Dimensions:</b>			
Length:	88 mm	88 mm	88 mm
Diameter:	26 mm	26 mm	26 mm
<b>Mass:</b>	0.27 kg	0.27 kg	0.27 kg
<b>Gear Ratio:</b>	1:1	1:1	1:1
<b>Ingress Protection (IP):</b>	IPX0		
<b>European Conformity:</b>			

## Keyed Chucks

<b>Model:</b>	1:1 1/4 inch (6.4 mm) Keyed Chuck	High Torque 1/4 inch (6.35 mm) Drill with Chuck 2:1 with Key	3.25:1 1/4 inch (6.4 mm) Keyed Chuck
<b>Catalog Number:</b>	4103-131-000	4103-180-000	4103-231-000
<b>Dimensions:</b>			
Length:	105.90 mm	111.25 mm	118.62 mm
Diameter:	28.60 mm	33.80 mm	28.60 mm
<b>Mass:</b>	0.313 kg	0.198 kg	0.182 kg
<b>Gear Ratio:</b>	1:1	2:1	3.25:1
<b>Ingress Protection (IP):</b>	IPX0		
<b>European Conformity:</b>	CE		

## Keyless Chucks

<b>Model:</b>	Keyless Chuck Adaptor	Keyless Drill Attachment
<b>Catalog Number:</b>	4103-082-131	4103-133-000
<b>Dimensions:</b>		
Length:	87.12 mm	126.50 mm
Diameter:	30.07 mm	31.95 mm
<b>Mass:</b>	0.185 kg	0.157 kg
<b>Gear Ratio:</b>	1:1	1:1
<b>Ingress Protection (IP):</b>	IPX0	
<b>European Conformity:</b>	CE	

## Drill Attachments

<b>Model:</b>	1:1 AO Small Drill	1:1 Hudson	1:1 Hudson/Modified Trinkle	1:1 Trinkle
<b>Catalog Number:</b>	4103-110-000	4103-113-000	4103-135-000	4103-160-000
<b>Dimensions:</b>				
Length:	89.50 mm	84.33 mm	78.25 mm	79.85 mm
Diameter:	26.04 mm	26.04 mm	26.04 mm	26.04 mm
<b>Mass:</b>	0.185 kg	0.183 kg	0.157 kg	0.154 kg
<b>Gear Ratio:</b>	1:1	1:1	1:1	1:1
<b>Ingress Protection (IP):</b>	IPX0			
<b>European Conformity:</b>	CE			

## Reamer Attachments

<b>Model:</b>	3.25:1 AO Large Reamer	3.25:1 Hudson	3.25:1 Hudson/Modified Trinkle	3.25:1 Trinkle	4:1 Hudson	4:1 Hudson/Modified Trinkle
<b>Catalog Number:</b>	4103-210-000	4103-213-000	4103-235-000	4103-260-000	4103-413-000	4103-435-000
<b>Dimensions:</b>						
Length:	95.45 mm	100.46 mm	94.00 mm	97.30 mm	100.46 mm	94.00 mm
Diameter:	27.95 mm	27.95 mm	27.95 mm	27.70 mm	32.40 mm	32.40 mm
<b>Mass:</b>	0.198 kg	0.256 kg	0.182 kg	0.180 kg	0.225 kg	0.208 kg
<b>Gear Ratio:</b>	3.25:1	3.25:1	3.25:1	3.25:1	4:1	4:1
<b>Ingress Protection (IP):</b>	IPX0					
<b>European Conformity:</b>						

## Transfer Shield

<b>Model:</b>	STRYKER SYSTEM G Transfer Shield
<b>Catalog Number:</b>	7304-120-000
<b>Dimensions:</b>	
Length:	150.9 mm
Width:	103.5 mm
Height:	37.0 mm
<b>Mass:</b>	77 g
<b>European Conformity:</b>	

## Product Safety Compliance

### International Electrotechnical Commission (IEC)

IEC 60601-1 A1: R:2012 (3.1 Edition): *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*; IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); Amendment 1 (2012)

IEC 60601-1-6: 2013: *USABILITY – Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability* – Edition 3.1 Consolidated Reprint.

IEC 62366: 2014: *USABILITY – Medical devices – Application of usability engineering to medical devices* - Edition 1.1 Consolidated Reprint.

IEC 60601-1-2: 2010: *EMC (Version 3.0) – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests* – Edition 3.0

IEC 62133: 2012/12/06 Edition: 2.0: *Secondary Cells and Batteries Containing Alkaline or other Non-Acid Electrolytes – Safety Requirements for Portable Sealed Secondary Cells and for Batteries Made from Them for use in Portable Applications*, Corrigendum 1:2013

### European Committee for Electrotechnical Standardization (CENELEC)

EN 60601-1:2006/A1:2013 (3.1 Edition): *Medical electrical equipment Part 1: General requirements for basic safety and essential performance* – Incorporates Amendment A12: 2014

## Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions		
The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration - electromagnetic immunity			
The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A N/A	N/A
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	N/A N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0,5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 seconds	N/A N/A N/A N/A	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

NOTE:  $U_T$  is the alternating current mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity			
The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment specified below. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance: $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ 150 kHz to 80 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz  Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:   (Non-ionizing electromagnetic radiation)
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	[ $E_1$ ] = 3 V/m	

NOTE 1: At 80 MHz and 800MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is used exceeds the applicable RF compliance level above, the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000).

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000)			
The STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the STRYKER SYSTEM G Rotary Drill (REF 7305-001-000) as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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IEC 60601-1-6: 2013: *USABILITY – Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability – Edition 3.1 Consolidated Reprint.*

IEC 62366: 2014: *USABILITY – Medical devices – Application of usability engineering to medical devices - Edition 1.1 Consolidated Reprint.*

IEC 60601-1-2: 2010: *EMC (Version 3.0) – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests – Edition 3.0*

IEC 62133: 2012/12/06 Edition: 2.0: *Secondary Cells and Batteries Containing Alkaline or other Non-Acid Electrolytes – Safety Requirements for Portable Sealed Secondary Cells and for Batteries Made from Them for use in Portable Applications, Corrigendum 1:2013*

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EN 60601-1:2006/A1:2013 (3.1 Edition): *Medical electrical equipment Part 1: General requirements for basic safety and essential performance – Incorporates Amendment A12: 2014*

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