**Cartiva® Synthetic Cartilage Implant (SCI) for CMC Instructions for Use**

**DESCRIPTION**

Cartiva® Synthetic Cartilage Implant (SCI) for CMC is intended to treat diseased or damaged articular surfaces of the 1st carpometacarpal (CMC) joint associated with joint pain or decreased range of motion. The implant, a cylindrical device made from an elastic biomaterial, may be used as a replacement for damaged cartilage and bone without requiring the destruction or removal of a patient's healthy tissue. It is intended for use during a single surgical procedure. The procedure is similar to that used for osteochondral autograft or allograft transplantation; a part is placed into a pre-drilled hole to resurface the damaged area of cartilage/bone.

Cartiva SCI for CMC is made from a proprietary biomaterial. The device, which is classified as a hydrated polymer, consists of water in similar proportion to human tissue. This organic polymer-based biomaterial is capable of withstanding repetitive loading typical of normal hand use, and its mechanical properties are similar to articular cartilage.

Cartiva SCI for CMC is supplied in one size. The device is supplied sterile and is packaged as a single unit.

**INDICATIONS FOR USE**

Cartiva SCI for CMC is indicated for use in treatment of diseased or damaged articular surface of the first carpometacarpal joint.

**CONTRAINDICATIONS**

Cartiva SCI for CMC is not designed, sold, or intended for use other than as indicated.

Use of the device is contraindicated for patients with the following:

- Any significant bone loss, avascular necrosis, or subchondral bone cyst > 8mm of the supporting bone structure;
- Inadequate cortical bone stock required to support placement of the implant;
- Inflammatory arthropathy and/or diagnosis of gout;
- Physical conditions that would tend to eliminate adequate implant support (e.g., insufficient quality or quantity of bone resulting from cancer, congenital dislocation, or osteoporosis), systemic and metabolic disorders leading to progressive deterioration of bone (e.g., cortisone therapies, immunosuppressive therapies, neuromuscular compromise, vascular deficiency in the affected limb, absence of musculoligamentous supporting structures, and joint neuropathy), tumors and/or cysts of the supporting bone structures, arthritis of the scaphotrapezotrapezoidal (STT) joint, or joints previously treated with trapeziectomy;
- Known or suspected allergic reaction to polyvinyl alcohol;
- Active infection of the surgery site.
- Patient is on chronic anticoagulation due to a bleeding disorder or has taken anticoagulants within 5 days prior to surgery.

**WARNINGS AND PRECAUTIONS**

It is the responsibility of each surgeon utilizing the Cartiva SCI for CMC device to consider the clinical and medical status of each patient and to be knowledgeable about all aspects of implant procedures and the potential complications that may occur. The benefits derived from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon. Surgeons must balance many considerations to achieve the best result in individual patients.

Clinical results depend on physician and technique, preoperative and postoperative care, the implant, patient pathology, and daily activity. It is suggested that surgeons use appropriate informed consent and discuss the potential for complications and the use of alternative techniques with each patient scheduled for surgery.

The Cartiva implant is a single-use only device. Each device should be used in one patient, during one procedure, and in only one implant site. Re-use of a previously implanted device is strictly prohibited. Material properties required for implant duration and longevity cannot be assured by the manufacturer if the device is re-used, and the potential cross-contamination between implant sites and/or patients poses a serious health risk.

Cartiva SCI for CMC has been sterilized. The implant is not compatible with gas or steam (autoclave) sterilization. **DO NOT RESTERILIZE THE IMPLANT.**

The implant is not compatible with storage or shipment temperatures in excess of 49°C (120°F). If the temperature-sensitive indicator on the container has turned dark gray to black, **DO NOT USE THE IMPLANT.**

Physical activity should be resumed according to the rehabilitation plan recommended by the physician. Vigorous activity may compromise the durability of clinical benefit from Cartiva® SCI for CMC.

**Product Limitations**

The long-term safety and effectiveness of cartilage replacement are unknown.

Safety and effectiveness of Cartiva SCI for CMC in the following populations have not been established:

- Pediatric population;
- Pregnant women;
- Patients with either systemic or local infections;
- Patients with a diagnosis of concomitant injury that the physician believes may interfere with healing;
- Patients with clinically significant (as defined by the surgeon) renal, hepatic, cardiac, endocrine, hematologic, autoimmune or any systemic disease which may make interpretation of the results difficult;
- Patients who have undergone systemic administration within 30 days prior to implantation of any type of corticosteroid in the thumb, antineoplastic, immunostimulating or immunosuppressive agents;
- Patients with evidence of osteonecrosis of the involved joint, or history of peripheral neuropathy, active or on-going neoplastic disease, or immunosuppression.

**Product Handling**

To avoid dehydration, maintain the Cartiva SCI for CMC in sterile saline until ready to implant.

**POTENTIAL COMPLICATIONS AND ADVERSE REACTIONS**

The following potential adverse effects may occur in association with Cartiva implant:

- An implant site may become infected, painful, swollen, or inflamed;
- Continued prolonged joint pain or recurrence of joint pain;
- Medial impingement between the first metacarpal and trapezium
- Risks associated with implantation of hemi-arthroplasty devices or Cartiva Synthetic Cartilage Implant include bony reaction, implant loosening, implant dislocation/dislodgement, cysts, bony erosion, and device wear.

Any of these adverse effects may require medical or surgical intervention.

Undesirable implant performance may be caused by:

- Inappropriate patient selection;
- Leaving an insufficient cortical bone stock around implant;
- Incorrect drill hole size for implant;
- Incorrect depth of implant placement;

- Uncorrected joint malalignment;
- Uncorrected joint instability;
- Unrecognized or uncorrected cartilage or bony pathology, including osteophytes;
- Excessive early post-operative load bearing.

Some preventive measures to consider in order to minimize the potential for complications:

- Follow accepted patient selection guidelines;
- Identify and correct co-existing pathology;
- Ensure appropriate implant placement;
- Sufficient joint capsule repair.

If complications develop, possible corrective procedures include:

- Replacement of implant(s);
- Implant removal or revision;
- Autograft or allografting of the implant site;
- Partial or complete trapeziectomy;
- Other subsequent secondary surgical intervention(s).

**PRIOR TO USE**

1. Before presentation to the operative field, inspect the package to ensure sterility has not been breached during transportation. If the package is obviously damaged or if there is a question of sterility for any other reason, **DO NOT USE THE IMPLANT.**
2. Before opening the carton, ensure that the temperature-sensitive indicator on the carton is light gray. The implant is not compatible with storage or shipment temperatures in excess of 49°C (120°F). If the temperature-sensitive indicator has turned dark gray to black, **DO NOT USE THE IMPLANT.**
3. The contents of the pouch, including the tray and implant, are sterile. The carton and exterior of the container are not sterile. **DO NOT PRESENT THE INTACT CARTON OR POUCH TO THE OPERATIVE FIELD.**
4. Inspect each Cartiva® SCI implant to ensure that it is not hard, brittle, torn, or otherwise defective. If any of these problems is detected, **DO NOT USE THE IMPLANT.**

**DIRECTIONS FOR USE**

Implantation of Cartiva SCI for CMC has been validated for use with surgical instrumentation distributed by Cartiva, Inc. Cartiva SCI Instrumentation for CMC is recommended for proper placement of the implant and is designed for use with a cannulated drill and 2 mm diameter non-threaded guide pin (minimum length 200 mm).

**General directions for carpometacarpal joint (CMC) implantation are as follows:**

The following table references optimal dimensions for successful implantation of Cartiva SCI for CMC implants slightly proud (~3mm) with the surrounding cartilage:

Implant Size	Catalog Number	Drill Part	Hole Diameter	Hole Depth
8 mm	CMC-08	DMC-08	7.7 mm	7.0 mm

1. Access the affected joint using standard surgical technique and expose the joint to gain access to the center of the base of the metacarpal. Ensure the incisions are conducive to capsular repair. During dissection, the abductor pollicis longus tendon should be protected and retracted radially while the metacarpal base is exposed.
  2. Confirm the fit of the 8 mm implant on the metacarpal using the concave end of the Placer on the base of the metacarpal.
- Warning: No more than one Cartiva SCI for CMC implant should be used in one articular area.**
3. Use the concave end of the Placer to create a perpendicular angle to the articular surface of the metacarpal base and to identify the target implantation site. Cartiva requires a minimum of 2mm of surrounding bone stock. Insert the Placement Guide Pin into the center of the metacarpal to a depth of approximately 3 cm ensuring it is perpendicular to the articular surface. The Drill Bit (DMC-08) will be used to drill a hole into the subchondral bone to the proper depth (recommendations can be found in the table above). Insert the Drill Bit over the Placement Guide Pin and advance the drill until the stop reaches the level of the adjacent tissue, to ensure appropriate depth is achieved. Care should be taken to advance only to the drill stop using light pressure and to irrigate while drilling. Note the Placement Guide Pins are single use.
  4. If necessary, remove any cartilage and/or bone debris from the recipient implant site. Advance the Placer over the Placement Guide Pin and confirm the depth of the cavity is approximately 7 mm in the area of the highest cartilage topography using the Placer's marking lines. Remove the Placement Guide Pin.
  5. Remove implant from packaging. Moisten the inner walls of the Introducer tube with sterile saline. Place Cartiva implant into the proximal (wide) end of the Introducer with the flat side of the implant facing distally (saddle-shaped side facing proximally) so that the flat side of the implant will be placed in the bottom of the recipient implant site. Align the highest points of the device with the markings on the Introducer. Insert the T-Handle into the proximal end of the Introducer matching its surface geometry with that of the implant taking care to maintain this alignment throughout the device insertion step. The markings on the Introducer and T-Handle should align. Rest the distal end of the Introducer on a flat, non-shedding, sterile surface and slowly advance the implant to the distal end of the Introducer using the T-Handle. Do not advance the flat side of the implant beyond the distal end of the Introducer until ready to deliver the implant to the recipient site. Place the distal end of the Introducer into the recipient implant site aligning the markings with the metacarpal topography. Advance the Cartiva® implant into the implant site using the T-Handle. Continue to apply steady pressure to the T-Handle until the flat side of the implant reaches the bottom of the recipient site. Remove the Introducer and T-Handle.
  6. Confirm that the implant is seated firmly in the base of the recipient site, is tight within it and the topography of the device is aligned with that of the surrounding cartilage. The implant should be slightly proud (~3 mm) in the implant site.
  7. Repair the joint capsule.

**Cartiva SCI Non-sterile Instrumentation for CMC Instructions for Use**

**INSTRUMENT DESCRIPTION**

This document pertains to the following instruments supplied by Cartiva, Inc.:

**Non-Sterile Instrumentation**

Part Description	Instrumentation REF	Classification
Drill Bit (Fabricated from 455 H900 or 17-4 H900 Stainless Steel) <b>Note: The Drill Bits (DMC-08) are compatible with drills having a chuck size of at least 0.25" or 6 mm.</b>	DMC-08	Reusable

Part Description	Instrumentation REF	Classification
Introducer (Fabricated from 17-4 H900 Stainless Steel)	IMC-08	Reusable
T-Handle (Fabricated from 17-4 H900 Stainless Steel)	TMC-08	Reusable
Placer (Fabricated from 17-4 H900 Stainless Steel)	PMC-08	Reusable
Non-threaded Placement Guide Pin, 2 mm (Fabricated from 316L Stainless Steel) <b>Note: Different Guide Pins may be provided with the instrumentation in the Australian market</b>	PNN-02	Single Use Only

#### INDICATIONS FOR USE

Instrumentation supplied by Cartiva, Inc is indicated for use with Cartiva for CMC implantable medical device products manufactured by Cartiva, Inc.

#### CONTRAINDICATIONS,

Instrumentation supplied by Cartiva, Inc. is not designed, sold or intended for use other than as indicated.

#### WARNINGS AND PRECAUTIONS

- Failure to properly clean instruments prior to sterilization may lead to inadequate sterilization.
- Surgical instruments are used with or on patients who may harbor both recognized and unrecognized infections. To prevent the spread of infection, all reusable instruments must be thoroughly cleaned and sterilized prior to initial use and after each patient use.
- Instruments may have sharp edges or features. Users and reprocessors must be cautious when handling instruments.

#### LIMITATIONS ON REPROCESSING

Repeated processing, according to these instructions, has minimal effect on and should not compromise the performance of reusable Cartiva CMC instruments. End of life is normally determined by wear and damage due to use.

In addition to the Cartiva for CMC Instrumentation that is labeled for re-use, Cartiva, Inc. provides single-use Placement Guide Pins in some markets for use during the Cartiva for CMC implantation procedure. Re-use of the single-use Placement Guide Pins is strictly prohibited. The material properties and reliability of these devices in a multi-use scenario have not been explicitly tested or demonstrated. Re-use of a single-use Placement Guide Pin could result in improper device placement (depth, alignment, etc.) and undesired clinical outcomes. Placement Guide Pins must be discarded after one use.

#### DAMAGE INSPECTION

- Inspect the instruments for damage, wear, and corrosion at all stages of handling.
- Cutting edges should be free of nicks and present a continuous edge.
- Check instruments with long slender features for distortion.
- If damage is detected, do not use instrument but consult Cartiva, Inc. for guidance.

#### SPECIAL PRECAUTIONS

##### Transmissible Spongiform Encephalopathy Agents

- It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy agents.
- The agents for transmission of Creutzfeldt-Jakob disease (CJD) are believed to be resistant to normal processing methods of disinfection and sterilization. Methods of decontamination and sterilization outlined below may not be appropriate where CJD transmission is a risk.
- Refer to the World Health Organization guidelines for a detailed listing of appropriate decontamination methods.

#### PACKAGING

The Cartiva SCI reusable sterilization tray and associated reusable surgical instruments are supplied non-sterile and must be cleaned and sterilized prior to use according to the instructions in this document. The reusable instruments and tray are shipped and stored in packaging that is labeled according to its contents. Store the sterilization tray in normal hospital environmental conditions. Store the instruments in the original packaging. Do not remove a reusable instrument from the packaging until it is ready to be placed in the sterilization tray.

#### CLEANING & STERILIZATION INSTRUCTIONS

The instructions provided below have been validated by Cartiva, Inc. as being capable of preparing non-sterile instruments for initial use or re-use. It is the responsibility of the reprocessor to ensure that reprocessing is performed using the appropriate equipment and materials, and that personnel in the reprocessing facility have been adequately trained in order to achieve the desired result. This normally requires validation and routine monitoring of the process. Any deviation by the reprocessor from these instructions should be properly evaluated for effectiveness to avoid potential adverse consequences.

#### Post-use

- Remove excess soil with disposable non-shedding wipe.
- Instruments should be covered with a damp cloth to prevent drying of soil prior to cleaning.

#### Containment and Transportation

- Observe universal precautions for handling contaminated/biohazardous materials.
- Instruments should be cleaned within 30 minutes of use to minimize the potential for drying prior to cleaning.

#### Cleaning

##### Preparation for Cleaning

- No assembly/disassembly of Cartiva instruments is required.
- Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and must not be used. Acidic cleaning agents should be avoided.
- Saline solution has a corrosive effect on stainless steel and should not be used to rinse, soak, or clean instruments.

##### Manual Cleaning Instructions

- For initial and subsequent uses, follow all cleaning and sterilization instructions.
- Prepare a neutral pH or nearly neutral pH enzymatic detergent at the use-dilution and temperature recommended by the agent's manufacturer.
- Submerge the instruments in enzymatic detergent and soak for 20 minutes.
- While submerged in enzymatic detergent, scrub each instrument with a soft-bristled brush, paying special attention to areas where debris might accumulate. Lumens and crevices should be cleaned with a long, narrow, soft-bristled brush. Avoid any harsh materials or cleaning motions that can scratch the surface of the instruments.
- Remove the instruments from the enzymatic detergent and rinse each instrument thoroughly in purified water (such as distilled or deionized water) for a minimum of 3 minutes. Thoroughly flush lumens and other difficult to reach areas.
- Sonicate instruments for a minimum of 10 minutes in an ultrasonic cleaner containing fresh enzymatic detergent, preferably at 45-50 kHz (according to the ultrasonic unit's directions).

- Remove the instruments from the enzymatic detergent and rinse each instrument thoroughly with purified water (such as distilled or deionized water) for at least 3 minutes and until there is no sign of soil in the rinse stream. Thoroughly flush lumens and other difficult to reach areas.

#### Verifying Cleaning

- Check instruments for visible soil. All exterior surfaces as well as inner lumens should be inspected.
- Repeat cleaning if soil is visible and re-inspect.

#### Disinfection

- Disinfection is only acceptable as an adjunct to full sterilization for reusable surgical instruments.
- Disinfectant solution Cidex® or equivalent may be used in accordance with solution supplier's instructions. Instruments should be thoroughly rinsed in purified water (such as distilled or deionized water) following disinfection.

#### Drying

- Instruments with inner lumens should be agitated or positioned so that liquid inside the lumens may drain.
- Dry the exterior of the instruments with a clean, disposable, non-shedding wipe.

#### Automated Cleaning

##### Manual Pre-Cleaning

- Prepare an enzymatic cleaner per manufacturer's directions. Soak all instruments in the enzymatic cleaning solution for a minimum of five (5) minutes. After five (5) minutes of soak time, remove gross soil using a soft bristled brush, paying special attention to areas where debris might accumulate. Lumens and crevices should be cleaned with a long, narrow, soft-bristled brush. Flush the inside of lumens using a syringe or pipette filled with the prepared cleaning solution for a minimum of three (3) times.
- Rinse each soiled instrument under running cold tap water for a minimum of one (1) minute. Thoroughly flush lumens and other difficult to reach areas.

##### Automated Process

- Load the instruments in a washer disinfectant so that the lumens can drain.
- Using a validated washer disinfectant and an alkaline cleaning agent (8-11 pH range) intended for use in an automated cleaning process, use the minimum cycle parameter set points below.

Cycle	Time (minutes:seconds)	Minimum Temperature	Type of Detergent/ Water
Pre-Cleaning	2:00	Cold	Tap
Cleaning	10:00	Heated 50°C - 60°C	Alkaline detergent (per the manufacturer's instructions) and tap water
Rinse 1	2:00	Cold	Tap
Rinse 2	1:00	Cold	Purified Water*
Dry	25:00	Heated (90°C)	N/A

\* Purified water for the final rinse may be distilled, deionized, or reverse osmosis water

- Visually inspect the instruments for remaining soil in a well-lit area; no visible soil should be left on any instrument surfaces.
- If instruments are still wet after automated cleaning cycle, thoroughly dry the instruments using a clean lint-free cloth. If needed, use filtered pressurized air to aid in drying.
- The instruments withstand steam exposure up to 137°C and, thus, may be thermally disinfected.

#### Cleaning of Reusable Tray

The instrument sterilization tray provided with reusable instrumentation should be cleaned, sterilized and inspected prior use in accordance with the tray's Instructions for Use.

#### Warnings

- Do not stack trays on top of one another. Be sure that ventilation holes are not obstructed, and that mats are correctly installed. For effective sterilization trays must have adequate steam circulation around all surfaces. They must also be placed upright on shelves for proper ventilation. Condensation can pool on non-absorbent surfaces. Do not place trays on their sides or at vertical angles in chamber, in order to ensure that proper drainage can occur during the cycle.
- Small baskets, trays, or other accessories with covers or lids should only be used in trays specifically designed and labeled for the purpose. Do not overload trays. Overloading may inhibit steam flow, cause excessive drying times, and make trays too heavy to safely handle. Load and sterilize instruments in trays in accordance with the instructions provided within this IFU.

#### STERILIZATION

##### Packaging for Sterilization

Instruments may be loaded into dedicated instrument trays or general-purpose sterilization trays. The maximum load configuration, regardless of instrument size, is as follows:

- 1 x Drill Bit (DMC-08)
- 1 x Introducer (IMC-08)
- 1 x T-Handle (TMC-08)
- 1 x Placer (PMC-08)
- 3 x Placement Guide Pins (PNN-02)

Use standard medical-grade steam sterilization wrap to double-wrap the tray.

##### Recommended Sterilization Parameters

- Steam sterilize using one of the three following steam cycles. Each has been found to demonstrate a sterilization assurance level (SAL) of 10<sup>-6</sup> for the maximum load configurations described above (AAMI TIR12:2010):

Cycle Number	1	2	3
Autoclave Type	Gravity	Pre-Vacuum	Pre-Vacuum
Sterilization Temperature	270°F / 132°C (+5°F / +3°C)	270°F / 132°C (+5°F / +3°C)	273°F / 134°C (+5°F / +3°C)
Exposure Time	25 minutes	4 minutes	3 minutes
Minimum Drying Time	30 minutes	20 minutes	20 minutes

- Sterilizers vary in design and performance characteristics, so cycle parameters should be verified against the sterilizer manufacturer's instructions for the specific sterilizer and load configuration being used.
- When sterilizing multiple instruments in one steam sterilization cycle, ensure that the sterilizer manufacturer's maximum load is not exceeded.
- Drying time may vary according to load size (larger loads require longer drying times).
- Instruments must be adequately cooled after removal from the sterilizer. Do not touch instruments during the cooling process.

## Storage

- Sterilized, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.
- Sterilized instrument packages should be examined closely prior to opening to ensure that there has been no loss of package integrity.

## Maintenance

- No lubrication is needed.
- Discard blunt, damaged, severely corroded or severely discolored instruments.

## Complaints

Any health care professional (e.g., customer or user of this system), who has a complaint or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Cartiva, Inc. USA. Further, if any of the implanted system ever "malfunctions" (i.e. does not meet any of its performance specifications or otherwise does not perform as intended) or may have caused or contributed to the death or serious injury of a patient, Cartiva, Inc. should be notified immediately by telephone, fax or written correspondence. When filing a complaint, please provide the device size, part number, lot number(s), your name and address, and the nature of the complaint. For US customers, complaints may also be reported directly to Medwatch at <http://www.fda.gov/medwatch>. For EU customers, any serious incident that has occurred in relation to the device should be reported to Cartiva, Inc. USA and the competent authority of the Member State in which the user and/or patient is established.

For the most current Cartiva® SCI product information, please visit our website at [www.cartiva.net](http://www.cartiva.net).

## REFERENCES

1. AAMI TIR12:2010. Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers
2. ISO 17664: 2004. Sterilization of medical devices — Information to be provided by the manufacturer for the processing of resterilizable medical devices
3. WHO/DCS/CSR/APH/2000.3, WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies, March 1999.

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Cartiva® is a registered trademark of Cartiva, Inc.

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MDSS GmbH
Tornier SAS – 161 rue Lavoisier, 38330
Stryker Australia Pty Ltd
Cartiva® est une marque déposée de Cartiva, Inc.