CE
0459

Stryker Spine
2 Pearl Ct. Allendale
NJ 07401-1677 USA

Manufactured by :
Stryker Spine SAS
Z.I Marticot 33610 CESTAS - FRANCE

Stryker European Operations B.V.
Herikerbergweg 110
Amsterdam 1101 CM Netherlands

NOLI115MV01 Rev.05 12/2016

stryker®

Spine

NOLI115MV01REV05

NOLI115MV01REV05

**IMPORTANT PRODUCT INFORMATION FOR
STRYKER SPINE VERTEBRAL BODY REPLACEMENT SYSTEMS
NON STERILE PRODUCT
and STERILE PRODUCT**

GB

Description

The Stryker Spine Vertebral Body Support System (VBOSS®) is intended for use as an aid in spinal fusion and consists of a hollow cylindrical tube. The sides of the cylinder are perforated by equally spaced round holes. The cylinder is segmented and the grooves can be used as cutting lines. The Stryker Spine Vertebral Body Support System implants are available in a variety of diameters from Ø10 to 25 mm and lengths from 8 to 120 mm. Ø10 mm is not available in the USA.

The end caps snap into each end of the Stryker Spine Vertebral Body Support System cage. The exterior side of the end cap features evenly spaced round spikes providing fixation. The end cap is available in round, oval, round angled and oval angled shapes depending on the diameter of the Stryker Spine Vertebral Body Support System. The end caps are also offered in different heights.

The VLIFT® Distractable In-situ (DIS) Vertebral Body Replacement system is intended for use as an aid in spinal fusion and consists of a single, pre-assembled implant, with a distractible center. The hollow cylindrical center may be adjusted continuously via an outer concentric ring. The hollow core of the cage allows for packing bone graft. The use of bone graft with VLIFT® is optional. As the implant is distracted via its outer concentric ring, additional, slotted openings appear. The VLIFT® implants are available in diameters Ø18 and Ø22 mm and cover a minimum range of heights from 20.5 to 60.5 mm

Extensions are available in diameters Ø18 and Ø22 mm and allow the implant height to increase 15mm or 30 mm, depending on whether one or two extensions are used respectively.

Similar to VBOSS®, modular end caps snap into each end of the VLIFT® implant. Spikes at the rim of the end cap attachments improve the anchorage of the implant to the end plates of the vertebral bodies. The end caps are round and are available in a variety of angles.

Material

All components of the VBOSS® system are manufactured out of commercially pure titanium (ASTM F67).

All components of the VLIFT® system are manufactured out of Titanium Alloy: Ti6Al4V according to ISO 5832-3 and ASTM F-136.

Indications

For USA only:

Stryker Spine Vertebral Body Replacement System (VBOSS®)

The Stryker Spine Vertebral Body Support System implant is a device intended to replace a vertebral body or an entire vertebra. It is for use in the thoracolumbar spine (T1-L5) to replace a collapsed, damaged, or unstable vertebral body or vertebra due to tumor or trauma (i.e., fracture). For both corpectomy and vertebrectomy procedures, the Stryker Spine Vertebral Body Support System is intended to be used with supplemental internal fixation systems. The use of bone graft with the Stryker Spine Vertebral Body Support System is optional.

Stryker Spine VLIFT® system

The VLIFT® system is a device intended to replace a vertebral body or an entire vertebra. It is for use in the thoracolumbar spine (T1-L5) to replace a collapsed, damaged, or unstable vertebral body or vertebra resected or excised during total and partial corpectomy and vertebrectomy procedures due to tumor or trauma (i.e., fracture). For both, corpectomy and vertebrectomy procedures, VLIFT® system is intended to be used with supplemental internal fixation systems. The supplemental internal fixation systems that may be used with VLIFT® system include, but are not limited to, Stryker Spine plate or rod systems (Xia® Spinal System, Spiral Radius 900® and Trio®). The use of bone graft with the VLIFT® system is optional.

Outside the USA:

Stryker Spine Vertebral Body Replacement System (VBOSS®)

The Stryker Spine Vertebral Body Support System implant is a device intended to replace a vertebral body or an entire vertebra. It is for use in all levels of spine (C1-L5) to replace a collapsed, damaged, or unstable vertebral body or vertebra due to tumor or trauma (i.e., fracture). For both corpectomy and vertebrectomy procedures, the Stryker Spine Vertebral Body Support System is intended to be used with supplemental internal fixation systems. The use of bone graft with the Stryker Spine Vertebral Body Support System is optional.

Stryker Spine VLIFT® system

The VLIFT® system is a device intended to replace a vertebral body or an entire vertebra. It is for use in the thoracolumbar spine (T1-L5) to replace a collapsed, damaged, or unstable vertebral body or vertebra resected or excised during total and partial corpectomy and vertebrectomy procedures due to tumor or trauma (i.e., fracture). For both, corpectomy and vertebrectomy procedures, VLIFT® system is intended to be used with supplemental internal fixation systems. The supplemental internal fixation systems that may be used with VLIFT® system include, but are not limited to, Stryker Spine plate or rod systems (Xia® Spinal System, Spiral Radius 90DTM and Trio®). The use of bone graft with the VLIFT® system is optional.

Warnings

The VLIFT® and VBOSS® implant systems have not been evaluated for safety and compatibility in the MR environment. The VLIFT® and VBOSS® implant systems have not been tested for heating or migration in the MR environment.

General conditions of use

Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and limitations of the spinal device. Knowledge of surgical techniques, proper reduction, selection and placement of implants, and pre- and post-operative patient management are considerations essential to a successful surgical outcome. Consult the medical literature for information regarding proper surgical techniques, precautions, and potential adverse effects associated with spinal fixation surgery.

The components of the system have been designed to work together. Do not substitute another manufacturer's device for any component of the System. Any such use will negate the responsibility of Stryker Spine for the performance of the resulting mixed component implant.

Do not mix metals (e.g. Titanium based devices with stainless steel items). Some corrosion occurs on all implanted metals and alloys. Contact of dissimilar metals, however, may accelerate corrosion. Corrosion may accelerate fatigue fracture of implants, and cause metal compounds to be released into the body.

Reuse

Never reuse or reimplant spinal surgical implants. These could become contaminated resulting in infection. In addition, even though the device appears undamaged, it may have small defects which could compromise structural integrity reducing its service life and/or leading to patient injury.

Surgeons must verify that the instruments are in good condition and operating order prior to use during surgery.

Contra-indications

- The VBOSS® and VLIFT® systems are not to be used for interbody fusion.
- Do not cut any VBOSS® cage to a segment that is less than 10mm in length.
- The VBOSS® and VLIFT® systems should not be implanted in patients with an active infection at the operative site.
- These devices are not intended for use except as indicated.
- Marked local inflammation.
- Any mental or neuromuscular disorder, which would create an unacceptable risk of fixation failure or complications in postoperative care.
- Bone stock compromised by disease, infection or prior implantation, which cannot provide adequate support and/or fixation to the devices.
- Open wounds.
- Rapid joint disease, bone absorption, osteopenia, osteomalacia, and/or osteoporosis. Osteoporosis or osteopenia are relative contraindications, since this condition may limit the degree of obtainable correction and/or the amount of mechanical fixation.
- Metal sensitivity, documented or suspected.
- Pregnancy.
- Anytime implant utilization would interfere with anatomical structures or physiological performance.
- Inadequate tissue coverage over the operative site.

Other medical or surgical conditions that could preclude the potential benefit of surgery, such as congenital abnormalities, immunosuppressive disease, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or marked left shift in the WBC differential count, must be carefully analyzed before surgery. These contra-indications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive. Surgeons must discuss all relative contra-indications including limited life time of the device when deemed appropriate.

Pre-operative precautions

The surgical indication and the choice of implants must take into account certain important criteria such as:

- Patients involved in an occupation or activity that applies excessive loading upon the implant (e.g., substantial walking, running, lifting, or muscle strain) may be at increased risk for failure of the fusion and/or the device.
- Patients should be instructed in detail about the limitations of the implants, including, but not limited to, the impact of excessive loading through patient weight or activity, and be taught to govern their activities accordingly. The procedure will not restore function to the level expected with a normal, healthy spine, and the patient should not have unrealistic functional expectations.
- A condition of senility, mental illness, chemical dependence or alcoholism. These conditions among others may cause the patients to ignore certain necessary limitations and precautions in the use of the implant, leading to failure and other complications.
- Foreign body sensitivity. Where material sensitivity is suspected appropriate tests should be made prior to material implantation.
- Patients who smoke have been shown to have an increased incidence of non-unions. Such patients should be advised of this fact and warned of the potential consequences.
- Care must be taken to protect the components from being marred, nicked, or notched as a result of contact with metal or abrasive objects. Alterations will produce defects in surface finish and internal stresses which may become the focal point for eventual breakage of the implant.

Surgeons must discuss these precautions pre-operatively with their patients when appropriate.

The choice of implants

- The choice of proper shape, size and design of the implant for each patient is crucial to the success of the surgery. The surgeon is responsible for this choice, which depends on each patient.
- The size and shape of the bone structures determine the size, shape and type of the implants. Once implanted, the implants are subjected to stresses and strains. These repeated stresses on the implants should be taken into consideration by the surgeon at the time of the choice of the implant, during implantation as well as in the post-operative follow-up period. Indeed, the stresses and strains on the implants may cause metal fatigue or fracture or deformation of the implants, before the bone graft has become completely consolidated. This may result in further side effects or necessitate the early removal of the osteosynthesis device.

Intra-operative precautions

- The insertion of the implants must be carried out using instruments designed and provided for this purpose and in accordance with the specific implantation instructions for each implant. Those detailed instructions are provided in the surgical technique brochure supplied by Stryker Spine.
- Discard all damaged or mishandled implants.
- Stryker Spine implants must not be reshaped, unless otherwise indicated in the surgical technique instructions. When VBOSS® implants need to be cut, the cutting must be carried out using the appropriate instruments, provided by Stryker Spine. The use of inappropriate instruments may result in scratches, notches, and sharp bending, causing the breakage of the implants. Improper seating of the implant may result in implant failure.
- Never reuse an implant, even though it may appear undamaged.
- Do not mix metals.

Post-operative precautions

- Physician instructions regarding full weight-bearing activities must be complied with until maturation of the fusion mass is confirmed. Failure to comply with physician instructions may result in failure of the implant, the fusion, or both.

Side effects

Include but are not limited to:

- Late bone fusion or no visible fusion mass and pseudarthrosis;
- Peripheral neuropathies, nerve damage, heterotopic bone formation and neurovascular compromise, including paralysis;
- Superficial or deep-set infection and inflammatory phenomena;
- Allergic reactions to the implanted materials, although uncommon, can occur;
- Metal sensitivity of allergic reactions to a foreign body have been reported, possibly leading to tumor formation;
- Decrease in bone density due to stress shielding;
- Neurological and spinal dura mater lesions from surgical trauma;
- Dural leak requiring surgical repair;
- Asymptomatic presence of microparticles may be observed around the implants as a result of interaction between the components as well as between the component and bone (i.e. wear).
- Cessation of growth of the fused portion of the spine;
- Loss of proper spinal curvature, correction, height and/or reduction;
- Pain, discomfort, or abnormal sensations due to the presence of the device;

- Early loosening may result from inadequate initial fixation, latent infection, premature loading of the device or trauma. Late loosening may result from trauma, infection, biological complications or mechanical problems, with the subsequent possibility of bone erosion, or pain.
- Serious complications may occur with any spinal surgery. These complications include, but are not limited to, genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; bursitis, hemorrhage, myocardial infarction, infection, paralysis or death.
- Inappropriate or improper surgical placement of this device may cause distraction or stress shielding of the graft or fusion mass. This may contribute to failure of an adequate fusion mass to form.
- Intraoperative fissure, fracture, or perforation of the spine can occur due to implantation of the components. Postoperative fracture of bone graft or the intervertebral body above or below the level of surgery can occur due to trauma, the presence of defects, or poor bone stock. Adverse effects may necessitate reoperation.
- Surgeons should discuss these side effects, including the finite life of the device, when appropriate.

Information for Patients

- The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion should be directed to the issues of premature weightbearing, activity levels, and the necessity for periodic medical follow-up.
The surgeon must warn the patient of the surgical risks and made aware of possible adverse effects. The surgeon must warn the patient that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) the surgeon must advise the patient that resultant forces can cause failure of the device. Patients who smoke have been shown to have an increased incidence of non-unions. Such patients should be advised of this fact and warned of the potential consequences. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief.

Removal


- Stryker Spine devices are indicated for treatment of fracture or stabilization of a surgical site even when fusion of the operated site is not sought.
- If fusion / bone graft growth occurs, the device will be deeply integrated into the bony tissues. As a result, the VBOSS® and VLIFT® systems are not intended to be removed unless the management of a complication or adverse event requires the removal.
- Standard instrument will be used to hold and disengage the device from the vertebrae.
- Any decision by a physician to remove the device should take into consideration such factors as the risk to the patient of the additional surgical procedure as well as the difficulty of removal.

Packaging

- The implants are single use devices, provided either non-sterile or sterile.
- Implants sold sterile are always presented in individual packaging and are clearly labelled as "sterile" on the package label. The sterilization method applied is indicated on the product label. Before utilization, it is recommended to verify the sterility expiration date. Stryker Spine cannot be held responsible for use of its products after their expiration date. It is recommended to verify the integrity of the original package before use. Sterility is ensured only if there is no trace of damage to the packaging. In case of damage to the packaging, or after opening of the packaging, re-sterilization of the implant is strictly forbidden, regardless of the method that might be employed.
- Implants sold non-sterile may be presented in individual packaging or implant containers as a complete set: implants and instruments arranged on trays and placed in specially designed storage boxes. The typical packaging used for implants sold non-sterile is clear plastic tubes and polyethylene bags. The packages must be intact at the time of receipt.
- The implants sold non-sterile must be removed entirely from their packaging prior to sterilization.

PRE-CLEANING / CLEANING AND STERILIZATION PROCEDURE RECOMMENDED FOR NON STERILE MEDICAL DEVICE

For safety reasons, non-sterile devices must be pre-cleaned, cleaned and sterilized prior to use. Moreover, for good maintenance, reusable instruments must be pre-cleaned, cleaned and sterilized immediately after surgery following the sequence of steps described in the following chart.

Sterilization procedure recommended for non-sterile medical devices including implants

Medical Devices should be sterilized in their container with water vapor in an autoclave in accordance with standard hospital procedure. The sterilization method suggested has been validated according to the AAMI TIR 12 in order to obtain a Sterility Assurance Level (SAL) of 10⁻⁶.

STERILIZATION CONDITIONS: 2 sets of low parameters have been validated on wrapped items:

- Prevacuum steam sterilization (Porous load autoclave): TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 4 minutes, DRY TIME: 45min.
- Gravity-displacement steam sterilization: TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 10 minutes, PRESSURE: 2.5 Bars/36-PSIG, DRY TIME: 45min.

CAUTION (U.S.A)

An FDA-cleared sterilization wrap is required when wrapping the sterilization containers.

The autoclave must be validated by the hospital and regularly checked to guarantee the recommended sterilization temperature is reached for the entire exposure time.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried and sterilization repeated.

Recommendations for implants provided sterile

Products delivered sterile have been exposed to a minimum of 25 kGy of gamma radiation process. Sterile products may be stored at room temperature and withstand the normal transportation conditions.

Products not labelled as sterile are not sterile.

Do not use if package is opened or damaged or after the "Use by" date on the label has expired.

Stryker shall not be responsible for the use of products presenting package deterioration or expiration of shelf life.

Re-sterilization of implants is strictly prohibited.

Care must be taken to prevent contamination of implant after opening of package.

Further Information

A surgical technique brochure is available on request through your Stryker agent or directly from Stryker Spine. Users with brochures that are over two years old at the time of surgery are advised to ask for an updated version.

CAUTION (U.S.A)

Based on the fatigue testing results, the physician/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

Federal law restricts this device to sale by or on the order of a licensed physician.

Complaints

Any health professional having a complaint or grounds for dissatisfaction relating to the quality of the product, its identity, its durability, its reliability, safety, effectiveness and / or its performance, should notify STRYKER Spine or its representative. Moreover, if a device has malfunctioned, or is suspected of having malfunctioned, STRYKER Spine or its representative must be advised immediately.

If a STRYKER Spine product has ever worked improperly and could have caused or contributed to the death of or serious injury to a patient, the distributor or STRYKER Spine must be informed as soon as possible by telephone, fax or in writing.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and an exhaustive description of the event to help STRYKER Spine understand the causes of the complaint.

For further information or complaints, please contact:

STRYKER SPINE
2 Pearl Court,
Allendale, NJ 07401-1677 USA
Tel. 201-760-8000
<http://www.stryker.com>

STRYKER SPINE
2 Pearl Ct
Allendale, N.J. 07401, USA
Tel: 201 760 8000
Fax. 201-760-8398 (Customer Service)

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