stryker[®]

Spine

NOLI135B01REV06

RR DONNELLEY

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Format strony X: 74.9999
Page Size X:

Format strony Y: 135
Page Size Y:

**IMPORTANT PRODUCT INFORMATION FOR
STRYKER SPINE ANTERIOR CERVICAL PLATING SYSTEMS
NON STERILE PRODUCT**

Description

The Stryker Spine Reflex ACP (Anterior Plate System) and Reflex Hybrid ACP are designed for use as an aid in cervical spinal fusion and are designed for unilateral fixation. The Stryker Spine ACP systems consist of bone screws and bone plates that are available in a variety of sizes in order to help accommodate individual patient physiology and pathology and to help facilitate anterior stabilization of the cervical spine.

The Reflex ACP and Reflex Hybrid ACP system components are only to be used within their designated system and must not be interchanged.

Material

The Reflex ACP plates are manufactured out of TiMZF alloy (ASTM F1813). The Reflex Hybrid plates and all other components for both systems, Reflex ACP and Reflex Hybrid ACP, are made out of Titanium alloy as defined in the ISO 5832-3 and ASTM F136 standards.

Indications

The ACP Systems are intended for anterior intervertebral screw fixation of the cervical spine from C2 – T1. These systems are indicated for temporary stabilization of the anterior spine during the development of cervical spine fusions in patients with the following indications:

- Degenerative disc disease (as defined by neck pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies)
- Decompression of the spinal cord following total or partial cervical vertebrectomy
- Trauma (including fractures)
- Tumors
- Deformities or curvatures (including kyphosis, lordosis or scoliosis)
- Pseudoarthrosis
- Failed previous fusions
- Spondylolisthesis
- Spinal stenosis

WARNING: These devices are not approved for screw attachment to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

Caution: Based on the fatigue testing results, the physician/ surgeon must consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

The Stryker Spine Reflex ACP (Anterior Plate System) and Reflex Hybrid ACP have not been evaluated for safety and compatibility in the MR environment. The Stryker Spine Reflex ACP (Anterior Plate System) and Reflex Hybrid ACP has not been tested for heating or migration in the MR environment

General conditions of use

Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and limitations of the spinal device. Knowledge of surgical techniques, proper reduction, selection and placement of implants, and pre- and post-operative patient management are considerations essential to a successful surgical outcome. Consult the medical literature for information regarding proper surgical techniques, precautions, and potential adverse effects associated with spinal fixation surgery

The components of each system have been designed to work together within each system; however, the components are not interchangeable between the two systems. Do not substitute another manufacturer's device for any component of the ACP Systems. Any such use will negate the responsibility of Stryker Spine for the performance of the resulting mixed component implant.

Do not mix metals (i.e. Titanium based devices with stainless steel items). Some corrosion occurs on all implanted metals and alloys. Contact of dissimilar metals, however, may accelerate corrosion. Corrosion may accelerate fatigue fracture of implants, and cause metal compounds to be released into the body.

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Format strony X: 75
Page Size X:
Format strony Y: 135
Page Size Y:

Anatomical Limitations

- The ACP Systems are intended for use in the cervical spine only. However, as with any orthopaedic implant, even when an implant's design does not expressly contraindicate its placement in a particular area, the surgeon may encounter certain patient physiologies which impose their own unique anatomic limitations.
- Anterior cervical plates are for use in the cervical region of the spine only and must not be used below T1.

Contra-indications

- Marked local inflammation.
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.
- Bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/or fixation to the devices.
- Bony abnormalities preventing safe screw fixation.
- Open wounds.
- Rapid joint disease, bone absorption, osteopenia, osteomalacia, and/or osteoporosis.
- Osteoporosis or osteopenia are relative contraindications, since this condition may limit the degree of obtainable correction and/or the amount of mechanical fixation.
- Metal sensitivity, documented or suspected.
- Pregnancy.
- Anytime implant utilization would interfere with anatomical structures or physiological performance.
- Inadequate tissue coverage over the operative site.

Other medical or surgical conditions which would preclude the potential benefit of surgery, such as congenital abnormalities, immunosuppressive disease, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or marked left shift in the WBC differential count.

These contra-indications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive.

Information for Patients

The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion should be directed to the issues of premature weightbearing, activity levels, and the necessity for periodic medical follow-up.

The surgeon must warn the patient of the surgical risks and make the patient aware of possible adverse effects. The surgeon must warn the patient that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) the surgeon must advise the patient that resultant forces can cause failure of the device. Patients who smoke have been shown to have an increased incidence of non-unions. Surgeons must advise of this fact and warned of the potential consequences. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief.

Pre-operative precautions

The surgical indication and the choice of implants must take into account certain important criteria such as:

- Patients involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) may be at increased risk for failure of the fusion and/or the device.
- Surgeons must instruct patients in detail about the limitations of the implants, including, but not limited to, the impact of excessive loading through patient weight or activity, and be taught to govern their activities accordingly. The procedure will not restore function to the level expected with a normal, healthy spine, and the patient should not have unrealistic functional expectations.
- A condition of senility, mental illness, chemical dependence or alcoholism. These conditions among others may cause the patients to ignore certain necessary limitations and precautions in the use of the implant, leading to failure and other complications.
- Foreign body sensitivity. Where material sensitivity is suspected appropriate tests must be made prior to material implantation.
- Surgeons must advise patients who smoke have been shown to have an increased incidence of non-unions. Surgeons must advise patients of this fact and warned of the potential consequences.

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Format strony X: 75
Page Size X:

Format strony Y: 135
Page Size Y:

- Care must be taken to protect the components from being marred, nicked, or notched as a result of contact with metal or abrasive objects. Alterations will produce defects in surface finish and internal stresses which may become the focal point for eventual breakage of the implant.

The choice of implants

- The choice of proper shape, size and design of the implant for each patient is crucial to the success of the surgery. The surgeon is responsible for this choice which depends on each patient.
- The size and shape of the bone structures determine the size, shape and type of the implants. Once implanted, the implants are subjected to stresses and strains. These repeated stresses on the implants must be taken into consideration by the surgeon at the time of the choice of the implant, during implantation as well as in the post-operative follow-up period. Indeed, the stresses and strains on the implants may cause metal fatigue or fracture or deformation of the implants, before the bone graft has become completely consolidated. This may result in further side effects or necessitate the early removal of the osteosynthesis device.

Intra-operative precautions

- The insertion of the implants must be carried out using instruments designed and provided for this purpose and in accordance with the specific implantation instructions for each implant. Those detailed instructions are provided in the surgical technique brochure supplied by Stryker Spine.
- Discard all damaged or mishandled implants.
- Stryker Spine implants must not be reshaped, unless otherwise indicated in the surgical technique instructions. When implants need to be bent, the bending must be carried out gradually using the appropriate instruments, provided by Stryker Spine. The use of inappropriate instruments may result in scratches, notches, and sharp bending, causing the breakage of the implants. Improper seating of the implant may result in implant failure.
- Never reuse an implant, even though it may appear undamaged.
- Do not mix metals.

Post-operative precautions

- Physician instructions regarding full weight-bearing activities must be complied with until maturation of the fusion mass is confirmed. Failure to comply with physician instructions may result in failure of the implant, the fusion, or both.

Side effects

Include but are not limited to:

- Late bone fusion or no visible fusion mass and pseudarthrosis;
- Peripheral neuropathies, nerve damage, heterotopic bone formation and neurovascular compromise, including paralysis;
- Superficial or deep-set infection and inflammatory phenomena;
- Allergic reactions to the implanted materials although uncommon can occur;
- Metal sensitivity of allergic reactions to a foreign body have been reported, possibly leading to tumor formation;
- Decrease in bone density due to stress shielding;
- Neurological and spinal dura mater lesions from surgical trauma;
- Dural leak requiring surgical repair;
- Asymptomatic presence of microparticles may be observed around the implants as a result of interaction between the components as well as between the component and bone (i.e. wear).
- Cessation of growth of the fused portion of the spine;
- Loss of proper spinal curvature, correction, height and/or reduction;
- Pain, discomfort, or abnormal sensations due to the presence of the device;
- Early loosening may result from inadequate initial fixation, latent infection, premature loading of the device or trauma.
- Late loosening may result from trauma, infection, biological complications or mechanical problems, with the subsequent possibility of bone erosion, or pain.
- Serious complications may occur with any spinal surgery. These complications include, but are not limited to, genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; bursitis, hemorrhage, myocardial infarction, infection, paralysis or death.
- Inappropriate or improper surgical placement of this device may cause distraction or stress shielding of the graft or fusion mass. This may contribute to failure of an adequate fusion mass to form.
- Intraoperative fissure, fracture, or perforation of the spine can occur due to implantation of the components. Postoperative fracture of bone graft or the intervertebral body above or below the level of surgery can occur due to trauma, the presence of defects, or poor bone stock.

Adverse effects may necessitate reoperation.

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Format strony Y: 135
Page Size Y:

Removal

- Stryker Spine devices are designed for treatment of fracture or stabilisation of a surgical site during the normal bone consolidation process. After this period, the presence of the device is no longer strictly required and its removal can be planned. Removal may also be necessary as a result of the above mentioned adverse effects.
- Removal of an ACP System may require special instruments to disengage the implant from the vertebrae. Appropriate recommendations are provided in the Surgical Technique brochure.
- Any decision by a physician to remove the internal fixation device must take into consideration such factors as the risk to the patient of the additional surgical procedure as well as the difficulty of removal.

Packaging

- The implants are delivered in packages; these must be intact at the time of receipt.
- The systems are sometimes supplied as a complete set: implants and instruments are arranged on trays and placed in specially designed storage boxes.

Reuse

An implant must never be reused. While it may appear undamaged, a used implant may have acquired blemishes or latent compromise of its integrity, which would reduce its service life. Surgeons must verify that the instruments are in good condition and operating order prior to use during surgery.

PRE-CLEANING / CLEANING AND STERILIZATION PROCEDURE RECOMMENDED FOR NON-STERILE MEDICAL DEVICE

For safety reasons, non-sterile devices must be pre-cleaned, cleaned and sterilized prior to use. Moreover, for maintenance, reusable instruments must be pre-cleaned, cleaned and sterilized immediately after surgery following the sequence of steps described in the following chart.

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
Format strony X: 75
 Page Size X:

Format strony Y: 135
 Page Size Y:

Sterilization procedure recommended for non-sterile medical devices including implants

Medical Devices should be sterilized in their container with water vapor in an autoclave in accordance with standard hospital procedure. The sterilization method suggested has been validated according to the AAMI TIR 12 in order to obtain a Sterility Assurance Level (SAL) of 10⁻⁶.

STERILIZATION CONDITIONS: Prevacuum steam sterilization (Porous load autoclave): TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 4 minutes, DRY TIME: 45min.

- Gravity-displacement steam sterilization: TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 10 minutes, PRESSURE: 2.5 Bars/36-PSIG, DRY TIME: 45min.

The autoclave must be validated by the hospital and regularly checked to guarantee the recommended sterilization temperature is reached for the entire exposure time.

If sterilization containers with paper filters are used, it is advisable to use a new filter for each sterilization. If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried and sterilization repeated.

Further Information

A surgical technique brochure is available on request through your Stryker agent or directly from Stryker Spine. Users with brochures that are over two years old at the time of surgery are advised to ask for an updated version.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Complaints

Any health professional, having a complaint or grounds for dissatisfaction relating to the quality of the product, its identity, its durability, its reliability, safety, effectiveness and / or its performance, should notify STRYKER Spine or its representative. Moreover, if a device has malfunctioned, or is suspected of having malfunctioned, STRYKER Spine or its representative must be advised immediately.

If a STRYKER Spine product has ever worked improperly and could have caused or contributed to the death of or serious injury to a patient, the distributor or STRYKER Spine must be informed as soon as possible by telephone, fax or in writing.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and an exhaustive description of the event to help STRYKER Spine understand the causes of the complaint.

For further information or complaints, please contact:

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B.V. Herikerbergweg 110
Stryker Spine SAS