stryker[®]

Spine

NOLI135C00REV04

RR DONNELLEY

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DESCRIPTION

Stryker Spine AVS® AL and AVS® ALign PEEK Spacers

The AVS® AL (Anterior Large) and AVS® ALign PEEK Spacers are intended for use as interbody fusion devices. They are offered in a variety of lengths, heights and lordotic angles. The hollow, ring shaped implant has serrations on the top and bottom for fixation.

Stryker Spine AVS® TL PEEK Spacers

The AVS® TL (Transforaminal Lumbar) PEEK Spacer is intended for use as an interbody fusion device. It is offered in both parallel (0°) and wedge (4°) shapes. The hollow implant has serrations on the top and bottom for fixation. It is offered in two medial/lateral widths, 25 & 30mm, and a variety of heights ranging from 7mm to 18mm.

Stryker Spine AVS® Navigator PEEK Spacers

The AVS® Navigator PEEK Spacer is intended for use as an interbody fusion device. It is offered in a variety of lengths, heights and lordotic angles. The hollow, oblong-shaped implant has serrations on the top and bottom for fixation. Radiopaque markers have been embedded within the implant to help allow for visualization in radiographic images.

MATERIAL

All components of the system are manufactured out of Polyetheretherketone (PEEK) (ASTM F2026) and Tantalum (ASTM F560).

INDICATIONS

Stryker Spine AVS® AL and AVS® ALign PEEK Spacers

The AVS® AL and AVS® ALign PEEK Spacers are indicated for the treatment of spondylolisthesis, degenerative spine disorders, discal and vertebral instability as well as in cases of spine revision surgery. It is recommended to pack bone graft material inside the implant.

The AVS® AL and AVS® ALign PEEK Spacer is intended for use with supplemental fixation. The supplemental fixation systems that may be used with the AVS® AL and AVS® ALign PEEK Spacers include, but are not limited to, STRYKER Spine plate or rod systems (XIA, Radius, Trio and Thor).

Stryker Spine AVS® TL PEEK Spacers

The AVS® TL PEEK Spacer is indicated for the treatment of spondylolisthesis, degenerative spine disorders, discal and vertebral instability as well as in cases of spine revision surgery. It is recommended to pack bone graft material inside the implant.

The AVS® TL PEEK Spacer is intended for use with supplemental fixation. The supplemental fixation systems that may be used with the AVS® TL PEEK Spacer include, but are not limited to: STRYKER Spine plate or rod systems (Xia, Spiral Radius 90D, and Trio).

Stryker Spine AVS® Navigator PEEK Spacers

The AVS® Navigator PEEK Spacers are indicated for the treatment of spondylolisthesis, degenerative spine disorders, discal and vertebral instability as well as in cases of spine revision surgery. It is recommended to pack bone graft material inside the implant.

The AVS® Navigator are to be implanted via posterior or posterolateral approach.

The AVS® Navigator PEEK Spacer is intended for use with supplemental fixation.

GENERAL CONDITIONS OF USE

The implantation of intervertebral body fusion devices must be performed only by experienced spinal surgeons having undergone the necessary specific training in the use of such systems because this is a technically demanding procedure presenting a risk of serious injury to the patient.

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Format strony X: 75
Page Size X:

Format strony Y: 135
Page Size Y:

The information contained in the Package Insert is necessary but not sufficient for the use of this device. This information is in no sense intended as a substitute for the professional judgment, skill and experience of the surgeon in careful patient selection, preoperative planning and device selection, knowledge of the anatomy and biomechanics of the spine, understanding of the materials and the mechanical characteristics of the implants used, training and skill in spinal surgery and the use of associated instruments for implantation, securing the patient's cooperation in following an appropriately defined post-operative management program and conducting scheduled post-operative follow-up examinations.

CAUTION

Based on the fatigue testing results, the physician/surgeon must consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

The implantation of the intervertebral body fusion device must be performed only by experienced spinal surgeons with specific training in the use of this device because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Potential risks identified with the use of this intervertebral body fusion device, which may require additional surgery, include: device component fracture, loss of fixation, pseudoarthrosis (i.e. non-union), fracture of the vertebrae, neurological injury, and vascular or visceral injury.

Instruments are provided by STRYKER Spine and must be used to assure accurate implantation of the intervertebral body fusion device. While rare, intraoperative fracture or breakage of instruments can occur, instruments, which have experienced extensive use or extensive force, are more susceptible to fracture depending on the operative precaution, number of procedures, and disposal attention. Instruments must be examined for wear or damage prior to surgery. Instruments for implantation of the AVS® PEEK Spacers are provided non-sterile and must be sterilized prior to use.

INFECTION

Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have been associated with transient bacteremia. To help prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.

INSTRUMENTS

Instruments are provided by STRYKER Spine and must be used to assure accurate implantation of the device. While rare, intraoperative fracture or breakage of instruments can occur. Instruments which have experienced extensive use or extensive force are more susceptible to fracture depending on the operative precaution, number of procedures, disposal attention. Instruments must be examined for wear or damage prior to surgery.

REUSE

Never reuse or reimplant spinal surgical implants. These could become contaminated resulting in infection. In addition, even though the device appears undamaged, it may have small defects which could compromise structural integrity reducing its service life and/or leading to patient injury. Surgeons must verify that the instruments are in good condition and operating order prior to use during surgery.

HANDLING

Correct handling of the implant is extremely important. The operating surgeon must avoid notching or scratching the device.

ALLERGY AND HYPERSENSITIVITY TO FOREIGN BODIES

When hypersensitivity is suspected or proven, it is highly recommended that the tolerance of the skin to the materials that make up the implants be checked before they are implanted

CONTRA-INDICATIONS FOR STRYKER SPINE AVS® AL AND AVS® ALIGN PEEK SPACERS

- The AVS® AL and AVS® ALign PEEK Spacers must not be implanted in patients with an active infection at the operative site.
- The AVS® AL and AVS® ALign PEEK Spacers are not intended for use except as indicated.
- Marked local inflammation.
- Any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis.
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.

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Format strony X: 75
Page Size X:

Format strony Y: 135
Page Size Y:

- Bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/ or fixation to the devices.
- Open wounds.
- Rapid joint disease, bone absorption, osteopenia, osteomalacia, and/or osteoporosis. Osteoporosis or osteopenia are relative contraindications, since this condition may limit the degree of obtainable correction and/or the amount of mechanical fixation.
- Pregnancy.
- Anytime implant utilization would interfere with anatomical structures or physiological performance.
- Inadequate tissue coverage over the operative site.
- Any neuromuscular deficit which places an unsafe load level on the device during the healing period.

Obesity. An overweight or obese patient can produce loads on the spinal system which can lead to failure of the fixation of the device or to failure of the device itself. Obesity is defined according to the W.H.O. standards.

A condition of senility, mental illness, or substance abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests must be made prior to material selection or implantation.

Other medical or surgical conditions that could preclude the potential benefit of surgery, such as congenital abnormalities, immunosuppressive disease, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or marked left shift in the WBC differential count, must be carefully analyzed before surgery.

These contra-indications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive. Surgeons must discuss the relative contraindications with the patients.

CONTRA-INDICATIONS FOR STRYKER SPINE AVS® TL PEEK SPACERS

The AVS® TL PEEK Spacer should not be implanted in patients with an active infection at the operative site.

The AVS® TL PEEK Spacers are not intended for use except as indicated.

Marked local inflammation.

Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.

Bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/ or fixation to the devices.

Open wounds.

Rapid joint disease, bone adsorption, osteopenia, osteomalacia, and/or osteoporosis. Osteoporosis or osteopenia are relative contraindications, since this condition may limit the degree of obtainable correction and/or the amount of mechanical fixation.

Pregnancy.

Anytime implant utilization would interfere with anatomical structures or physiological performance

Inadequate tissue coverage over the operative site.

Other medical or surgical conditions that could preclude the potential benefit of surgery, such as congenital abnormalities, immunosuppressive disease, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or marked left shift in the WBC differential count, must be carefully analyzed before surgery.

These contraindications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive.

CONTRA-INDICATIONS FOR STRYKER SPINE AVS® NAVIGATOR PEEK SPACERS

• The AVS® Navigator PEEK Spacers must not be implanted in patients with an active infection at the operative site.

• The AVS® Navigator PEEK Spacers are not intended for use except as indicated.

• Marked local inflammation.

• Any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis.

• Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.

• Open wounds.

• Pregnancy.

• Inadequate tissue coverage over the operative site.

• Any neuromuscular deficit which places an unsafe load level on the device during the healing period.

• Obesity. An overweight or obese patient can produce loads on the spinal system which can lead to failure of the fixation of the device or to failure of the device itself. Obesity is defined according to the W.H.O. standards.

RR DONNELLEY

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Product (language version):

NOLI135C00REV04

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Product ID & Issue:

cNOLI135C00REV04

Plik:

File Name:

UG_NOLI135C00REV04.p3.pdf

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Creation Date:

16-04-11

Strona w egzemplarzu: c3
Page position in booklet:

Format strony X: 75
Page Size X:

Format strony Y: 135
Page Size Y:

A condition of senility, mental illness, or substance abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests must be made prior to material selection or implantation.

Other medical or surgical conditions that could preclude the potential benefit of surgery, such as congenital abnormalities, immunosuppressive disease, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or marked left shift in the WBC differential count, must be carefully analyzed before surgery. Prior fusion at the levels to be treated.

These contra-indications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive. Surgeons must discuss the relative contraindications with the patients.

INFORMATION FOR PATIENTS

The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion should be directed to the issues of premature weightbearing, activity levels, and the necessity for periodic medical follow-up.

The surgeon must warn the patient of the surgical risks and made aware of possible adverse effects. The surgeon must warn the patient that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) the surgeon must advise the patient that resultant forces can cause failure of the device. Patients who smoke have been shown to have an increased incidence of non-unions. Such patients should be advised of this fact and warned of the potential consequences. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief. Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.

PRE-OPERATIVE PRECAUTIONS

The surgical indication and the choice of implants must take into account certain important criteria such as:

- Patients involved in an occupation or activity that applies excessive loading upon the implant (e.g., substantial walking, running, lifting, or muscle strain) may be at increased risk for failure of the fusion and/or the device.
- Surgeons must instruct patients in detail about the limitations of the implants, including, but not limited to, the impact of excessive loading through patient weight or activity, and be taught to govern their activities accordingly. The procedure will not restore function to the level expected with a normal, healthy spine, and the patient should not have unrealistic functional expectations.
- A condition of senility, mental illness, chemical dependence or alcoholism. These conditions among others may cause the patients to ignore certain necessary limitations and precautions in the use of the implant, leading to failure and other complications.
- Foreign body sensitivity. Where material sensitivity is suspected appropriate tests should be made prior to material implantation.
- Patients who smoke have been shown to have an increased incidence of non-unions. Such patients must be advised of this fact and warned of the potential consequences.
- Care must be taken to protect the components from being marred, nicked, or notched as a result of contact with metal or abrasive objects.

THE CHOICE OF IMPLANTS

The choice of proper shape, size and design of the implant for each patient is crucial to the success of the surgery. The surgeon is responsible for this choice, which depends on each patient.

Patients who are overweight may add additional stresses and strains on the device which can lead to implant fatigue and/or lead to deformation or failure of the implants.

The size and shape of the bone structures determine the size, shape and type of the implants. Once implanted, the implants are subjected to stresses and strains. These repeated stresses on the implants should be taken into consideration by the surgeon at the time of the choice of the implant, during implantation as well as in the post-operative follow-up period. Indeed, the stresses and strains on the implants may cause fatigue, fracture or deformation of the implants, before the bone graft has become completely consolidated. This may result in further side effects or necessitate the early removal of the device.

RR DONNELLEY

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Product (language version):

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Product ID & Issue:

NOLI135C00REV04

cNOLI135C00REV04

Plik:

File Name:

UG_NOLI135C00REV04.p4.pdf

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Creation Date:

16-04-11

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Page position in booklet:

Format strony X: 75

Page Size X:

Format strony Y: 135

Page Size Y:

INTRA-OPERATIVE PRECAUTIONS

The insertion of the implants must be carried out using instruments designed and provided for this purpose and in accordance with the specific implantation instructions for each implant. Those detailed instructions are provided in the surgical technique brochure supplied by STRYKER Spine.

Discard all damaged or mishandled implants.

Never reuse an implant, even though it may appear undamaged.

POST-OPERATIVE PRECAUTIONS

Physician instructions regarding full weight-bearing activities must be complied with until maturation of the fusion mass is confirmed. Failure to comply with physician instructions may result in failure of the implant, the fusion, or both.

PATIENT CARE FOLLOWING TREATMENT

Prior to adequate maturation of the fusion mass, implanted spinal instrumentation may need additional help to accommodate full load bearing. External support may be recommended by the physician from two to four months postoperatively or until x-rays or other procedures confirm adequate maturation of the fusion mass; external immobilization by bracing or casting be employed. Surgeons must instruct patients regarding appropriate and restricted activities during consolidation and maturation for the fusion mass in order to prevent placing excessive stress on the implants which may lead to fixation or implant failure and accompanying clinical problems. Surgeons must instruct patients to report any unusual changes of the operative site to his/her physician. The physician must closely monitor the patient if a change at the site has been detected.

ADVERSE EFFECTS

Include but are not limited to:

- Late bone fusion or no visible fusion mass and pseudarthrosis;
- While the expected life of spinal implant components is difficult to estimate, it is finite. These components are made of foreign materials which are placed within the body for the potential fusion of the spine and reduction of pain. However, due to the many biological, mechanical and physicochemical factors which affect these devices but cannot be evaluated in vivo, the components cannot be expected to indefinitely withstand the activity level and loads of normal healthy bone;
- Superficial or deep-set infection and inflammatory phenomena;
- Allergic reactions to the implanted materials although uncommon can occur;
- Decrease in bone density due to stress shielding;
- Dural leak requiring surgical repair;
- Peripheral neuropathies, nerve damage, heterotopic bone formation and neurovascular compromise, including paralysis;
- Cessation of growth of the fused portion of the spine;
- Loss of proper spinal curvature, correction, height and/or reduction;
- Pain, discomfort, or abnormal sensations due to the presence of the device;
- Delayed Union or Nonunion: Internal fixation appliances are load sharing devices which are used to obtain alignment until normal healing occurs. In the event that healing is delayed, does not occur, or failure to immobilize the delayed/nonunion results, the implant will be subject to excessive and repeated stresses which can eventually cause loosening, bending or fatigue fracture. The degree or success of union, loads produced by weight bearing, and activity levels will, among other conditions, dictate the longevity of the implant. If a nonunion develops or if the implants loosen, bend or break, the device(s) must be revised or removed immediately before serious injury occurs;
- Neurological and spinal dura mater lesions from surgical trauma;
- Early loosening may result from inadequate initial fixation, latent infection, premature loading of the device or trauma. Late loosening may result from trauma, infection, biological complications or mechanical problems, with the subsequent possibility of bone erosion, or pain.
- Serious complications may occur with any spinal surgery. These complications include, but are not limited to, genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; bursitis, hemorrhage, myocardial infarction, infection, paralysis or death.
- Inappropriate or improper surgical placement of this device may cause distraction or stress shielding of the graft or fusion mass. This may contribute to failure of an adequate fusion mass to form.
- Intraoperative fissure, fracture, or perforation of the spine can occur due to implantation of the components. Postoperative fracture of bone graft or the intervertebral body above or below the level of surgery can occur due to trauma, the presence of defects, or poor bone stock. Adverse effects may necessitate reoperation.

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Format strony Y: 135
Page Size Y:

REMOVAL

If fusion / bone graft growth occurs, the device will be deeply integrated into the bony tissues. As a result, the AVS® PEEK Spacer is not intended to be removed unless the management of a complication or adverse event requires the removal. Any decision by a physician to remove the device should take into consideration such factors as:

- The risk to the patient of the additional surgical procedure as well as the difficulty of removal.
- Migration of the implant, with subsequent pain and/or neurological, articular or soft tissue lesions
- Pain or abnormal sensations due to the presence of the implants
- Infection or inflammatory reactions
- Reduction in bone density due to the different distribution of mechanical and physiological stresses and strains.

PACKAGING

- The implants are single use devices, provided either non-sterile or sterile.
- Implants sold sterile are always presented in individual packaging and are clearly labelled as "sterile" on the package label. The sterilization method applied is indicated on the product label. Before utilization, it is recommended to verify the sterility expiration date. Stryker Spine cannot be held responsible for use of its products after their expiration date. It is recommended to verify the integrity of the original package before use. Sterility is ensured only if there is no trace of damage to the packaging. In case of damage to the packaging, or after opening of the packaging, re-sterilization of the implant is strictly forbidden, regardless of the method that might be employed.
- Implants sold non-sterile may be presented in individual packaging or implant containers as a complete set: implants and instruments arranged on trays and placed in specially designed storage boxes. The typical packaging used for implants sold non-sterile is clear plastic tubes and polyethylene bags. The packages must be intact at the time of receipt.
- The implants sold non-sterile must be removed entirely from their packaging prior to sterilization.

PRE-CLEANING / CLEANING AND STERILIZATION PROCEDURE RECOMMENDED FOR NON STERILE MEDICAL DEVICE

For safety reasons, non-sterile devices must be pre-cleaned, cleaned and sterilized prior to use. Moreover, for good maintenance, reusable instruments must be pre-cleaned, cleaned and sterilized immediately after surgery following the sequence of steps described in the following chart.

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
Format strony X: 75
 Page Size X:

Format strony Y: 135
 Page Size Y:

Sterilization procedure recommended for non-sterile medical devices including implants

Medical Devices should be sterilized in their container with water vapor in an autoclave in accordance with standard hospital procedure. The sterilization method suggested has been validated according to the AAMI TIR 12 in order to obtain a Sterility Assurance Level (SAL) of 10⁻⁶.

STERILIZATION CONDITIONS: 2 sets of low parameters have been validated on wrapped items:

- Vacuum steam sterilization (Porous load autoclave): TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 4 minutes, DRY TIME: 45min.
- Gravity-displacement steam sterilization: TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 10 minutes, PRESSURE: 2.5 Bars/36-PSIG, DRY TIME: 45min.

The autoclave must be validated by the hospital and regularly checked to guarantee the recommended sterilization temperature is reached for the entire exposure time.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried and sterilization repeated.

RECOMMENDATIONS FOR IMPLANTS PROVIDED STERILE

- Products delivered sterile have been exposed to a minimum of 25 kGy of gamma radiation process. Sterile products may be stored at room temperature and withstand the normal transportation conditions.
- Products not labelled as sterile are not sterile.
- Do not use if package is opened or damaged or after the "Use by" date on the label has expired.
- Stryker shall not be responsible for the use of products presenting package deterioration or expiration of shelf life.
- Re-sterilization of implants is strictly prohibited.
- Care must be taken to prevent contamination of implant after opening of package.

FURTHER INFORMATION

A surgical technique brochure is available on request through your STRYKER representative or directly from STRYKER Spine. Users with brochures that are over two years old at the time of surgery are advised to ask for an updated version.

COMPLAINTS

Any health professional having a complaint or grounds for dissatisfaction relating to the identity, quality, durability, reliability, safety, effectiveness or performance of a device should notify STRYKER Spine or its representative. Moreover, if a device has malfunctioned, or is suspected of having malfunctioned, STRYKER Spine or its representative must be advised immediately.

If a STRYKER Spine product has ever worked improperly and could have caused or contributed to the death of or serious injury to a patient, the distributor or STRYKER Spine must be informed as soon as possible by telephone, fax or in writing.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and a complete description of the event to help STRYKER Spine understand the causes of the complaint.

For further information or complaints, please contact:

STRYKER SPINE
2 Pearl Court,
Allendale, NJ 07401-1677 USA
Tel. 201-760-8000
<http://www.stryker.com>

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B.V. Herikerbergweg 110
Stryker Spine SAS