stryker[®]

Spine

NOLI135P00REV05

RR DONNELLEY

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Format strony X: 74.9999
Page Size X:

Format strony Y: 135
Page Size Y:

**IMPORTANT PRODUCT INFORMATION FOR
SOLIS® CERVICAL Cage and SOLIS® AS PEEK Spacer**

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DESCRIPTION

The SOLIS® CERVICAL Cage & SOLIS® AS PEEK Spacers are intended for use as interbody fusion devices. They are offered in various lengths, heights and lordotic angles to adapt to a variety of patient anatomies. These hollow, ring shaped implants have serrations on the top and bottom for fixation.

MATERIAL

The SOLIS® CERVICAL Cage is manufactured out of Polyetheretherketone (PEEK) (ASTM F2026), Titanium Alloy Ti6Al4V (ASTM F136; ISO 5832-3) and Commercially Pure Titanium (ASTM F67, ISO 5832-2).

The SOLIS® AS PEEK Spacer is manufactured out of Polyetheretherketone (PEEK) (ASTM F2026) and Tantalum (ASTM F560).

INDICATIONS

The SOLIS® CERVICAL Cage & SOLIS® AS PEEK Spacer are indicated for the treatment of cervical degenerative spine disorders, discal and vertebral instability as well as in cases of cervical spine revision surgery. It is recommended to pack bone graft material inside the implant.

The SOLIS® CERVICAL Cage incorporates spikes for fixation and is suitable for stand-alone use on one level fusion only. An anterior fixation is recommended while using cages in more than one disc level. This decision is the surgeon's responsibility according to the situation of each patient (bone quality, pathology, stability of the spine). When an anterior plating fixation is combined with the use of the cage, the surgeon must refer to the specific recommendations for this type of implant.

The SOLIS® AS PEEK Spacer is intended for use with supplemental fixation.

GENERAL CONDITIONS OF USE

Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and limitations of the spinal device. Knowledge of surgical techniques, proper reduction, selection and placement of implants, and pre- and post-operative patient management are other considerations essential to a successful surgical outcome. Consult the medical literature for information regarding proper surgical techniques, precautions, and potential adverse effects associated with spinal surgery.

INFORMATION FOR PATIENTS

The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion should be directed to the issues of premature weightbearing, activity levels, and the necessity for periodic medical follow-up.

The surgeon must warn patients of the surgical risks and make them aware of possible adverse effects. The patient must be warned that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) the surgeon must warn the patient that resultant forces can cause the resulting adverse effects after surgery. Patients who smoke have been shown to have an increased incidence of non-unions. Such patients should be advised of this fact and warned of the potential consequences, including device breakage. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief.

REUSE

Never reuse or reimplant spinal surgical implants. These could become contaminated resulting in infection. In addition, even though the device appears undamaged, it may have small defects which could compromise structural integrity reducing its service life and/or leading to patient injury.

Surgeons must verify that the instruments are in good condition and operating order prior to use during surgery.

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Format strony Y:
Page Size Y: 135

CAUTION

Based on the fatigue testing results, the physician/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

CONTRA-INDICATIONS

- The SOLIS® CERVICAL Cage & SOLIS® AS PEEK Spacer should not be implanted in patients with an active infection at the operative site.
- The SOLIS® CERVICAL Cage & SOLIS® AS PEEK Spacer are not intended for use except as indicated.
- Marked local inflammation.
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.
- Bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/or fixation to the devices.
- Open wounds.
- Rapid joint disease, bone absorption, osteopenia, osteomalacia, and/or osteoporosis. Osteoporosis or osteopenia are relative contraindications, since this condition may limit the degree of obtainable correction and/or the amount of mechanical fixation.
- Pregnancy.
- Anytime implant utilization would interfere with anatomical structures or physiological performance.
- Inadequate tissue coverage over the operative site.

Other medical or surgical conditions that could preclude the potential benefit of surgery, such as congenital abnormalities, immunosuppressive disease, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or marked left shift in the WBC differential count, must be carefully analyzed before surgery. These contra-indications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive.

PRE-OPERATIVE PRECAUTIONS

The surgical indication and the choice of implants must take into account certain important criteria such as:

- Patients involved in an occupation or activity that applies excessive loading upon the implant (e.g., substantial walking, running, lifting, or muscle strain) may be at increased risk for failure of the fusion and/or the device.
- Surgeons must instruct patients in detail about the limitations of the implants, including, but not limited to, the impact of excessive loading through patient weight or activity, and be taught to govern their activities accordingly. The procedure will not restore function to the level expected with a normal, healthy spine, and the patient should not have unrealistic functional expectations.
- A condition of senility, mental illness, chemical dependence or alcoholism. These conditions among others may cause the patients to ignore certain necessary limitations and precautions in the use of the implant, leading to failure and other complications.
- Foreign body sensitivity. Where material sensitivity is suspected appropriate tests should be made prior to material implantation.
- Surgeons must advise patients who smoke have been shown to have an increased incidence of non-unions. Such patients should be advised of this fact and warned of the potential consequences.
- Care must be taken to protect the components from being marred, nicked, or notched as a result of contact with metal or abrasive objects.
- Surgeons must discuss these precautions pre-operatively with their patients when appropriate

THE CHOICE OF IMPLANTS

- The choice of proper shape, size and design of the implant for each patient is crucial to the success of the surgery. The surgeon is responsible for this choice, which depends on each patient.
- The size and shape of the bone structures determine the size, shape and type of the implants. Once implanted, the implants are subjected to stresses and strains. These repeated stresses on the implants should be taken into consideration by the surgeon at the time of the choice of the implant, during implantation as well as in the post-operative follow-up period. Indeed, the stresses and strains on the implants may cause fatigue, fracture or deformation of the implants, before the bone graft has become completely consolidated. This may result in further side effects or necessitate the early removal of the device.

INTRA-OPERATIVE PRECAUTIONS

- The insertion of the implants must be carried out using instruments designed and provided for this purpose and in accordance with the specific implantation instructions for each implant. Those detailed instructions are provided in the surgical technique brochure supplied by STRYKER Spine.
- Discard all damaged or mishandled implants.
- Never reuse an implant, even though it may appear undamaged.

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Format strony X:
Page Size X: 74.9997

Format strony Y:
Page Size Y: 135

- Ensure the disc space is adequately distracted while the implant is inserted. Trials must be used to evaluate the correct size. To avoid risk of breakage, the SOLIS® CERVICAL cage & SOLIS® AS PEEK Spacer must not be impacted.

POST-OPERATIVE PRECAUTIONS

- Physician instructions regarding full weight-bearing activities must be complied with until maturation of the fusion mass is confirmed. Failure to comply with physician instructions may result in failure of the implant, the fusion, or both.

SIDE EFFECTS

Include but are not limited to:

- Late bone fusion or no visible fusion mass and pseudarthrosis;
- Peripheral neuropathies, nerve damage, heterotopic bone formation and neurovascular compromise, including paralysis;
- Superficial or deep-set infection and inflammatory phenomena;
- Allergic reactions to the implanted materials although uncommon can occur;
- Decrease in bone density due to stress shielding;
- Neurological and spinal dura mater lesions from surgical trauma;
- Dural leak requiring surgical repair;
- Cessation of growth of the fused portion of the spine;
- Loss of proper spinal curvature, correction, height and/or reduction;
- Pain, discomfort, or abnormal sensations due to the presence of the device;
- Early loosening may result from inadequate initial fixation, latent infection, premature loading of the device or trauma. Late loosening may result from trauma, infection, biological complications or mechanical problems, with the subsequent possibility of bone erosion, or pain.
- Serious complications may occur with any spinal surgery. These complications include, but are not limited to, genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; bursitis, hemorrhage, myocardial infarction, infection, paralysis or death.
- Inappropriate or improper surgical placement of this device may cause distraction or stress shielding of the graft or fusion mass. This may contribute to failure of an adequate fusion mass to form.
- Intraoperative fissure, fracture, or perforation of the spine can occur due to implantation of the components. Postoperative fracture of bone graft or the intervertebral body above or below the level of surgery can occur due to trauma, the presence of defects, or poor bone stock. Adverse effects may necessitate reoperation.
- Surgeons should discuss these side effects, including the finite life of the device, when appropriate.

REMOVAL

- If fusion / bone graft growth occurs, the device will be deeply integrated into the bony tissues. As a result, the SOLIS® CERVICAL Cage & SOLIS® AS PEEK Spacer are not intended to be removed unless the management of a complication or adverse event requires the removal.
- The SOLIS® or SOLIS® AS Inserter will be used to hold and disengage the device from the vertebrae. Additional general instruments may be required to assist in the removal of the device.
- Any decision by a physician to remove the device should take into consideration such factors as the risk to the patient of the additional surgical procedure as well as the difficulty of removal.

PACKAGING

- The implants are single use devices, provided either non-sterile or sterile.
- Implants sold sterile are always presented in individual packaging and are clearly labelled as "sterile" on the package label. The sterilization method applied is indicated on the product label. Before utilization, it is recommended to verify the sterility expiration date. Stryker Spine cannot be held responsible for use of its products after their expiration date. It is recommended to verify the integrity of the original package before use. Sterility is ensured only if there is no trace of damage to the packaging. In case of damage to the packaging, or after opening of the packaging, re-sterilization of the implant is strictly forbidden, regardless of the method that might be employed.
- Implants sold non-sterile may be presented in individual packaging or implant containers as a complete set: implants and instruments arranged on trays and placed in specially designed storage boxes. The typical packaging used for implants sold non-sterile is clear plastic tubes and polyethylene bags. The packages must be intact at the time of receipt.
- The implants sold non-sterile must be removed entirely from their packaging prior to sterilization.

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
Format strony X: 74.9997
Page Size X:

Format strony Y: 135
Page Size Y:

PRE-CLEANING / CLEANING AND STERILIZATION PROCEDURE RECOMMENDED FOR NON STERILE MEDICAL DEVICE

For safety reasons, non-sterile devices must be pre-cleaned, cleaned and sterilized prior to use. Moreover, for maintenance, reusable instruments must be pre-cleaned, cleaned and sterilized immediately after surgery following the sequence of steps described in the following chart.

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Format strony X: **74.9997**
 Page Size X:

Format strony Y: **135**
 Page Size Y:

Sterilization procedure recommended for non-sterile medical devices including implants

Medical Devices should be sterilized in their container with water vapor in an autoclave in accordance with standard hospital procedure. The sterilization method suggested has been validated according to the AAMI TIR 12 in order to obtain a Sterility Assurance Level (SAL) of 10⁻⁶.

STERILIZATION CONDITIONS: 2 sets of low parameters have been validated on wrapped items:

- Precavacuum steam sterilization (Porous load autoclave): TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 4 minutes, DRY TIME: 45min.
- Gravity-displacement steam sterilization: TEMPERATURE: 132°C (270°F), EXPOSURE TIME: 10 minutes, PRESSURE: 2.5 Bars/36-PSIG, DRY TIME: 45min.

The autoclave must be validated by the hospital and regularly checked to guarantee the recommended sterilization temperature is reached for the entire exposure time.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried and sterilization repeated.

RECOMMENDATIONS FOR IMPLANTS PROVIDED STERILE

- Products delivered sterile have been exposed to a minimum of 25 kGy of gamma radiation process. Sterile products may be stored at room temperature and withstand the normal transportation conditions.
- Products not labelled as sterile are not sterile.
- Do not use if package is opened or damaged or after the "Use by" date on the label has expired.
- Stryker shall not be responsible for the use of products presenting package deterioration or expiration of shelf life.
- Re-sterilization of implants is strictly prohibited.
- Care must be taken to prevent contamination of implant after opening of package.

FURTHER INFORMATION

A surgical technique brochure is available on request through your STRYKER representative or directly from STRYKER Spine. Users with brochures that are over two years old at the time of surgery are advised to ask for an updated version.

COMPLAINTS

Any health professional having a complaint or grounds for dissatisfaction relating to the identity, quality, durability, reliability, safety, effectiveness or performance of a device should notify STRYKER Spine or its representative. Moreover, if a device has malfunctioned, or is suspected of having malfunctioned, STRYKER Spine or its representative must be advised immediately.

If a STRYKER Spine product has ever worked improperly and could have caused or contributed to the death of or serious injury to a patient, the distributor or STRYKER Spine must be informed as soon as possible by telephone, fax or in writing.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and a complete description of the event to help STRYKER Spine understand the causes of the complaint.

For further information or complaints, please contact:

STRYKER SPINE
2 Pearl Court,
Allendale, NJ 07401-1677 USA
Tel. 201-760-8000
<http://www.stryker.com>

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B.V. Herikerbergweg 110
Stryker Spine SAS