stryker
Howmedica
OSTEONICS

Non-Sterile Instruments

Howmedica Osteonics Corp.
325 Corporate Drive
Mahwah, NJ 07430, USA
A subsidiary of Stryker Corporation

Telephone #: +1 201-831-5000

CE 0086

© 2017 Howmedica Osteonics Corp.

QIN 4333 Rev. AB

QIN 4333, Ver. AB EON Release Date: Jun 02, 2017

Print Date: Feb 08, 2018 07:11:49 GMT

Labeling Symbols

Attention, See Instructions for Use

Do Not Reuse

Sterilized using Irradiation

Sterilized using Hydrogen Peroxide

Sterilized using Ethylene Oxide

Use by Date

Date of Manufacture

Legal Manufacturer

Authorized Representative in the European Community

Catalog Number

Batch Code

Serial Number

Non-Sterile Instruments.....	1
Instruments non stériles	4
Nicht sterile Instrumente.....	7
Strumenti non sterili.....	10
Instrumentos no estériles	13
Instrumentos não esterilizados.....	16
lcke-sterila instrument	19
Steriloimattomat instrumentit	22
Usterile instrumenter	25
Niet-steriele instrumenten	28
Ikke-sterile instrumenter.....	31
Oprzrządowanie niesterylne	34
Μη αποστειρωμένα εργαλεία	37
非消毒器械.....	40
비멸균 기구.....	43
非滅菌性手術器具.....	47

English

NON-STERILE INSTRUMENTS

HOWMEDICA OSTEONICS Instrumentation consists of manual surgical instruments intended for use during hip, knee, shoulder, elbow or wrist arthroplasty, spinal fusion or trauma surgeries.

Utilization

- Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- HOWMEDICA OSTEONICS Trial Components should be used for size determination, preparation evaluation, trial reduction and range of motion evaluation, thus helping to preserve the integrity of the actual implants and their sterile packaging.
- Radiographic templates are available to help assist in the preoperative prediction of component size and style.
- The Howmedica Osteonics Corp. Surgical Protocols provide additional procedural information.

Warnings

- Care must be taken not to cut through surgical gloves when handling any sharp-edged orthopaedic device.
- Due to different manufacturers employing differing design parameters, varying tolerances, different materials and manufacturing specifications, HOWMEDICA OSTEONICS instrumentation should not be used to implant any other manufacturer's components. Any such use will negate the responsibility of Howmedica Osteonics Corp. for the performance of the resulting implant.
- Instruments made of non-metallic material(s) and fragments thereof may not be visible using certain forms of external imaging (e.g. x-ray) unless otherwise specified, such as radiopaque femoral head trials that are visible.
- Unless otherwise specified on the instrument, markings on instruments with a measuring function have a linear accuracy of ± 0.01 inches (± 0.25 mm) or $\pm 5^\circ$ for angular measurements.

Adverse Effects

- Malfunctions of instruments, including, but not limited to, weld fracture, thread galling and material fatigue, have occurred in a small percentage of cases.
- Serious complications may be associated with any joint replacement, spinal or trauma surgery. These complications include, but are not limited to: genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Peripheral neuropathies, nerve damage, circulatory compromise and heterotopic bone formation may occur.
- Intraoperative fissure, fracture, or perforation of the bone can occur due to numerous factors including the presence of defects, poor bone stock, the use of the surgical instruments and the impaction of the component into the preparation.

Cleaning and Sterilization

- This reusable instrument is **not** sterile.
- Remove all instruments from their packaging prior to cleaning and sterilization.
- All parts must be cleaned and sterilized before use. Cleaning of the reusable instrument prior to sterilization should include: presoaking, manual cleaning, ultrasonication (e.g. for complex parts), and washing/disinfecting. Enzymatic cleaners, manual detergents, and neutral cleaners should be used.

- During and after cleaning, inspect hard-to-reach areas such as lumens within the instrument to help ensure that entrapped debris and/or soil are completely removed. If debris and/or soil are evident, then reclean the instrument until there is no evidence of debris and/or soil.
- Rinsing should be performed after each cleaning step to help ensure cleaners and/or debris/soil are removed.
- Howmedica Osteonics reusable instruments are recommended to be sterilized by moist heat sterilization.
- Instruments with articulating surfaces must be tested for movement. A moist heat compatible, medical grade lubricant should be applied to all articulating joints prior to sterilization.
- After cleaning and prior to sterilization, the reusable instrument should be double-wrapped or packaged in CSR sterilization wraps or pouches. Wrapping should be performed using the appropriate wrapping method (e.g. AAMI CSR wrapping technique).
- Select Stryker Orthopaedics trays and instruments are compatible with Aesculap SterilContainer System (JN441 or JN442 perforated bottom with JK48X series of aluminum 2000 lids). Please refer to LSTPI-B for a list of compatible parts, and detailed instructions for use.
- Biological indicators and/or steam integrators should be used within the sterilizer load to help ensure the proper sterilization conditions of time, temperature, and saturated steam have been met.
- The parameters for sterilization are:

Method:	Moist-Heat Sterilization
Cycle:	Pre-Vacuum (Pre-Vac)
Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Pressure:	2-15 PSIA
Dry-Time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Howmedica Osteonics Corp. has validated the above recommended sterilization cycle based on AAMI/ANSI/ISO guidelines and recommendations. Other sterilization methods and cycles may also be suitable. However, individuals or hospitals are advised to validate whichever method they deem appropriate at their institution.

- EtO sterilization and cold sterilization techniques are not recommended.

Select Stryker Orthopaedics trays and instruments may be sterilized with the following parameters **outside of the United States**:

Method:	Moist-Heat Sterilization
Cycle:	Pre-Vacuum
Temperature:	273° F - 279° F (134° C - 137° C)
Exposure Time:	3 minutes (minimum)
Dry-time:	30 minutes (minimum, in chamber)
Cool-Time:	60 minutes (minimum, at room temperature)

Howmedica Osteonics Corp. has validated the above recommended OUS sterilization cycle based on ISO guidelines and recommendations. Please refer to LSTPI-B for a list of instruments/trays that may be sterilized using this cycle.

Detailed manual and automatic cleaning instructions, including times and temperatures, are provided in Stryker Orthopaedics Literature LSTPI-B available at www.stryker.com/orthopaedics/cleaning. Disassembly and assembly instructions applicable to all Stryker Orthopaedics re-usable instruments are all also provided in this brochure.

Refer to Aesculap Instructions for Use for Care & Handling of Aesculap SterilContainer Systems (JN441 or JN442).

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a licensed physician.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademark(s) or service mark(s): Howmedica, Osteonics, Stryker. All other trademarks are trademarks of their respective owners or holders.

Refer to product label for CE Mark Status and Legal Manufacturer. The CE mark is only valid if also found on the product label.

The following table contains a list of abbreviations that are used on Howmedica Osteonics Corp. product labeling:

Term	Abbreviation	Term	Abbreviation
Alpha Code	ALPH CDE	Neck	NK
Angle	ANG	Offset	OFFST
Degree	DEG or °	Outer Diameter	OD
Diameter	DIA	Right	RT ▶
Extra Deep	XDP	Screw Holes	SCR HLS
Extra Large	XLGE	Side	SDE
Extra Small	XSM	Size	SZE
Head	HD	Small	SM
Height	HT	Standard	STD
Inner Diameter	ID	Taper	TPR
Insert	INSR	Thickness	THKNS
Large	LGE	Type	TYP
Left	◀ LFT	With	W/
Length	LNTH	Without	W/O
Medium	MED		

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