DOCUMENT COVER SHEET

TITLE : TD6002190713 Drill
Guide Set IFU Multilingual

Name: TD6002190713 Drill Guide
Set IFU Multilingual

Number: 0000063011

Created By: Hassler Nadine

Created: 2017-09-29 13:04:46.0 GMT

Updated By: Hassler Nadine

Last Updated: 2017-09-29 13:04:46.0 GMT

Version: C.0

Signature History Report

Version	Task	User Name	Event Date	Role
C	Give your signature information and complete the task to approve the attached review document. Upon completion a composite PDF will be created using the signature given here and the primary content of the attached document.	Brunello Paola Maria	2017-09-29 19:06:14 CEST	Technical Documentation

Instructions for Use





Navigated Drill Guide Set




Rx Only

Navigated Short Drill Guide	REF 6002-140-000
Navigated Long Drill Guide	REF 6002-190-000
Navigated Short Drill Guide Calibrator	REF 6002-140-001
Navigated Long Drill Guide Calibrator	REF 6002-190-001
Navigated Short Ø 2.5 mm Drill Bit	REF 6002-140-025
Navigated Short Ø 3.0 mm Drill Bit	REF 6002-140-030
Navigated Short Ø 3.5 mm Drill Bit	REF 6002-140-035
Navigated Long Ø 3.0 mm Drill Bit	REF 6002-190-030
Navigated Long Ø 3.5 mm Drill Bit	REF 6002-190-035
Navigated Long Ø 4.0 mm Drill Bit	REF 6002-190-040
Navigated Long Ø 4.5 mm Drill Bit	REF 6002-190-045
Navigated Long Ø 5.0 mm Drill Bit	REF 6002-190-050
Navigated Long Ø 5.5 mm Drill Bit	REF 6002-190-055
Navigated Drill Bit Stop	REF 6002-140-100
Navigated Drill Guide Set Container	REF 6002-140-200

English	1
Español.....	10
Deutsch.....	20
Français.....	30
Italiano	40
Nederlands	50
Svenska.....	59
Dansk.....	68
Norsk.....	77
Suomi.....	86
Português	95
Polski.....	105
Ελληνικά	115
Русский.....	125
한국어	136
简体中文.....	145
日本語	154
Türkçe	163
Română.....	172

1 Introduction

This document is the most comprehensive source of information for the safe and effective use of the product. It is directed at medical staff.

Read this document carefully. Familiarization with the user documentation for the components of the Stryker navigation system prior to use is important.

Pay special attention to safety information. Keep this document accessible to users.

If you need additional information, assistance or training, contact your Stryker sales representative (USA) or your nearest Stryker subsidiary (outside USA).

The healthcare provider performing any procedure is responsible for the adequate training of the medical staff, for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.

2 Safety Information

WARNING:

Use

- The product is intended for use by healthcare professionals only.
- Do not use the product outside of its indications for use. Failure to comply may compromise navigation accuracy.
- Handle the product with care to avoid injury to yourself or to the patient. The drill guides and the drill bits have sharp tips.

Damage or Loose Components

- Prior to each use, check the product for loose components, damage and bent or deformed parts. Do not use if these conditions exist. Failure to comply may lead to patient or medical staff injury. If service is required, contact your Stryker sales representative.

Accessories

- Use only Stryker approved parts and accessories unless otherwise specified. Failure to comply may lead to patient or medical staff injury.

Single-Use Parts

- Do not re-use the drill bits. The drill bits are single use only. Re-use may cause mechanical failure, such as bending or fracture.

Cleaning and Sterilization

- Clean and sterilize all parts for the first use, and clean and sterilize reusable parts for every subsequent use.

Calibrator

A bent calibrator tip or otherwise changed calibrator geometry may lead to incorrect calibration of the drill guide. Incorrect calibration leads to navigation inaccuracy, which may lead to patient injury.

- Do not use the calibrator tip as a lever or bend the tip.
- Use the calibrator only for calibration. Do not use the calibrator as an obturator or awl. Hammering impact may damage the calibrator such that all subsequent calibrations with the damaged calibrator become inaccurate.

3 Indications for Use

The Navigated Drill Guide Set is intended to be used as an accessory to the Stryker Spine Navigation System. The Navigated Drill Guide Set consists of manual instruments that are intended to be used in spine surgery by providing guidance during drilling.

The Navigated Drill Guide Set may be used as part of the Stryker Spine Navigation System, which is indicated for any medical condition in which the use of computer assisted surgery may be appropriate. The system can be used for intraoperative guidance where a reference to a rigid anatomical structure can be identified.

4 Product Description

Drill Guide Set

The drill guide set consists of navigated drill guides, calibrators, drill bits, a drill bit stop and a container. The drill guides are available in two lengths. Refer to Chapter 10 "Technical Specifications".

The set serves to drill holes of predefined depths during navigated surgery. The user marks the desired depth on the drill bit with the drill bit stop. The user then directs the drill bit to the target anatomy through the drill guide. The position of the drill guide relative to the patient is tracked and displayed in the patient images on the navigation screen.

Drill Guide Set Container

The container holds up to 30 instruments or instrument parts for sterilization and storage.

5 For Use With

The drill guide set is intended for exclusive use with the nGenius Universal Tracker using the Stryker Spine Navigation System.

6 Overview

Figure ❶ - Drill guide and nGenius Universal Tracker	A Tracker interface
Figure ❷ - Drill guide and calibrator	B Rotatable handle
Figure ❸ - Rotatable handle	C Insertion point (for calibrator and drill bits)
Figure ❹ - Container	D Calibrator tip
Figure ❺ - Drill bit and drill bit stop	E Calibrator
Figure ❻ - Setting the depth mark on the drill bit	F Calibrator knob
Figure ❼ - Lubrication (Maintenance)	G Interface button
	H Auxiliary space
	I Release button
	J Drill bit stop
	K Drill bit
	L Depth mark

7 Use

- ✓ Choose the correct drill guide size for your patient (short or long) and the corresponding calibrator size (short or long). See Chapter 10 "Technical

Specifications". The long drill guide is only intended to be used on the thoracic and lumbar spine.

Attach Tracker to Drill Guide (1)

1. Press and hold the interface button (G) on the tracker.
2. Slide the tracker onto the drill guide's tracker interface (A) and release the interface button (G).
3. Verify that the tracker is locked into place.

Calibrate Drill Guide (2)

▲ WARNING: Calibrate the drill guide only with the calibrator fully inserted. Incorrect calibration leads to navigation inaccuracy and may lead to patient injury.

1. Screw the calibrator (E) into the drill guide until the calibrator is fully tightened.
2. Turn on the tracker and follow the calibration workflow on the navigation screen. For detailed instructions, refer to the user manual for the software application.
3. After calibration and validation, remove the calibrator from the drill guide.

Adjust Tracker Position (Optional) (3)

For navigation, the LEDs on the tracker must be visible to the navigation camera. If required, you can adjust the position of the tracker by rotating the handle of the drill guide:

1. Hold the handle (B) and place the thumb of the same hand on the insertion point (C) for the calibrator and drill bits.
2. Push your thumb against the insertion point to pull back the handle.
3. With the handle pulled back, rotate the handle into the desired position.
4. Verify that the handle is locked into place.

Use Drill Bit Stop To Mark Drilling Depth on Drill Bit (4/5)

✓ Select the drill bit of appropriate size based on the intended screw diameter and length.

1. Press and hold the release button on the drill bit stop (J).
2. Slide the drill bit stop onto the drill bit (K). **Note:** The drill bit stop is bidirectional and may be slid on in either orientation.
3. Place the drill bit stop at the intended depth mark (L) as read on the distal end of the drill bit. **▲ WARNING:** Ensure that the drill bit stop clicks into a slot on the drill bit and verify that it is locked into place at the intended depth mark. Otherwise the intended drilling depth may be surpassed and the patient may be injured.
4. Attach a Quick Release Handle.

Insert Drill Guide and Drill Bit To Drilling Location

▲ WARNING:

- Insert drill bits only through the drill guide. If you do not use the drill guide, navigation is not possible and the drill bit may breach sensitive anatomy.
- Be careful when you insert the drill guide. The drill guide has a sharp serrated tip.
- Dock the drill guide tip firmly into place on the bone to prevent the drill guide from moving off the bone during drilling.

- Do not use excessive force on the drill bits to avoid bending or fracture. Broken pieces that come in contact with the patient may damage sensitive anatomy or cause adverse reactions.

Note: The navigation system tracks and visualizes the position and orientation of the drill guide, not of the drill bit. Plan and adjust the drilling depth with the drill bit stop before drilling. For additional information on navigating the drill guide, refer to the user manual for the software application.

1. Insert the drill guide through the prepared incision site to the intended drilling location.
2. Navigate the drill guide to establish the final screw trajectory using the navigation system.
3. Guide the drill bit through the drill guide to the intended drilling location.
4. During drilling, confirm the final drilling depth via visual and tactile feedback when the drill bit stop contacts the insertion point of the drill bit.

8 Cleaning, Sterilization and Storage

When to Clean and Sterilize

⚠ WARNING: Clean and sterilize all parts for the first use, and clean and sterilize reusable parts for every subsequent use.


As part of routine maintenance, preclean, clean and sterilize the parts immediately after surgery.

Instructions

Follow the sequence of steps in Chart 1 and the additional recommendations for sterilization. Consider the following product-specific instructions:

- Preparation for cleaning: Remove the tracker from the drill guide.
- Sterilization and storage: Place the parts in the dedicated location in the drill guide set container.

Chart 1: Cleaning and Sterilization Steps

Sterilization Procedure Recommended for Non-Sterile Medical Devices

Medical Devices should be sterilized in their container by means of moist heat steam sterilization in an autoclave in accordance with standard hospital procedure. The recommended sterilization methods have been validated according to the ANSI/AAMI ST79 in order to obtain a Sterility Assurance Level (SAL) of 10^{-6} .

The autoclave must be validated by the hospital and regularly checked to guarantee the recommended sterilization temperature is reached for the entire exposure time.

If sterilization containers with paper filters are used, it is advisable to use a new filter for each sterilization.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried and sterilization repeated.

The sterilization parameters have been validated on containers in fully-loaded autoclaves.

Sterilization Conditions - USA

The product can be sterilized individually using FDA cleared sterilization wrap or in an appropriate instrument tray or container wrapped in FDA cleared sterilization wrap or in an FDA cleared rigid sterilization container. See Table 1 for recommended sterilization parameters. Stryker container lids must be removed prior to use with the Aesculap reusable, rigid sterilization container. Sterilization using the Aesculap container configuration in Table 2 with the pre-vacuum sterilization parameters in Table 1 is recommended.

Table 1: Sterilization Parameters - USA

	Pre-vacuum	Gravity
Temperature	132°C (270°F)	132°C (270°F)
Exposure time	4 minutes	15 minutes
Drying time	30 minutes	30 minutes

Table 2: Rigid Container Configuration - USA

Aesculap Part Number	Description
JN442	Full size container, perforated bottom with retention plates (6-inch height)
JK488	Full-size container lid, gold
JK489	Full-size container lid, silver
US751	Single use 7.5-inch diameter filter with indicator dot
US994	Single use 7.5-inch diameter filter without indicator dot

Sterilization Conditions – Outside of USA

The product can be sterilized individually using cleared sterilization wrap or in the appropriate instrument tray or container wrapped in cleared sterilization wrap or in a cleared rigid sterilization container. See Table 3 for recommended sterilization parameters. Stryker container lids must be removed prior to use with the Aesculap reusable, rigid sterilization container.

Sterilization using the Aesculap container configuration in Table 4 with the pre-vacuum sterilization parameters in Table 3 is recommended.

Table 3: Sterilization Parameters - Outside of USA

	Pre-vacuum
Temperature	134°C
Exposure time	3 minutes
Drying time	30 minutes

Table 4: Rigid Container Configuration - Outside of USA

Aesculap Part Number	Description
JK442	Full-size container, solid bottom (6-inch height)
JK488	Full-size container lid, gold
JK489	Full-size container lid, silver
JK095	Single use filter

9 Maintenance

⚠ WARNING: To ensure correct function, do not service the product other than lubrication. If service is required, contact your Stryker sales representative.

Lubrication Instructions

If moveable parts of the product grind, make noise or become difficult to move, lubricate the product. Repeat lubrication as required.

1. Lubricate moveable parts with a drop of light mineral oil REF 1605-010-000 as shown in Figure 7.
2. Wipe off excessive oil.
3. Sterilize the product.

10 Technical Specifications

Table 5: Technical Specifications

	Size Height/Width/Length	Material
Navigated Short Drill Guide	143 mm/22 mm/220 mm	Stainless steel
Navigated Long Drill Guide	143 mm/22 mm/270 mm	Stainless steel
Navigated Short Drill Guide Calibrator	18 mm/18 mm/160 mm	Stainless steel
Navigated Long Drill Guide Calibrator	18 mm/18 mm/210 mm	Stainless steel
Navigated Short Drill Bits	6 mm/6 mm/226 mm	Stainless steel

	Size Height/Width/Length	Material
Navigated Long Drill Bits	6 mm/6 mm/307 mm	Stainless steel
Navigated Drill Bit Stop	27 mm/25 mm/14 mm	Stainless steel

11 Definition of Symbols

The following table defines the symbols used in this document, on the product and on the product label.

Table 6: Definition of Symbols

ANSI Z535.6 American National Standard for Product Safety Information in Product Manuals, Instructions, and Other Collateral Materials

Symbol	Name: Definition
	Safety alert symbol: Alerts the user to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.

EN ISO 15223-1 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied - Part 1 General requirements

Symbol/ number	Name: Definition
 5.1.1	Manufacturer: Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
 5.1.3	Date of manufacture: Indicates the date when the medical device was manufactured.
 5.1.6	Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.
 5.4.3	Consult instructions for use: Indicates the need for the use to consult the instructions for use.
 5.1.5	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
 5.2.7	Non-Sterile: Indicates a medical device that has not been subjected to a sterilization process.
 5.4.2	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

Device-Specific Symbols

Symbol	Name: Definition
QTY	Quantity: Indicates the number of medical devices in the packaging.
S	Stainless steel alloys: Indicates that the device contains stainless steel alloys.

21 Code of Federal Regulations (CFR), section 801.109(b)(1)

Symbol	Name: Definition
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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for the following trademarks or service
All other trademarks are trademarks of
Copyright © 2017 Stryker
Stryker Corporation or its divisions or other
Stryker Leibinger GmbH & Co. KG