## INSTRUCTIONS FOR USE

## Reinforced Anesthesia Needle

REF RAN-027-2  
RAN-027-5

0459

Entellus Medical Inc.  
3600 Holly Lane North, Suite 40  
Plymouth, MN 55447 USA  
www.ent.stryker.com

EC REP  
ICON (LR) Limited  
South County Business Park  
Leopardstown, Dublin 18,  
D18 X5R3, Ireland  
patent http://www.ent.stryker.com/patents  
Entellus Medical, Inc., a Stryker Company

2021-07

5117-001-r8

## EN

INSTRUCTIONS FOR USE  
Reinforced Anesthesia Needle

Read all Instructions prior to use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Sterility: Provided Sterile, Gamma Irradiation Sterilization

Single Use: Disposable, For Single Patient Use Only, Do Not

Resterilize and/or Reuse

Storage: Store in a cool, dry place.

Symbol Definitions: Refer to Table 1.

## Indication for Use

For use in injecting local anesthetics into a patient to provide regional anesthesia.

## Device Description

The Reinforced Anesthesia Needle is a 3.5" long disposable medical device that allows the user to transnally administer anesthetic solutions. The Reinforced Anesthesia Needle device is a 27-gauge needle with a reinforcing sleeve that supports the needle to reduce needle flex during administration of anesthetic solutions. A slight bend at the distal end of the needle is designed to improve access to nasal anatomy. The Quincke tip extends approximately 2 mm past the reinforcing sleeve creating a distal stop. The needle hub features a standard luer-lock connector. The needle may be reshaped along its length. Refer to Figure 1.

## Contents

The Reinforced Anesthesia Needle consists of:

- 1 needle
- 1 needle protector sleeve

## Contraindications

None known.

## Intended Users

Trained ear, nose and throat (ENT) physicians.

## Warnings

- Visually inspect all sterile barrier systems immediately prior to use in order to determine if breaches in sterile barrier system integrity are evident. Do not use breached or damaged packages, since the sterility and functionality of the device may be compromised.
- Single use only. Do not resterilize or reuse, as it may result in compromised device performance and risk improper sterilization and cross contamination.
- Read all warnings and precautions in pharmaceutical literature prior to administration of anesthetic solutions; administer per state and local regulations.
- Do not over-insert the needle when injecting anesthetic – this may increase risk of delivering anesthetic solutions directly to blood stream.
- If the needle tip is bent during procedure, remove needle and replace with new device.
- When injecting anesthetic, if working channel is occluded, remove needle and replace with new device.
- To avoid infection, dispose of needle safely after single use, according to federal, state, and/or local regulations and procedures for disposal of sharps and/or biomedical waste. Biomedical waste should be considered infectious and requires special management/treatment.

## Precautions

- Maintain aseptic conditions during preparation of the needle and throughout the procedure.
- Handle all used needles according to current blood borne pathogens procedures, such as Standard Precautions.
- Do not attempt to reinsert needle into needle protector sleeve.
- Do not reshape needle beyond a 25 degree bend angle as this may cause device damage.
- Administration of anesthetic solutions may cause elevated heart rate – monitor patient accordingly.

## Adverse Effect

Possible adverse effects include, but are not limited to, the following:

- Complications from anesthetic solutions
- Infection
- Tissue inflammation, irritation, swelling or trauma
- Bleeding
- Needle stick injury
- Orbital damage

## Supplies

The following supplies are not provided with the Reinforced Anesthesia Needle and should be available for the procedure: medical syringes, injectable anesthetic solutions, sharps disposal container, and gloves.

## Instructions for Use

1. Fill a medical syringe with anesthetic solution.
2. Remove the Reinforced Anesthesia Needle device from packaging.
3. Connect the needle hub onto filled syringe and discard clear needle protector.
4. Carefully insert needle transnally into tissue at desired location.
  - The needle may be reshaped up to a 25 degree bend to access various anatomy
  - To avoid unnecessary tissue trauma, do not insert the needle into the tissue beyond the distal stop, which is 1.5 – 2 mm from the needle distal tip.

5. Administer anesthetic solution into desired location.
  - The same needle may be used multiple times within the same patient – repeat injections at additional sites within anatomy as necessary.
6. Remove the Reinforced Anesthesia Needle.
7. Discard used needle in accordance with federal, state or local regulations and procedures for disposal of sharps and/or biomedical waste. Biomedical waste should be considered infectious and requires special management/treatment

## Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to Entellus Medical by calling Customer Service at (866) 620-7615. If dialing outside the United States, precede the number with your country's exit code and the US country code (+1). EU Member States should also notify the competent authority of the Member State in which the incident occurred.

## Limited Warranty

Refer to Entellus Medical, Inc. Standard Terms and Conditions.

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