

# Augment<sup>®</sup>

## Regenerative Solutions



Augment  
Bone Graft



Augment  
Injectable

# Choosing the appropriate graft for Hindfoot and Ankle Fusion Surgery

With the myriad of bone graft options available and continuing to grow, it has become increasingly difficult for surgeons and hospitals to avoid confusion and maintain a clear understanding of the clinical value, specific features, benefits, and appropriate uses of available graft choices.

This guide is intended to offer some perspective on how to better understand the three basic categories of commercially available bone graft products, and to offer a method of applying four simple, overarching questions that will allow you to make better informed purchasing decisions:

- 1 What is the evidence?
- 2 How was it approved?
- 3 How should it be used?
- 4 How does it work?

	Categorization	Characteristics	Regulatory pathway and clinical burden of proof	Examples
<b>Recombinant growth factors</b>	Proven alternative to autograft <sup>1</sup>	True biologic: Combination products with consistent, highly active signal proteins that drive bone regeneration <sup>4</sup>	Pre-market approval (PMA) Safety and efficacy in large pivotal clinical trial <sup>5</sup>	Augment Regenerative Solutions (rhPDGF-BB/β-TCP)  InFuse Bone Graft (rhBMP-2/ACS)
<b>Allograft tissue</b>	Bone void fillers <sup>2</sup>	Osteoconductive and weakly osteoinductive putties with variable handling  Cell-containing products include various claims of cell viability at point of use	Human tissue products (361 HCT/Ps) – No proof needed <sup>6</sup>  Devices containing human tissue (351 HCT/Ps—more than minimally manipulated)—510(k) <sup>11</sup>	BIO4, Allomatrix, DBX, Trinity Elite, Accell, Osteocel Plus, Grafton, Osteosponge, Fusionflex, Allograft chips, femoral head allograft, Allopure wedges
<b>Synthetic scaffolds</b>	Bone void fillers <sup>3</sup>	Passive osteoconduction and fills space	Synthetics—510(k) <sup>7</sup>	Pro-Dense, Osteoset, Vitoss BA, Hydroset XT, Norian Drillable Inject

1. Combination Growth Factor products composed of biologically active signals that promote chemotaxis, mitogenesis, angiogenesis, and/or osteoinductivity, rigorously reviewed by FDA and proven to be non-inferior to the Gold Standard, autograft in specific approved indications.
2. Void fillers with mineralized or demineralized bone, with or without cryopreserved cells from same donor intended for treatment of musculoskeletal defects.
3. Physical scaffolds composed of synthetic materials (e.g. calcium phosphate) intended to be used to fill bone voids.
4. Platelet-derived growth factor (rhPDGF-BB) and bone morphogenetic protein-2 (rhBMP-2) have been tested, reviewed and established as alternatives to autograft in multiple clinical studies for specific indications in foot & ankle (rhPDGF-BB only), spine/orthopaedic trauma (rhBMP-2 only) and dental (rhPDGF-BB & rhBMP-2).
5. Multiple clinical trials culminating in a large pivotal trial are typically required to prove that the combination device is both safe and efficacious in the specified indications.
6. Human tissues designated as 361 HCT/Ps are not regulated as medical devices and do not require a submission and/or review for commercialization. Tissue processors are required to register with FDA and follow Good Tissue Practices (GTP) per 21 CFR 1271.
7. A 510(k) clearance demonstrates that the device is substantially equivalent to a legally marketed device.

## 1 What is the evidence?

**Augment is the first and only** proven alternative to autograft in hindfoot and ankle arthrodesis.

- The role of PDGF in bone repair and regeneration is documented in over 20 peer-reviewed publications including three large-scale randomized, controlled clinical trials involving rhPDGF-BB.
- Level-I clinical data is considered the most reliable quality of data based upon standards for peer-reviewed evidence.



## 2 How was it approved?

**Augment** was approved via the rigorous PMA pathway by the U.S. FDA as a Class III combination medical device/drug product.

- The PMA application was supported by two pilot clinical trials, the largest prospective, randomized, controlled clinical trial in foot and ankle history, and numerous preclinical studies. The PMA supplement for Augment Injectable was supported by two additional randomized, controlled clinical trials.




● Trial sites

## 3 How should it be used?

Among the criteria for which patients undergoing hindfoot and ankle arthrodesis should receive **Augment**, the following factors should be considered:

- In the randomized, controlled pivotal trial conducted to support U.S. FDA approval of Augment Bone Graft, 75% of the treated patients had one or more risk factors for non-union. Augment Bone Graft treated patients were found to have equivalent improvements in clinical outcomes and a better overall safety profile versus autograft (due to the elimination of harvest site pain and morbidity).<sup>1,\*\*</sup>
- A natural deficiency in PDGF-BB has been found to be correlated with non-union and poor bone formation in clinical and pre-clinical studies involving subjects with diabetes, osteoporosis and smokers.<sup>2,3,4</sup>
- Unlike growth factors such as the bone morphogenetic proteins (BMPs), which lead solely to osteoblastic differentiation of cells at the implantation site, there are no reported incidences of ectopic bone formation.<sup>5</sup>

### Augment® vs. autograft results



	Equivalent to Autograft	Compared to Autograft
Clinical healing	✓	
CT Fusion Rates*	✓	
Therapeutic Failure Rate	✓	
No Graft Harvest Site Pain		✓
2x odds of Fusion success ≥65 years old <sup>2</sup>		✓

\*FDA did not base its approval of Augment Bone Graft on radiologic findings from the pivotal study, but instead relied on clinical outcomes.

- DiGiovanni CW, et al., JBJS (2013)
- Verma, et al., Curr Orthop Pract (2011)
- Hollinger, et al., JOR (2008)
- Al-Zube, et al., J Orthop Res (2009)
- Carragee EJ, et al., Spine J (2011)

\*\* A majority of patients studied in the two pivotal studies evaluating Augment Injectable also presented with one or more risk factors for non-union.

# 4 How does it work?

## Implantation

Following preparation of the bony surfaces, Augment is applied to the fusion site. The rhPDGF-BB releases from the carrier, forming a concentration gradient as it migrates throughout the local environment.




## Completion of the Bone Formation Process

These mature osteoblasts will then lay down new bone to create a continuous scaffold, fusing the bony surfaces.

## Chemotaxis

(cell movement)


Mesenchymal Stem Cells (MSCs) are attracted to the fusion site by the increase concentration of rhPDGF-BB from bleeding bone, muscle and the periosteum.



## Mitogenesis

(cell division)


MSCs are stimulated to divide and proliferate in the presence of the higher concentration of rhPDGF-BB in the graft site.



## Morphogenesis

(Cell Transformation)


Once the colony of MSCs has divided repeatedly, native bone morphogenetic proteins (BMPs) induce MSCs to mature into osteoblasts.



## Angiogenesis

(blood vessel formation)

In parallel with the effects of rhPDGF-BB in the bone formation cascade, the protein also promotes angiogenesis by increasing vascular endothelial cell, pericyte, and smooth muscle responses. Pericytes then synthesize VEGF, thereby enhancing the neovascular drive.




These newly formed blood vessels support the formation of bone by supplying oxygen and nutrients, carrying additional cells and signals to the healing environment, and eliminating local waste.

# Competitive summary

	Augment	Infuse Bone Graft	DBM w/ allogenic cells	DBM	Allograft chips + BMA
Level I evidence demonstrating safe and effective use as alternative to autograft	✓	✓			
FDA approved for hindfoot and ankle fusions	✓				
MOA: Works early in healing cascade to impact multiple tissues involved in clinical healing (more than bone)	✓				
Consistent and reliable concentration, composition and activity	✓	✓			
Easy to administer/apply to surgical site	✓	✓		✓	
Off-the-shelf (minimal preparation required)	✓	✓		✓	
Regulatory classification	<b>Class III – PMA</b> Required safety and efficacy evidence required	<b>Class III – PMA</b> Required safety and efficacy evidence required	HCT/P (human tissue) for homologous use No safety or efficacy evidence required	<b>510k clearance</b> No clinical trials HCT/P (Human Tissue) for homologous use No safety or efficacy evidence required	HCT/P (human tissue) for homologous use No safety or efficacy evidence required

# The same powerful protein **rhPDGF-BB**

Augment<sup>®</sup>  
Injectable



- Carrier  $\beta$ -TCP is 70-95% smaller than  $\beta$ -TCP in Augment<sup>®</sup> Bone Graft
- $\beta$ -TCP is paired with collagen matrix to make flowable
- Resorbs and replaced with bone as site heals

Augment<sup>®</sup>  
Bone Graft



- Carrier  $\beta$ -TCP is an osteomimetic scaffold with long history of orthopedic use
- Resorbs and replaced with bone as site heals

# Augment

## Peer-reviewed Bibliography

### Peer-reviewed clinical evidence

- Loveland JD, McMillen RL, Cala MA. A Multicenter, Retrospective, Case Series of Patients With Charcot Neuroarthropathy Deformities Undergoing Arthrodesis Utilizing Recombinant Human Platelet-derived Growth Factor With Beta-Tricalcium Phosphate. *The Journal of Foot and Ankle Surgery*. 2021 Jan 1;60(1):74-9.
- Berlet GC, Baumhauer JE, Glazebrook M, Haddad SL, Younger A, Quiton JD, Fitch DA, Daniels TR, DiGiovanni CW. The impact of patient age on foot and ankle arthrodesis supplemented with autograft or an autograft alternative (rhPDGF-BB/ $\beta$ -TCP). *JBJS Open Access*. 2020 Oct;5(4).
- Loveland JD, Basile P, Collier BN, Manning ES. Revisional hindfoot and ankle arthrodesis using recombinant human platelet-derived growth factor and beta-tricalcium phosphate. *Clin Res Foot Ankle* 2020; 8(4):302-306.
- Baumhauer JE, Glazebrook M, Younger A, Quiton JD, Fitch DA, Daniels TR, DiGiovanni CW. Long-term autograft harvest site pain after ankle and hindfoot arthrodesis. *Foot & Ankle International*. 2020 Aug;41(8):911-5.
- Daniels TR, Anderson J, Swords MP, Maislin G, Donahue R, Pinsker E, Quiton JD. Recombinant Human Platelet-Derived Growth Factor BB in Combination With a Beta-Tricalcium Phosphate (rhPDGF-BB/ $\beta$ -TCP)-Collagen Matrix as an Alternative to Autograft. *Foot & ankle international*. 2019 Sep;40(9):1068-78.
- Scott RT, McAlister JE, Rigby RB. Allograft Bone\*: What Is the Role of Platelet-Derived Growth Factor in Hindfoot and Ankle Fusions. *Clin Podiatr Med Surg*. 2018 Jan;35(1):37-52. doi: 10.1016. PMID: 29156166
- DiGiovanni CW, Lin SS, Daniels TR, Glazebrook M, Evangelista P, Donahue R, Beasley W, Baumhauer JE. The Importance of Sufficient Graft Material in Achieving Foot or Ankle Fusion. *J Bone Joint Surg Am*. 2016 Aug 3;98(15):1260-7. doi: 10.2106/JBJS.15.00879. PMID: 27489316
- Krell ES, Digiovanni CW. The Efficacy of Platelet-Derived Growth Factor as a Bone-Stimulating Agent. *FootAnkle Clin*. 2016;21(4):763-770. doi: 10.1016/j.fcl.2016.07.002. PMID: 27871409
- Younger A, Penner M, Montijo HE. Vancouver Experience of Recombinant Human Platelet-Derived Growth Factor. *Foot Ankle Clin*. 2016;21(4):771-776. doi: 10.1016/j.fcl.2016.07.009. PMID: 27871410
- Dodd A, Daniels TR. Injectable Recombinant Human Platelet-derived Growth Factor in Collagen Carrier for Hindfoot Fusion. *Foot Ankle Clin*. 2016;21(4):777-791. doi: 10.1016/j.fcl.2016.07.013. PMID: 27871411
- Daniels TR, Younger AS, Penner MJ, Wing KJ, Le IL, Russell IS, Lalonde KA, Evangelista PT, Quiton JD, Glazebrook M, DiGiovanni CW. Prospective Randomized Controlled Trial of Hindfoot and Ankle Fusions Treated With rhPDGF-BB in Combination With a  $\beta$ -TCP-Collagen Matrix. *Foot Ankle Int*. 2015 Jul;36(7):739-48. doi: 10.1177/1071100715576370. PMID: 25848134
- Glazebrook M, Beasley W, Daniels T, Evangelista PT, Donahue R, Younger A, Pinzur MS, Baumhauer JE, DiGiovanni CW. Establishing the Relationship Between Clinical Outcome and Extend of Osseus Bridging Between Computed Tomography Assessment in Isolated Hindfoot and Ankle Fusions. *Foot Ankle Int*. 2013 Dec;34(12):1612-8. doi: 10.1177/1071100713504746. Epub 2013 Sep 16. PMID: 24043351
- Baumhauer J, Pinzur MS, Donahue R, Beasley W, DiGiovanni C. Site Selection and Pain Outcome After Autologous Bone Graft Harvest. *Foot Ankle Int*. 2014 Feb;35(2):104-7. doi: 10.1177/1071100713511434. doi: 10.1177/1071100713511434. PMID: 24227683
- DiGiovanni CW, Lin S, Baumhauer JE, Daniels T, Younger A, Glazebrook M, Anderson J, Evangelista P, Lynch SE, et al. Recombinant Human Platelet-Derived Growth Factor-BB and Beta-Tricalcium Phosphate (rhPDGF-BB/ $\beta$ -TCP): An Alternative to Autogenous Bone Graft. *J Bone Joint Surg Am*. 2013 Jul 3;95(13):1184-92. doi: 10.2106/JBJS.K.01422. PMID: 23824386
- Note: FDA did not base its approval of Augment Bone Graft on radiologic findings from the pivotal study, but instead relied on clinical outcomes.**
- Friedlaender GE, Lin S, Solchaga LA, Snel LB, Lynch SE. The Role of Recombinant Human Platelet Derived Growth Factor BB(rhPDGF-BB) in Orthopaedic Bone Repair and Regeneration. *Curr Pharm Des*. 2013;19(19):3384-90. PMID: 23432673
- Solchaga LA, Daniels T, Roach S, Beasley W, Snel LB. Effect of implantation of Augment Bone Graft on serum concentrations of platelet-derived growth factors: a pharmacokinetic study. *Clin Drug Investig*. 2013 Feb;33(2):143-9. doi: 10.1007/s40261-013-0053-5. PMID: 23334906
- Glazebrook M; Daniels T; Abidi N; Younger A; DiGiovanni CW. Role of Platelet-derived Growth Factor in Hindfoot Fusion. *Techniques in Foot & Ankle Surgery*. 11(1):34-38, March 2012. doi: 10.1097/BTF.0b013e3182469626
- DiGiovanni CW, Lin S, Pinzur M. Recombinant human PDGF-BB in foot and ankle fusion. *Expert Rev Med Devices*. 2012 Mar;9(2):111-22. doi: 10.1586/erd.11.76. Review. PMID: 22404772
- DiGiovanni CW, Baumhauer J, Lin SS, Berberian WS, Flemister AS, Enna MJ, Evangelista P, Newman J. Prospective, Randomized, Multi-center Feasibility Trial of rhPDGF-BB Versus Autologous Bone Graft in a Foot and Ankle Fusion Model. *Foot Ankle Int*. 2011 Apr;32(4):344-54. doi: 10.3113/FAL.2011.0344. PMID: 21733435
- DiGiovanni CW, Petricek JM. The Evolution of rhPDGF-BB in Musculoskeletal Repair and Its Role in Foot and Ankle Fusion Surgery. *Foot Ankle Clin*. 2010 Dec;15(4):621-40. doi: 10.1016/j.fcl.2010.07.001. Review. PMID: 21056861
- Daniels T, DiGiovanni C, Lau JT, Wing K, Younger A. Prospective Clinical Pilot Trial in a Single Cohort Group of rhPDGF in Foot Arthrodeses. *Foot Ankle Int*. 2010 Jun;31(6):473-9. doi: 10.3113/FAL.2010.0473. PMID: 20557811
- Hollinger JO, Hart CE, Hirsch SN, Lynch S, Friedlaender GE. Recombinant Human Platelet-derived Growth Factor: Biology and Clinical Applications. *J Bone Joint Surg Am*. 2008 Feb;90 Suppl 1:48-54. doi: 10.2106/JBJS.G.01231. Review. PMID: 18292357
- Nevins M, Giannobile WV, McGuire MK, Kao RT, Mellonig JT, Hinrichs JE, McAllister BS, Murphy KS, McClain PK, Nevins ML, Paquette DW, Han TJ, Reddy MS, Lavin PT, Genco RJ, Lynch SE. Platelet-derived Growth Factor Stimulates Bone Fill and Rate of Attachment Level Gain: Results of a Large Multicenter Randomized Controlled Trial. *J Periodontol*. 2005 Dec;76(12):2205-15. PMID: 16332231

### Peer-reviewed pre-clinical evidence

- Gadomski B, C. Labus, K. M., Puttlitz, C. M., McGillvray, K. C., Regan, D. P, Nelson, B., Seim, H. B. III, & Easley, J. T. (2021). Evaluation of lumbar spinal fusion utilizing recombinant human platelet derived growth factor-B chain homodimer (rhPDGF-BB) combined with a bovine collagen/ $\beta$ -tricalcium phosphate ( $\beta$ -TCP) matrix in a bovine model. *JOR Spine*, 4(3), e1166.
- Yamaguchi, J. T., Weiner, J. A., Minardi, S., Greene, A. C., Ellenbogen, D. J., Hallman, M. J., Shah, V. P., Weisz, K. M., Jeong, S., Nandurkar, T., Yun, C., Hsu, W. K., & Hsu, E. L. (2021). Characterizing the host response to rhPDGF-B in a rat spinal arthrodesis model. *JOR Spine*, e1173.
- Perrien DS, Young CS, Alvarez-urena PE, Dean DD, Lynch SE, Hollinger JO. Percutaneous Injection of Augment Injectable Bone Graft (rhPDGF-BB and  $\beta$ -tricalcium Phosphate [ $\beta$ -TCP]/Bovine Type I Collagen Matrix) Increases Vertebral Bone Mineral Density in Geriatric Female Baboons. *Spine J*. 2013;13(5):580-6. doi: 10.1016/j.spinee.2013.02.029. PMID: 23541446
- Solchaga LA, Hee CK, Roach S, Snel LB. Safety of Recombinant Human Platelet derived Growth Factor-BB in Augment Bone Graft. *J Tissue Eng*. 2012;3(1):2041731412442668. doi: 10.1177/2041731412442668. PMID: 22511993; PMCID: PMC3324841
- Caplan AI, Correa D. PDGF in Bone Formation and Regeneration: New Insights Into a Novel Mechanism Involving MSCs. *J Orthop Res*. 2011 Dec;29(12):1795-803. doi: 10.1002/jor.21462. Epub 2011 May 25. PMID: 21618276
- Verma R, Koerner J, Breitbart E, Paglia D, Vaidya S, Pinzur M, Sheldon L. Correlation of Growth Factor Levels at the Fusion Site of Diabetic Patients Undergoing Hindfoot Arthrodesis and Clinical Outcome. *Current Ortho Practice* May/June 2011; 22(3):251-56. doi: 10.1097/BCO.0b013e318212584c
- Young CS, Bradica G, Hart CE, Karunanidhi A, Street RM, Schutte L, Hollinger JO. Preclinical Toxicology Studies of Recombinant Human Platelet-Derived Growth Factor-BB Either Alone or in Combination with Beta-Tricalcium Phosphate and Type I Collagen. *J Tissue Eng*. 2011 Jan 10;2010:246215. doi: 10.4061/2010/246215. PMID: 21350649
- Pountos I, Georgouli T, Henshaw K, Bird H, Jones E, Giannoudis PV. The Effect of Bone Morphogenetic Protein-2, Bone Morphogenetic Protein-7, Parathyroid Hormone, and Platelet-derived Growth Factor on the Proliferation and Osteogenic Differentiation of Mesenchymal Stem Cells Derived from Osteoporotic Bone. *J Orthop Trauma*. 2010 Sep;24(9):552-6. doi: 10.1097/BOT.0b013e3181efa8fe. PMID: 20736793
- Moore DC, Ehrlich MG, McAllister SC, Machan JT, Hart CE, Voigt C, Lesieur-Brooks AM, Weber EW. Recombinant Human Platelet-derived Growth Factor-BB Augmentation of New-bone Formation in a Rat Model of Distraction Osteogenesis. *J Bone Joint Surg Am*. 2009 Aug;91(8):1973-84. doi: 10.2106/JBJS.H.00540. PMID: 19651957
- Al-Zube I, Breitbart EA, O'Connor JP, Parsons JR, Bradica G, Hart CE, Lin SS. Recombinant Human Platelet-derived Growth Factor BB (rhPDGF-BB) and Betatricalcium Phosphate/Collagen Matrix Enhance Fracture Healing in a Diabetic Rat Model. *J Orthop Res*. 2009 Aug;27(8):1074-81. doi: 10.1002/jor.20842. PMID: 19170096
- McCarthy HS, Williams JH, Davie MW, Marshall MJ. Platelet-derived Growth Factor Stimulates Osteoprotegerin Production in Osteoblastic Cells. *J Cell Physiol*. 2009 Feb;218(2):350-4. doi: 10.1002/jcp.21600. PMID: 8814141
- Tokunaga A, Oya T, Ishii Y, Motomura H, Nakamura C, Ishizawa S, Fujimori T, Nabeshima Y, Umezawa A, Kanamori M, Kimura T, Sasahara M. PDGF Receptor Beta Is a Potent Regulator of Mesenchymal Stromal Cell Function. *J Bone Miner Res*. 2008 Sep;23(9):1519-28. doi: 10.1359/jbmr.080409. PMID: 18410236
- Hollinger JO, Onikepe AO, MacKrell J, Einhorn T, Bradica G, Lynch S, Hart CE. Accelerated Fracture Healing in the Geriatric, Osteoporotic Rat with Recombinant Human Platelet-derived Growth Factor-BB and an Injectable Beta-tricalcium Phosphate/Collagen Matrix. *J Orthop Res*. 2008 Jan;26(1):83-90. PMID: 17676626
- Sanchez-Fernandez MA, Gallois A, Riedl T, Jurdic P, Hoffack B. Osteoclasts Control Osteoblast Chemotaxis via PDGF-BB/PDGF Receptor Beta Signaling. *PLoS One*. 2008;3(10):e3537. doi: 10.1371/journal.pone.0003537. PMID: 18953417; PMCID: PMC2569415
- Ozaki Y, Nishimura M, Sekiya K, Suehiro F, Kanawa M, Nikawa H, Hamada T, Kato Y. Comprehensive Analysis of Chemotactic Factors for Bone Marrow Mesenchymal Stem Cells. *Stem Cells Dev*. 2007 Feb;16(1):119-29. PMID: 15660418
- Kilian O, Alt V, Heiss C, Jonuleit T, Dingeldein E, Flesch I, Fidorra U, Wenisch S, Schnettler R. New Blood Vessel Formation and Expression of VEGF Receptors after Implantation of Platelet Growth Factor-enriched Biodegradable Nanocrystalline Hydroxyapatite. *Growth Factors*. 2005 Jun;23(2):125-33. PMID: 16019434
- Fiedler J, Etzel N, Brenner RE. To Go or Not To Go: Migration of Human Mesenchymal Progenitor Cells Stimulated by Isoforms of PDGF. *J Cell Biochem*. 2004 Nov 15;93(5):990-8. PMID: 15389881
- Mehrotra M, Krane SM, Walters K, Pilbeam C. Differential Regulation of Platelet-derived Growth Factor Stimulated Migration and Proliferation in Osteoblastic cells. *J Cell Biochem*. 2004 Nov 1;93(4):741-52. PMID: 15660418
- Kilian O, Flesch I, Wenisch S, Taborski B, Jork A, Schnettler R, Jonuleit T. Effects of Platelet Growth Factors on Human Mesenchymal Stem Cells and Human Endothelial Cells In Vitro. *Eur J Med Res*. 2004 Jul 30;9(7):337-44. PMID: 15337634
- Bouletreau PJ, Warren SM, Spector JA, Steinbrech DS, Mehrara BJ, Longaker MT. Factors in the Fracture Microenvironment Induce Primary Osteoblast Angiogenic Cytokine Production. *Plast Reconstr Surg*. 2002 Jul;110(1):139-48. PMID: 12087245
- Fiedler J, Roderer G, Günther KP, Brenner RE. BMP-2, BMP-4, and PDGF-BB Stimulate Chemotactic Migration of Primary Human Mesenchymal Progenitor Cells. *J Cell Biochem*. 2002;87(3):305-12. PMID: 12397612
- Zhang Z, Chen J, Jin D. Platelet-derived Growth Factor (PDGF)-BB Stimulates Osteoclastic Bone Resorption Directly: The Role of Receptor Beta. *Biochem Biophys Res Commun*. 1998 Oct 9;251(1):190-4. PMID: 9790928

# Trust Augment®

## Proven

- Level 1 evidence from largest F&A clinical trial ever conducted
- Equivalent improvements in clinical outcomes, compared to the gold standard autograft<sup>1</sup>

## Labeled

- Class III combination product labeled for ankle and hindfoot arthrodesis

## Unique

- Bioengineered human PDGF-BB stimulates multiple aspects of healing in response to injury
- rhPDGF-BB is highly purified with consistent biological potency; allograft and autograft are highly variable in quality and potency<sup>2,3,4,5</sup>

## Safe

- In commercial use since 2009 (Canada)
- Eliminates risks, morbidities, and costs associated with autograft harvest

1. DiGiovanni, et al. JBJS (2013).  
2. Fiedler, et al. J Cell Biochem (2002).  
3. Ozaki, et al. J Stem Cells & Dev (2007).  
4. Bouletreau, et al. Plast Reconstr Surg (2002).  
5. Hollinger, et al. JBJS (2008).



Part no.	Description	Volume
K30003010	Augment Injectable	3.0cc
K30001510	Augment Injectable	1.5cc
K20003010	Augment Bone Graft	3.0cc
K20001510	Augment Bone Graft	1.5cc

This document is intended solely for healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Augment, Stryker. All other trademarks are trademarks of their respective owners or holders.

Content ID: AP-017085  
Copyright © 2023 Stryker

**Manufactured by:**

BioMimetic Therapeutics, LLC  
389 Nichol Mill Lane  
Franklin, TN 37067

[stryker.com](http://stryker.com)