



DARCO™

3.2mm Headless Compression Screw



SURGICAL TECHNIQUE

---

**DARCO™**

**3.2mm Compression  
Screws**

**SURGICAL TECHNIQUE**

---

**Surgical Technique as described by:**

Robert Anderson, MD  
Bruce Cohen, MD  
W. Hodges Davis, MD

---

# Contents

<b>Chapter 1</b>	<b>4</b>	<b>Introduction</b>
		Device Description
		Indications
<b>Chapter 2</b>	<b>5</b>	<b>Surgical Technique</b>
<b>Appendix A</b>	<b>6</b>	<b>Explant Information</b>
<b>Appendix B</b>	<b>7</b>	<b>Ordering Information</b>

Proper surgical procedures and techniques are the responsibility of the medical professional. The following guidelines are furnished for information purposes only. Each surgeon must evaluate the appropriateness of the procedures based on his or her personal medical training and experience. Prior to use of the system, the surgeon should refer to the product package insert for complete warnings, precautions, indications, contraindications and adverse effects. Package inserts are also available by contacting Wright Medical Technology, Inc. Contact information can be found on the back of this surgical technique and the package insert is available on the website listed.

Please contact your local Wright representative for product availability.

# Introduction



## Device Description

The DARCO™ 3.2mm Headless Screw is a self drilling, self tapping and self countersinking 3.2mm hex drive cannulated screw with a reverse cutting nib designed into the thread pattern. The cannulated feature allows for the use of a drill guide for precise placement while the smooth shank between the threaded portions of the screw allows the bone surfaces to be compressed to facilitate healing. The screws are made of Ti 6-Al 4-V biocompatible titanium alloy and coated with an anodized finish.

## Indications

The DARCO™ 3.2mm Headless Screw is to be used on indications that are common for currently marketed compression screws. The primary indication for use is the fixation and stabilization of fractures and non-unions of small bones and small bone arthrodeses including but not limited to intra-articular fractures of the tarsals, metatarsals, carpals and metacarpals, bunionectomies and osteotomies, and arthrodeses of small joints (i.e. phalanges).

## Contraindications

No product specific contraindications.

Prior to use of the system, the surgeon should refer to the product package insert for complete warnings, precautions, indications, contraindications and adverse effects. Package inserts are also available by contacting the manufacturer. Contact information can be found on the back of this surgical technique and the package insert is available on the website listed.

## K-wire Placement

Use a powered driver to place a 1mm K-Wire (P/N N02228210) through the osteotomy or fusion site.



## Screw Selection

Using the ruler (P/N MT1009), measure the length of the K-wire exposed from the bone. Subtract this length from the total length of the K-wire (100mm) to obtain the screw length needed.



**Note:** The DARCO™ 3.2mm Headless Compression Screw is designed to be self-tapping and self-drilling. However, in dense cortical bone, drilling and countersinking may be required. In these cases, use the provided drill/countersink instruments over the K-wire prior to screw insertion.

## Screw Placement

Load the 2mm Cannulated Hex Driver (P/N 44112001) into the AO driver handle (P/N 44112009). Use the driver to advance the appropriate screw until the head is completely countersunk within the bone.



# Explant Information

If the removal of the implant is required due to revision or failure of the device, the surgeon should contact the manufacturer using the contact information located on the back cover of this surgical technique to receive instructions for returning the explanted device to the manufacturer for investigation.

## Postoperative Care

Postoperative care is the responsibility of the medical professional

# Ordering Information



## DARCO™

### 3.2mm Screw

NON-STERILE #	DESCRIPTION
DC4110	3.2mm X 10mm
DC4112	3.2mm X 12mm
DC4114	3.2mm X 14mm
DC4116	3.2mm X 16mm
DC4118	3.2mm X 18mm
DC4120	3.2mm X 20mm
DC4122	3.2mm X 22mm
DC4124	3.2mm X 24mm
DC4126	3.2mm X 26mm
DC4128	3.2mm X 28mm
DC4130	3.2mm X 30mm
DC4132	3.2mm X 32mm

## DARCO™

### SNAP-Off Screws

PART NUMBER	DESCRIPTION
45110001	2.0mm X 11mm
45110002	2.0mm X 12mm
45112013	2.0mm X 13mm
45110003	2.0mm X 14mm
45112015	2.0mm X 15mm
45112711	2.7mm x 11mm
45112712	2.7mm x 12mm
45110004	2.7mm x 13mm
45112714	2.7mm x 14mm
45110005	2.7mm x 15mm

## DARCO™

### Instrumentation

PART NUMBER	DESCRIPTION
NO2228210	K-wire 1.0mm
MT10010	Drill Bit 33mm (GOLD)
MT10020	Drill Bit 24mm (BLUE)
MT10160	Drill Bit 15mm (RED)
DC2602016	Scarf Fixation Forceps
MT1009	Ruler
MT1010	Tweezers
2225000	Gauge for screws & k-wire
MT1015	Sterilization Box
44112001	2.0mm Hex
41112017	AO Quick Connect
44112009	AO Quick Connect Handle
45112004	Snap-Off Screwdriver AO-QC





1023 Cherry Road  
Memphis, TN 38117  
800 238 7117  
901 867 9971  
[www.wright.com](http://www.wright.com)

62 Quai Charles de Gaulle  
69006 Lyon  
France  
+33 (0)4 72 84 10 30  
[www.tornier.com](http://www.tornier.com)

18 Amor Way  
Letchworth Garden City  
Hertfordshire SG6 1UG  
United Kingdom  
+44 (0)845 833 4435

Proper surgical procedures and techniques are the responsibility of the medical professional. This material is furnished for information purposes only. Each surgeon must evaluate the appropriateness of the material based on his or her personal medical training and experience. Prior to use of any Tornier implant system, the surgeon should refer to the product package insert for complete warnings, precautions, indications, contraindications, and adverse effects. Package inserts are also available by contacting Wright. Contact information can be found in this document and the package insert.

™ and ® denote Trademarks and Registered Trademarks of Wright Medical Group N.V. or its affiliates.  
©2016 Wright Medical Group N.V. or its affiliates. All Rights Reserved.