
## **Endoscope Tracker, Rotatable**

REF 8000-060-004

Instructions for Use

Rx Only

English	EN	1
Français	FR	18
日本語	JA	34

# 1. About this document

This document is the comprehensive source of information for the safe and effective use of the product. Read this document carefully. Familiarization with the user documentation for the components of the guidance system before use is important. Pay special attention to safety information. Keep this document accessible to users.

The following conventions are used in this document:

- The signal word **WARNING** indicates a hazardous situation that, if not avoided, could result in death or serious injury.
- The signal word **CAUTION** indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
- The symbol  indicates a note that supplements or clarifies information.

## 1.1. Other applicable documents

In addition to this document, the following documents are provided:

- User manual for the software application
- Safety information for the software application
- Guide for Cleaning, Disinfection, and Steam-Based Sterilization (TD6000005750)

## 1.2. Definition of symbols

The following tables define the symbols used in this document, on the product, and on the product label.

### EN ISO 7010: Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol/number      Name: definition

**General warning sign:** To signify a general warning.

### EN ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labeling, and information to be supplied - Part 1 General requirements

Symbol/number      Name: definition

**Manufacturer:** Indicates the medical device manufacturer as defined in the European Union harmonization legislation.

**Date of manufacture:** Indicates the date when the medical device was manufactured.

Symbol/number      Name: definition

 5.1.5	<b>Batch code:</b> Indicates the manufacturer's batch code so that the batch or lot can be identified.
 5.1.6	<b>Catalog number:</b> Indicates the manufacturer's catalog number so that the medical device can be identified.
 5.2.7	<b>Non-Sterile:</b> Indicates a medical device that has not been subjected to a sterilization process.
 5.3.2	<b>Keep away from sunlight:</b> Indicates a medical device that needs protection from light sources.
 5.3.4	<b>Keep dry:</b> Indicates a medical device that needs to be protected from moisture.
 5.3.7	<b>Temperature limit:</b> Indicates the temperature limits to which the medical device can be safely exposed.

Symbol/number      Name: definition

5.3.8

**Humidity limitation:** Indicates the range of humidity to which the medical device can be safely exposed.

5.3.9

**Atmospheric pressure limitation:** Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

5.4.3

**Consult instructions for use:** Indicates the need to consult the instructions for use.

### **ASTM F2503: Standard practice for marking medical devices and other items for safety in the MR environment**

Symbol      Name: definition

**MR Unsafe:** Indicates a product that is known to pose hazards in all MRI environments.

## Product-specific symbols

Symbol                      Name: definition

**QTY**

**Quantity:** Indicates the number of medical devices in the packaging.

**GTIN**

Global Trade Item Number

## 21 CFR 801.109

Symbol                      Name: definition

**Rx Only**

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

## Regulatory marks and logos

Symbol                      Name: definition

CSA certified for Canada and USA

## 2. Safety information

### WARNING

- Before each use, check the product for loose components, damage, malfunction, and bent or deformed parts. Do not use if these conditions exist. Failure to comply may lead to patient or medical staff injury. If service is required, contact your Stryker sales representative.
- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, it must be reprocessed according to a validated procedure. Refer to *Section "Reprocessing."*
- In case of known pre-existing conditions of variants of Creutzfeldt-Jakob-Disease (CJD) of the patient, ensure that the product is quarantined and not being used further. Follow the recommendation of the national authorities (e. g. WHO, RKI or CDC) concerning Creutzfeldt-Jakob-Disease.
- Avoid serious strains, for instance heavy impacts. After any excessive strain, the affected instrument must be checked for defects and recalibrated (if possible) or replaced.
- Do not expose the product to a high magnetic field such as from a magnetic resonance imaging (MRI) device.

### CAUTION

- Unauthorized modifications of the product are forbidden for safety reasons.
- The product may only be used for its intended purpose.

The user or patient must report any serious product-related incident to both the manufacturer and the national competent authority where the user or patient is established.

---

## 2.1. User group

Healthcare professionals (surgeon, resident, nurse, professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product. To request additional in-service instruction, contact Stryker.

## 2.2. Indications for use and contraindications

### **United States and Canada\***

Refer to the user manual supplied with the Cranial Guidance Software for system indications and contraindications.

### **Rest of World\***

The Endoscope Tracker Rotatable is indicated for use as an accessory to the Cranial Guidance Software. It is intended to be used to connect and navigate conventional rigid endoscopes used during cranial neurosurgical procedures.

The system is indicated as an aid for locating anatomical structures in either open or percutaneous neurosurgical procedures. The system is indicated for any medical condition in which the use of computer-assisted surgery may be appropriate and where reference to a rigid anatomical structure can be identified.

**Contraindications:** None known.

\*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

### 3. Product description

The Endoscope Tracker is attached to the optical unit of a rigid endoscope and allows the tip position of the endoscope to be calibrated with a calibration device and tracked during navigated surgery.


Figure 1: Endoscope Tracker, Rotatable. **A** - Shaft, **B** - Tracker, **C** - Fastening Nut.

### 3.1. For use with

#### **WARNING**

Use only the Stryker-approved products specified in this section.

---

For primary related products, consult the following table:

<b>Product</b>	<b>Catalog number</b>
Navigation Spheres (100 Pack)	8000-100-004
Cranial Guidance Software	6000-670-000
Stryker Q Guidance System	8900-100-000
Passive Trackers Insert Tray	8000-810-000

Refer to the software user manual for a complete list of compatible products on a system level.

---

## 4. Using the product

### 4.1. Preparation for use

1. Mount three navigation spheres onto the mounting pins of the tracker. The spheres produce a click sound when they are properly fastened in position.
2. Pass the shaft (A) through the ring of the tracker (B) and fix the shaft in place with the fastening nut (C). Ensure that three small knobs on the ring of the tracker match with concavities on the shaft. Ensure that a black mark on the shaft matches with an incision on the ring of the tracker.


Figure 2: Preassembly

3. Slide the tracker onto your endoscope and fix the tracker with the bayonet catch. Refer to *Figures 3 and 4*.


Figure 3


Figure 4

The Endoscope Tracker, Rotatable is suited for endoscopes with 4 mm diameter and 120 mm length.

---

4. If it is necessary to adjust the position of the tracker on the endoscope, loosen the fastening nut, bring the tracker into the desired position, and fix it again with the fastening nut. The black mark on the shaft and the incision on the tracker do not have to match anymore.

#### **CAUTION**

- The tracker must be attached to a compatible endoscope in a way that allows the endoscope to be used in an unrestricted and usual manner.
  - Make sure that the endoscope tracker is positioned in the direction of the 3D camera during intraoperative use. A direct line of sight between the camera and the tracker must not be blocked.
- 

#### **WARNING**

- Use the navigation spheres only if they have been taken directly out of intact and sterile packaging. Check the navigation spheres for damage before use.
  - Liquids on navigation spheres reduce the visibility of the navigated instruments. When possible, clean up the navigation spheres using a sterile cloth or replace the affected navigation spheres. Damaged or liquid-covered navigation spheres may compromise the navigation accuracy and lead to a serious damage of the nervous system.
-

- Make sure that the navigation spheres are properly attached to the tracker. If you suspect that a navigation sphere has moved, perform instrument tip validation again to verify the navigation accuracy. Failure to attach the navigation spheres correctly may compromise the navigation accuracy and lead to a serious damage of the nervous system.
- 

## 4.2. Instructions for application

### **WARNING**

Do not exert force when calibrating or validating.

---

## 5. Reprocessing

This section contains only product-specific information. For general reprocessing instructions, safety directives, and reprocessing equipment, refer to the Guide for Cleaning, Disinfection, and Steam-Based Sterilization (GCDS) (TD6000005750).

The **cleaning group** of this product according to the GCDS is: IV.

Remove the navigation spheres from the tracker before reprocessing. Do not reprocess the navigation spheres.

Unscrew the fastening nut and remove the shaft from the tracker before reprocessing. Refer to *Figure 2* for a picture of the tracker in disassembled state.

This device can be reprocessed in the Passive Trackers Insert Tray (REF 8000-810-000). Refer to the tray instructions for use for handling of the tray and placement of the instruments.

---

## **6. Transport, storage, and disposal**

### **6.1. Transport and storage**

Store the product within the specified environmental conditions throughout its shelf life. For more information, refer to *Section "Technical specifications."*

### **6.2. Disposal**

Products that have been in contact with material of human origin may be infectious. Dispose of with the necessary precautionary measures in accordance with local regulations. Ensure that infected products are decontaminated before recycling.

## 7. Technical specifications

### 7.1. Environmental conditions

Environmental limitations	Operation	Storage and transportation
Temperature	between 10°C and 30°C	between -10°C and +50°C
Relative humidity	30 % to 75 %	keep dry
Atmospheric air pressure	between 80 kPa and 106 kPa	-

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This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

U.S. Patents: [stryker.com/patents](https://stryker.com/patents)

Stryker Corporation or its divisions or other affiliated entities own, use, or have applied for the following trademarks or service marks: Stryker, Leibinger, Q Guidance. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2024 Stryker

Stryker Leibinger GmbH & Co. KG  
Bötzingen Straße 41  
79111 Freiburg (Germany)  
t: +49 761 4512 0 (Germany)  
t: +1 269 323 7700 (USA)