Instrument Tracker 1	REF 8000-060-014
Instrument Tracker 2	REF 8000-060-015
Instrument Tracker 3	REF 8000-060-016
Instrument Tracker L	REF 8000-060-017
Instrument Tracker XL	REF 8000-060-018
Elbow, 45°	REF 8000-060-019

Rx Only

Instructions for Use

English	EN	1
Français	FR	23
日本語	JA	45

1. About this document

This document is the comprehensive source of information for the safe and effective use of the product. Read this document carefully. Familiarization with the user documentation for the components of the guidance system before use is important. Pay special attention to safety information. Keep this document accessible to users.

The following conventions are used in this document:

- The signal word **WARNING** indicates a hazardous situation that, if not avoided, could result in death or serious injury.
- The signal word **CAUTION** indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
- The symbol  indicates a note that supplements or clarifies information.

1.1. Other applicable documents

In addition to this document, the following documents are provided:

- User manual for the software application
- Safety information for the software application
- Guide for Cleaning, Disinfection, and Steam-Based Sterilization (TD6000005750)

1.2. Definition of symbols

The following tables define the symbols used in this document, on the product, and on the product label.

EN ISO 7010: Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol/number Name: definition

General warning sign: To signify a general warning.

EN ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labeling, and information to be supplied - Part 1 General requirements

Symbol/number Name: definition

Manufacturer: Indicates the medical device manufacturer as defined in the European Union harmonization legislation.

Symbol/number Name: definition

5.1.3

Date of manufacture: Indicates the date when the medical device was manufactured.

LOT

5.1.5

Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.

REF

5.1.6

Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.

5.2.7

Non-Sterile: Indicates a medical device that has not been subjected to a sterilization process.

5.3.2

Keep away from sunlight: Indicates a medical device that needs protection from light sources.

5.3.4

Keep dry: Indicates a medical device that needs to be protected from moisture.

Symbol/number Name: definition

5.3.7

Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

5.3.8

Humidity limitation: Indicates the range of humidity to which the medical device can be safely exposed.

5.3.9

Atmospheric pressure limitation: Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

5.4.3

Consult instructions for use: Indicates the need to consult the instructions for use.


ASTM F2503: Standard practice for marking medical devices and other items for safety in the MR environment

Symbol Name: definition

MR Unsafe: Indicates a product that is known to pose hazards in all MRI environments.

Product-specific symbols

Symbol	Name: definition
	Quantity: Indicates the number of medical devices in the packaging.
GTIN	Global Trade Item Number

21 CFR 801.109

Symbol	Name: definition
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Regulatory marks and logos

Symbol	Name: definition
	CSA certified for Canada and USA

2. Safety information

WARNING

- Before each use, check the product for loose components, damage, malfunction, and bent or deformed parts. Do not use if these conditions exist. Failure to comply may lead to patient or medical staff injury. If service is required, contact your Stryker sales representative.
- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, it must be reprocessed according to a validated procedure. Refer to *Section "Reprocessing."*
- In case of known pre-existing conditions of variants of Creutzfeldt-Jakob-Disease (CJD) of the patient, ensure that the product is quarantined and not being used further. Follow the recommendation of the national authorities (e. g. WHO, RKI or CDC) concerning Creutzfeldt-Jakob-Disease.
- Avoid serious strains, for instance heavy impacts. After excessive strain, the affected instrument must be checked for defects and recalibrated (if possible) or replaced.
- Do not expose the product to a high magnetic field such as from a magnetic resonance imaging (MRI) device.

CAUTION

- Unauthorized modifications of the product are forbidden for safety reasons.
- The product may only be used for its intended purpose.

The user or patient must report any serious product-related incident to both the manufacturer and the national competent authority where the user or patient is established.

2.1. User group

Healthcare professionals (surgeon, resident, nurse, professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product. To request additional in-service instruction, contact Stryker.

2.2. Indications for use and contraindications

United States and Canada*

Refer to the user manual supplied with the Cranial Guidance Software for system indications and contraindications.

Rest of World*

– Instrument Trackers

The Instrument Tracker 1, Instrument Tracker 2, Instrument Tracker 3, Instrument Tracker L, and Instrument Tracker XL are indicated for use as accessories to the Cranial Guidance Software. They are intended to be attached to the Stryker Instrument Clamps to track surgical instruments during cranial neurosurgical procedures.

– Elbow, 45°

The Elbow, 45° is indicated for use as an accessory to the Stryker Instrument Clamps. When used with the Cranial Guidance Software, the Elbow, 45° adjusts the Instrument Tracker's alignment relative to the surgical instrument allowing for an unrestricted view of the tracker.

The system is indicated as an aid for locating anatomical structures in either open or percutaneous neurosurgical procedures. The system is indicated for any medical condition in which the use of computer-assisted surgery may be appropriate and where reference to a rigid anatomical structure can be identified.

Contraindications: None known.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3. Product description


Instrument trackers

There is a choice of five different instrument trackers, which can be freely combined with the instrument clamps. Thus it is possible for five different surgical instruments to be navigated during one surgical procedure.


Elbow, 45°

The elbow serves to adjust the chosen instrument tracker's alignment relative to the surgical instrument, so that the camera has an unrestricted view of the instrument tracker at all times. The elbow can be used to widen the gap between the instrument and the instrument tracker, thereby changing their ergonomic features.


Product name	Catalog number	Figure
Instrument Tracker 1	8000-060-014	1a
Instrument Tracker 2	8000-060-015	1b
Instrument Tracker 3	8000-060-016	1c
Instrument Tracker L	8000-060-017	1d
Instrument Tracker XL	8000-060-018	1e
Elbow, 45°	8000-060-019	1f


1a: Instrument Tracker 1


1b: Instrument Tracker 2

1c: Instrument Tracker 3

1d: Instrument Tracker L

1e: Instrument Tracker XL

1f: Elbow, 45°

Figure 1: Instrument trackers and Elbow, 45°

3.1. For use with

WARNING

Use only the Stryker-approved products specified in this section.

For primary related products, consult the following table:

Product	Catalog number
Instrument Clamp, 2-6 mm	8000-060-011
Instrument Clamp, 6-10 mm	8000-060-012
Instrument Clamp, 10-16 mm	8000-060-013
Instrument Clamp, Forceps	8000-060-010
Navigation Spheres (100 Pack)	8000-100-004
Cranial Guidance Software	6000-670-000
Stryker Q Guidance System	8900-100-000
Passive Trackers Insert Tray	8000-810-000

Refer to the software user manual for a complete list of compatible products on a system level.

4. Using the product

4.1. Suitability of surgical instruments for navigation

Numerous surgical instruments can be navigated if attached to the instrument trackers by means of the instrument clamps. Not all surgical instruments are suited for navigation. Soft or flexible instruments as well as instruments with vibrating parts are unsuitable for navigation.

Only those points on a surgical instrument that remain rigidly stable in relation to the instrument tracker are suitable for calibration. If the position of the calibrated point relative to the instrument tracker changes, the navigation becomes inaccurate. The closer the instrument tracker is placed to the navigated point, the more precise the navigation information will be. Refer to *Figure 2*.

CAUTION

- Make sure that the attachment of the instrument clamp and the instrument tracker does not restrict the function and usability of the surgical instrument to be navigated. Otherwise this instrument must not be navigated via the instrument tracker.
- If it is not possible to fix the clamp on a surgical instrument without damaging it, the instrument is not suitable for navigation. Do not use any tools for tightening the fastening nut.
- Verify that the instrument clamp is rigidly attached to the surgical instrument and cannot slip during the procedure.
- Make sure that the instrument tracker is positioned in the direction of the 3D camera during the intraoperative use. The direct line of sight between the camera and the instrument tracker must not be blocked. Check if this is the case before calibrating

the instrument and, if necessary, adjust the instrument clamp position on the instrument and the positioning of the instrument tracker.


Figure 2: Instrument with low precision vs. instrument with high precision

4.2. Preparation for use

4.2.1. Attaching the instrument clamp to a surgical instrument

1. Unscrew the fastening nut of the instrument clamp by rotating it counterclockwise.
2. Press the clamp together and position it on the surgical instrument.
3. Rotate the fastening nut clockwise until the clamp is rigidly attached to the instrument.


Figure 3: Clamp attachment

4.2.2. Instrument tracker preassembly

1. Mount the navigation spheres onto the mounting pins of the instrument tracker. The spheres produce a click sound when they are properly fastened in position.
2. Rotate the knurled screw through the tracker until the knurled screw is loosely connected with the tracker.


Figure 4: Tracker preassembly

WARNING

- Use the navigation spheres only if they have been taken directly out of intact and sterile packaging. Check the navigation spheres for damage before use.
- Liquids on the navigation spheres reduce the visibility of navigated instruments. When possible, clean up the navigation spheres using a sterile cloth or replace the affected navigation spheres. Damaged or liquid-covered navigation spheres may compromise navigation accuracy and lead to serious damage of the nervous system.
- Make sure that the navigation spheres are properly attached to the tracker. If you suspect that a navigation sphere has moved, perform the instrument tip validation

to verify the navigation accuracy. Failure to attach the navigation spheres correctly may compromise the navigation accuracy and lead to a serious damage of the nervous system.

4.2.3. Tracker to clamp attachment

1. Mount the preassembled instrument tracker onto the instrument clamp.
2. Rotate the knurled screw clockwise until the tracker is rigidly attached to the clamp.


Figure 5: Tracker to clamp attachment

4.2.4. Tracker to clamp attachment by means of the Elbow, 45°

1. Mount the preassembled instrument tracker onto the elbow and rotate the knurled screw clockwise until the tracker is rigidly attached to the elbow. Refer to *Figure 6 A*.
2. Rotate the knurled screw through the elbow until the screw is loosely connected with the elbow. Refer to *Figure 6 B*.
3. Mount the elbow onto the instrument clamp.
4. Rotate the knurled screw clockwise until the elbow is rigidly attached to the clamp. Refer to *Figure 7*.


Figure 6: Tracker to elbow attachment


Figure 7: Elbow to clamp attachment

4.3. Instructions for application

WARNING

Do not exert force when calibrating or validating.

CAUTION

Make sure that the instrument tracker is in the direct line of sight with the camera.

5. Reprocessing

This section contains only product-specific information. For general reprocessing instructions, safety directives, and reprocessing equipment, refer to the Guide for Cleaning, Disinfection, and Steam-Based Sterilization (GCDS) (TD6000005750).

The **cleaning group** of this product according to the GCDS is: IV.

Remove the navigation spheres from the instrument trackers before reprocessing. Do not reprocess the navigation spheres.

Remove the knurled screws out of the trackers and out of the elbow before reprocessing.

i The instrument trackers and the Elbow, 45° can be reprocessed in the Passive Trackers Insert Tray (REF 8000-810-000). Refer to the tray instructions for use for handling of the tray and placement of the devices.

6. Transport, storage, and disposal

6.1. Transport and storage

Store the product within the specified environmental conditions throughout its shelf life. For more information, refer to *Section "Technical specifications."*

6.2. Disposal

Products that have been in contact with material of human origin may be infectious. Dispose of with the necessary precautionary measures in accordance with local regulations. Ensure that infected products are decontaminated before recycling.

7. Technical specifications

7.1. Environmental conditions

Environmental limitations	Operation	Storage and transportation
Temperature	between 10°C and 30°C	between -10°C and +50°C
Relative humidity	30 % to 75 %	keep dry
Atmospheric air pressure	between 80 kPa and 106 kPa	-

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This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

U.S. Patents: stryker.com/patents

Stryker Corporation or its divisions or other affiliated entities own, use, or have applied for the following trademarks or service marks: Stryker, Leibinger, Q Guidance. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2024 Stryker

Stryker Leibinger GmbH & Co. KG
Bötzingen Straße 41
79111 Freiburg (Germany)
t: +49 761 4512 0 (Germany)
t: +1 269 323 7700 (USA)