Patient Tracker Electromagnetic -10 uses

REF 8000-040-002

Instructions For Use

Rx Only

EN Instructions for use	1
DA Brugsanvisning	19
DE Gebrauchsanweisung	37
EL Οδηγίες χρήσης	55
ES Instrucciones de uso	73
FI Käyttöohjeet	91
FR Notice d'utilisation	109
IT Istruzioni per l'uso	127
JA 使用説明書	145
KO 사용 설명서	163
NL Gebruiksaanwijzing	181
NO Bruksanvisning	199
PL Instrukcja użycia	217
PT Instruções de utilização	235
RO Instrucțiuni de utilizare	253
RU Инструкции по применению	271
SV Bruksanvisning	289
TR Kullanma Talimatı	307
ZH 使用说明	325

1. How to Use this Document

1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the respective system user manual before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient and medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

 Supplements or clarifies information.

1.2. Symbol Definitions

EN ISO 7010: Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol	Definition
 W001	General warning sign: To signify a general warning
 M002	Refer to instruction manual/booklet: To signify that the instruction manual/ booklet must be read.

ISO 15223-1: Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements

Symbol	Definition
 5.1.1	Manufacturer: Indicates the medical device manufacturer.

Symbol	Definition
 5.1.3	Date of manufacture: Indicates the date when the medical device was manufactured.
 5.1.5	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
 5.1.6	Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.
 5.1.7	Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.
 5.2.7	Non-sterile: Indicates a medical device that has not been subjected to a sterilization process.
 5.3.2	Keep away from sunlight: Indicates a medical device that needs protection from light sources.

Symbol	Definition
 5.3.4	Keep dry: Indicates a medical device that needs to be protected from moisture.
 5.3.7	Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.
 5.3.8	Humidity limitation: Indicates the range of humidity to which the medical device can be safely exposed.
 5.3.9	Atmospheric pressure limitation: Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
 5.4.3	Consult instructions for use or consult electronic instructions for use: Indicates the need for the user to consult the instructions for use.
 5.7.7	Medical device: Indicates the item is a medical device.

Product-specific symbols

Symbol	Definition
	Quantity: Indicates the number of products in the packaging.
	Note: Used to supplement or clarify information.
GTIN	Global Trade Item Number.

IEC 60417: Graphical symbols for use on equipment

Symbol	Definition
	Type BF applied part (5333): To identify a type BF applied part complying with IEC 60601-1.

ASTM F2503: Standard practice for marking medical devices and other items for safety in the MR environment

MR Unsafe: Indicates an item that poses unacceptable risks to the patient, medical staff, or other persons within the MR environment.

21 CFR 801.109

Symbol	Definition
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol	Definition
	Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

Regulatory marks and logos

Symbol	Definition
	Indicates conformity to Annex I of Medical Device Regulation (EU) 2017/745 without involvement of the notified body

2. Safety Information

WARNING

- The product is delivered non-sterile. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments.(e.g. MRI)

- Do not use the device in the environment of a magnetic resonance imaging (MRI) scanner
 - The product is a high-precision instrument. Avoid subjecting it to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
 - The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
-

i The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Stryker ENT Navigation System

3.1.1. United States / Rest of World*

Intended Use

The Patient Tracker Electromagnetic – 10 Uses is an accessory to the Electromagnetic Navigation Unit and is intended for locating the patient.

Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures:
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

Contraindications

None known.

3.2. Cranial Guidance System

3.2.1. United States and Canada*

Refer to the user manual supplied with the Cranial Guidance Software (REF 6000-670-000) for system indications and contraindications.

3.2.2. Rest of World*

The Patient Tracker Electromagnetic – 10 Uses is indicated for use as an accessory with adult and pediatric patients when using the Cranial Guidance Software using

electromagnetic navigation. It is intended to be used as a patient tracker with electromagnetic navigation.

The system is indicated as an aid for locating anatomical structures in either open or percutaneous neurosurgical procedures where the use of stereotactic surgery may be appropriate. The system can be used for intraoperative guidance where reference to a rigid anatomical structure, such as the skull, can be identified.

Contraindications

None known.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3.3. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

4. Product Overview

Figure 1: Patient Tracker Electromagnetic-10 uses

4.1. Instructions for Application

The patient tracker provides a reference for the patient during surgery. For this purpose, the tracker is securely fixed onto the patient with appropriate accessories as near as possible to the operation area.

WARNING

- Prior to and during the surgical procedure, regularly verify that the patient tracker is firmly attached to the patient. The patient tracker must not move with respect to the patient. If the patient tracker moves, the system displays inaccurate positional information of the instruments. If you suspect that the patient tracker has moved, check landmarks to verify navigation accuracy.
-

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

5. For Use With

 **WARNING**

- Use only Stryker-approved products, unless otherwise specified.
-

For information related to compatible software applications, refer to the user manual supplied with the respective software application. For information related to product-specific compatibility, refer to the table below.

Description	REF
Patient Tracker Tabs	8000-100-001

6. Disassembly

1. Carefully remove the patient tracker from the patient skin.
2. Remove the Patient Tracker Tab from the patient tracker, and dispose of the Patient Tracker Tab.

7. Product Lifetime

The instrument is validated for ten reprocessing cycles. Thus, the instrument contains a usage counter that disables the instrument after running the tenth usage cycle in surgery. The usage counter is updated as soon as the instrument is connected during surgery and shown in the software.

To maintain accurate use count, connect the instrument to the electromagnetic port at least once during surgery. This also applies to replacement instruments if they have been taken from sterile packaging.

8. Reprocessing

If you use the dedicated EM Instruments Insert Tray REF 8000-820-000* for cleaning, disinfection and sterilization, refer to the Guide for Cleaning, Disinfection and Steam-Based Sterilization (TD6000005750) for reprocessing instructions, safety directives, and reprocessing equipment. The Patient Tracker Electromagnetic– 10 uses is designed for machine cleaning and disinfection like instruments without electronics. The cleaning group of the Patient Tracker Electromagnetic– 10 uses is: V.

If you use sterilization pouches for sterilization, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706) for reprocessing instructions, safety directives, and reprocessing equipment.

9. Maintenance

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition.

10. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

- After using the electromagnetic system, make sure that instruments that cannot be reprocessed again are disposed of.
- In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
- The plugs used in the electronic products contain the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH). Handling Instructions: No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
- The silicone glue used in the products contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH).

11. Technical Specifications

Environmental Limitations	Operation	Storage and Transportation
Temperature	 A diagram showing a thermometer with a scale. The minimum temperature is marked as 10°C and the maximum as 30°C.	 A diagram showing a thermometer with a scale. The minimum temperature is marked as -10°C and the maximum as 50°C.
Relative humidity	 A circular icon containing a percentage sign and wavy lines representing water. The range is indicated as 30% to 75%.	 An icon of an umbrella with rain falling on it, representing weather conditions.
Atmospheric air pressure	 A circular icon with arrows pointing inward and outward, representing air pressure. The range is indicated as 80 kPa to 106 kPa.	-

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11. Tekniske specifikationer

Miljømæssige begrænsninger	Drift	Opbevaring og transport
Temperatur	 A line graph showing a temperature range from 10 °C to 30 °C. The line starts at 10 °C, rises to a peak at 30 °C, and then drops. A small black circle is positioned at the peak of the line.	 A line graph showing a temperature range from -10 °C to 50 °C. The line starts at -10 °C, rises to a peak at 50 °C, and then drops. A small black circle is positioned at the peak of the line.
Relativ luftfugtighed	 A circular icon containing a percentage sign and wavy lines representing water. Two lines point to the values 30 % and 75 %.	 An icon of an umbrella with rain falling on it.
Atmosfærisk lufttryk	 A circular icon containing a central dot and two horizontal arrows pointing outwards. Two lines point to the values 80 kPa and 106 kPa.	-

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Einschränkungen in Bezug auf die Umgebungsbedingungen	Betrieb	Lagerung und Transport
Temperatur	 A line graph showing a temperature range from 10 °C to 30 °C. The line starts at 10 °C, rises to a peak at 30 °C, and then drops. A small black circle is positioned at the peak of the graph.	 A line graph showing a temperature range from -10 °C to 50 °C. The line starts at -10 °C, rises to a peak at 50 °C, and then drops. A small black circle is positioned at the peak of the graph.
Relative Luftfeuchte	 A circular icon containing a percentage sign and wavy lines representing water. Two lines point to the values 30 % and 75 %.	 An icon of an open umbrella with several raindrops falling from it.
Luftdruck	 A circular icon containing a central dot and two horizontal arrows pointing outwards. Two lines point to the values 80 kPa and 106 kPa.	-

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Limitaciones ambientales	Funcionamiento	Almacenamiento y transporte
Temperatura	 A diagram showing a temperature range for operation. It consists of a horizontal line with a vertical tick mark at the left end labeled "10 °C" and another vertical tick mark at the right end labeled "30 °C". A diagonal line rises from the right end of the horizontal line to a point above the "30 °C" label, and then a vertical line drops down to a small black circle.	
Humedad relativa	 A diagram showing a relative humidity range for operation. It features a circle containing a percentage sign (%) and a wavy line representing water. Two lines extend from the circle: one to the top right labeled "75 %" and one to the bottom left labeled "30 %".	
Presión atmosférica	 A diagram showing an atmospheric pressure range for operation. It features a circle with a central dot and two horizontal arrows pointing towards each other. Two lines extend from the circle: one to the top right labeled "106 kPa" and one to the bottom left labeled "80 kPa".	
Almacenamiento y transporte	 A diagram showing a temperature range for storage and transport. It consists of a horizontal line with a vertical tick mark at the left end labeled "-10 °C" and another vertical tick mark at the right end labeled "50 °C". A diagonal line rises from the right end of the horizontal line to a point above the "50 °C" label, and then a vertical line drops down to a small black circle.	
	 An icon of an umbrella with several raindrops falling from it, representing weather conditions during storage and transport.	

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Ограничения по условиям окружающей среды	Эксплуатация	Хранение и транспортировка
Температура воздуха	 A line graph showing a temperature range from 10 °C to 30 °C. The line starts at 10 °C, rises to a peak at 30 °C, and then falls. A thermometer icon is placed at the peak.	 A line graph showing a temperature range from -10 °C to 50 °C. The line starts at -10 °C, rises to a peak at 50 °C, and then falls. A thermometer icon is placed at the peak.
Относительная влажность	 A circular icon containing a percentage sign and wavy lines representing humidity. Two lines point to the values 30 % and 75 %.	 An icon of an open umbrella with rain falling from it.
Атмосферное давление	 A circular icon containing a pressure gauge symbol with two arrows pointing towards a central dot. Two lines point to the values 80 kPa and 106 kPa.	—

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11. Tekniska specifikationer

Gränser avseende miljöfaktorer	Drift	Förvaring och transport
Temperatur	 A line graph showing a temperature range from 10 °C to 30 °C. The line starts at 10 °C, rises to a peak at 30 °C, and then falls. A thermometer icon is positioned at the peak.	 A line graph showing a temperature range from -10 °C to 50 °C. The line starts at -10 °C, rises to a peak at 50 °C, and then falls. A thermometer icon is positioned at the peak.
Relativ luftfuktighet	 A circular icon containing a percentage sign and wavy lines representing humidity. Two lines point to the values 30 % and 75 %.	 An icon of an umbrella with rain falling on it.
Atmosfärstryck	 A circular icon containing a central dot and two horizontal arrows pointing outwards. Two lines point to the values 80 kPa and 106 kPa.	-

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Leibinger, Stryker.
All other trademarks are trademarks of their respective owners or holders.
Copyright © 2024 Stryker
Stryker Leibinger GmbH & Co. KG