Precision Pointer Electromagnetic

REF 8000-050-001

Pointer Electromagnetic

REF 8000-050-002

Instructions for Use

Rx Only

EN Instructions for use	1
DA Brugsanvisning	21
DE Gebrauchsanweisung	41
EL Οδηγίες χρήσης	61
ES Instrucciones de uso	81
FI Käyttöohjeet	101
FR Notice d'utilisation	121
IT Istruzioni per l'uso	141
JA 使用説明書	161
KO 사용 설명서	181
NL Gebruiksaanwijzing	201
NO Bruksanvisning	221
PL Instrukcja użycia	241
PT Instruções de utilização	261
RO Instrucțiuni de utilizare	281
RU Инструкции по применению	301
SV Bruksanvisning	321
TR Kullanım talimatları	341
ZH 使用说明	361

1. How to Use this Document

1.1. About this Document

This manual is the most comprehensive source of information for the safe, effective, and compliant use and/or maintenance of the product. Read and understand this manual as well as the instructions for use supplied with the software application before using the product or any component compatible with the product. When combined with other medical devices, the instructions for use of these devices are to be considered as well. This manual is a permanent part of the product. Keep this manual for future reference.

The following conventions are used in this document:

The signal word **WARNING** highlights a safety-related issue. Comply with this information to prevent patient or medical staff injury.

The signal word **CAUTION** highlights a product reliability issue. Comply with this information to prevent product damage.

Supplements or clarifies information.

1.2. Symbol Definitions

EN ISO 7010: Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol	Definition
	General warning sign (W001): To signify a general warning.
	Refer to instruction manual/booklet (M002): To signify that the user instruction manual/booklet must be read.

ISO 15223-1: Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements

Symbol	Definition
	Manufacturer (5.1.1): Indicates the medical device manufacturer
	Date of manufacture (5.1.3): Indicates the date when the medical device was manufactured.
	Batch code (5.1.5): Indicates the manufacturer's batch code so that the batch or lot can be identified.

Symbol	Definition
	Catalog number (5.1.6): Indicates the manufacturer's catalog number so that the medical device can be identified.
	Serial number (5.1.7): Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Non-Sterile (5.2.7): Indicates a medical device that has not been subjected to a sterilization process.
	Keep away from sunlight (5.3.2): Indicates a medical device that needs protection from light sources.
	Keep dry (5.3.4): Indicates a medical device that needs to be protected from moisture.
	Temperature limit (5.3.7): Indicates the temperature limits to which the medical device can be safely exposed.
	Humidity limitation (5.3.8): Indicates the range of humidity to which the medical device can be safely exposed.

Symbol	Definition
	Atmospheric pressure limitation (5.3.9): Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
	Consult instructions for use or consult electronic instructions for use (5.4.3): Indicates the need for the user to consult the instructions for use.
	Contains hazardous substances (5.4.10): Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties.
	Medical device (5.7.7): Indicates the item is a medical device.

IEC 60417 Graphical symbols for use on equipment

Symbol

Definition

Type BF applied part (5333): To identify a type BF applied part complying with IEC 60601-1.

Product-Specific Symbols

Symbol

Definition

Quantity: Indicates the number of products in the packaging

Note: Used to supplement or clarify information

GTIN

Global Trade Item Number

21 CFR 801.109

Symbol

Definition

Rx Only

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Directive 2012/19/EU on waste electrical and electronic equipment (WEEE)

Symbol

Definition

Indicates that the product must be collected separately and must not be disposed of as unsorted municipal waste.

Regulatory marks and logos

Symbol

Definition

CE 0197

Indicates conformity to Annex I of Medical Device Regulation (EU) 2017/745 or Medical Device Directive 93/42 EEC with involvement of the notified body.

2. Safety Information

WARNING

- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, the product must be processed according to a validated procedure.
- Unauthorized modifications of the product are forbidden for safety reasons.
- Due to its composition, the product is not to be used together with other magnetically sensitive medical products, devices or instruments (e.g. MRI).
- Before using a medical product, assure yourself of its functional safety and proper condition via a visual inspection. Do not use the product if any defects, such as breaks, cracks, deformations or wear, are detected. Particularly important parts such as tips, notches and all moveable components should be checked with extra care.
- The product is a high precision instrument. Avoid subjecting the product to serious strains, such as heavy impacts. The product cannot be used if there are visible defects. After a heavy impact, the product must be checked for defects.
- The product has a sharp tip. Please pay special attention to sharp areas in order to avoid injuries to medical personnel or patients.

- The healthcare provider performing any procedure is responsible for determining the appropriateness of using the product and for the specific technique for each patient. Stryker, as a manufacturer, does not recommend a specific surgical procedure.
 - Do not use the system during the discharge of a defibrillator. All navigated instruments must be removed from the patient and the operating area beforehand.
-

The user and/or patient should report any serious product-related incident to both the manufacturer and the national competent authority where the user and/or patient is established.

3. Product Information

WARNING

- The product may only be used for its intended purpose and in accordance with this document and all current versions of relevant system and software application documentation.
 - This document is part of the product and must be accessible to personnel at all times. It must be provided to subsequent owners or users.
-

3.1. Stryker ENT Navigation System

3.1.1. United States / Rest of World*

Precision Pointer Electromagnetic / Pointer Electromagnetic Intended Use

The Precision Pointer Electromagnetic / Pointer Electromagnetic is an accessory to the Electromagnetic Navigation Unit and is intended for pointing at anatomical structures.

Indications for Use

The Stryker ENT Navigation System is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure in the field of ENT surgery, such as the paranasal sinuses, mastoid anatomy, can be identified relative to a CT- or MR-based model of the anatomy.

Example procedures include, but are not limited to the following ENT procedures:

- Transsphenoidal access procedures
- Intranasal procedures
- Sinus procedures, such as maxillary antrostomies, ethmoidectomies, sphenoidotomies/ sphenoid explorations, turbinate resections, and frontal sinusotomies
- ENT-related anterior skull base procedures

Contraindications

None known.

3.2. Cranial Guidance System

3.2.1. United States and Canada*

Precision Pointer Electromagnetic

Refer to the user manual supplied with the Cranial Guidance Software (REF 6000-670-000) for system indications and contraindications.

3.2.2. Rest of World*

The Precision Pointer Electromagnetic is indicated for use as an accessory with adult and pediatric patients when using the Cranial Guidance Software using electromagnetic navigation. It is intended to be used for precisely locating anatomical structures.

The system is indicated as an aid for locating anatomical structures in either open or percutaneous neurosurgical procedures where the use of stereotactic surgery may be appropriate. The system can be used for intraoperative guidance where reference to a rigid anatomical structure, such as the skull, can be identified.

Contraindications: None known.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3.3. User Group

Healthcare professionals (surgeon/resident, nurse/professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product.

To request an additional in-service instruction, contact Stryker.

4. Product Overview

Precision Pointer Electromagnetic (REF 8000-050-001)

The Precision Pointer Electromagnetic is a pointing device with a tracked tip. The steel tube (Figure 1, B) may be reshaped once per procedure. The distal 20 mm up to the black ring marking of the instrument (Figure 1, A) must remain unbent. For detailed information about permissible bending refer to Table 1.


Figure 1: Precision Pointer Electromagnetic

- A Black ring marking
- B Flexible steel tube
- C Non-flexible tube

Pointer Electromagnetic (REF 8000-050-002)

The Pointer Electromagnetic is a pointing device with a tracked tip and with a long design for use. The steel tube (Figure 2, B) may be reshaped once per procedure. The distal 20 mm up to the black ring marking of the instrument (Figure 2, A) and the proximal 20 mm must remain unbent. For detailed information about permissible bending refer to Table 1.


Figure 2: Pointer Electromagnetic

- A Black ring marking
- B Steel tube

Table 1: Product specifications

	Precision Pointer Electromagnetic	Pointer Electromagnetic
Tube Length	120 mm	180 mm
Diameter	1.5 mm (from tip – 80 mm) 2.5 mm (80 mm – 120 mm)	2.7 mm
Flexibility	<ul style="list-style-type: none">• Flexible section length: 60 mm (Figure 1, B); min. bending radius: 15 mm• Not flexible parts:<ul style="list-style-type: none">- first 20 mm of the tip (Figure 1, A)- proximal tube (Figure 1, C)	<ul style="list-style-type: none">• Flexible section length: 140 mm; min. bending radius: 30 mm• Not flexible parts: tip and end of the tube (each 20 mm)

4.1. Instructions for Application

WARNING

Do not bend the distal 20 mm of the Precision Pointer Electromagnetic and the Pointer Electromagnetic.

For instructions on how to use the product for a surgical procedure, refer to the instructions for use supplied with the respective software application.

5. For Use With

WARNING

Use only Stryker-approved products, unless otherwise specified.

For information related to compatible software applications, refer to the instructions for use supplied with the software application.

6. Product Lifetime

The instrument is validated for ten reprocessing cycles. Thus, the instrument contains a usage counter that disables the instrument after running the tenth usage cycle in surgery. The usage counter is updated as soon as the instrument is connected during surgery and shown in the software.

To maintain accurate use count, connect the instrument to the electromagnetic port at least once during surgery. This also applies to replacement instruments if they have been taken from sterile packaging.

7. Reprocessing

7.1. Precision Pointer Electromagnetic REF 8000-050-001

If you use the dedicated EM Instruments Insert Tray REF 8000-820-000* for cleaning, disinfection and sterilization, refer to the Guide for Cleaning, Disinfection and Steam-Based Sterilization (TD6000005750) for reprocessing instructions, safety directives, and reprocessing equipment. The Precision Pointer Electromagnetic is designed for machine cleaning and disinfection like instruments without electronics. The cleaning group of the Precision Pointer Electromagnetic is: V.

If you use sterilization pouches for sterilization, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706) for reprocessing instructions, safety directives, and reprocessing equipment.

7.2. Pointer Electromagnetic REF 8000-050-002

For reprocessing instructions, safety directives, and reprocessing equipment, refer to the Stryker ENT Navigation System Reprocessing Instructions (TD8000010706).

8. Maintenance

Maintenance and repair may only be conducted by the manufacturer or authorized partners. The product components may only be sent back to the manufacturer in a cleaned, disinfected, and sterilized condition. Sharp or pointed components need to be sent back in a protected state.

9. Disposal

Risk of injury or infection can be avoided by safe disposal of the product components. Sharp and pointed components need to be collected and locked in a tight and break-proof container. They must be stored in such a way that they are protected from unauthorized use. Contaminated products are to be supplied to a hazardous waste site and handled in a way that contamination of third parties is excluded.

In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.


10. Technical Specifications

Environmental Limitations

Operation

Storage and
Transportation

Temperature

Relative humidity

Atmospheric air pressure

—

Materials

- This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0
Current scientific evidence supports the position that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- The plugs used in the electronic products contain the following substance: Lead, CAS No. 7439-92-1 (1907/2006 REACH).
No specific precautions are required for handling items manufactured from alloys containing lead in the supplied condition.
- The silicone glue used in the products contains the following substances: Decamethylcyclopentasiloxane, CAS No. 541-02-6 and Dodecamethylcyclohexasiloxane, CAS No. 540-97-6 (1907/2006 REACH)

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Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks:
All other trademarks are trademarks of their respective owners or holders.
Copyright © 2024 Stryker
Stryker Leibinger GmbH & Co. KG