
Universal Base Skull

REF 8000-040-017

Instructions for Use

Rx Only

English	EN	1
Français	FR	19
日本語	JA	37

1. About this document

This document is the comprehensive source of information for the safe and effective use of the product. Read this document carefully. Familiarization with the user documentation for the components of the guidance system before use is important. Pay special attention to safety information. Keep this document accessible to users.

The following conventions are used in this document:

- The signal word **WARNING** indicates a hazardous situation that, if not avoided, could result in death or serious injury.
- The signal word **CAUTION** indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
- The symbol  indicates a note that supplements or clarifies information.

1.1. Other applicable documents

In addition to this document, the following documents are provided:

- User manual for the software application
- Safety information for the software application
- Guide for Cleaning, Disinfection, and Steam-Based Sterilization (TD6000005750)

1.2. Definition of symbols

The following tables define the symbols used in this document, on the product, and on the product label.

EN ISO 7010: Graphical symbols – Safety colors and safety signs – Registered safety signs

Symbol/number Name: definition

 W001	General warning sign: To signify a general warning.
---	--

EN ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labeling, and information to be supplied - Part 1 General requirements

Symbol/number Name: definition

 5.1.1	Manufacturer: Indicates the medical device manufacturer as defined in the European Union harmonization legislation.
 5.1.3	Date of manufacture: Indicates the date when the medical device was manufactured.

Symbol/number Name: definition

 5.1.5	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
 5.1.6	Catalog number: Indicates the manufacturer's catalog number so that the medical device can be identified.
 5.2.7	Non-sterile: Indicates a medical device that has not been subjected to a sterilization process.
 5.3.2	Keep away from sunlight: Indicates a medical device that needs protection from light sources.
 5.3.4	Keep dry: Indicates a medical device that needs to be protected from moisture.
 5.3.7	Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

Symbol/number Name: definition

5.3.8

Humidity limitation: Indicates the range of humidity to which the medical device can be safely exposed.

5.3.9

Atmospheric pressure limitation: Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

5.4.3

Consult instructions for use: Indicates the need to consult the instructions for use.

ASTM F2503: Standard practice for marking medical devices and other items for safety in the MR environment

Symbol Name: definition

MR Unsafe: Indicates a product that is known to pose hazards in all MRI environments.

Product-specific symbols

Symbol Name: definition

QTY

Quantity: Indicates the number of medical devices in the packaging.

GTIN

Global Trade Item Number

21 CFR 801.109

Symbol Name: definition

Rx Only

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Regulatory marks and logos

Symbol Name: definition

CSA certified for Canada and USA

2. Safety information

WARNING

- Before each use, check the product for loose components, damage, malfunction, and bent or deformed parts. Do not use if these conditions exist. Failure to comply may lead to patient or medical staff injury. If service is required, contact your Stryker sales representative.
- The product is delivered in a non-sterile condition. Before first use, as well as before each following use, it must be reprocessed according to a validated procedure. Refer to *Section "Reprocessing."*
- In case of known pre-existing conditions of variants of Creutzfeldt-Jakob-Disease (CJD) of the patient, ensure that the product is quarantined and not being used further. Follow the recommendation of the national authorities (e. g. WHO, RKI or CDC) concerning Creutzfeldt-Jakob-Disease.
- Avoid serious strains, for instance heavy impacts. After an excessive strain, the affected instrument must be checked for defects and recalibrated (if possible) or replaced.
- Do not expose the product to a high magnetic field such as from a magnetic resonance imaging (MRI) device.

CAUTION

- The product has three sharp tips at the base. Pay special attention to these areas in order to avoid injury to yourself or others.
- Unauthorized modifications of the product are forbidden for safety reasons.
- The product may only be used for its intended purpose.

The user or patient must report any serious product-related incident to both the manufacturer and the national competent authority where the user or patient is established.

2.1. User group

Healthcare professionals (surgeon, resident, nurse, professional caregiver) educated in computer-assisted surgery and thoroughly familiar with the instructions for use and with the operation of this product. To request additional in-service instruction, contact Stryker.

2.2. Indications for use and contraindications

United States and Canada*

Refer to the user manual supplied with the Cranial Guidance Software for system indications and contraindications.

Rest of World*

The Universal Base Skull is indicated for use as an accessory to the Cranial Guidance Software. It is intended for fixation of the Patient Tracker Advanced to the patient. The system is indicated as an aid for locating anatomical structures in either open or percutaneous neurosurgical procedures. The system is indicated for any medical condition in which the use of computer-assisted surgery may be appropriate and where reference to a rigid anatomical structure can be identified.

Contraindications: None known.

*Note: Product may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Contact your Stryker representative for product availability.

3. Product description

The Universal Base Skull is an invasive patient tracker adapter. It allows the Patient Tracker Advanced (8000-040-018) to be fixed directly on the patient's head. This allows the system to track the position of the patient during the surgery even if the patient needs to be moved.


Figure 1: Universal Base Skull

3.1. For use with

WARNING

Use only the Stryker-approved products specified in this section.

For primary related products, consult the following table:

Product	Catalog number
Patient Tracker Advanced	8000-040-018
Cranial Guidance Software	6000-670-000
Stryker Q Guidance System	8900-100-000
Passive Trackers Insert Tray	8000-810-000

Refer to the software user manual for a complete list of compatible products on a system level.

4. Using the product

4.1. Preparation for use

You need:

- Titanium bone screw 1.5 mm x 6 mm (or 1.5 mm x 8 mm)
- Surgical screwdriver incl. screwdriver blade compatible with the bone screw

WARNING

Use only a titanium single use bone screw.

4.1.1. Base preassembly (Figure 2)

1. Insert the tube (a) into the tracker foot (b) from below and fix it with the adjusting screw (c) (with one or two rotations).
2. Attach the connecting piece (d) and secure it with the knurled nut (e).

Figure 2:


a - Tube

b - Tracker Foot

c - Adjusting Screw

d - Connecting Piece

e - Knurled Nut

4.1.2. Fixing the base on the patient's head (Figure 3)

1. Select a suitable location on the patient's head.
2. Carry out an incision of approximately 7 mm to 8 mm.
3. Insert the tube into the incision and fix it with a bone screw at the neurocranium.

4. Loosen the bone screw by rotating it half a turn counterclockwise so that the base does not tilt when the adjusting screw is tightened.
5. Lightly press the pins of the base through the skin until they lie on the neurocranium.
6. Tighten the adjusting screw until the base is rigidly attached to the neurocranium.

WARNING

Do not pull the adjusting screw too tight, as otherwise the bone screw could be pulled out.


Figure 3

4.1.3. Attaching the Patient Tracker Advanced to the base

1. Preassemble the Patient Tracker Advanced and attach it to the Universal Base Skull. Refer to the Patient Tracker Advanced Instructions for Use.

WARNING

- Fix the navigation spheres to the Patient Tracker Advanced before attaching it to the Universal Base Skull.
 - Hold the Universal Base Skull firmly while attaching the Patient Tracker Advanced to the base so that the base neither twists nor tilts on the patient.
-

4.2. Instructions for application

WARNING

Do not exert force when using the patient tracker for instrument validation.

4.3. Disassembly

1. Remove the Patient Tracker Advanced from the Universal Base Skull by loosening the knurled screw.
2. Unscrew the bone screw with a screwdriver and remove the base from the patient's head.
3. If the bone screw does not come unscrewed, completely loosen the adjusting screw and pull the tracker foot out of the tube. Now, remove the tube by pulling the bone screw through the opening of the tube. Then remove the bone screw.

WARNING

When completely loosening the adjusting screw, ensure that it does not fall uncontrollably into the operating area.

5. Reprocessing

This section contains only product-specific information. For general reprocessing instructions, safety directives, and reprocessing equipment, refer to the Guide for Cleaning, Disinfection, and Steam-Based Sterilization (GCDS) (TD6000005750).

The **cleaning group** of this product according to the GCDS is: IV.

Disassemble the device entirely before reprocessing. Refer to *Figure 2* for a picture of the device in disassembled state.

 This device can be reprocessed in the Passive Trackers Insert Tray (REF 8000-810-000). Refer to the tray instructions for use for handling of the tray and placement of the instruments.

6. Transport, storage, and disposal

6.1. Transport and storage

Store the product within the specified environmental conditions throughout its shelf life. For more information, refer to *Section "Technical specifications."*

6.2. Disposal

Products that have been in contact with material of human origin may be infectious. Dispose of with the necessary precautionary measures in accordance with local regulations. Ensure that infected products are decontaminated before recycling.

7. Technical specifications

7.1. Environmental conditions

Environmental limitations	Operation	Storage and transportation
Temperature	between 10°C and 30°C	between -10°C and +50°C
Relative humidity	30 % to 75 %	keep dry
Atmospheric air pressure	between 80 kPa and 106 kPa	-

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This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

U.S. Patents: [stryker.com/patents](https://www.stryker.com/patents)

Stryker Corporation or its divisions or other affiliated entities own, use, or have applied for the following trademarks or service marks: Stryker, Leibinger, Q Guidance. All other trademarks are trademarks of their respective owners or holders.

Copyright © 2024 Stryker

Stryker Leibinger GmbH & Co. KG
Bötzingen Straße 41
79111 Freiburg (Germany)
t: +49 761 4512 0 (Germany)
t: +1 269 323 7700 (USA)