AIS TAYORI

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- Goyal M, Menon BK, van Zwam WH, Dippel DW, Mitchell PJ, Demchuk AM, Davalos A, Majoie CB, van der Lugt A, de Miquel MA, et al. Endovascular thrombectomy after large-vessel ischaemic stroke: a meta-analysis of individual patient data from five randomised trials. *Lancet* (London, England). 2016;387:1723-1731. doi: 10.1016/s0140-6736(16)00163-x
- Duloquin G, Graber M, Garnier L, Crespy V, Comby PO, Baptiste L, Mohr S, Delpont B, Guénat J, Blanc-Labarre C, et al. Incidence of Acute Ischemic Stroke With Visible Arterial Occlusion: A Population-Based Study (Dijon Stroke Registry). *Stroke*. 2020;51:2122-2130. doi: 10.1161/strokeaha.120.029949
- Ospel JM, Menon BK, Demchuk AM, Almekhlafi MA, Kashani N, Mayank A, Fainardi E, Rubiera M, Khaw A, Shankar JJ, et al. Clinical Course of Acute Ischemic Stroke Due to Medium Vessel Occlusion With and Without Intravenous Alteplase Treatment. *Stroke*. 2020;51:3232-3240. doi: 10.1161/strokeaha.120.030227
- Bilgin C, Hardy N, Hutchison K, Pederson JM, Mebane A, Olaniran P, Kobeissi H, Kallmes KM, Fiorella D, Kallmes DF, et al. First-line thrombectomy strategy for distal and medium vessel occlusions: a systematic review. *J Neurointerv Surg*. 2023;15:539-546. doi: 10.1136/jnis-2022-019344
- Nogueira RG, Doheim MF, Al-Bayati AR, Lee JS, Haussen DC, Mohammaden M, Lang M, Starr M, Rocha M, da Câmara CP, et al. Distal Medium Vessel Occlusion Strokes: Understanding the Present and Paving the Way for a Better Future. *J Stroke*. 2024;26:190-202. doi: 10.5853/jos.2023.02649
- Kuwahara S, Uchida K, Sakai N, Yamagami H, Imamura H, Takeuchi M, Shirakawa M, Sakakibara F, Haraguchi K, Kimura N, et al. Technical and clinical outcomes of thrombectomy in patients with acute medium vessel occlusion and large vessel occlusion: sub-analyses of Japan Trevo registry. *J Neurol Sci*. 2024;459:122956. doi: 10.1016/j.jns.2024.122956
- Haussen DC, Rebello LC, Nogueira RG. Optimizing Clot Retrieval in Acute Stroke: The Push and Fluff Technique for Closed-Cell Stent Retrievers. *Stroke*. 2015;46:2838-2842. doi: 10.1161/strokeaha.115.010044
- Pinheiro AC, Nogueira RG, Grandfield RM, Lin SP, Majidhoo AQ, Aghaebrahim AN, Abraham MG, Mazaris P, Singer JA, Al-Bayati AR, et al. Push and Fluff technique for optimization of clot integration with stent-retriever: An in vitro model. *Interventional neuroradiology : journal of peritherapeutic neuroradiology, surgical procedures and related neurosciences*. 2023;15910199231175348. doi: 10.1177/15910199231175348
- Kyselyova AA, Fiehler J, Leischner H, Flottmann F, Buhk JH, Frölich AM. Vessel diameter and catheter-to-vessel ratio affect the success rate of clot aspiration. *J Neurointerv Surg*. 2021;13:605-608. doi: 10.1136/neurintsurg-2020-016459
- Halama D, Merkel H, Werdehausen R, Gaber K, Schob S, Quäsching U, Ziganshyna S, Hoffmann KT, Lindner D, Richter C. Reference Values of Cerebral Artery Diameters of the Anterior Circulation by Digital Subtraction Angiography: A Retrospective Study. *Diagnostics* (Basel). 2022;12. doi: 10.3390/diagnostics12102471

All Photographs taken by Kindai University Hospital.

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

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Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: AXS Catalyst, AXS Universal, AXS Vecta, Stryker, Synchro SELECT, Trevo, Trevo NXT, Trevo Trak, TriGlide. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

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