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Bibliografia

1. Jackson J.D., Pagnano M.W. (2012) Cement fixation for total knee arthroplasty. In: *The Knee Joint*. Springer, Paris. pagg. 759-764.
2. Dalury, DF. Cementless total knee arthroplasty. *Bone Joint J* 2016;98-B:867-73.
3. Webb, JCJ and Spencer, RF. The role of polymethylmethacrylate bone cement in modern orthopaedic surgery. *BJJ*. VOL. 89-B, n. 7, luglio 2007
4. Barrack, RL. Early Failure of Modern Cemented Stems. *The Journal Of Arthroplasty*. Vol. 15, n. 8, ago 2000.
5. Lettich, T. et al., Primary Total Hip Arthroplasty with an Uncemented Femoral Component: Two- to Seven-Year Results. *The Journal of Arthroplasty* Vol. 22, n. 7, suppl. 3 2007
6. Fleischman, AN et al. Reduced Incidence of Intraoperative Femur Fracture With a Second-Generation Tapered Wedge Stem. *The Journal of Arthroplasty* 32 (2017), pagg. 3457-3461
7. National Joint Registry for England, Wales and Northern Ireland. 15th annual report, 2018. <http://www.njrreports.org.uk/Portals/0/PDFdownloads/NJR%2011th%20Annual%20Report%202014.pdf> (ultimo accesso: 11 novembre 2018).
8. Australian Orthopaedic Association National Joint Replacement Registry: Annual Report, 2018 (ultimo accesso: 11 novembre 2018).
9. American Joint Replacement Registry: Annual Report, 2017
10. Memtsoudis SG, Della Valle AG, Besculides MC, Gaber L, Laskin R. (2009) Trends in Demographics, Comorbidity Profiles, InHospital Complications and Mortality Associated with Primary Knee Arthroplasty. *Journal of Arthroplasty*, Jun; 24(4):51827
11. Schreurs, BW, Hannik, G. Total Joint Arthroplasty in Younger Patients: Heading for Trouble? *The Lancet* (2017), vol. 389, n. 10077, pagg. 1374-1375.
12. Nam, D., et al. Perioperative and Early Postoperative Comparison of a Modern Cemented and Cementless Total Knee Arthroplasty of the Same Design. *The Journal of Arthroplasty* 32 (2017), pagg. 2151-2155
13. Nilsson K G, Henricson A, Norgren B, Dalen T. Uncemented HA-coated implant is the optimum fixation for TKA in the young patient. *Clin Orthop Relat Res* 2006; 448: 129-39.
14. Meneghini RM, Hanssen A. Cementless Fixation in Total Knee Arthroplasty: Past, Present, and Future. *J Knee Surg*. 2008; 21:307-314.
15. Bayley JC, Scott RD, Ewald FC, Holmes GB Jr. Failure of the metal-backed patellar component after total knee replacement. *J Bone Joint Surg Am*. 1988 Jun; 70(5):668-74.
16. Bircher, A. et al. Contact Dermatitis. Allergic complications from orthopaedic joint implants: the role of delayed hypersensitivity to benzoyl peroxide in bone cement. *Contact Dermatitis*. 2011 August, 66, pagg. 20-26
17. Granchi D., Cenni E., Tigani D., Trisolino G., Baldini N., Giunti A. (2008). Sensitivity to implant materials in patients with total knee arthroplasties. *Science Direct*, 29, 1494-1500. doi:10.1016/j.biomaterials.2007.11.038
18. Vega, F Aseptic loosening of a total knee prosthesis caused by delayed hypersensitivity to bone cement. *Annals of Allergy, Asthma & Immunology*. Volume 117, Issue 1, July 2016, pagg. 89-91
19. Keene, R. et al. Occupational Hazards to the Pregnant Orthopaedic Surgeon. *J Bone Joint Surg Am*. 2011;93:e141(1-5)
20. Pacheco, Karin A. Allergy to Surgical Implants. *The Journal of Allergy and Clinical Immunology: In Practice*, Volume 3, Issue 5, pagg. 683-695
21. Nam, D. et al. Cemented Versus Cementless Total Knee Arthroplasty of the Same Modern Design. *J Bone Joint Surg Am*. 2019; 101:1185-92
22. Sharpe, K., Robinson, K., Cohen, R., Barnett TM., Rastogi, A., Masini, M. Does Implant Fixation Affect Early Return to Function Following Primary Total Knee Arthroplasty. *ISTA 2018*. Londra, Regno Unito.
23. Sinicrope BJ, Feher AW, Bhimani SJ, Smith LS, Harwin SF, Yakkanti MR, Malkani AL, Increased Aseptic Failures in Cemented versus Cementless Total Knee Arthroplasty in Morbidly Obese Patients. Minimal 5 Year Follow-up, *The Journal of Arthroplasty* (2018), doi: <https://doi.org/10.1016/j.arth.2018.10.016>.
24. Miller, AJ, et al., Results of Cemented vs Cementless Primary Total Knee Arthroplasty Using the Same Implant Design. *J Arthroplasty*. 2018 Apr;33(4): 1089-1093. doi: 10.1016/j.arth.2017.11.048. Epub 2017 Dec 2.
25. Cohen et al. "Early Clinical Outcomes of a New Cementless Total Knee Arthroplasty Design. *Orthopedics*. 2018.

Bibliografia

26. Buzhardt, P. et al. "Clinical and radiographic results of a highly porous titanium cementless tibial baseplate in TKA." *Bone Joint J 99-B.SUPP 3*(2017): 56. Web. 13 mar 2017."
27. Harwin, SF. et al. Outcomes of Newer Generation Cementless Total Knee Arthroplasty: Beaded Periapatite-Coated vs Highly Porous Titanium-Coated Implant. *J Arthroplasty*. 2017 Jul;32(7):2156-2160. doi: 10.1016/j.arth.2017.01.044. Pubbl. elettronica 3 febbraio 2013.
28. Bhowmik-Stoker et al. "Clinical Performance of a Novel 3D Printed Cementless Titanium Tibial Baseplate, 2-4-year Follow-up" *World Arthroplasty Congress*, 4/19-21, 2018, Rome.
29. King, J.D., et al. The Hidden Cost of Commercial Antibiotic-Loaded Bone Cement: A Systematic Review of Clinical Results and Cost Implications Following Total Knee Arthroplasty. *The Journal of Arthroplasty* 33 92018, pagg. 3789-3792.
30. Swedish Knee Arthroplasty Register. Annual Report 2017. <http://www.myknee.se/en/> (ultimo accesso: 11 novembre 2018).
31. Sultan, A.A., et al. Routine use of commercial antibiotic-loaded bone cement in primary total joint arthroplasty: a critical analysis of the current evidence. *Ann Transl Med* 2019;7(4):73
32. McCalden et al "Comparison of Outcomes and Survivorship Between Patients of Different Age Groups Following TKA." *The Journal of Arthroplasty*, Churchill Livingstone, 24 July 2013, www.sciencedirect.com/science/article/pii/S0883540313004774.
33. Meehan, JP et al. Younger Age Is Associated with a Higher Risk of Early Periprosthetic Joint Infection and Aseptic Mechanical Failure After Total Knee Arthroplasty. *J Bone Joint Surg Am*. 2014; 96:529-35.
34. Mont MA, Gwam C, Newman JM, Chughtai M, Khlopas A, Ramkumar PN, Harwin SF. Outcomes of a newer-generation cementless total knee arthroplasty design in patients less than 50 years of age. *Ann Transl Med* 2017; 5 (Suppl 3):S24. doi: 10.21037/atm.2017.08.20
35. Newman, J. et al. Cementless Total Knee Arthroplasty in Patients Older Than 75 Years. *The Journal of Knee Surgery*. Vol. 30, n. 9/2017, pagg. 930-934.
36. Martin, JR, et al. Morbid Obesity and Total Knee Arthroplasty: A Growing Problem. *Journal of the American Academy of Orthopaedic Surgeons*: March 2017 - Volume 25 - Issue 3 - pagg. 188-194. doi: 10.5435/JAAOS-D-15-00684
37. Derman, PB, Fabricant, PD, David, G. The Role of Overweight and Obesity in Relation to the More Rapid Growth of Total Knee Arthroplasty Volume Compared with Total Hip Arthroplasty Volume. *J Bone Joint Surg Am*. 2014 Jun 4; 96 (11): 922-928. Pubbl. elettronica 4 giu 2014.
38. Abdel, MA. et al. Increased aseptic tibial failures in patients with a BMI ≥ 35 and well-aligned total knee arthroplasties. *J Arthroplasty*, 30 (2015), pagg. 2181-2184
39. Sharpe, K., Robinson, K., Cohen, R., Barnett TM., Cohen, R., Masini, M. Prospective Assessment of the Impact of Obesity on Early Postoperative Outcomes of Cementless Total Knee Arthroplasty AAHKS 2017
40. Sinicropo BJ, Feher AW, Bhimani SJ, Smith LS, Harwin SF, Yakkanti MR, Malkani AL, Increased Aseptic Failures in Cemented versus Cementless Total Knee Arthroplasty in Morbidly Obese Patients. Minimal 5 Year Follow-up, *The Journal of Arthroplasty* (2018), doi: <https://doi.org/10.1016/j.arth.2018.10.016>.
41. Bhowmik-Stoker et al. "Clinical Performance of a Novel 3D Printed Cementless Titanium Tibial Baseplate, 2-4-year Follow-up" *World Arthroplasty Congress*, 4/19-21, 2018, Rome.
42. Firestein GS. Evolving concepts of rheumatoid arthritis. *Nature*. 2003;423(6937):356-361. doi:10.1038/nature01661
43. Helmick CG, et al. Estimates of the Prevalence of Arthritis and Other Rheumatic Conditions in the United States. Part I. *Arthritis Rheumatism*. 2008;58(1):15-25
44. Patel, N. et al. Outcomes of Cementless Total Knee Arthroplasty in Patients With Rheumatoid Arthritis. *Othopedics*, vol. 41, n. 2, 2018, pagg.103-106
45. Ryd et al. 1995 Ryd L, Albrektsson B E, Carlsson L, Dansgard F, Herberts P, Lindstrand A, Regner L, Toksvig-Larsen S. Roentgen stereophotogrammetric analysis as a predictor of mechanical loosening of knee prostheses. *J Bone Joint Surg Br* 1995; 77(3): 377-83.

46. Nilsson K et al. Evaluation of Micromotion in Cemented vs Uncemented Knee Arthroplasty in Osteoarthritis and Rheumatoid Arthritis. *Journal of Arthroplasty*. Vol 6. N. 3. Settembre, 1991. 265-278.
47. Bhimji et al. The effect of fixation design on micromotion of cementless tibial baseplates. ORS 2012. Poster 1977.
48. Stryker Test Report RD-13-107.
49. Stryker Test Protocol 92911.
50. Nevelos, J. et al. Design, Migration and Early Clinical Results of the First Mass Produced 3D Printed Cementless Total Knee Implants. Scientific Exhibit presented at AAOS. Las Vegas, NV, 12-16 Mar 2019.
51. Pijls, BG, Plevier J, Nelissen, R. RSA migration of total knee replacements: A systematic review and meta-analysis. *Acta Orthopaedica* 2018; 89 (3): 320-328
52. Dunbar, MJ et al. Stable migration of peri-apatite-coated uncemented tibial components in a multicentre study. *Bone Joint J* 2017;99-B:1596-1602.
53. Robinson JC, Pozen A, Tseng S, Bozic KJ. Variability in costs associated with total hip and knee replacement implants. *J Bone Joint Surg Am*. 2012 Sep 19;94(18):1693-8.
54. Macario, A. What does one minute of operating room time cost? *Journal of Clinical Anesthesia* (2010) 22, 233–236
55. Volpin A, Khan O, Haddad FS. Theater Cost Is £16/Minute So What Are You Doing Just Standing There? *J Arthroplasty*. 2016 Jan;31(1):22-6. doi: 10.1016/j.arth.2015.08.008. Pubbl. elettronica 14 ago 2015.
56. Chen, A. et al. Cementless versus Cemented TKA: Shorter Operative Time and Similar Complications. Eastern Orthopedic Association 2016. New Orleans, LA.
57. Lawrie, C. et al. The Cost of Implanting a Cemented vs. Cementless Total Knee Arthroplasty. Poster n. 72. AAHKS 2018. Dallas, TX
58. Hansen, EN. et al. Routine use of antibiotic laden bone cement for primary total knee arthroplasty: impact on infecting microbial patterns and resistance profiles. *J Arthroplasty*. 2014 Jun;29(6):1123-7. doi: 10.1016/j.arth.2013.12.004. Epub 2013 Dec 10.
59. Inabathula, A. et al. Extended Oral Antibiotic Prophylaxis in High-Risk Patients Substantially Reduces Primary Total Hip and Knee Arthroplasty 90-Day Infection Rate. *J Bone Joint Surg Am*. 2018 Dec 19;100(24):2103-2109.
60. Kallala, R., Anderson, P., Morris, S., Haddad, FS. The cost analysis of cemented versus cementless total hip replacement operations on the NHS. *Bone Joint J* 2013;95-B:874-6.

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Copyright © 2020 Stryker