Trevo NXT™ ProVue Retriever

Instructions for Use	2
Instrucciones de uso	9
Mode d'emploi	17
Gebrauchsanweisung	25
Istruzioni per l'uso	33
Gebruiksaanwijzing	41
Instruções de Utilização	49

Trevo NXT™ ProVue Retriever

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Stryker Neurovascular representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or reesterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or reesterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Retriever consists of a flexible, tapered core wire with a shaped section at the distal end. Platinum markers at the distal end allow fluoroscopic visualization. In addition, the shaped section is also radiopaque. Retriever dimensions are indicated on product label. The Retriever delivery wire has a hydrophilic coating on the distal 101cm to reduce friction during use. The Retriever has a shaft marker to indicate proximity of Retriever tip relative to Microcatheter tip. A torque device is provided with the Retriever to facilitate manipulation and retrieval. The Retriever comes preloaded in an insertion tool to introduce the Retriever into a Microcatheter.

User Information

The Trevo NXT ProVue Retriever should only be used by physicians trained in interventional endovascular procedures.

Contents

One (1) Trevo NXT ProVue Retriever

One (1) Torque Device

One (1) Insertion Tool

Compatibility

Retriever Size	Trevo Pro14 Microcatheter Inner Diameter 0.017 inches (0.43 mm)	Trevo Trak 21 Microcatheter Inner Diameter 0.021 inches (0.53 mm)	Trevo™ Pro18 Microcatheter Inner Diameter 0.021 inches (0.53 mm)	Excelsior® XT-27® Microcatheter (REF XT275081) Inner diameter 0.027 inches (0.69 mm)	Recommended Minimum Vessel ID (mm)
Trevo NXT 3x32	✓	✓	✓	✓	1.4
Trevo NXT 4x28		✓	✓	✓	
Trevo NXT 4x41		✓	✓	✓	
Trevo NXT 6x37		✓	✓	✓	2.5

Compatibility of the Retriever with other microcatheters has not been established. Performance of the Retriever device may be impacted if a different microcatheter is used.

Balloon Guide Catheters and Aspiration Catheters (commercially available aspiration catheters with minimum inner diameter 0.058 inches (1.47 mm)) are recommended for use during thrombus removal procedures.

INTENDED USE/INDICATIONS FOR USE

- The Trevo Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion, and smaller core infarcts who have first received intravenous tissue plasminogen activator (IV t-PA). Endovascular therapy with the device should start within 6 hours of symptom onset.

- The Trevo Retriever is intended to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.
- The Trevo Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion of the internal carotid artery (ICA) or middle cerebral artery (MCA)-M1 segments with smaller core infarcts (0-50cc for age <80 years, 0-20cc for age ≥80 years). Endovascular therapy with the device should start within 6-24 hours of time last seen well in patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy.

CONTRAINDICATIONS

None known.

WARNINGS

SPECIFIC WARNINGS FOR INDICATION 1

- The safety and effectiveness of the Trevo Retrievers in reducing disability has not been established in patients with large core infarcts (i.e., ASPECTS ≤ 7). There may be increased risks, such as intracerebral hemorrhage, in these patients.
- The safety and effectiveness of the Trevo Retrievers in reducing disability has not been established or evaluated in patients with occlusions in the posterior circulation (e.g., basilar or vertebral arteries) or for more distal occlusions in the anterior circulation.

SPECIFIC WARNINGS FOR INDICATION 2

- To reduce risk of vessel damage, take care to appropriately size Retriever to vessel diameter at intended site of deployment.

SPECIFIC WARNINGS FOR INDICATION 3

- The safety and effectiveness of the Trevo Retrievers in reducing disability has not been established in patients with large core infarcts (i.e., ASPECTS ≤ 7). There may be increased risks, such as intracerebral hemorrhage, in these patients.
- The safety and effectiveness of the Trevo Retrievers in reducing disability has not been established or evaluated in patients with occlusions in the posterior circulation (e.g., basilar or vertebral arteries) or for more distal occlusions in the anterior circulation.
- Users should validate their imaging software analysis techniques to ensure robust and consistent results for assessing core infarct size.

GENERAL WARNINGS APPLIED TO ALL INDICATIONS

- Administration of IV t-PA should be within the currently approved window
- To reduce risk of vessel damage, adhere to the following recommendations:
 - Do not perform more than six (6) retrieval attempts in same vessel using Retriever devices.
 - Maintain Retriever position in vessel when removing or exchanging Microcatheter.
- To reduce risk of kinking/fracture, adhere to the following recommendations:
 - Immediately after unsheathing Retriever, position Microcatheter or Aspiration Catheter tip marker over the proximal section of the Retriever. Maintain this position during manipulation and withdrawal.
 - Do not rotate or torque Retriever.
 - Use caution when passing Retriever through stented arteries.
- The Retriever is a delicate instrument and should be handled carefully. Before use and when possible during procedure, inspect device carefully for damage. Do not use a device that shows signs of damage. Damage may prevent device from functioning and may cause complications.
- Do not advance or withdraw Retriever against resistance or significant vasospasm. Moving or torquing device against resistance or significant vasospasm may result in damage to vessel or device. Assess cause of resistance using fluoroscopy and if needed resheath the device to withdraw.
- If Retriever is difficult to withdraw from the vessel, do not torque Retriever. Advance Microcatheter or Aspiration Catheter over the Retriever and remove devices as a unit. If undue resistance is met when withdrawing the Retriever into the Microcatheter, consider exchanging for a larger diameter Aspiration Catheter. Gently withdraw the Retriever and larger diameter catheter as a unit.
- Administer anti-coagulation and anti-platelet medications per standard institutional guidelines.
- Do not use open or damaged packages.
- Do not expose Retriever to organic solvents.

PRECAUTIONS

- Store in cool, dry, dark place
- Use by "Use By" date.
- Exposure to temperatures above 54°C (130°F) may damage device and accessories. Do not autoclave.
- Use Retriever in conjunction with fluoroscopic visualization and proper anti-coagulation agents.
- To prevent thrombus formation and contrast media crystal formation, maintain a constant infusion of appropriate flush solution through all catheter lumens.
- Users should take all necessary precautions to limit X-ray radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors whenever possible.

ADVERSE EVENTS

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with possible complications which may occur during or after the procedure. Possible complications include, but are not limited to, the following:

- air embolism
- hematoma or hemorrhage at puncture site
- infection
- distal embolization
- pain/headache
- vessel spasm
- thrombosis
- dissection
- perforation
- emboli
- acute occlusion
- ischemia
- intracranial hemorrhage
- false aneurysm formation
- neurological deficits including stroke
- death

Use of device requires fluoroscopy which presents potential risks to physicians and patients associated with x-ray exposure. Possible risks include, but are not limited to, the following:

- alopecia,
- burns ranging in severity from skin reddening to ulcers,
- cataracts,
- delayed neoplasia

Adverse Event Reporting

Please notify your Stryker Neurovascular representative immediately if a device malfunctions or patient complication or injury is experienced or suspected. Please make every attempt to retain any suspect device, its associated components and their packaging for return to Stryker Neurovascular.

HOW SUPPLIED

Stryker Neurovascular products are sterile and non-pyrogenic in unopened packaging that is designed to maintain sterility unless the primary product pouch has been opened or damaged.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Recommended Preparation Procedure

1. Flush Retriever package hoop with saline. Hydrate for 2 minutes minimum before removing Retriever from hoop. Keep Retriever hydrated.
2. Prepare and place Balloon Guide Catheter or Guide Catheter into internal carotid artery, common carotid artery or subclavian artery as specified in guide catheter labeling.
3. Connect rotating hemostasis valve to Microcatheter hub.
4. Using conventional catheterization techniques, place Microcatheter into target vessel using a standard neurovascular guidewire. Anatomy permitting, position Microcatheter tip distal to thrombus.

5. Remove guidewire from Microcatheter. Inject contrast media through Microcatheter to visualize distal vasculature. Flush Microcatheter.
6. Remove insertion tool and preloaded Retriever as a unit from package hoop. Do not allow Retriever to exit insertion tool tip or to retract further into insertion tool.
7. Important: Introduce insertion tool half way into the rotating hemostasis valve and use syringe or infusion line to flush insertion tool until saline exits the proximal end of the insertion tool. If insertion tool is not properly flushed, it may be difficult to advance the Retriever through the insertion tool.
8. Seat insertion tool into Microcatheter hub and close rotating hemostasis valve tightly to secure insertion tool in place.
9. Advance Retriever until half of its length has been inserted into Microcatheter. Remove insertion tool.

Recommended Retrieval Procedure

1. Continue to advance the Retriever into the Microcatheter. Once the fluoro-saver markers on the proximal Retriever wire reach the Microcatheter hub or proximal end of the Rotating Hemostasis Valve, advance the Retriever under fluoroscopic guidance.
2. Advance Retriever until the distal tip aligns with the distal tip of the Microcatheter.
3. Retract Microcatheter while applying gentle strut forward force to Retriever to deploy shaped section of Retriever within clot. Position Microcatheter tip marker over the proximal section of the Retriever.
4. If using an Aspiration Catheter, position Aspiration Catheter just proximal to the Microcatheter tip, then remove the Microcatheter.

WARNING: To reduce risk of kinking/fracture, maintain Microcatheter or Aspiration Catheter tip marker over the proximal section of the Retriever during manipulation and withdrawal of the Retriever.

5. After deploying Retriever, visualize strut expansion and allow sufficient time for clot to integrate into the Retriever (approximately 5 minutes).
6. If using a Balloon Guide Catheter, inflate balloon to occlude vessel as specified in Balloon Guide Catheter labeling.
7. If using an Aspiration Catheter, position Aspiration Catheter tip marker over the proximal section of the Retriever while applying aspiration using a 60 mL syringe or an aspiration pump.
8. Position and lock torque device onto Retriever wire at Microcatheter or Aspiration Catheter rotating hemostasis valve.
9. Continue aspirating through Aspiration Catheter while the Retriever and Aspiration Catheter are withdrawn as a unit.
10. Slowly withdraw Retriever with Microcatheter or Aspiration Catheter as a unit to Balloon Guide Catheter tip while applying aspiration to Balloon Guide Catheter using 60 mL syringe.

Note: Aspiration through the balloon guide catheter with an aspiration pump is not recommended.

11. Continue vigorous aspiration until Retriever with Microcatheter or Aspiration Catheter are withdrawn to guide catheter hub.
-
- Note:** If withdrawal into Balloon Guide Catheter or Guide Catheter is difficult, deflate Balloon Guide Catheter balloon and simultaneously withdraw all devices as a unit through sheath. Remove sheath if necessary.
-

12. Disconnect Balloon Guide Catheter, Guide Catheter rotating hemostasis valve and fully remove Retriever with Microcatheter or Aspiration Catheter as a unit from guide catheter.
13. Attach 60 mL syringe to Balloon Guide Catheter hub and aspirate.
14. Deflate Balloon Guide Catheter balloon.
15. Clean the device with saline. Inspect Retriever for damage. Do not reuse Retriever if core wire, shaped section or platinum weave appears damaged. If not damaged, the Retriever may be used for up to three (3) retrieval attempts. A retrieval attempt is one (1) advancement and complete withdrawal cycle.

CLINICAL STUDIES

Safety and Effectiveness Information

DAWN™ Study to support Indication 3

DAWN (DWI and CTP Assessment with Clinical Mismatch in the Triage of Wake-Up and Late Presenting Strokes Undergoing Neurointervention) was a prospective, randomized, multi-center, controlled trial designed to evaluate the hypothesis that Trevo™ Retriever thrombectomy plus medical management leads to superior clinical outcomes at 90 days as compared to medical management alone in appropriately selected subjects experiencing an acute ischemic stroke when treatment is initiated within 6-24 hours after last seen well. A total of 206 subjects were enrolled across 26 investigative sites: 107 subjects (51.9%) were randomized to the Treatment Arm and 99 subjects (48.1%) were randomized to the Control Arm. The mean UW-mRS values at 90 days were 5.5 in the Treatment arm vs 3.4 in the Control arm; the core infarct adjusted group difference in utility-weighted degree of disability was 2.0 (95% Bayesian Credible Interval, 1.1 to 3.0, posterior probability

of superiority >0.9999). Treatment with Trevo™ Retriever led to higher rates of functional independence (mRS 0-2) at 90 days, 48.6% vs. 13.1%; core infarct adjusted group difference 32.8% (95% Credible Interval, 21.1% to 44.1%, posterior probability of superiority >0.9999). The rate of symptomatic intracranial hemorrhage at 24 hours was higher in patients undergoing endovascular therapy. Patients and/or their legally authorized representatives should be advised of this risk during the treatment decision-making process. There is greater level of uncertainty in the results shown for subjects in the 18-24 hour time from last seen well.

Key Inclusion Criteria

Acute ischemic stroke with confirmed occlusion of intracranial ICA and/or M1, that failed or were contraindicated for IV t-PA, treatment within 6-24 hours after time last known well, baseline NIHSS ≥10, informed consent given and age 18 or over; In patients younger than 80, a stroke score equal to or greater than 10 must be associated with a core volume less than 51ml. In patients 80 or older, the stroke score must be greater than 10 and volume less than 21ml. Clinical-core mismatch which is defined as the mismatch between baseline infarct volume (or core) on CT or MRI imaging and the extent of total brain tissue at risk. Infarct volume was assessed with the use of diffusion-weighted magnetic resonance imaging (MRI) or perfusion CT and was measured with one of the commercially available automated software packages (RAPID, IschemaView).

Key Exclusion criteria

Rapid improvement in neurological status to a NIHSS <10 or evidence of vessel recanalization prior to randomization; Arterial blood pressure > 185/110 mmHg; laboratory evidence of electrolyte imbalance (i.e. sodium <130 mmol/L, potassium <3 mEq/L or >6 mEq/L); laboratory evidence of renal failure (i.e. serum creatinine >3.0 mg/dL (264 μmol/L); laboratory evidence of coagulation abnormalities (i.e., platelet count <50,000/uL, APTT>3 times normal or INR>3.0, if given factor Xa inhibitor 24-48 hours ago must have normal PTT); laboratory evidence of bleeding (i.e. hemoglobin <7mmol/L).

The primary endpoint was 90-day disability assessed by modified Rankin scale (mRS). The secondary objective is to provide evidence that endovascular reperfusion with Trevo Retriever is associated with a significant reduction in median infarct size compared to the control group at 24 (-6/+24) hours post randomization. Secondary endpoints include proportion of subjects with good functional outcome at 90 days (mRS of 0-2); proportion of subjects with NIHSS drop of ≥10 points from baseline or NIHSS score 0 or 1; difference between treatment and control arm in all cause mortality; difference in median final infarct size at 24 hours; difference in revascularization rates at 24 hours; and rate of reperfusion to TIC1 >2b.

The safety outcomes include incidence of stroke-related mortality at 90 days; incidence of SICH, neurological deterioration, procedure-related and device-related serious adverse events.

90 Day Modified Rankin Scale ITT Population

Modified Rankin Scale (mRS)	Treatment Arm N=107	Control Arm N=99
0 - No symptoms / UW = 10	9.3% (10/107)	4.0% (4/99)
1 - No significant disability / UW = 9.1	22.4% (24/107)	5.1% (5/99)
2 - Slight disability / UW = 7.6	16.8% (18/107)	4.0% (4/99)
3 - Moderate disability / UW = 6.5	13.1% (14/107)	16.2% (16/99)
4 - Moderately severe disability / UW = 3.3	13.1% (14/107)	34.3% (34/99)
5 - Severe disability / UW = 0	6.5% (7/107)	18.2% (18/99)
6 - Dead / UW = 0	18.7% (20/107)	18.2% (18/99)
90 Day mRS [0-2]	48.6% (52/107)	13.1% (13/99)
90 Day Weighted mRS*		
Mean ± SD (N)	5.5 ± 3.8 (107)	3.4 ± 3.1 (99)
Median (Q1, Q3)	6.5 (0.0, 9.1)	3.3 (0.0, 6.5)
Range (min, max)	(0.0, 10.0)	(0.0, 10.0)
[95% Conf. Interval]†	[4.8, 6.3]	[2.7, 4.0]

†By normal approximation
**mRS = (0, 1, 2, 3, 4, 5, 6) was assigned a corresponding numerical value = (10, 9.1, 7.6, 6.5, 3.3, 0.0), which represents its clinical utility*

Co-Primary outcomes ITT population

Outcome	Intervention (N=107)	Control (N=99)	Mean Absolute Difference (95% CI)	Posterior mean benefit, core-adjusted (95% credible interval) *	Probability of superiority
Utility-Weighted mRS-mean (SD)	5.5 (3.8)	3.4 (3.1)	2.1 (1.2, 3.1)	2.0 (1.1, 3.0)	>0.9999 ¹
Functional Independence (mRS 0-2)	48.6%	13.1%	35.5% (23.9%, 47.0%)	32.8% (21.1%, 44.1%)	>0.9999

**Estimated by Bayesian general linear model adjusting for Core Infarct*

†The overall probability of a device benefit is then the average of the benefit probabilities for the imputed data sets, weighted by the imputation probabilities of those data sets. This probability is 0.999986.

The second secondary effectiveness endpoint was the NIHSS drop of ≥10 points from baseline or an NIHSS score of 0 or 1.

Early Neurological Response at Day 5-7/Discharge

Outcome	Intervention N=107	Control N=99	Mean Absolute Difference (95% CI)	Risk Ratio (95% CI)	p-value
Early Response†	47.7%	19.2%	28.5% (16.2%, 40.7%)	2.5 (1.6, 3.9)	<0.001

†Defined as NIHSS drop of ≥ 10 points from baseline or NIHSS 0 or 1 at day 5-7 or discharge (whichever was earlier).The p-value was not adjusted for multiplicity.

The secondary effectiveness outcome was the revascularization rates at 24 hours defined as the presence of partial or complete recanalization.

Revascularization rates at 24 hours

	Total N=206	Treatment Arm N=107	Control Arm N=99	Difference [95% CI]	p-value
Difference in revascularization rates at 24 (-6/+24) hours from randomization	58.7% (121/206)	76.6% (82/107)	39.4% (39/99)	37.2% [24.7%, 49.8%]	< 0.0001*

Revascularization at 24 hours was defined as the presence of partial or complete recanalization.

a. t-test. The p-value was not adjusted for multiplicity.

The primary safety endpoint was the incidence of stroke-related mortality at 90 days. The secondary safety outcomes for both Treatment and Control Arms are incidence of SICH, by ECASS definition, within 24 (-6/+24) hours post randomization (time zero), and incidence of neurological deterioration from baseline NIHSS score through Day 5-7/discharge (whichever is earlier) post randomization. Neurological deterioration was defined as ≥ 4 point increase in the NIHSS score from the baseline score.

Stroke-Related Mortality at 90 Days

	Total N=206	Treatment Arm N=107	Control Arm N=99	Difference [95% CI]†	p-value
Incidence of all stroke-related mortality at 90 days	17.0% (35/206)	15.9% (17/107)	18.2% (18/99)	-2.3% [-12.6%, 8.0%]	0.7126*

†By normal approximation

a. t-test. The p-value was not adjusted for multiplicity.

The secondary safety outcome was the Incidence of SICH, by ECASS III definition, within 24 (-6/+24) hours post randomization (time zero).

Incidence of SICH by CEC Adjudication

	Total N=206	Treatment Arm N=107	Control Arm N=99	Difference [95% CI]†	p-value
SICH at 24 (-6/+24) hrs post randomization	4.4% (9/206)	5.6% (6/107)	3.0% (3/99)	2.6% [-2.9%, 8.1%]	0.5011*

†By normal approximation

a. Fisher's exact test. The p-value was not adjusted for multiplicity.

Incidence of neurological deterioration was defined as a ≥ 4 point increase in the NIHSS score from baseline through Day 5-7/Discharge post randomization (time zero).

Incidence of Neurological Deterioration by CEC Adjudication

	Total N=206	Treatment Arm N=107	Control Arm N=99	Difference [95% CI] ¹	p-value
Incidence of neurological deterioration between baseline and Day 5 to 7/ Discharge	19.9% (41/206)	14.0% (15/107)	26.3% (26/99)	-12.2% [-23.1%, -1.4%]	0.0358 ^a


¹By normal approximation
a. t-test. The p-value was not adjusted for multiplicity.

Functional independence (90 Day mRS 0-2) Sub-grouped by Time

TLSW Intervals	Treatment	Control	Difference	p-value
6-9 hrs	50.0% (9/18)	14.3% (3/21)	35.7% [8.2%, 63.2%]	0.0346 ^a
9-12 hrs	58.1% (18/31)	25.0% (6/24)	33.1% [8.5%, 57.6%]	0.0273 ^a
6-12 hrs Total	55.1% (27/49)	20.0% (9/45)	35.1% [16.9%, 53.3%]	0.0006^a
12-15 hrs	37.5% (6/16)	9.1% (2/22)	28.4% [1.8%, 55.0%]	0.0498 ^a
15-18 hrs	50.0% (13/26)	5.3% (1/19)	44.7% [23.1%, 66.4%]	0.0025 ^a
18-21 hrs	36.4% (4/11)	0.0% (0/8)	36.4% [7.9%, 64.8%]	0.1032 ^a
21-24 hrs	40.0% (2/5)	20.0% (1/5)	20.0% [-35.4%, 75.4%]	1.0000 ^a
12-24 Total*	43.1% (25/58)	7.4% (4/54)	35.7% [21.2%, 50.2%]	< 0.0001^a

a. Fisher's exact test

* Greater level of uncertainty in the results shown for subjects in the 18-24 hour time from last seen well
The p-value was not adjusted for multiplicity.

The following table is a summary of all serious and non-serious adverse events for all Control Arm and Treatment Arm subjects by MedDRA Dictionary System Organ Classes (SOC).

MedDRA* System Organ Class	Treatment arm N=107						Control Arm N=99					
	SAEs	Subjects with SAEs	Non SAEs	Subjects with Non SAEs	Total AEs	Subjects with Total AEs	SAEs	Subjects with SAEs	Non SAEs	Subjects with Non SAEs	Total AEs	Subjects with Total AEs
Any Adverse Event (AE)	65	41 (38.3%)	476	102 (95.3%)	541	103 (96.3%)	83	47 (47.5%)	473	91 (91.9%)	556	91 (91.9%)
Blood and lymphatic system disorders	0	0	11	10 (9.3%)	11	10 (9.3%)	1	1 (1.0%)	8	8 (8.1%)	9	8 (8.1%)
Cardiac disorders	8	7 (6.5%)	37	29 (27.1%)	45	33 (30.8%)	6	5 (5.1%)	29	24 (24.2%)	35	28 (28.3%)
Endocrine disorders	0	0	1	1 (0.9%)	1	1 (0.9%)	0	0	0	0	0	0
Eye disorders	1	1 (0.9%)	0	0	1	1 (0.9%)	0	0	2	2 (2.0%)	2	2 (2.0%)
Gastrointestinal disorders	7	6 (5.6%)	42	30 (28.0%)	49	34 (31.8%)	12	11 (11.1%)	38	28 (28.3%)	50	34 (34.3%)
General disorders and administration site conditions	0	0	22	18 (16.8%)	22	18 (16.8%)	0	0	17	15 (15.2%)	17	15 (15.2%)
Hepatobiliary disorders	0	0	0	0	0	0	0	0	1	1 (1.0%)	1	1 (1.0%)
Immune system disorders	0	0	0	0	0	0	0	0	1	1 (1.0%)	1	1 (1.0%)
Infections and infestations	6	6 (5.6%)	46	34 (31.8%)	52	38 (35.5%)	7	7 (7.1%)	54	41 (41.4%)	61	43 (43.4%)
Injury, poisoning and procedural complications	2	2 (1.9%)	8	7 (6.5%)	10	9 (8.4%)	2	2 (2.0%)	8	6 (6.1%)	10	7 (7.1%)
Investigations	0	0	5	5 (4.7%)	5	5 (4.7%)	0	0	10	10 (10.1%)	10	10 (10.1%)
Metabolism and nutrition disorders	0	0	76	47 (43.9%)	76	47 (43.9%)	0	0	62	34 (34.3%)	62	34 (34.3%)
Musculoskeletal and connective tissue disorders	0	0	11	10 (9.3%)	11	10 (9.3%)	2	2 (2.0%)	22	15 (15.2%)	24	16 (16.2%)
Neoplasms benign, malignant and unspecified	0	0	1	1 (0.9%)	1	1 (0.9%)	0	0	0	0	0	0
Nervous system disorders	21	21 (19.6%)	95	63 (58.9%)	116	71 (66.4%)	34	30 (30.3%)	88	53 (53.5%)	122	66 (66.7%)
Psychiatric disorders	1	1 (0.9%)	26	24 (22.4%)	27	24 (22.4%)	2	2 (2.0%)	30	27 (27.3%)	32	28 (28.3%)
Renal and urinary disorders	0	0	18	15 (14.0%)	18	15 (14.0%)	1	1 (1.0%)	19	14 (14.1%)	20	15 (15.2%)
Reproductive system and breast disorders	0	0	1	1 (0.9%)	1	1 (0.9%)	0	0	2	2 (2.0%)	2	2 (2.0%)
Respiratory, thoracic and mediastinal disorders	14	11 (10.3%)	33	28 (26.2%)	47	33 (30.8%)	11	9 (9.1%)	42	34 (34.3%)	53	37 (37.4%)
Skin and subcutaneous tissue disorders	0	0	3	3 (2.8%)	3	3 (2.8%)	0	0	10	7 (7.1%)	10	7 (7.1%)
Surgical and medical procedures	1	1 (0.9%)	1	1 (0.9%)	2	2 (1.9%)	1	1 (1.0%)	2	2 (2.0%)	3	3 (3.0%)
Vascular disorders	4	3 (2.8%)	39	33 (30.8%)	43	34 (31.8%)	4	4 (4.0%)	28	21 (21.2%)	32	23 (23.2%)

*MedDRA v17.0 was used.

MR CLEAN study to support Indication 1

MR CLEAN was a large, prospective, randomized, open label, controlled, multicenter trial in which every endovascular hospital center in the Netherlands participated. Intra-arterial treatment plus usual care (which could include intravenous administration of IV t-PA) was compared with usual care alone (control group) in patients with acute ischemic stroke and a proximal intracranial arterial occlusion of the anterior circulation that was confirmed on vessel imaging. MR CLEAN randomized 500 patients (233 treatment ["IAT"], 267 control) at 16 medical centers. Retrievable stents (including the Trevo™ Retriever product family) were used in 190 of the 233 patients (81.5%) assigned to the IAT arm, of which 120/190 subjects received first line intra-arterial treatment with a Trevo Retriever and 104/120 received IV t-PA within 3 hours from symptoms onset compared to the entire MR CLEAN Medical Management (MM) control group, of which 224 of 249 were treated with IV t-PA in accordance with FDA labeling. Of the 120 subjects who received first line intra-arterial treatment with a Trevo Retriever, 24 subjects had concomitant carotid stenting. These subjects are excluded from the efficacy summary below resulting in a primary efficacy cohort of 96 Trevo Retriever treated subjects. Subjects with IA Lytic usage or use of other mechanical thrombectomy devices after Trevo Retriever were counted as failures in efficacy summaries below.

Key Inclusion Criteria for MR CLEAN

Acute ischemic stroke with confirmed occlusion of ICA, M1, M2, A1 or A2, treatment within 6 hours of symptom onset, and a NIHSS >2, CT or MRI scan ruling out intracranial hemorrhage; informed consent given; and age 18 or over.

Key Inclusion Criteria for Trevo Retriever Subset

Additional key inclusion criteria for Trevo Retriever subset (FDA Cohort) analysis from MR CLEAN: first line treatment with Trevo Retriever and IV t-PA treatment within 3 hours of symptoms onset.

Key Exclusion criteria for MR CLEAN

Arterial blood pressure > 185/110 mmHg; blood glucose <2.7 or > 22.2 mmol/L; intravenous treatment with thrombolytic therapy in a dose exceeding 0.9 mg/kg alteplase or 90 mg; intravenous treatment with thrombolytic therapy despite contra-indications (i.e., major surgery, gastrointestinal bleeding or urinary tract bleeding within the previous 2 weeks, or arterial puncture at a non-compressible site within the previous 7 days; cerebral infarction in the distribution of the relevant occluded artery in the previous 6 weeks; laboratory evidence of coagulation abnormalities (i.e., platelet count <40x10⁹/L, APTT>50 sec or INR>3.0).

The primary effectiveness endpoint was to demonstrate that Trevo Retriever thrombectomy plus MM leads to superior functional independence (mRS <3) at 90 days as compared to MM alone in eligible subjects experiencing an acute ischemic stroke. Based on the analysis, the IAT Trevo Retriever treatment resulted in a highly significant treatment effect when adjusted for sites, IV t-PA use, and baseline NIHSS.

Primary Effectiveness Endpoint

Variable	Ratio for Binary Variables	Odds Ratio [95% CI]	p-value
Trevo Excluding Concomitant Carotid Stenting	Trevo vs MM	1.88 [1.07, 3.29]	0.014 ¹
IV t-PA	Yes vs No	1.45 [0.60, 3.52]	0.4134 ²
NIHSS	>14 or ≤14	0.41 [0.24, 0.69]	0.0010 ²

¹One-sided t-test.

²Two-sided t-test.

Distribution of 90-day mRS*

mRS	Trevo Excluding Concomitant Carotid Stenting (N=96)	90 DAY Mr Clean/Control (N=249)
0	4.2% (4/96)	0.4% (1/249)
1	8.3% (8/96)	5.2% (13/249)
2	17.7% (17/96)	13.7% (34/249)
3	15.6% (15/96)	16.5% (41/249)
4	32.3% (31/96)	31.3% (78/249)
5	7.3% (7/96)	10.4% (26/249)
6	14.6% (14/96)	22.5% (56/249)
mRS [0, <3] with Any IA Lytic use Considered Treatment Failure*	29.2% (28/96)	19.3% (48/249)

*1 subject with IA lytic use was counted as a failure.

*Subjects in both arms received IV tPA within 3 hours (81 out of 96 subjects in the Trevo Retriever cohort and 224 out of 249 subjects in the medical management arm).

The primary safety objective analysis was defined as the all-cause mortality rate between Trevo Retriever plus medical management in comparison to the all-cause mortality rate in the medical management alone control group.

	Trevo FDA Cohort % (x/n) (LCL, UCL) ¹	MM % (x/n) (LCL, UCL) ¹
Mortality by 90 Days	13.33 (16/120) (7.82, 20.75)	22.89 (57/249) (17.82, 28.62)

¹Two-sided 95% Exact Clopper Pearson confidence intervals.

The first secondary effectiveness endpoint was the percentage of patients with no intra-cranial occlusion assessed by a consensus review by up to three readers after 24 hours. The endpoint was assessed by CTA MRA using the Arterial Occlusive Lesion (AOL) scale.

	Trevo Excluding Concomitant Carotid Stenting % (x/N ¹) (LCL, UCL) ²	MM % (x/N ¹) (LCL, UCL) ²
Occlusion Free (24 Hours)	77.5% (62/80) (66.8, 86.1)	33.51 (65/194) (26.90, 40.62)

¹Subjects without 24 hour CTA MRA were excluded from this analysis

²Two-sided 95% Exact Clopper Pearson confidence intervals.

The second secondary effectiveness endpoint is the percentage of Trevo Retriever patients that achieved recanalization (TICI 2a or better). This endpoint was also assessed by a consensus review by at least three readers.

	Trevo Excluding Concomitant Carotid Stenting % (x/N) (LCL, UCL) ¹
Revascularization Success (TICI ≥ 2a)	81.3% (78/96) (72.0, 88.5)

¹Two-sided 95% Exact Clopper Pearson confidence intervals.

The third effectiveness endpoint was defined as the difference in neurological outcomes assessed by NIHSS at 24 hours and at 5-7 days post-randomization between Trevo Retriever plus medical management in comparison to the MM control group.

NIHSS	Trevo Excluding Concomitant Carotid Stenting Mean ± SD (N) Median (Min, Max) (LCL, UCL) ²	MM Mean ± SD (N) Median (Min, Max) (LCL, UCL) ²
NIHSS (24 Hours)	14.23 ± 9.23 (96) 14.00 (1, 42) (12.36, 16.10)	16.19 ± 7.77 (240) 16 (0, 42) (15.20, 17.18)
NIHSS (5-7 Days) ¹	13.02 ± 12.16 (94) 12.00 (0, 42) (10.53, 15.51)	15.76 ± 11.49 (237) 15 (0, 42) (14.29, 17.23)

¹Subjects with missing NIHSS were excluded from this analysis. Subjects who expired prior to discharge were assigned a NIHSS value of 42.

²Two-sided 95% confidence limits by normal approximation.

The following table summarizes all adverse events through 90 days post randomization, for the Trevo™ Retriever cohort of subjects, utilizing the MR CLEAN study categorization. For all safety analyses, subjects exposed to Trevo Retriever with or without carotid stenting are included.

MR CLEAN Class	Trevo FDA Cohort	MM
	120 Subjects w/189 Events % (x/N) [LCL, UCL] ¹ Number of Events	249 Subjects w/331 Events % (x/N) [LCL, UCL] ¹ Number of Events
Ischemic stroke	2.5% (3/120) [0.5%, 7.1%] 3	0.8% (2/249) [0.1%, 2.9%] 2
Symptomatic intracranial hemorrhage	7.5% (9/120) [3.5%, 13.8%] 9	6.8% (17/249) [4.0%, 10.7%] 17
Extracranial hemorrhage	2.5% (3/120) [0.5%, 7.1%] 3	4.4% (11/249) [2.2%, 7.8%] 11
Cardiac Ischemia	0	2.0% (5/249) [0.7%, 4.6%] 5
Allergic reaction	1.7% (2/120) [0.2%, 5.9%] 2	0.4% (1/249) [0.0%, 2.2%] 1
Pneumonia	18.3% (22/120) [11.9%, 26.4%] 24	19.3% (48/249) [14.6%, 24.7%] 52
Other infection	23.3% (28/120) [16.1%, 31.9%] 34	13.7% (34/249) [9.6%, 18.6%] 41
Other complication	47.5% (57/120) [38.3%, 56.8%] 89	38.6% (96/249) [32.5%, 44.9%] 154
Progression of stroke	20.8% (25/120) [14.0%, 29.2%] 25	18.9% (47/249) [14.2%, 24.3%] 48

¹Two-sided 95% Exact Clopper Pearson confidence intervals

The following table summarizes all MR CLEAN adverse events through 90 days post randomization, for the FDA cohort, by MedRA coding.

MR CLEAN MedRA v18.0 Coded Adverse Events for FDA cohort over 1% Frequency

MedDRA Preferred Term	Trevo FDA Cohort (N=120) Patients with Events (%)	Control (N=249) Patients with Events (%)
Stroke in evolution	19 (15.8%)	41 (16.5%)
Haemorrhage intracranial	13 (10.8%)	20 (8.0%)
Urinary tract infection	11 (9.2%)	17 (6.8%)
Atrial fibrillation	10 (8.3%)	10 (4.0%)
Pyrexia	6 (5.0%)	6 (2.4%)
Neurological decompensation	6 (5.0%)	6 (2.4%)
Pneumonia aspiration	5 (4.2%)	8 (3.2%)
Delirium	4 (3.3%)	6 (2.4%)
Bone graft	4 (3.3%)	1 (0.4%)
Cardiac failure	3 (2.5%)	7 (2.8%)
Diarrhoea	3 (2.5%)	4 (1.6%)
Urosepsis	3 (2.5%)	4 (1.6%)
Fall	3 (2.5%)	4 (1.6%)
Headache	3 (2.5%)	1 (0.4%)
Ischaemic stroke	3 (2.5%)	1 (0.4%)
Cardiac arrest	2 (1.7%)	0
Hypersensitivity	2 (1.7%)	1 (0.4%)
Clostridium difficile infection	2 (1.7%)	0
Gout	2 (1.7%)	3 (1.2%)
Carotid artery dissection	2 (1.7%)	0
Epilepsy	2 (1.7%)	3 (1.2%)
Depression	2 (1.7%)	4 (1.6%)
Urinary retention	2 (1.7%)	2 (0.8%)
Pulmonary embolism	2 (1.7%)	6 (2.4%)
Haematoma	2 (1.7%)	2 (0.8%)
Anaemia	1 (0.8%)	3 (1.2%)
Bradycardia	1 (0.8%)	3 (1.2%)
Seizure	1 (0.8%)	3 (1.2%)
Respiratory failure	1 (0.8%)	3 (1.2%)
Gastrointestinal tube insertion	1 (0.8%)	5 (2.0%)
Hypotension	1 (0.8%)	3 (1.2%)
Phlebitis	1 (0.8%)	3 (1.2%)
Vasospasm	1 (0.8%)	0
Arrhythmia	0	3 (1.2%)
Myocardial ischaemia	0	3 (1.2%)
Renal failure	0	3 (1.2%)
Deep vein thrombosis	0	3 (1.2%)

TREVO 2 study to support Indication 2

The TREVO 2 (Thrombectomy Revascularization of Large Vessel Occlusions in acute ischemic stroke) clinical trial assessed the safety and efficacy of the Trevo™ Retriever device in comparison to the Merci™ Retriever. One hundred and seventy-eight (178) patients were enrolled and randomized to treatment with either the Trevo Retriever (N=88) or the Merci Retriever (N=90). Key inclusion criteria were: patients with a large-vessel occlusion who failed or were ineligible for intravenous tissue plasminogen activator (IV t-PA) and who could be treated within 8 hours of stroke symptom onset; age 18-85; 8 ≤ NIHSS ≤ 29; and angiographic confirmation of a persistent large vessel occlusion in the internal carotid, middle cerebral (M1 and/or M2 segments), basilar and/or vertebral arteries. Key exclusion criteria were: stenosis in a proximal vessel that requires treatment or that prevents access to the thrombus with the assigned study device.

The primary effectiveness endpoint was revascularization defined as at least TIC1 2 flow after use of the assigned device, as measured by a Core Lab. Subjects with a baseline TIC1 2a by Core Lab and subjects in whom the device was never attempted were excluded from the analyses shown in the table below.

The primary safety endpoint was incidence of procedure-related serious adverse events (PRSAEs) through 24 hours post procedure defined as the following: vascular perforation or intramural arterial dissection, symptomatic ICH, embolization to a previously uninvolved territory, access site complication requiring surgical repair or blood transfusion, mortality within 24 hours, device failure (in-vivo breakage), or any other complications judged by the Clinical Events Committee to be related to the procedure.

TREVO 2 Study Endpoints

(Patients with Baseline TIC1 0 or 1, Study Device Attempted)

Endpoint	Trevo N=79 pts	Merci N=81 pts	Difference [95% CI] ^a	p-value
Primary Effectiveness Endpoint				
Post-Device Revascularization Success (TIC1 ≥ 2a)				
Any Adjunctive Therapy Use Considered Treatment Failure	79.7% (63/79)	49.4% (40/81)	30.4%	<0.0001 ^c
[95% Conf. Interval] ^a	[69.2%, 88.0%]	[38.1%, 60.7%]	[15.0%, 44.1%]	<0.0001 ^d
Post-Device Revascularization Success (TIC1 ≥ 2a)				
Any IA Lytic use Considered Treatment Failure	87.3% (69/79)	58.0% (47/81)	29.3%	<0.0001 ^c
[95% Conf. Interval] ^a	[78.0%, 93.8%]	[46.5%, 68.9%]	[15.0%, 42.4%]	<0.0001 ^d
Primary Safety Endpoint				
Composite Events	13.9% (11/79)	23.5% (19/81)	-9.5% [-22.1%, 2.8%]	0.1567 ^a
[95% Conf. Interval] ^a	[7.2%, 23.6%]	[14.8%, 34.2%]		
Vessel Perforation	0.0% (0/79)	9.9% (8/81)	-9.9% [-18.5%, -3.9%]	
Intramural Arterial Dissection	0.0% (0/79)	1.2% (1/81)	-1.2% [-6.7%, 3.5%]	
Symptomatic ICH	5.1% (4/79)	9.9% (8/81)	-4.8% [-14.1%, 3.8%]	
Embolization to Previously Uninvolved Territory	7.6% (6/79)	4.9% (4/81)	2.7% [-5.6%, 11.4%]	
Access Site Complication Requiring Surgical Repair or Blood Transfusion	1.3% (1/79)	0.0% (0/81)	1.3% [-3.3%, 6.9%]	
Mortality within 24 hrs	1.3% (1/79)	0.0% (0/81)	1.3% [-3.3%, 6.9%]	
In-vivo Device Failure	0.0% (0/79)	0.0% (0/81)	0.0% [-4.6%, 4.6%]	
Other PR-SAE	0.0% (0/79)	0.0% (0/81)	0.0% [-4.6%, 4.6%]	
mRS 0-2 at 90 days	38.2% (29/76)	17.9% (14/78)	20.2% [6.4%, 34.1%]	0.0068 ^e

a: Exact Clopper Pearson confidence intervals on individual proportions; b: Exact confidence intervals on differences in proportions computed with StatXact Version 8; c: Non-inferiority hypothesis using Blackwelder's method with non-inferiority margin of 10%; d: One-sided Fisher's exact test of superiority; e: Two-sided Fisher's exact test.

WARRANTY

Stryker Neurovascular warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Stryker Neurovascular's control directly affect the instrument and the results obtained from its use. Stryker Neurovascular's obligation under this warranty is limited to the repair or replacement of this instrument and Stryker Neurovascular shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Stryker Neurovascular neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **Stryker Neurovascular assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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Copyright © 2020 Stryker
Stryker Australia Pty Ltd